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- Shimojo N, Naka K, Nakajima C, Yoshikawa C, Okuda K, Okada K. Test-strip method for measuring lactate in whole blood. Clin Chem. 1989;35(9):1992-1994.
- Shimojo N, Fujino K, Kitahashi S, Nakao M, Naka K, Okuda K. Lactate analyzer with continuous blood sampling for monitoring blood lactate during physical exercise. Clin Chem. 1991;37(11):1978-1980
- Lehninger AL. Principles of Biochemistry. New York, NY: Worth Publishers; 1993:416.
- Mascini M. Lactate and pyruvate electrochemical biosensors for whole blood in extracorporeal experiments with an endocrine artificial pancreas. Clin Chem. 1987;33(4):591-593.
- Jacobs DS, Kasten BL, DeMott WR, eds. Laboratory Test Handbook. 2nd ed. Stow, OH: LexiComp Inc; 1990:245.
- Sakushima K, Hayashino Y, Kawaguchi T, Jackson JL, Fukuhara S. Diagnostic accuracy of cerebrospinal fluid lactate for differentiating bacterial meningitis from aseptic meningitis: A meta-analysis. J Infect. 2011;62: 255-262.
- Clinical and Laboratory Standards Institute. Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Seventh Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2017. CLSI Document GP41-A7.
- Clinical and Laboratory Standards Institute. Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Clinical and Laboratory Standards Institute. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Tietz NW. Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, PA: WB Saunders Company; 1995:520-521.
- Data on file at Siemens Healthcare Diagnostics.
- Guder WG, Narayanan S, Wisser H, Zwata B. List of Analytes Preanalytical Variables. German Society of Clinical Chemistry. 1st

ed. 1996.

- Young DS. Effects of Drugs on Clinical Laboratory Tests. 3rd ed. Washington, DC: AACC Press; 1990.
- Tintu A, Rouwet E, Russcher H.: Interference of ethylene glycol with L-lactate measurement is assay-dependent. Annals of Clinical Biochemistry 2013; 50:70-72.
- Clinical and Laboratory Standards Institute. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
- Leen WG, Willemsen MA, Wevers RA, Verbeek MM. Cerebrospinal Fluid Glucose and Lactate: Age-Specific Reference Values and Implications for Clinical Practice. PLOS ONE. 2012;7(8): e42745. <https://doi.org/10.1371/journal.pone.0042745>.
- Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
- Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
- Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP09c.
- Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP07-ED3.

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Siemens Healthcare Diagnostics Inc.