**IMMULITE® 2000**

Systems

**EBV-EBNA IgG**

# IMMULITE® 2000 EBV-EBNA IgG

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the qualitative detection of IgG antibodies to the Epstein-Barr virus nuclear antigen (EBNA) in human serum or plasma (heparinized or EDTA). This assay may also be used as an aid in the diagnosis of EBV-associated diseases when used in conjunction with the IMMULITE 2000 EBV-VCA IgG and IMMULITE 2000 EBV-VCA IgM immunoassays.

Catalog Number: **L2KEB2** (200 tests)

Test Code: **EBN** Color: **Dark Pink**

### Summary and Explanation

Epstein-Barr virus (EBV) is a member of the herpesvirus family and a worldwide infectious agent that is the classic cause of infectious mononucleosis (IM).<sup>1</sup> IM is often diagnosed based on the presentation of characteristic symptoms of sore throat, lymphadenopathy and fever, along with a positive heterophile antibody test and a characteristic hematological picture of lymphocytosis (> 50%) with 10% or more of the lymphocytes identified as atypical.<sup>2</sup>

EBV is generally acquired by oral transmission of infected saliva. Following multiplication in epithelial cells of the oropharynx, the virus selectively infects B lymphocytes in the peripheral blood and other reticuloendothelial tissues.

EBV has been intimately associated with several human tumors, including African Burkitt lymphoma, nasopharyngeal carcinoma, Hodgkin's disease, and B-cell lymphomas and Leiomyosarcomas in organ transplant recipients and patients with immunodeficiency disorders.<sup>3</sup>

Serological methods have been an important aid in the diagnosis of EBV infections. Several distinct EBV-associated antigen systems and their corresponding antibodies have been characterized and developed into important diagnostic tools. The EBV

antigen systems are classified by the phase of the viral replicative cycle during which they are expressed. The latent-phase antigens include EBNA proteins, of which EBNA-1 is the principal constituent. Anti-EBNA antibodies typically appear much later, gradually emerging over several (2 to 4) months following the onset of illness and remaining for life. Some immunocompromised patients demonstrate diminishing or undetectable levels of anti-EBNA antibodies. A past EBV infection is characterized by the concurrent presence of moderate but stable titers of IgG antibodies to VCA and EBNA and the absence of IgM antibody to VCA.<sup>4</sup>

### Principle of the Procedure

IMMULITE 2000 EBV-EBNA IgG is a solid-phase, two-step chemiluminescent enzyme immunoassay.

The solid phase (bead) is coated with EBNA-1 antigen. The liquid phase consists of two reagents: 1) protein based buffer and 2) alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-human IgG.

In the first cycle, the on-board diluted patient sample (1-in-20) and the protein-based buffer are incubated together with the coated bead for 30 minutes. During this time, IgG specific to EBNA in the sample binds to the EBNA-1 antigen on the bead. Unbound sample is then removed by centrifugal washes.

In the second cycle, the enzyme-conjugated monoclonal murine anti-human IgG antibody is added to the original reaction tube for an additional 30-minute incubation. The enzyme-conjugated monoclonal murine anti-human IgG antibody binds to immobilized IgG to form the antibody sandwich complex. The unbound enzyme conjugates are removed by centrifugal washes. Finally, chemiluminescent substrate is added to the reaction tube containing the bead, and the signal is generated in proportion to the bound enzyme.

**Incubation Cycles:** 2 × 30 minutes

**Time to first result:** 65 minutes

## Specimen Collection

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

The use of an ultracentrifuge is recommended to clear lipemic samples. Samples which are cloudy or have particulate material should be clarified by low-speed centrifugation.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 EBV-EBNA IgG has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 10 µL serum or plasma (heparinized or EDTA). After automatic predilution, 50 µL of diluted sample is used for the immunoassay.

**Storage:** 3 days at 2–8°C or 6 months at –20°C.<sup>5</sup>

**Automatic Predilution Factor:** 20

## Warnings and Precautions

For *in vitro* diagnostic use.

For Professional Use.

Safety data sheets (SDS) available on [siemens-healthineers.com](http://siemens-healthineers.com).

### WARNING! POTENTIAL BIOHAZARD

Contains human source material.

No known test method can ensure that products derived from human source materials will not transmit infection. These materials should be handled using good laboratory practices and universal precautions.<sup>14-16</sup>

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

### Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

**H317  
P280,  
P302+P352,  
P333+P313,  
P362+P364**

**Contains:** 2-methyl-2H-isothiazol-3-one; EBV-EBNA IgG Reagent Wedge

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1%, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing. Disposal into drain systems must be in compliance with prevailing regulatory requirements.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

Dispose of hazardous or biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with prevailing regulatory requirements.

## Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

### **EBV-EBNA IgG Bead Pack (L2EB12)**

With barcode. 200 beads, coated with EBNA antigen. Stable at 2–8°C until expiration date.

**L2KEB2:** 1 pack

### **EBV-EBNA IgG Reagent Wedge (L2EBA2)**

With barcode. Two reagents: 11.5 mL buffer and 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-human IgG antibody in a protein-based buffer, with preservative. Stable at 2–8°C until expiration date.

**L2KEB2:** 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

### **EBV-EBNA IgG Adjustor (LEBR)**

4 mL human serum with IgG reactive to EBV-EBNA, in a protein-based buffer, with preservative. (No predilution is required.) Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

**L2KEB2:** 1 vial

### **EBV-EBNA IgG Controls (LEBC1, LEBC2)**

**Negative Control (LEBC1):** One vial containing 2 mL human serum nonreactive to EBV-EBNA, with preservative.

**Positive Control (LEBC2):** One vial containing 2 mL human serum with IgG reactive to EBV-EBNA, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

**L2KEB2:** 1 set

The IMMULITE 2000 software performs automatic on-board dilution of control samples, and the results will be tracked in the QC database. Enter controls as controls.

For the current control ratio ranges, refer to the Control insert.

Before running adjustors or controls, place the appropriate Aliquot Labels (supplied with the kit) on test tubes, so that the barcodes can be read by the on-board reader.

### **IgG/IgM Sample Diluent (L2IGZ2)**

For the on-board dilution of patient samples and controls. 55 mL concentrated (ready-to-use) nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or 6 months (aliquotted) at –20°C.

**L2KEB2:** 1 vial

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2KEB2:** 3 labels

## Materials Required but Not Provided

### **IgG/IgM Sample Diluent (L2IGZ2)**

For the on-board dilution of patient samples and controls. 55 mL concentrated (ready-to-use) nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or 6 months (aliquotted) at –20°C.

**L2IGZ2:** 1 vial

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2IGZ2:** 5 labels

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**L2ZT:** 250 Sample Diluent Test Tubes (16 × 100 mm)

**L2ZC:** 250 Sample Diluent Tube Caps

Also Required

Distilled or deionized water; test tubes

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

### Recommended Adjustment Interval:

2 weeks

**Quality Control Samples:** The EBV-EBNA IgG controls supplied with the kit should be used as quality control material to monitor the performance of the assay.

The Controls should be processed at (or near) the beginning of every run containing patient samples to be tested for EBV-EBNA IgG, and also when readjusting.

It is also recommended that known reactive and nonreactive specimens be run periodically to assure pipetting accuracy for the dilution step.

Users of the IMMULITE system are advised to consult the NCCLS document C24-A, Internal Quality Control Testing: Principles and Definitions, for additional guidance on quality control for the basic principles and definitions in dealing with internal quality control testing.

Additional quality control material can be used at the discretion of the laboratory. Use the quality control material in accordance with the quality control instructions for use.

In addition, perform quality control:

- Following a valid calibration
- With use of a new lot of reagent
- When troubleshooting test results that do not match clinical conditions or symptoms

Follow government regulations or accreditation requirements for quality control frequency. Individual laboratory quality control programs and procedures may require more frequent quality control testing.

Acceptable performance is achieved when the analyte values obtained are within the expected control interval for the system, as indicated by the manufacturer of the control material or within the interval determined by an internal laboratory quality control procedure.

Follow your laboratory's quality control procedures if the results obtained do not fall within the acceptable limits. For information about entering quality control definitions, refer to the system online help.

### Calculation of Cutoff and S/CO Ratio:

The Master Cutoff of the assay was determined from representative samples to achieve optimal sensitivity and specificity for the assay.

The cutoff is set equal to the average counts per second (mean cps) of the Adjustor (from the most recent adjustment) multiplied by Curve Parameter 1. (See the "Low Adjustor CPS" and "Curve Parameter 1" fields in the IMMULITE 2000 Kit Information screen, which can be accessed from the menu via Data Entry: Kit Entry.)

Calculation of a signal/cutoff (s/co) ratio is done by using the following formula:

$$\text{S/CO Ratio} = \frac{\text{Sample or Control cps}}{\text{Mean Adjustor cps} \times \text{P1}}$$

Calculation and reporting of qualitative (reactive/nonreactive/indeterminate) and s/co ratio results are handled automatically by the IMMULITE 2000.

The result for a sample is reported as "Indeterminate" if the counts per second for that sample fall within  $\pm 10\%$  of the cutoff. The result is reported as "Reactive" if the sample's counts are *above* the indeterminate range, and "Nonreactive" if *below* this range.

Additional controls may be tested in accordance with guidelines or requirements of local, state and/or federal regulations or accrediting organizations.

It is also recommended that known reactive and nonreactive specimens be run periodically to assure pipetting accuracy for the dilution step.

## Interpretation of Results

The cutoff of the IMMULITE 2000 EBV-EBNA IgG assay was determined with reactive and nonreactive patient samples by a ROC analysis, with a balanced consideration of sensitivity and specificity.

A result of “**Reactive**” (ratio of  $\geq 1.1$ ) indicates that EBV-EBNA IgG antibodies were detected in the patient sample.

A result of “**Nonreactive**” (ratio of  $< 0.9$ ) indicates that EBV-EBNA IgG antibodies were not detected in the patient sample.

Any result of “**Indeterminate**” (ratio between 0.9 and  $< 1.1$ ) should be retested. Samples which still test as “Indeterminate” should be tested by an alternate method, or a second sample should be taken — if possible — within a reasonable period of time (e.g., one week).

The presence of IgG antibodies to EBV-EBNA is an indication of previous exposure to the virus. A single specimen can only be used to determine the serological status of the individual.

The magnitude of the measured results (cps) above the Cutoff is not indicative of the total amount of antibodies detected.

Reports by the laboratory to the physician should include “The following results were obtained with the IMMULITE 2000 EBV-EBNA IgG EIA. Values obtained from other manufacturers' assay methods may not be used interchangeably.”

An interpretation of the results is not recommended if the positive or negative control falls outside the range specified in the EBV-EBNA IgG Control Insert.

Clinical interpretation requires knowledge of the patient's medical history and clinical condition (signs and symptoms), as well as other diagnostic results including presence or absence of antibodies against VCA IgM and VCA IgG utilizing the IMMULITE 2000 EBV-VCA IgM and IMMULITE 2000 EBV-VCA IgG immunoassays.

Although EBV serodiagnosis requires measurement of more than one analyte, universal agreement of a serological profile does not exist. Nonreactive results do not rule out the diagnosis of infectious mononucleosis (IM). The specimen may

have been drawn before appearance of detectable antibodies. Nonreactive results in suspected early IM infection should be retested in 4–5 weeks. A serological profile used for diagnostic interpretation is shown below.<sup>12,13</sup>

### Interpretation Criteria:

Interpretation	VCA IgG	EBNA IgG	VCA IgM
Past Infection	+	+	–
Current / Acute	+/-	–	+
No Infection	–	–	–

## Expected Values

Antibodies to EBV have been detected in all population groups studied. However, the prevalence may vary with the population group under review (e.g. age, race, geographic location, socioeconomic status, etc.).<sup>6</sup> Approximately 50% of the early childhood ( $< 5$  years) population in the United States and Great Britain demonstrate EBV seroconversion. Further seroconversion occurs in adolescence and young adulthood, reaching greater than 90% seroconversion by the 4<sup>th</sup> decade of life.<sup>7, 8, 9</sup> The incidence of infectious mononucleosis is 30 times higher for whites than for blacks.<sup>9</sup> Persons in tropical and developing nations, as well as lower socioeconomic persons in industrialized nations, acquire EBV in childhood, while the acquisition of EBV by economically advantaged persons in industrialized nations does not occur until adolescence or early adulthood.<sup>10</sup> After EBV infection, VCA IgG antibodies persist indefinitely.<sup>11</sup>

## Limitations

This assay is not intended for use in screening blood or plasma donors.

The results of the test must be taken within the context of the patient's clinical history, symptomology and other laboratory findings.

The presence of IgG antibodies in a single specimen is not sufficient to distinguish between active or past infection.

Diagnosis should not be made on the basis of anti-EBV-EBNA IgG results alone.

Test results should be interpreted in conjunction with the results of IMMULITE 2000 EBV-VCA IgM and IMMULITE 2000 EBV-VCA IgG. Caution should be used if interpreting a profile using more than one manufacturer's kits or different assay methods.

Samples containing antinuclear antibodies or other anticell antibodies may give false positive results.

The performance characteristics of this assay have not been established for use with specimens from neonates, cord blood, or pretransplant patients.

The results in HIV patients, in patients undergoing immunosuppressive therapy, or in patients with other disorders leading to immunosuppression, should be interpreted with caution.

The performance of IMMULITE 2000 EBV-EBNA IgG has not been established for sample matrices other than human serum and plasma.

IMMULITE 2000 EBV-EBNA IgG ELISA assays are not intended to replace virus isolation and/or identification.

Heterophilic antibodies in human serum and plasma can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

The drug asfotase alfa (STRENSIQ®), a recombinant form of alkaline phosphatase, is expected to interfere with *in vitro* diagnostic assays utilizing an alkaline phosphatase detection system. Test samples from patients taking asfotase alfa with a non-alkaline phosphatase methodology.<sup>17</sup>

## Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed as a signal-to-cutoff ratio. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Precision:** Samples were processed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Based on internal testing, the overall reproducibility is estimated to be  $\leq 15\%$  CV for samples tested  $\geq 0.9$  S/CO and  $\leq 20\%$  CV for samples tested  $0.6$  to  $< 0.9$  S/CO and includes multiple reagent lots, instruments, days, operators, and replicates. Performance of the assay at individual laboratories may vary.

**Crossreactivity:** A study was conducted on 57 samples to evaluate whether the measurement of EBV-EBNA IgG antibody is affected by closely related microorganisms. Sera containing antibodies to the microorganisms listed in the table below were tested by IMMULITE 2000 EBV-EBNA IgG and by a commercially available ImmunoWell™ EBV-EBNA IgG.

Disease state	Total <i>n</i>	IMMULITE 2000 nonreactive/ tested	ImmunoWell™ nonreactive/ tested
ANA	10	0/10	0/10
Anti-HIV	10	0/10	0/10
Rheumatoid Factor	10	0/10	0/10
Toxoplasma IgG	9	1/9	1/9
Anti-VZV	3	0/3	0/3
CMV IgG	5	0/5	0/5
HSV IgG	5	0/5	0/5
Rubella IgG	5	0/5	0/5

Due to the high prevalence of EBV in a normal population, fifty-six out of 57 samples yielded reactive results. The 56 reactive samples were further tested by a commercially available ImmunoWell™ EBV-EBNA IgG assay and yielded similar results.

**Alternate Sample Types:** To assess the effect of alternate sample types, blood was collected from 52 volunteers into plain, sodium heparin, EDTA and Becton Dickinson SST® vacutainer plastic tubes and assayed by the IMMULITE 2000 EBV-EBNA IgG procedure. Results are expressed as a signal-to-cutoff ratio. By linear regression:

$$(NaHeparin) = 0.98 (\text{Serum}) + 0.21$$

$$r = 0.98$$

$$(EDTA) = 0.98 (\text{Serum}) + 0.48$$

$$r = 0.98$$

$$(SST) = 1.00 (\text{Serum}) + 0.12$$

$$r = 0.98$$

Means:  
27.0 (Serum)  
26.8 (NaHeparin)  
27.2 (EDTA)  
27.0 (SST)

**Bilirubin:** Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Presence of hemoglobin in concentrations up to 537 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

**Method Comparison 1:** The assay was compared to a commercially available DiaSorin LIAISON EBNA IgG kit across 200 samples.

IMMULITE 2000 EBV-EBNA IgG

DiaSorin LIAISON EBNA IgG	Re-active	Ind	Non-reactive	Relative Sensitivity	Relative Specificity
Positive	99	0	1		
Equivocal	12	0	1	99%	58%
Negative	30	7	50		

Total Agreement: 74%

Forty-nine discordant samples were resolved by western blot testing. The table below shows the resolved relative sensitivity and specificity after immunoblot resolution testing.

IMMULITE 2000 EBV-EBNA IgG

DiaSorin LIAISON EBNA IgG	Re-active	Ind	Non-reactive	Relative Sensitivity	Relative Specificity
Positive	130	0	1		
Equivocal	4	2	0	99%	88%
Negative	6	1	50		

Total Agreement: 94%

## References

- 1) Jenson HB, Baltimore RS. Infectious mononucleosis. Pediatric infectious diseases: principles and practice. Norwalk, Conn. Appleton & Lange, 1995: 565-76. 2) Henle W, Henle G. Epstein-Barr virus and infectious mononucleosis. In: Human herpesvirus infections: clinical aspects. New York: Marcel Dekker, 1982. 3) Gaffey MJ, Weiss LM. Association of Epstein-Barr virus with human neoplasia. Pathol Ann 1992: 27:55-74. 4) Sumaya CV, Ench Y. Heterophile antibody and viral-specific responses. Pediatrics 1985: 75:1011-19. 5) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995: 358. 6) Schooley RT, Dolin R. Epstein-Barr virus (infectious mononucleosis). In: Mandell G, editor. Principles and practice of infectious diseases. 3rd ed. New York: Churchill Livingstone, 1990: 1172-85. 7) Pereira MS, et al. EB virus antibody at different ages. Br. Med J 1969; 4:526. 8) Porter DD, et al. Prevalence of antibodies to EB virus and other herpesviruses. JAMA 1969; 208:1675. 9) Tamir D. Infectious mononucleosis and EB virus in childhood. Pediatrics 1974; 53:330. 10) De-The G. Epidemiology of Epstein-Barr virus and associated diseases in man. In: Roizman B. editor. The herpes viruses. Vol.1. New York: Plenum Press, 1982: 25-103. 11) Pochedly C. Laboratory testing for infectious mononucleosis. Postgrad Med 1987; 81:335. 12) Lennette ET. Epstein-Barr Virus (EBV). In: Lennette EH, Lennette DA, Lennette ET. Diagnostic Procedures for Viral, Rickettsial, and Chlamydial Infections, 7th Ed. American Public Health Association, Washington D.C., 1995. 13) Lennette ET. Epstein-Barr Virus. In: Murray PR, et al. Manual of Clinical Microbiology, 6th Ed. ASM Press, Washington D.C., 1995. 14) US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: US Government Printing Office; December 2009. 15) World Health Organization. Laboratory Biosafety Manual. 3rd ed. Geneva: World Health Organization; 2004. 16) Clinical and Laboratory Standards Institute. Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4. 17) STRENSIQ® (asfotase alfa) injection [package insert]. Boston, MA: Alexion

## Technical Assistance

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

Available outside the United States only.  
For technical assistance, contact your National Distributor.

siemens-healthineers.com

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Shading indicates technical content that differs from the previous version.

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IMMULITE is a trademark of Siemens Healthineers.
All other trademarks and brands are the property of their respective owners.
All rights reserved.
O Sistema da Qualidade da Siemens Healthcare Diagnostics Products Ltd. está registado sob a norma ISO 13485.