MAGNETOM Skyra

Operator Manual - TimTX TrueShape

syngo MR D13



www.siemens.com/healthcare



Manufacturer's notes:

This product bears a CE marking in accordance with the provisions of regulation 93/42/EEC of June 14, 1993 for medical products.

The CE marking applies only to medico-technical products/ medical products introduced in connection with the above-mentioned comprehensive EC regulation.

Overview of Contents

TimTX TrueShape

А



A TimTX TrueShape

- A.1 Safety information
- A.2 Overview
- A.3 B1 shimming
- A.4 Zoomed imaging



Introduction

This operator manual is an addendum to the MR system operator manual.

This manual may include descriptions covering standard as well as optional hardware and software. Contact your Siemens Sales Organization with respect to the hardware and software available for your system. The description of an option does not infer a legal requirement to provide it.

The graphics, figures, and medical images used in this operator manual are examples only. The actual display and design of these may be slightly different on your system.

Male and female patients are referred to as "the patient" for the sake of simplicity.

References to "Siemens Service" include service personnel authorized by Siemens.

Important icons

For readability, certain contents are highlighted. In the following sections, you will find the symbols and their contents used:

- Prerequisites for the operating steps to follow
- Request for action
- Item in list

Notes for optimal use of the system.

Remarks that facilitate working with the system.



Problem

Description of possible sources of errors

Requests for action to solve problems

A.1	Safety information	A.1-1
	General	A.1-1
	Patient padding	A.1-2
A.2	Overview	A.2-1
A.3	B1 shimming	A.3-1
	Mitigating B1 inhomogeneities	A.3-1
	Performing B1 shim	A.3-2
	Performing B1 shim manually	A.3-5
A.4	Zoomed imaging	A.4-1
	syngo ZOOMit applications	A.4-1
	Zoomed SPACE	A.4-3
	Zoomed BOLD imaging in the brain (Zoomed fMRI)	A.4-4
	Zoomed diffusion-weighted imaging	A.4-5
	Performing zoomed imaging	A.4-8

Safety information

General

All safety information of the valid operator manuals must be adhered to without fail. It is merely supplemented by the notes contained herein.

All personnel should be instructed in the proper connection, use and handling of the system.

Always follow proper safety, operating and maintenance procedures to prevent exposure to electrical or mechanical hazards that may cause injury.

Patient padding

Dangerous current loops may be generated when parts of the patient's body touch. For information on how to position the patient appropriately, see (> Operator Manual MR system, chapter Safety)

Furthermore, when using the TimTX TrueShape option, an adequate patient padding is necessary to keep a proper distance between the body parts (for example, legs and arms).

Overview

A.2

Overview

With syngo MR D13C, the TimTX TrueShape option is introduced for MAGNETOM Skyra, a new architecture for parallel transmit (pTx) technology. TimTX TrueShape provides multi-channel TX array imaging, resulting in higher image quality and shorter scan times.

The interaction of independent transmit channels in TimTX TrueShape enables the selective excitation of a specific body area. Besides a focused and more homogeneous *B1 Shimming*, TimTX TrueShape brings entirely new applications into the clinical world:

syngo ZOOMit allows Zoomed Imaging for highlighting regions and organs, with optimized protocols for neurology, oncology, and orthopedics applications.

B1 shimmin

A.3

B1 shimming

Mitigating B1 inhomogeneities

Depending on the patient, B1 inhomogeneities can present a problem for imaging in specific areas of the body, such as the liver, breast, the trunk (including hip and lumbar spine), and in the contrast-enhanced imaging of blood lumen, in particular in the pelvis. B1 field inhomogeneities are induced by the interaction of the patient with the RF energy introduced into a strong main magnetic field $B0 \ge 3$ T. They may result in image artifacts such as shading and a general inhomogeneous impression in certain regions of the image.

Performing B1 shim

The B1 shim allows you to correct inhomogeneities of the RF excitation (B1 field). The adjustment is performed in two steps:

- Measurement of a 3-dimensional B1 map for each Tx channel
- Adaptations of RF pulses for each Tx channel

High B1 field applied over a relatively long time heats the coil housing of the Spine 32. In order to avoid heating make sure that the coil cushion is mounted before performing B1 shim.



Selecting the B1 shim mode

• Open the **System pTx Volumes** parameter card.

Coils	Mis	cellaneous	Adjustme	ents Adju	st Volume	pTx Volu	imes	Tx/Rx
B1 SI	him mode	Volume-se	lective 💌		рТх	Volume 2	- + -	
	Excital	TrueForm Patient-spe	Form ent-specific		Vol. Prop	erty B1 S		Y
		vulume-sei	ecuve		Posi	tion L31.5	5 P17.7 H0.0	
					Orienta	tion Trans	sversal	_
					i	Rotation C).00	📑 deg
						A >> P 🗧	30	🕂 mm
						R >> L 🕄	30	🕂 mm
						F >> H 🗧	30	🕂 mm
					Vol. Visib	oility On		_
					ana ar 17			100 A 01
Program	Routine	Contrast	Resolution	Geometry	System	Physio	BOLD	Sequence

- Set the **B1 Shim mode**.
- TrueForm: The B1 field is optimized for body regions, so the optimization is anatomy-specific.
- Patient-specific: The RF excitation is parameterized on basis of B1 maps of the subject. The optimization aims at the whole volume covered by the slice group.
- Volume-selective: The B1 field is optimized in a specified volume according to the chosen pTx volumes.

Selecting second B1 shim volume

If anatomically necessary you can add a second B1 shim volume.

• Increase the **pTx Volume** parameter to **2**.

The B1 shim volumes are displayed in the GSP segments.

Set the Vol. Property of the selected volume(s) to B1 Shim Vol.

Performing B1 shim manually

The following description applies to *expert users* only.

With the B1 Shim you can correct inhomogeneities of the RF excitation (B1 field).

• Select **Options > Adjustments** from the main menu.

The Manual Adjustments dialog window is opened.

• Select the **B1 Shim** subtask card.

÷

Ma	Manual Adjustments							
	NO TX 0 1 0.86	()[°] T× 1 \$70 0.71	(°) Tx 2 (°) 7777	Tx 3 [][⁰]	Tx 4 [][°]	Tx 5 ([[*] Tx 6 ([[*] Tx 7 [[[*] Converged Yes	
Т×	Channel	Tx 0 💌				B1 Map (sys	5) 22	
Г					Max 0	B1 Map (temp) 22	Apply
					livino 1	Amplitude [\	/] 356.6	<u>M</u> easure
						B1 Mappin	g Standard 🔽	
						B1 Shimmin	g 📃 🗾	
				10		Temporary	System	
	1	An and				Tx Scale 0 0.86 - 0.00	0.86 0.00	Apply
					-	Tx Scale 2		<u>C</u> alculate
	\bigcirc				•	Tx Scale 3		Load Circular
						Tx Scale 5		Load E <u>l</u> liptical
		۲				Tx Scale 6		Load System
						mag pha	mag pha	
								Reset
	Magnit	ude	Phase		Mask			
	Freque	ncy	Transmitter		3D Shim	Inter. Shim	B1 Shim	Show
S	uccessfully	applied adj	ustment parameter	S.			Close	Help

Entering the amplitude

 Enter the transmit amplitude for the adjustment to be performed in the Amplitude [V] input field.

Starting the measurement

• Start the measurement with **Measure**.

After the measurement, the B1 field map is displayed in the image area of the **Manual Adjustments** dialog window.

Selecting the image type	 Display the magnitude images by selecting the Magnitude card. 				
	– or –				
	 Display the phase images by selecting the Phase card. 				
	– or –				
	 Display the mask images by selecting the Mask card. 				
Windowing images	You set the window values for the image, by dragging the mouse over the image with the center button pressed.				
	 Adjust the contrast by moving the mouse vertically. 				
	 Adjust the width by moving the mouse horizontally. 				
Selecting the transmit chan- nel	 Load the signals of specific transmit channels by selecting them from the Tx Channel list. 				
Applying the B1 map	 Transfer the current temporary B1 map to the MR system with Apply (to the right of the B1 Map (temp)). 				

Calculating the Tx scale factors

Loading the Tx scale factors

 Calculate the Tx scale factors from the current B1 map with Calculate.

The Tx scale factors are recalculated. They are displayed in the list of numeric results and used as the basis for the next B1 map measurement.

The numeric results are listed below the running number of the iteration. The following parameters are displayed.

Parameters	Comment
No.	Running number of the adjustment cal- culation
Tx0 Tx7	Tx scale factors in amplitude/phase dis- play
Converged	Results of adjustment

Transferring the Tx scale fac-
torsTransfer the current temporary Tx scale factors to the MR
system with Apply (to the right of System).

- Display the circular Tx scale factors as temporary parameters with Load Circular.
 - Display the elliptical Tx scale factors as temporary parameters with Load Elliptical.
 - Display the current Tx scale factors of the system with Load System.



Zoomed imaging

syngo ZOOMit applications

Zoomed MR imaging can be performed using *syngo* ZOOMit. Selecting a small stripe (reduced FoV) instead of a large excitation plane allows you to perform a targeted excitation by highlighting regions, organs, or even features of an organ. The reduced FoV can be imaged in a shorter scan time.

The excitation with ZOOMit is integrated in BOLD and diffusion EPI sequences (ep2_bold, ep2_diff) as well as SPACE. Using ZOOMit for EPI has the additional advantage that a shorter echo train for imaging is used which reduces distortions and in addition can shorten the echo time TE.



Example: Zoomed EPI in the brain, using a reduced FoV in phase-encoding direction, compared to a standard full FoV acquisition.

- (1) Reduced FoV using syngo ZOOMit
- (2) Full FoV

Explanation:

- PE: phase-encoding direction
- SS: slice-selection direction
- RO: readout direction

A.4

Zoomed SPACE

The SPACE sequence allows inner volume excitation. The resulting spin echo provides only signal from this inner volume.

The patient's anatomy does not need to be fully covered in order to avoid aliasing artefacts. In practice, you can reduce the FoV in all directions to the volume of interest, as a result, a smaller cuboid volume is defined.



The zoomed SPACE sequence allows you to reduce the FoV in phase-encoding direction resulting in shorter scan time.

Zoomed BOLD imaging in the brain (Zoomed fMRI)

Zoomed imaging is especially useful for BOLD imaging in the brain. Compared to a full FoV excitation, zoomed images show visibly less distortions without aliasing artifacts.



(1) Full FoV

(2) Reduced FoV using syngo ZOOMit





- (1) Inferior slice
- (2) Superior slice
- (3) Full FoV
- (4) Reduced FoV with aliasing artifacts using conventional techniques
- (5) Reduced FoV without aliasing artifacts using syngo ZOOMit

Zoomed diffusion-weighted imaging

Zoomed diffusion allows imaging of smaller organs or smaller parts of an organ with a reduced FoV, eliminating aliasing artifacts. A selected organ can be imaged without signals from surrounding tissue (for example fat).



Diffusion-weighted MRI: Full FoV (left side), reduced FoV using *syngo* ZOOMit (right side)

- (1) DWI with b-value 0
- (2) DWI with b-value 1000
- (3) ADC map

Images of the prostate and the pancreas in particular will be improved. These organs cannot be imaged well with small FoVs using conventional means, because they are surrounded by tissue of several times their size.



Trace-weighted images of the prostate (left side), ADC maps (right side)

- (1) DWI using conventional techniques
- (2) DWI using *syngo* ZOOMit

÷.

Performing zoomed imaging

• Open the System pTx Volumes parameter card.

• Switch Excitation to ZOOMit.

Coils	Misc	ellaneous	Adjustme	ents Adju	st Volume	pTx Volun	nes	Tx/Rx
B1 Sł	nim mode	TrueForm	-		рТ×	Volume 1	• + -	
	Excitati	on ZOOM	t 🗾		Vol. Prop	erty Optimi	ization Vol	~
					Posi	tion Isocen	iter	-
					Orienta	tion Transv	versal	
					I	Rotation 🛛	ł	🚽 deg
						R >> L 19	32 -	mm
						A >> P 25	50	- mm
						F >> Н 📴)	mm
					Vol. Visit	oility On		-
		_						
Program	Routine	Contrast	Resolution	Geometry	System	Physio	BOLD	Sequence

The pTx volume is shown in the graphical slice positioning (GSP) segments. The optimization pTx volume can be varied in the phase-encoding direction. The visibility can be switched on/ off.

When switching on the **Excitation** to **ZOOMit** in EPI or SPACE while reducing the FoV phase, it is recommended also to switch off the **PAT mode** to **None** (**Resolution iPAT** parameter card).

Preparing the zoom

A.4

Zooming the FoV

The zoomed FoV (reduced FoV in the phase-encoding direction) determines the region of excitation.

The FoV is displayed active (yellow) in the image. An arrow indicates the phase-encoding direction.



 Position the zoomed FoV according to your desired volume of interest.

Adapting the suppression volume

For zoomed EPI, in addition to positioning the zoomed FoV, you define the Field of Excitation (FoE) outer border. The FoE is an optimization volume used for suppressing undesired signals from outside the FoV.

Select the pTx volume with **Vol. Property Optimization**. ٠

The FoE is displayed in the image.

• Adapt the FoE outer border in the phase-encoding direction so that the expected object size of the patient's anatomy is completely covered.

The FoE outer border influences the total duration of the selective RF pulse and as a result the minimally achievable echo time TE. Therefore, it is recommended not to extend the FoE too far outside of the patient's anatomy.



pTx volume parameters	Comment
Position	Position of the selected pTx volume
Orientation	Orientation of the selected pTx volume
Rotation	In-plane rotation of the selected pTx vol- ume
R >> L	Extension of the volume in the right-left direction
A >> P	Extension of the volume in the ante- rior-posterior direction
F >> H	Extension of the volume in the feet-head direction
Vol. Visibility	Specify if the selected volume shall be visible in the GSP segments or not

Coupling the adjustment volume to the pTx volume

The system's adjustment volume is normally calculated automatically, and coupled to the acquisition volume. Alternatively, a sequence, e.g. EPI in ZOOMit mode, may permit to couple the automatic calculation of the adjustment volume to the FoE/ pTx volume.

- Open the **Adjust Volume** parameter card.
- Set the **Couple to** parameter to **pTx Volume**.

© Siemens AG 2012 Order No. MR-06015.640.03.01.02 11/2012

Global Business Unit

Siemens AG Medical Solutions Magnetic Resonance Henkestr. 127 DE-91052 Erlangen Germany Phone: +49 9131 84-0 www.siemens.com/healthcare

Global Siemens Headquarters

Siemens AG Wittelsbacherplatz 2 80333 Muenchen Germany

Global Siemens Healthcare Headquarters

Siemens AG Healthcare Sector Henkestrasse 127 91052 Erlangen Germany Phone: +49 9131 84-0 www.siemens.com/healthcare

Legal Manufacturer

Siemens AG Wittelsbacherplatz 2 DE-80333 Muenchen Germany