SIEMENS

CLINITEK Atlas® Positive and Negative Control Strips for Urinalysis

For use with CLINITEK Atlas Automated Urine Chemistry Analyzer and CLINITEK Novus® Automated Urine Chemistry Analyzer

SUMMARY AND EXPLANATION / INTENDED USE: When used with the CLINITEK Atlas® Automated Urine Chemistry Analyzer and CLINITEK Atlas 10 Reagent Pak or the CLINITEK Novus® Automated Urine Chemistry Analyzer and CLINITEK Novus 10 Urinalysis Cassette, CLINITEK Atlas Positive and Negative Control Strips provide defined results for color, clarity, glucose, bilirubin, ketone (acetoacetic acid), specific gravity, blood, pH, protein, urobilinogen, nitrite, and leukocytes.

Each CLINITEK Atlas Positive Control (☞ Control + ☞) Strip is a firm plastic strip on which 7 separate analyte areas are affixed. These each contain 1 or more natural or synthetic ingredients that, when dissolved out of the analyte areas in a measured quantity of distilled or deionized* water, produce a Positive Control Solution. Each CLINITEK Atlas Negative Control (☞ Control - ☞) Strip contains 6 separate analyte areas. The ingredients are dissolved out of the analyte areas to produce a Negative Control Solution.

Control solutions should be used on a regular basis to verify proper functioning of the CLINITEK Atlas Automated Urine Chemistry Analyzer and CLINITEK Atlas 10 Reagent Pak or the CLINITEK Novus Automated Urine Chemistry Analyzer and CLINITEK Novus 10 Urinalysis Cassette. Read the insert carefully before use (☞ Insert ☞).

This product is for professional, *in vitro* diagnostic use (☞ Intended Use ☞).

PROCEDURE:

Materials required, but not provided:

- URIN-TEK® Specimen Tubes (available in the USA) or equivalent (refer to your analyzer's operating manual for specific tube requirements)

- Distilled or deionized water*

- Directions:**
- Place 12 mL of distilled or deionized water in an appropriately labeled specimen tube. Do not use tap water.
 - Remove a CLINITEK Atlas Control Strip from the bottle and replace the cap immediately and tightly.
 - Place the strip into the tube. Cap tightly.
 - Repeat Steps 1–3 if using a second control.
 - Gently invert the tube(s) back and forth for 2 minutes.
 - Allow the tube(s) to stand for 30 minutes at room temperature.
 - Invert one more time, then remove and discard the strip(s), according to your standard laboratory procedures.
 - For the CLINITEK Atlas Automated Urine Chemistry Analyzer:** If controls are to be assayed at the beginning or end of a run, place the tubes in Positions 48–49 on the Carousel Tray (or Positions 1–2 of a specimen rack, if using a Rack Sample Handler). If a third control is being used, place it in Position 50 (Carousel Tray) or Position 3 (Rack). Alternatively, controls may be assayed within a run by placing them in any desired position(s) among the specimens.

TEST	EXPECTED RESULTS WITH CLINITEK ATLAS OR CLINITEK NOVUS ANALYZER	
	POSITIVE CONTROL	NEGATIVE CONTROL
Color ☞ Color ☞ (instrumentally determined)	Yellow–Orange	Yellow, Green, or Other
Clarity ☞ Clarity ☞	Clear	Clear
Glucose ☞ Glucose ☞	100–500 mg/dL (5.5–28 mmol/L)	Negative
Bilirubin ☞ Bilirubin ☞	Positive	Negative
Ketone ☞ Ketone ☞	Positive	Negative
Specific Gravity ☞ SG ☞	≤ 1.006	1.004–1.014
Occult Blood ☞ BLO ☞	Large	Negative
pH ☞ pH ☞	≥ 8.0	6.0–7.0
Protein ☞ PRO ☞	Trace–100 mg/dL (Trace–1.0 g/L)	Negative
Urobilinogen ☞ URO ☞	1–2 E.U./dL (mg/dL) (16–33 µmol/L)	0.2–1 E.U./dL (mg/dL) (3.2–16 µmol/L)
Nitrite ☞ NIT ☞	Positive	Negative
Leukocytes ☞ LEU ☞	Small–Large	Negative

If test results are questionable or inconsistent with expected findings, refer to the reagent package insert and your analyzer's operating manual for complete information on testing procedures and troubleshooting.

* If using deionized water, it must meet the minimum quality specified in CLSI Guideline C03-A4 *Preparation and Testing of Reagent Water in the Clinical Laboratory*.

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TRADEMARKS

CLINITEK Atlas, CLINITEK Novus, and URIN-TEK are trademarks of Siemens Healthcare Diagnostics. All other trademarks are the property of their respective owners.

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SIEMENS

ARTWORK DETAILS	SPECIFICATIONS	Notes/Comments:
File Name: 19005/E_10929095/B_Atlas Controls_IFU_draft.indd	Languages: EN, DE, FR, IT, ES, PT, EL	Updated to had the document # and it's title block and file name - 10/31/2013 - AD
Project #: POC_100	Dielne #: N/A	Revised copy per redline - AD - 07/19/2013.
Product Brand: CLINITEK Atlas®	Print Process: N/A	Received the translations for "The reconstituted solution" statement - AD - 10/01/2013
Description: Positive and Negative Control Strips for Urinalysis	Size: 13.5 x 13.5	Added the color and clarity boxes - AD - 10/3/2013
SMN #: 10929095 Rev. B	COLORS:	09/04/2014 -WW
REF #: 5029ADS	BLACK ONLY 0/0/100	10/14/2014 - AD Made edits per the redline.
Package Type: IFU (Instructions)		
Document #: 19005 Rev. E		
Artwork Part#: 10929095 Rev. B		

10929095 Rev. B EN/DE/FR/IT/ES/PT/EL

19005 Rev. E, 2014-09

© 2008–2014 Siemens Healthcare Diagnostics. All rights reserved.

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Siemens Healthcare Diagnostics Inc.
Siemens Healthcare Diagnostics Ltd.