Diff-Quik™

English

Rapid Staining Set for differential hemogram

IVD

Intended use and principle

The Diff-Quik™ Staining Set contains solutions for the rapid staining of blood smears, yielding results comparable to the Pappenheim technique (Giemsa-May-Grünwald). By consecutive brief immersions in the Diff-Quik™ Solutions, the previously air-dried smears are fixed and stained in only 15 seconds.

By varying the number of dips in the appropriate staining Solution, different degrees of shading and intensity are easily obtained.

Reagents

Diff-Quik I: Stain Solution I

Contains Eosin Y (1.22 g/l) in phosphate buffer (pH 6.6) and sodium azide (<0.1%) as preservative.

Diff-Quik II: Stain Solution II

Contains Thiazine Dye (1.1 g/l) in phosphate buffer (pH 6.6).

Diff-Quik Fix: Fixative Solution

Contains Fast Green (0.002 g/l) in Methanol.

Warning

1. Stain Solution I contains sodium azide. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. If discarded into sink, flush with a large amount of water to prevent azide build-up.

2. Fixative Solution contains Methanol >50%.

Danger: Highly flammable liquid and vapour. Toxic if inhaled. Causes damage to organs. Wear protective gloves/protective clothing/eye protection/face protection. Keep away from heat, sparks, open flames and hot surfaces. – No smoking. Do not breathe vapour. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower. Keep cool. Dispose of contents and container in accordance with all local regulations.

Stability

If stored appropriately at 15-30°C all 3 components of the staining set are stable after the first opening until the indicated expiration date. The stability is monitored throughout the dating period.

Sample

Capillary blood; freshly collected venous blood (EDTA).

Procedure

1. Dispense Diff-Quik™ Solutions into staining jars with lid or any other dipstaining device.
2. Prepare the slides in the same way as for a Pappenheim technique (smear blood on a degreased slide and let air dry).
3. Dip slide 5 x 1 second in Fixative Solution. Allow excess to drain after each dip.
4. Dip slide 5 x 1 second in Stain Solution I. Allow excess to drain after each dip.
5. Dip slide 5 x 1 second in Stain Solution II. Allow excess to drain after each dip.
6. Rinse slide with distilled water or Weise's buffer pH 7.2 (Merck).
7. Allow to air dry and proceed with differentiation.

Remark

The slides must be dipped several times in and out of the solutions. Continuous immersion during 5 seconds results in unsatisfactory staining.

Limitations of procedure

1. The intensity of the staining can be varied by increasing or decreasing the number of dips, but never go below 3 x 1 seconds.
2. To increase eosinophilic (red) staining, increase the number of dips in Stain Solution I. To increase basophilic (blue) staining, increase the number of dips in Stain Solution II.

3. If the standard procedure results in a too intensive eosinophilic or basophilic staining, the pH of the distilled water should be checked. Normally best staining is obtained at an approximatively neutral pH. Weise's buffer pH 7.2 (Merck) may be used.
4. Diff-Quik™ Solutions usually need to be renewed only once a week. However the solutions should also be renewed if staining differs from normal.

Warranty

This product is warranted to perform as described in its labeling and in the product literature, and Medion Diagnostics AG disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event shall Medion Diagnostics AG be liable for any consequential damages arising out of aforesaid express warranty.

Medion Grifols Diagnostics AG
Bonnstrasse 9, CH-3186 Düringen / Switzerland

© Copyright 2006, Medion Diagnostics AG
All rights reserved.

Siemens Healthcare Headquarters
Siemens Healthcare GmbH
Henkestraße 127, 91052 Erlangen / Germany
☎ +49 9131 84-0

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