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1. Hayes DF, et al. Elevated circulating c-neu oncogene product in patients with breast cancer. Abstract presented at the 12th Annual San Antonio Breast Cancer Symposium (December 8-9. 1989).

2. Carney WP, et al. Detection and quantitation of the human neu oncoprotein. *J Tumor Marker Oncol* 1991;6:53-72.

3. Zabrecky JR, et al. The extracellular domain of p185/neu is released from the surface of human breast carcinoma cells, SK-BR-3. *J Biol Chem* 1991; 266:1716-1720.

4. Kath R, et al. The neu-oncogene product in serum and tissue of patients with breast carcinoma. *Ann Oncol* 1993; 4:585-590.

5. Isola JJ, et al. Elevated erbB-2 oncoprotein levels in preoperative and follow-up serum samples define an aggressive disease course in patients with breast cancer. *Cancer* 1994; 73:652-658.

6. Kandl H, et al. Soluble c-erbB-2 fragment in serum correlates with disease stage and predicts for shortened survival in patients with early-stage and advanced breast cancer. *Brit J Cancer* 1994; 70:739-742.

7. Leitzel K, et al. Elevated soluble c-erbB-2 antigen levels in serum and effusions of a proportion of breast cancer patients. *J Clin Oncol* 1992; 10:1436-1443.

8. Brandt-Rauf PW, et al. The c-erbB-2 protein in oncogenesis: Molecular structure to molecular epidemiology. *Crit Rev Oncol* 1994; 5:313-329.

9. Watanabe N, et al. Serum c-erbB-2 in breast cancer patients. *Acta Oncol* 1994; 33:901-904.

10. Muss HB, et al. c-erbB-2 expression and response to adjuvant therapy in women with node-positive early

breast cancer. *N Eng J Med* 1994; 330:1260-1266.

11. Anderson TI, et al. Detection of c-erbB-2 related protein in sera from breast cancer patients. Relationship to ERBB2 gene amplification and c-erbB-2 protein overexpression in tumour. *Acta Oncol* 1995; 34:499-504.

12. Molina R, et al. Utility of c-erbB-2 in tissue and in serum in the early diagnosis of recurrence in breast cancer patients: Comparison with carcinoembryonic antigen and CA 15.3. *Br J Cancer* 1996; 74:1126-1131.

13. Brandt-Rauf PW, et al. Detection of increased amounts of the extracellular domain of the c-erbB-2 oncoprotein in serum during pulmonary carcinogenesis in humans. *Int J Cancer* 1994; 56:383-386.

14. Vogel W, et al. The neu oncogene product in serum and tissue of patients with metastatic gastrointestinal carcinomas. *J Cancer Res Clin Oncol* 1996; 122:118-121.

15. Okada N, et al. Elevated serum c-erbB-2 protein levels in patients with pancreatic cancer: Correlation to metastasis and shorter survival. *Oncol* 1995; 52:392-396.

16. Meden H, et al. Overexpression of the oncogene c-erbB-2 (HER-2/neu) in ovarian cancer: a new prognostic factor. *Euro J Obst, Gyn and Repro Biol* 1997; 71:173-179.

17. Luo JC, et al. Serum c-erbB-2 oncopeptide in hepatocellular carcinogenesis. *Med Sci Res* 1993; 21:305-307.

18. Myers RB, et al. Serum levels of erbB-2 protein in patients with prostatic adenocarcinoma. *Proceedings AACR* 1995; 36:645.

19. Wu JT, et al. Detection of the extracellular domain of c-erbB-2 oncoprotein in sera from patients with various carcinomas: Correlation with tumor markers. *J Clin Lab Anal* 1993; 7:31-40.

20. Chariyalertsak S, et al. Comparison of c-erbB-2 oncoprotein expression in tissue and serum of patients with stomach cancer. *Tumor Biol* 1994; 15:294-303.

21. Fajac A, et al. c-erbB-2 gene amplification and protein expression in ovarian epithelial tumors: Evaluation of their respective prognostic significance by multivariate analysis. *Int J Cancer* 1995; 64:146-151.

22. Nakano T, et al. Correlation of cervical carcinoma c-erbB-2 oncogene with cell proliferation parameters in patients treated with radiation therapy for cervical carcinoma. *Cancer* 1997; 79:513-520.

23. Reagent Water Technical Bulletin. Siemens Healthcare Diagnostics, 107060.

24. National Committee for Clinical Laboratory Standards. Procedures for Handling and Processing of Blood Specimens; approved guideline. NCCLS Document H18-A2. Wayne, PA: NCCLS; 1999 Oct. 40p.

25. Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988; 34:27-33.

26. Inloes R, Clark D, Drobnies A. Interference of fluorescein,

- used in retinal angiography with certain clinical laboratory tests. Clin Chem 1987; 33:2126-2127.
27. Meden H, Mielke S, Schauer A, Kuhn W. Serum levels of the c-erbB-2 (HER-2/neu) encoded oncoprotein fragment p105 in normal pregnancies. In Vivo 1997; 11(1):51-4.
28. National Committee for Clinical Laboratory Standards. How to define, determine, and utilize reference intervals in the clinical laboratory; approved guideline. NCCLS Document C28-A. Wayne, PA: NCCLS; 1995 June 59p.
29. National Committee for Clinical Laboratory Standards. Interference testing in clinical chemistry. Proposed Guideline. NCCLS Document EP7-P. Wayne, PA: NCCLS; 1986.

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