

EMIT® II Plus Buprenorphine Assay
Alliance Application Sheet
**Emit® II Plus Buprenorphine Assay
Application Sheet**
For the Atellica® CH Analyzer

Refer to the appropriate Instructions for Use for information regarding these reagents. Also refer to the online help for instructions on filling and loading empty CH reagent packs.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

The parameters defined in this application sheet have been developed by Siemens Healthcare Diagnostics to optimize product performance. Any modification to these parameters may affect performance of this and other assays in use on your system and the resulting assay values. It is the responsibility of the user to validate any modifications and their impact on all assay results.

Materials

These reagents are qualified for use with the Calibrators/Controls listed below only. Other material may be used for quality control purposes.

These reagents are for human urine use only.

<u>Emit® II Plus Buprenorphine Assay</u>	<u>REF</u>	<u>Catalog Number</u>
28 mL	10720048	9S039UL
EMPTY reagent pack* (4 x Empty P1s and 4 x Empty P2s)	11097534	N/A**
EMPTY1 reagent pack* (8 x Empty P1s)	11538114	N/A**
EMPTY2 reagent pack* (8 x Empty P2s)	11538115	N/A**

*Materials required but not provided.

**Not Applicable.

<u>Emit® Calibrator/Control</u>	<u>REF</u>	<u>Catalog Number</u>
Level 0	10445406	9A509UL
<u>Emit® II Plus Specialty Drug Calibrators/Controls</u>	<u>REF</u>	<u>Catalog Number</u>
Level 1 (2.5 ng/mL)	10720049	9S529UL
Level 2 (5 ng/mL)	10720050	9S549UL
Level 3 (15 ng/mL)	10720051	9S569UL
Level 4 (25 ng/mL)	10720052	9S589UL
<u>Emit® II Plus Specialty Drug Control</u>	<u>REF</u>	<u>Catalog Number</u>
Negative (3 ng/mL)	10718700	9S609UL
Positive (7 ng/mL)	10718701	9S619UL

**Syva Reagent Preparation, Storage, and
Stability**

Note For reagent preparation, storage information, and full product details, refer to the Syva Instructions for Use.

**Atellica Reagent Preparation, Storage, and
Stability**
Preparation Instructions

Assay reagents are provided as liquid and are ready to use. See reagent IFU for storage instructions of original reagent containers when not in use.

Transfer the appropriate reagent volume to Empty reagent packs (refer to Table 1).

Manually fill the Empty reagent packs with R1 and R2 (refer to Table 1).

EMIT® II Plus Buprenorphine Assay

Alliance Application Sheet

Table 1

Contents	Number of Tests
Pack 1 Well 1 (W1) 14.0 mL of Bupn* Well 2 (W2) Empty	2 x 100
Pack 2 Well 1 (W1) 7.0 mL of Bupn* Well 2 (W2) Empty	

*Do not exceed required fill volume.

Refer to the online help for instructions on filling and loading empty CH reagent packs.

Assay parameters have been predetermined in the instrument software. The reagent pack must be manually loaded on the Atellica® CH Analyzer.

Important Information

For performance characteristics, intended use, limitations and a detailed description how to perform this method, refer to the Emit® II Plus Buprenorphine Assay Instructions for Use.

Calibration and Stability

Table 2

Pack Calibration	Onboard Stability
8 days	28 days

Qualitative

Results

Results are reported as positive or negative.

Semi-Quantitative

Calibration

Prepare a calibration curve by running the calibrators/controls listed below. Recalibrate whenever a new lot of reagents are used or as indicated by control results.

Table 3

Assay	Level 0 ng/mL	Level 1 ng/mL	Level 2 ng/mL	Level 3 ng/mL	Level 4 ng/mL
Bupn	0	2.5	5	15	25

An automatic repeat condition for this assay extends the measuring interval to 100 ng/mL for urine. You may configure the system to trigger an automatic repeat. Automatic repeat results will be flagged **Autorepeat**.

Results

Results are reported as ng/mL.

Table 4

Assay	TDEF Code	Time to First Result
Bupn	8049	6 minutes

EMIT® II Plus Buprenorphine Assay

Alliance Application Sheet

Syva®, Emit®, and Atellica® are trademarks of Siemens Healthcare Diagnostics.

For technical assistance, call Siemens Healthcare

Diagnostics:

1-877-229-3711 in the USA and Canada

In other countries, please contact your local representative.

© 2017-2020 Siemens Healthcare Diagnostics

All rights reserved.

Symbols Key

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
LOT	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
REF	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante
EC REP	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos
IVD	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
CONTENTS	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
LEVEL	Level / Konzentration / Niveau / Livello / Nivel
RxOnly	Prescription Device (US Only) / Verschreibungspflichtiges Medizinprodukt (nur USA) / Dispositif soumis à prescription (États-Unis uniquement) / Dispositivo su prescrizione (solo USA) / Dispositivo de prescripción médica (solo EE. UU.)

2015-07_EFIGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA
siemens.com/healthineers

Siemens Healthineers Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthineers



2020-11-24_11111583_03_US