



Artis one


Basic Planning Information

| Table of contents | |
|-------------------------|----|
| Planning Example | 3 |
| Room Dimensioning | 5 |
| Statics and Transport | 7 |
| Environment | 9 |
| Electrical Installation | 10 |
| General Information | 11 |

| Legend | |
|--|--|
|  Motion area / Swivelling range / Minimal room size / Safety distance | |
|  Service area |  Wall mounted |
|  Floor mounted |  Additional equipment |
|  Ceiling mounted |  Demolition |

| Dimensioning | |
|--|--|
| All installation measurements apply to finished wall/floor/ceiling and are to be checked prior to assembling the unit. | |
|  | |
|  Orientation point = reference point of the Siemens Healthineers unit for planning and installation | |
| Please note: The drawing parts in this document are not to scale! | |

Planning Example



| Artis one – Equipment Legend | | | | |
|------------------------------|---|--|------|------------------------|
| Pos. | Description | Weight (kg), Heat dissipation to the air (W) | | |
| | | kg | W | Remark |
| 1.01 | Floor stand (incl. Mounting plate 45 kg) | 570 | 200 | |
| 1.02 | Patient table (incl mounting plate 25 kg) | 255 | 200 | |
| 1.03 | Control console and ECC I on patient table | 4 | | |
| 1.04 | DCS-EAS 2x 30" Display | 117 | 270 | |
| 1.05 | Upper body radiation shield with OR light LED, moveable | 85 | | optional |
| 1.05a | Transformer for LED-light | | | optional |
| 1.06 | Display with keyboard and mouse | 10 | 175 | |
| 1.07 | MEDRAD Injector Pedestal | 97 | | optional |
| 1.08 | Axis smart image system | 75 | 570 | 88 kg / 960 W with UPS |
| 1.09 | Generator POLYDOROS A100 Plus | 300 | 1000 | |
| 1.10 | Cooling unit - SMC one4all | 28 | 2400 | |
| 1.11 | System control cabinet | 297 | 1600 | |
| 1.12 | Cable cabinet | 120 | | optional |
| 1.13b | 30" Color Display with Mouse and Keyboard | 16 | 170 | |
| 1.13b | Cockpit Controller MLC | 15 | 400 | |


Room Dimensioning

Room dimensioning

The indicated room dimensions have to be checked on site. The planning department has to be informed about possible deviations. Otherwise we cannot assume any guarantee for the accurate implementation of the dimensions indicated in the planning documents.

Room height

not to scale



Room height (RH) : 2400 to 3100 mm

Measured from the highest point of the finished floor (incl. floor covering) to the lowest point of the ceiling substructure.

Required minimum height (MH) 150 mm in suspended ceiling for corrugated hose holder and cable routing.

Positioning of the DCS depends on the room height.

With a room height of < 2800 mm, the upper body radiation shield cannot be secured to the DCS.

The moveable upper-body radiation shield is required!

Rigid collisions of the C-arm with the DCS in case of room height equal 2700 and lower 2800 mm ($> = 2700\text{mm} < 2800\text{mm}$) may lead to damages to the DCS and the Stand. Danger of personal injury due to falling parts ! The distance of 2590 mm between the floor stand's orientation point and DCS mounting point shall not be reduced.

Safety Areas

not to scale

All components that do not belong to the Siemens Healthineers system are labeled as external system components.


When planning non-system components in the examination room, pay attention to the movement ranges of the system components (risk of collision!).

The max. rotation range of the floor stand is limited by the configuration of the room limits during system startup.

Safety Areas for ceiling-mounted Components in the motorized movement Range of the Floor Stand


Area **A** - In this area, no components that are not part of the system may be installed on the ceiling.



Area **B** - Non-system components that are installed in this area must meet the following requirements:
 minimum clear space between the bottom edge of the fixed part (e.g. support arm column) of the component and the finished floor: 2318 mm.
 Moving parts (e.g. support arm) of the components that are below the minimum height of 2318 mm must be able to perform evasion to outside of the max. movement range (including the safety distance) of the floor stand. Otherwise there is a risk of crushing for personnel or a risk of equipment damage.



Area **C** - represents the safety distance of 500 mm between motor-driven parts and personnel as defined by EN 60601.

Statics and Transport

| Statics | not to scale |
|--|---|
| <p>Floor Stand</p> <p>Mounting holes 1 to 7 for installation directly on the floor. The other holes are spare mounting holes.</p> <p>Max. tension forces 8.13 kN at fastening points 1 to 7.</p> <p>This maximum value does not appear simultaneously on all mounting points. Tensile forces depend on the operating position and on the unit movement.</p> <p>Patient Table</p> <p>Load on each mounting point at the table foot-end: max. tensile force equaling 8.55 kN</p> <p>The max. admissible patient weight is 250 kg.</p> <p>DCS-Stand</p> <p>The installation of the DCS depends on the ceiling type.</p> <p>For installation direct on a concrete ceiling the "DCS Artis One ceiling column" is necessary.</p> <p>Max. tensile force of 5.96 kN per mounting point (assuming that only one anchor is supporting the entire load).</p> <p>For installation on a suspended ceiling (heavy load ceiling) the "DCS Artis One adapter plate" is necessary.</p> <p>Max. 2.9 kN pull force per mounting point.</p> |  |


| Transport | | not to scale |
|--|--|---|
| The transport route (doors and hallway) needs sufficient dimensions for the following parts: | | |
| Floor stand with packaging | 250 x 95 x 218 cm (W x D x H), Weight 911 kg | |
| Floor stand with transport carriage | Weight 675 kg | |
| Patient table with packaging | 262,5 x 86 x 197,6 cm (W x D x H), Weight 605 kg | |
| Patient table with transport carriage | Weight 354 kg | |
| The door must have a final clearance of 125 cm if bed entrance is requested. | | |
| Comply with the following transport conditions to prevent damages to the system: | | |
| Maximum acceleration spectral density according to DIN EN 60721-3-2, class 2M2 | | |
| - 1 m2/s3 (0,01 g2/Hz) for oscillation frequencies from 10 to 200 Hz | | |
| Maximum acceleration according to DIN EN 60721-3-2, class 2M2 | | |
| Peak acceleration permitted \hat{a} (full shock response spectrum type I): | | |
| - horizontal shock (longitudinal / transversal):40 m/s2 (corresponds to the shunting of a train car) | | |
| - Vertical shock: 100 m/s2 | | |
| Note: | | |
| All AX components are to be transported in the transport packaging provided by the supplier. | | |
| Only vehicles with air suspension may be used for road transport of AX systems. | | |
| Only rail cars with shock absorbing bumpers ma be used for rail transport of AX systems. | | |
| Floor stand with transport carriage | | Patient table with transport carriage |
|  | |  |


Environment

| Environmental conditions for transport / storage | |
|---|--|
| Temperature range Relative humidity Barometric pressure | -20° to 70 °C 10 to 95 % without condensation 50 to 106 kPa |
| Acceleration | There must be absolute compliance with the following transport conditions to prevent damage to the system: Maximum acceleration spectral density according to DIN EN 60721-3-2 class 2M2 1 m ² /s ³ (0.01 g ² /Hz) for oscillation frequencies from 10 to 200 Hz Maximum acceleration according to DIN EN 60721-3-2 class 2M2 Permissible peak acceleration \hat{a} (full shock response spectrum type I): Horizontal shock (longitudinal / transverse): 40 m/s ² (corresponds to the shunting of a train car) Vertical shock: 100 m/s ² |

| Environment | | |
|--|--|--|
| Examination room Control room Equipment room | Temperature range Relative humidity Max. temperature gradient Barometric pressure | 15 bis 30 °C (recommended 22 °C) 20 bis 75 % non-condensing 5 °C / h 70 bis 104 kPa |
| Image system | Air flow rate Max. noise generation | 156 m ³ / h 53 dB(A) |
| Generator | Air flow rate Max. noise generation | 160 m ³ / h 55 dB(A) |
| System control cabinet | Air flow rate Max. noise generation | 295 m ³ / h 55 dB(A) |
| Cooling unit | Cooling air Flow rate Max. noise generation | 15 to 30 °C, frost-free room 4,2 l / min. 55 dB(A), 57 dB(A) at altitude > 3000m |
| Floor stand | Schocks Vibrations | max. 10 g / 16 ms max. 0,1 g / 10 bis 200 Hz |
| UPS 15 kVA | Temperature range Relative humidity Air change Max. noise generation | 0 to 40 °C 20 to 25 °C recommended 5 to 95 % non-condensing 0,18 m ³ / h 60 dB(A) |
| UPS 40 kVA | Temperature range Relative humidity Air change Max. noise generation | 0 to 40 °C 20 to 25 °C recommended 5 to 95 % non-condensing 0,35 m ³ / h 60 dB(A) |

Electrical Installation

| Power requirements for POLYDOROS A100 Plus  | | | | |
|--|--|----------------------------|--|------------------|
| Power Line: | | 3/N/PE AC, 50/60 Hz ± 1 Hz | Connection Value: | 34,6 kVA |
| Line Voltage: | | 400 V ± 10 % | Power Consumption: Fluoroscopy Radiography | 8 kVA 160 kVA |
| Line Impedance: | | ≤ 135 mΩ | | |
| Cable cross section to be determined by calculation, max 70 mm² | | | | |

| Power requirements for system control cabinet | |  | |
|--|----------------------------|---|---------|
| Power Line: | 3/N/PE AC, 50/60 Hz ± 1 Hz | Connection Value: | 5,5 kVA |
| Line Voltage: | 400 V ± 10 % | Power Consumption: | 7,5 kVA |
| Cable cross section to be determined by calculation, min. 16 mm² | | | |

| Room lighting | |
|---|--|
| <p>Ambient lighting in rooms with diagnostics or with workstations must comply with the respective local and national regulations.</p> <p>General requirements like the needed intensity of illumination - adjustable, reproducible, flicker-free or a limitation of dazzlings and reflections etc. have to be observed (EN 12464-1, DIN 5035-7).</p> | |

General Information

Display screen workstations

For setting up display screen workstations, take account of the guidelines in the Display Screen Workstation directive as well as any national regulations (e.g. EN ISO 9241-5).

Smart Remote Services (SRS)

Smart Remote Services (SRS) is used for remote diagnostics as well as remote service to provide highest system availability.

Requirements:

- Broadband connection (minimum 4 MBit/s down- and 768 kBit/s upstream, optimum 30 MBit/s down- and 2 MBit/s upstream) without time or volume limitations
- Router (for exclusive use with SRS)

Data protection and security is defined in the Smart Remote Services security concept.

Network Integration

The Siemens Healthineers components are using TCP/IP Protocol, a 100/1000 Mbit/s switched Ethernet network and static IP addresses.

The required network cabling (min. CAT 5 TP) has to be provided on site. Media converters, which are needed for using fibre optic cabling, are not in scope of delivery.

To prepare the implementation of the new system into the existing network environment, the availability of the needed network data at least two weeks before starting the installation is mandatory.

This is the only way to ensure a seamless integration of the new system into the workflow of the department.

Notes on preparations for installation

Contracts for performing and supervising on-site installation preparations should be concluded with technically competent companies by the customer. The customer is responsible for timely and proper completion and supervision of all preparations for installation at the construction site in observance of all applicable legal regulations (e.g. X-ray regulations, radiation protection regulations) and all applicable general recognized rules of technology (e.g. VDE regulations, DIN standards).

Execution and supervision of installation preparations at the construction site and later observance of the standard operating conditions are not included in our duties. The customer is responsible for checking the static calculations and, where applicable, the air conditioning in the building to be equipped.

Safety distances

Distances from moving parts of the medical device to walls, furniture and other equipment have to be kept to avoid injuries by crushing in compliance with local regulations, e.g. a minimum distance of 50 cm according to DIN EN ISO 13854.

It is the customer's responsibility to ensure the above requirements are followed. This is to avoid the risk of injury.

Radiation protection

The structural radiation protection depends on the location of the unit and the function of the surrounding rooms. By order, the planning departments of Siemens Healthineers prepare radiation protection calculation and radiation protection plan.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice. Some/all of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features which do not always have to be present in individual cases.

The statements by Siemens' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g. hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results. The customers cited are employed by an institution that might provide Siemens product reference services, R&D collaboration or other relationship for compensation pursuant to a written agreement.

Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Not all features shown in this brochure are necessarily standard and available in all countries.

Published by
Siemens Healthineers AG
SHS ES FD

Siemens Healthineers Headquarters
Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
Phone: +49 9191 180
siemens-healthineers.com