

EMIT® 2000 Phenytoin Assay

Application Sheet

Shading indicates technical content that differs from the previous version.

Emit® 2000 Phenytoin Assay Application Sheet

For the DxC 700 AU® Clinical Chemistry System

Refer to the appropriate Instructions for Use for information regarding these reagents. Also refer to the Operator's Guide for additional instructions.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

The parameters defined in this application sheet have been developed by Siemens Healthineers to optimize product performance. Any modification to these parameters may affect performance of this and other assays in use on your system and the resulting assay values. It is the responsibility of the user to validate any modifications and their impact on all assay results.

Reagents

These reagents are qualified for use with the Calibrator(s) listed below only.

Emit® 2000 Phenytoin Assay

Syva® Kit

Beckman Coulter Kit

Catalog Number

4A019UL

OSR4A229

Calibrators

Emit® 2000 Phenytoin Calibrators

Catalog Number

4A109UL

Storage

Reagents which are in use may be stored on board the analyzer for up to 9 weeks or as long as quality control results are within acceptable limits.

Important Information

For performance characteristics, intended use, limitations and a detailed description how to perform the method, refer to the Emit® 2000 Phenytoin Assay Instructions for Use.

Instrument

Calibration

Prepare a calibration curve whenever a new lot of reagent is used or as indicated by control results. Calibrate by placing the appropriate calibrators in the assigned positions in the calibration rack (yellow rack). Run a reagent blank (blue rack) with each calibration curve.

Table 1

Calibrator			PHY1G-1	PHY1G-2	PHY1G-3	PHY1G-4	PHY1G-5
Phenytoin	µg/mL	0	2.5	5.0	10	20	40
Concentration	µmol/L	0	10	20	40	79	158

Results

Results are reported in µg/mL (µmol/L).

EMIT® 2000 Phenytoin Assay

Application Sheet

Instrument Settings

General Screen

Reagent ID: 506

General		LIH		ISE		Calculated Test		Range	
Test Name: PHY1G ▼		Type: Serum ▼		Operation: Yes ▼					
Sample Volume	1.6 µL	Dilution	0 µL	OD Limit					
Pre-Dilution Rate	1 ▼			Min. OD	-2.0000	Max. OD	2.5000		
Reagent Volume	R1 (R1-1) 72 µL	Dilution	0 µL	Reagent OD Limit					
	R1-2	Dilution		1st	Low -2.0000	High	2.5000		
	R2 (R2-1) 52 µL	Dilution	10 µL	Last	Low -2.0000	High	2.5000		
Common Reagent	Type None	Name	None	Analytical Measuring Range	Low 0.5†	High	40.0†		
Wavelength	Pri 340 nm	Sec	410 nm	Correlation Factor	A 1	B	0		
Method	RATE ▼			Manufacturer Factor	A 1	B	0		
Reaction Slope	+ ▼			Onboard Stability Period	63 Day	0 Hour			
Measuring Point-1	1st 19	Last	27	LIH Influence Check	Yes ▼				
Measuring Point-2	1st	Last		Lipemia	+++++ ▼				
Linearity Limit	100 %			Icterus	++++ ▼				
Lag Time Check	No ▼			Hemolysis	+++++ ▼				

† Values set for working in µg/mL. To work in µmol/L, multiply by 3.96.

Range Screen

General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="PHY1G"/> Type: <input type="text" value="Serum"/>				
Value/Flag <input type="text" value="Value"/> Level Low <input type="text" value="-999999.9"/> High <input type="text" value="999999.9"/>				
Specific Ranges				
	From	To		
	Sex	Year	Month	Year
	Year	Month	Other Type	Low
	High			High
<input type="checkbox"/> 1:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 2:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 3:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 4:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 5:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 6:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
7:	Standard demographics			<input type="text" value="#"/>
8:	Not within expected values			<input type="text" value="#"/>
				<input type="text" value="#"/>
Critical Limits	Low	<input type="text" value="#"/>	High	<input type="text" value="#"/>
			Unit	<input type="text" value="µg/mL"/>
			Select	<input type="text" value="1"/>
			Decimal places	<input type="text" value="1"/>

User Defined

EMIT® 2000 Phenytoin Assay

Application Sheet

Calibration Specific Screen

Calibrators	General	ISE																																																		
Test Name: PHY1G ▼ Type: Serum ▼																																																				
Calibration Type 5AB ▼ Formula EIA Type 1 ▼ Counts 2																																																				
< Calibrator Parameters > <input type="checkbox"/> Use Serum Cal.																																																				
	<table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th rowspan="2"></th> <th rowspan="2">Calibrator</th> <th rowspan="2">OD</th> <th rowspan="2">Conc</th> <th colspan="2">Range</th> </tr> <tr> <th>Low</th> <th>High</th> </tr> </thead> <tbody> <tr> <td>Point-1</td> <td>PHY1G-1*</td> <td></td> <td>2.5††</td> <td>-2.0000</td> <td>2.5000</td> </tr> <tr> <td>Point-2</td> <td>PHY1G-2*</td> <td></td> <td>5.0††</td> <td>-2.0000</td> <td>2.5000</td> </tr> <tr> <td>Point-3</td> <td>PHY1G-3*</td> <td></td> <td>10††</td> <td>-2.0000</td> <td>2.5000</td> </tr> <tr> <td>Point-4</td> <td>PHY1G-4*</td> <td></td> <td>20††</td> <td>-2.0000</td> <td>2.5000</td> </tr> <tr> <td>Point-5</td> <td>PHY1G-5*</td> <td></td> <td>40††</td> <td>-2.0000</td> <td>2.5000</td> </tr> <tr> <td>Point-6</td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Point-7</td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>		Calibrator	OD	Conc	Range		Low	High	Point-1	PHY1G-1*		2.5††	-2.0000	2.5000	Point-2	PHY1G-2*		5.0††	-2.0000	2.5000	Point-3	PHY1G-3*		10††	-2.0000	2.5000	Point-4	PHY1G-4*		20††	-2.0000	2.5000	Point-5	PHY1G-5*		40††	-2.0000	2.5000	Point-6						Point-7						Slope Check + ▼ Allowable Range Check <input type="checkbox"/> Reagent Blank <input type="checkbox"/> Calibration Advanced Calibration Operation Yes ▼ Interval (RB) Lot ▼ Interval (ACAL) Lot ▼ Stability Reagent Blank # Day Hour Calibration # Day Hour
	Calibrator					OD	Conc	Range																																												
		Low	High																																																	
Point-1	PHY1G-1*		2.5††	-2.0000	2.5000																																															
Point-2	PHY1G-2*		5.0††	-2.0000	2.5000																																															
Point-3	PHY1G-3*		10††	-2.0000	2.5000																																															
Point-4	PHY1G-4*		20††	-2.0000	2.5000																																															
Point-5	PHY1G-5*		40††	-2.0000	2.5000																																															
Point-6																																																				
Point-7																																																				
MB Type Factor 1-Point Calibration Point None ▼ <input type="checkbox"/> with Conc-0																																																				

* Syva® Calibrator

†† Calibrator concentrations in µg/mL. Refer to Table 1 for µmol/L.

User Defined

Emit® and Syva® are trademarks of Siemens Healthineers.

DxC 700 AU® Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For technical assistance:

**Beckman Coulter customers, contact the
Customer Technical Support Center at
1-800-854-3633 (USA & Canada)**

**In other countries, please contact your local
Beckman Coulter representative.**

**Siemens Healthineers customers,
contact the Technical Solutions Center at
1-800-227-8994 in the USA**

Technical Assistance

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

[siemens-healthineers.com](https://www.siemens-healthineers.com)

© 2017–2024 Siemens Healthineers

All rights reserved.