IMMULITE®

IL-10

For use on the IMMULITE®
and IMMULITE® 1000 systems

SIEMENS

IMMULITE® / IMMULITE® 1000 IL-10

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® and IMMULITE 1000 Analyzers — for the quantitative measurement of interleukin 10 (IL-10) in serum or heparinized plasma, as an aid in the study of inflammatory diseases.

Catalog Number: **LKXP1** (100 tests)

Test Code: **I10** Color: **Red**

Summary and Explanation

The cytokine interleukin-10 (IL-10) was discovered in 1989. It is a 35-40 kDa protein and consists of 160 aminoacids. TH0 and TH2 cells, which represent different subpopulations of T helper cells, monocytes/macrophages and B cells are capable of IL-10 synthesis.

It has been shown that IL-10 inhibits cytokine synthesis by TH1 cells. These cells secrete cytokines like interferon- γ , TNF- β and IL-2. They tend to activate macrophages and respond to the antigens presented by these cells. Furthermore, IL-10 inhibits the stimulation of TH1 cells by macrophages and monocytes. Additionally, IL-10 has numerous other effects on the monocyte/macrophage system. It suppresses the production of the so-called inflammatory cytokines IL-1 α , IL-6, IL-8 and GM-CSF, and endogenous IL-10 inhibits the synthesis of IL-10 by monocytes/macrophages in a negative feedback mechanism. Besides these inhibitory effects, IL-10 has no effect on the production of TGF- β by monocytes and stimulates the expression of IL-1RA, another anti-inflammatory protein. In contrast to the effects mentioned above, IL-10 shows immunostimulatory effects, too. In the presence of IL-3 and IL-4 it initiates the growth of mast cells and their progenitor cells. Furthermore, IL-10 drives the proliferation and differentiation of B cells into antibody-producing cells.

Due to the inhibition of the production of inflammatory cytokines and at the same time the upregulation of anti-inflammatory proteins, like IL-1RA, IL-10 may play an important role in the regulation of a number of different diseases, including bacterial sepsis, rheumatoid arthritis and psoriasis.¹⁻⁹

Principle of the Procedure

IMMULITE/IMMULITE 1000 IL-10 is a solid-phase enzyme-labeled chemiluminescent immunometric assay.

Incubation Cycles: 1 × 60 minutes

Specimen Collection

EDTA plasma causes a depression in apparent IL-10 concentrations and is not recommended as a sample type for IL-10 determination.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Lipemic, hemolyzed or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 IL-10 has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 100 μ L serum or heparinized plasma. (Sample cup must contain at least 250 μ L more than the total volume required.)

Storage: 6 hours at 2–8°C or up to 6 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412	Harmful to aquatic life with long lasting effects.
P273, P501	Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations. Contains: 2-methyl-2H-isothiazol-3-one; IL-10 Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

IL-10 Test Units (LXP1)

Each barcode-labeled unit contains one bead coated with a murine monoclonal anti-IL-10 antibody. Stable at 2–8°C until expiration date.

LKXP1: 100 units

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

IL-10 Reagent Wedge (LXP2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to a murine monoclonal anti-IL-10 antibody in buffer. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is

within 30 days after opening when stored as indicated.

LKXP1: 1 wedge

IL-10 Adjustors (LXPL, LXPH)

Two vials (Low and High) of lyophilized IL-10 in a nonhuman serum matrix. Reconstitute each vial with **4.0 mL** distilled or deionized water. Mix by gentle, intermittent swirling. Aliquot and freeze any unused material at –20°C within 3 hours of reconstitution. Stable at –20°C for 6 months after reconstitution.

LKXP1: 1 set

Kit Components Supplied Separately

IL-8 Sample Diluent (L8PZ)

For the manual dilution of patient samples. One vial (25 mL) of an IL-8 and IL-10-free nonhuman serum matrix, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

LXPCM: IMMULITE IL-10 Control Module (Bi-level, synthetic matrix control).

Also Required

Sample transfer pipets, distilled or deionized water, controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

The sample cup must contain at least 250 µL more than the total volume required. Each sample cup holder can be followed by up to four test units.

Recommended Adjustment Interval:
2 weeks

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of IL-10.

Expected Values

A reference range study performed on 55 apparently healthy adults yielded a median of 1.5 pg/mL and a 95th percentile of 9.1 pg/mL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without anticoagulants, gel barriers or clot-promoting additives.)

Calibration Range: Up to 1,000 pg/mL

Analytical Sensitivity: 1 pg/mL

High-dose Hook Effect:

None up to 102,300 pg/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Specificity: The assay is highly specific for IL-10. (see "Specificity" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three IL-10 solutions (550, 5,500 and 11,000 pg/mL) were assayed. (See "Recovery" table for representative data.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples

Hemolysis: Grossly hemolyzed samples may show slight elevations in apparent IL-10 concentration. (See "Hemolysis" table.)

Lipemia: Presence of triglycerides in concentrations up to 2,000 mg/dL has no effect on results, within the precision of the assay. Higher concentrations may cause a depression in apparent IL-10 concentration. (See "Lipemia" table.)

Method Comparison 1: The IMMULITE IL-10 procedure was compared to a commercially available assay for IL-10 (Kit A) on 167 patient samples. (Concentration range: approximately 5 to 300 pg/mL. See graph 1.) By linear regression:

(IML) = 2.12 (Kit A) + 0.4 pg/mL
r = 0.988

Means:
61 pg/mL (IMMULITE)
29 pg/mL (Kit A)

Method Comparison 2: The calibration material used to construct the IMMULITE IL-10 master curves was compared to the IL-10 WHO reference material (WHO 93/722). (See graph 2.) By linear regression:

(IML Calibration) = 0.86 (WHO 93/722) + 0.7 pg/mL
r = 1.00

Means:
198 pg/mL (IMMULITE Calibration)
229 pg/mL (WHO 93/722)

Alternate Sample Type: Samples ($n = 26$) were collected into plain, heparinized, EDTA and SST® vacutainer tubes. Matched samples were supplemented with IL-10 and then assayed by the IMMULITE IL-10 procedure. By linear regression:

(Heparin) = 1.02 (Serum) + 2.2 pg/mL
 $r = 0.998$

(EDTA) = 0.49 (Serum) – 8.1 pg/mL
 $r = 0.982$

(SST) = 0.94 (Serum) – 1.2 pg/mL
 $r = 0.999$

Means:

279 pg/mL (Plain Serum)

286 pg/mL (Heparin)

128 pg/mL (EDTA)

260 pg/mL (SST)

EDTA plasma is not recommended as a sample type in the IMMULITE IL-10 assay.

References

- 1) Dummer W, Becker JC, Schwaaf A, Leverkus, M, Moll T, Brcker EB. Elevated serum levels of interleukin-10 in patients with metastatic malignant melanoma. Short Communication: Melanoma Res 1995;5:67-8.
- 2) Belardelli F. Role of interferons and other cytokines in the regulation of the immune response. APMIS 1995; 103:161-79.
- 3) Doherty TM. T-cell regulation of macrophage function. Curr Opin Immunol 1995;7:400-4.
- 4) Gomez-Jimenez J, Martin MC, Sauri R, et al. Interleukin-10 and the monocyte/macrophage-induced inflammatory response in septic shock. J Infect Dis 1995;171:472-5.
- 5) Stasi R, Zinzani L, Galieni P, et al. Clinical implications of cytokine and soluble receptor measurements in patients with newly-diagnosed aggressive non-Hodgkin's lymphoma. Eur J Haematol 1995;54:9-17.
- 6) Cush JJ, Splawski JB, Thomas R, et al. Elevated interleukin-10 levels in patients with rheumatoid arthritis. Arthritis Rheum 1995;38:96-104.
- 7) Ohmen JD, Hanifin JM, Nickoloff BJ, et al. Overexpression of IL-10 in atopic dermatitis. Contrasting cytokine patterns with delayed-type hypersensitivity reactions. J Immunol 1995;154:1956-63.
- 8) Krams SM, Cao S, Hayashi M, et al. Elevations in IFN- γ , IL-5, and IL-10 in patients with the autoimmune disease primary biliary cirrhosis: association with autoantibodies and soluble CD30. Clin Immunol Immunopathol 1996;80:311-20.
- 9) Kim DS, Lee HK, Noh GW, et al. Increased serum interleukin-10 level in Kawasaki disease. Yonsei Med J 1996;37:125-30.
- 10) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (pg/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	18.2	0.52	2.9%	1.8	9.9%
2	46.0	1.4	3.0%	2.2	4.8%
3	177	6.1	3.4%	8.0	4.5%
4	258	8.0	3.1%	12	4.7%
5	358	10	2.8%	17	4.7%
6	617	19	3.1%	26	4.2%

Linearity (pg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	145	—	—
	4 in 8	69.4	72.5	96%
	2 in 8	36.0	36.3	99%
	1 in 8	17.3	18.1	96%
2	8 in 8	275	—	—
	4 in 8	141	138	102%
	2 in 8	73.5	68.8	107%
	1 in 8	39.1	34.4	114%
3	8 in 8	455	—	—
	4 in 8	233	228	102%
	2 in 8	116	114	102%
	1 in 8	57.0	56.9	100%
4	8 in 8	648	—	—
	4 in 8	316	324	98%
	2 in 8	161	162	99%
	1 in 8	83.8	81.0	103%
5	8 in 8	900	—	—
	4 in 8	439	450	98%
	2 in 8	232	225	103%
	1 in 8	118	113	104%

Recovery (pg/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	ND	—	—
	A	30.3	27.5	110%
	B	276	275	100%
	C	534	550	97%
2	—	27.9	—	—
	A	55.2	54.0	102%
	B	306	302	101%
	C	637	577	110%
3	—	188	—	—
	A	213	206	103%
	B	426	454	94%
	C	690	729	95%
4	—	339	—	—
	A	353	350	101%
	B	581	597	97%
	C	819	872	94%

ND: not detectable⁵

Specificity



Compound ¹	pg/mL Added ²	Apparent pg/mL ³	% Cross-reactivity ⁴
IL-1 α	10,000	ND	ND
IL-1 β	10,000	ND	ND
IL-2	10,000	ND	ND
IL-4	10,000	ND	ND
IL-6	10,000	ND	ND
IL-8	10,000	ND	ND
H-IFN γ	10,000	ND	ND
TNF α	10,000	ND	ND

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IMMULITE is a trademark of
All rights reserved.
CONTROL AB En: Control Antibody