Total IgE

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Total IgE

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of immunoglobulin type E (IgE) in serum.

Catalog Number: **L2KIE2** (200 tests)
L2KIE6 (600 tests)

Test Code: **TIE** Color: **Violet**

Summary and Explanation

Modern laboratory methods for allergy testing are based on the comparatively recent discovery that many allergies are mediated by immunoglobulins of the IgE class acting as points of contact between the allergen and specialized cells.^{4,5,7,9,10}

The IgE molecules, which have a molecular mass of approximately 200,000, bind to the surface of mast cells and basophilic granulocytes.^{2,3} Subsequent binding of allergens to cell-bound IgE causes these cells to release histamines and other vasoactive substances, thereby initiating the events which we recognize as an allergic reaction.¹⁰

In deciding on a course of therapy, it is important to distinguish between IgE-mediated and non-IgE-mediated reactions.^{4,7} Measurement of the total circulating IgE level, in conjunction with other supporting diagnostic information, can aid in making this diagnosis. The supporting information should include appropriate tests for allergen-specific IgE. Measurement of the total circulating IgE level may also be of value in the early detection of allergy in infants, and as a means for predicting future atopic manifestations.^{4,6,7}

IgE levels normally show a slow increase during childhood, reaching adult levels in the second decade of life.^{1,8} In general, the total IgE level increases with the number of allergies which a person has and with the amount of exposure to relevant allergens.^{5,7}

Significant elevations may be encountered not only in sensitized individuals, but also

in cases of IgE myeloma, pulmonary aspergillosis, and during the active stage of parasitic infestations.^{4,5,7,9}

Principle of the Procedure

IMMULITE 2000 Total IgE is a solid-phase, chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes

Time to First Result: 35 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Total IgE has not been tested with all possible variations of tube types.

Volume Required: 5 µL serum

Storage: 3 days at 2–8°C,^{13,15} or 6 months at –20°C.¹⁵

Warnings and Precautions

For *in vitro* diagnostic use.

**CAUTION! POTENTIAL BIOHAZARD**

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.¹⁶⁻¹⁸

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Total IgE Bead Pack (L2IE12)

With barcode. 200 beads, coated with monoclonal mouse anti-IgE. Stable at 2–8°C until expiration date.

L2KIE2: 1 pack **L2KIE6:** 3 packs

Total IgE Reagent Wedge (L2IEA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-IgE antibody in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KIE2: 1 wedge **L2KIE6:** 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Total IgE Adjustors (LIEL, LIEH)

Two vials (Low and High), 2.5 mL each, of IgE in a nonhuman serum/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KIE2: 1 set **L2KIE6:** 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately**Multi-Diluent 2 (L2M2Z, L2M2Z4)**

For on-board dilution of samples. One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C. Dispose of in accordance with applicable laws.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M2Z: 3 labels **L2M2Z4:** 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

IECM: Bi-level, human serum-based control module containing IgE

Also Required

Distilled or deionized water; test tubes; controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks

Quality Control Samples: Use controls or serum pools with at least two levels (low and high) of total IgE.

Expected Values

Based on its relationship to IMMULITE Total IgE (see Method Comparison), the assay can be expected to have essentially the same reference ranges for nonatopic individuals.

Age (years)	Median (IU/mL)	95%ile (IU/mL)
0–1	6.6	29
1–2	10.1	49
2–3	12.9	45
3–9	14.4	52
Adult	20.4	87

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Serum total IgE levels for the majority of individuals with IgE-mediated disease can be expected to be elevated compared to the reference range for healthy adults. However, not all allergic patients exhibit elevated serum total IgE levels.¹¹

Since not all atopic reactions are IgE-mediated, a total IgE result in the reference range should always be interpreted in light of other clinical observations.

Tests for allergen-specific IgE may clarify the finding of an elevated total IgE result.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing

interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in IU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

$\text{IU/mL} \times 2.4 \rightarrow \text{ng/mL}$

$\text{IU/mL} \times 1 \rightarrow \text{kIU/L}$

Calibration Range: up to 2000 IU/mL. Standardized in reference to WHO 2nd IRP 75/502.

Analytical Sensitivity: 1.0 IU/mL

High-dose Hook: No effect up to 13,000 IU/mL

Precision: Within-Run (Intraassay) precision data was generated from 20 replicates in a single run. Total (Interassay) precision data was generated from 2 replicates per run across 29 different runs. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. The results show that the IMMULITE 2000 Total IgE kit maintains good linearity throughout its calibration range; but due to the heterogeneity of circulating IgE with different affinities, there may be a loss of dilutional parallelism in certain patients. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three total IgE solutions (1750, 3500 and 7000 IU/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The assay is highly specific for human IgE and exhibits no crossreactivity to other human immunoglobulin classes.

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of hemoglobin in concentrations up to 384 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 5000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison: The assay was compared to IMMULITE Total IgE on 175 samples. (Concentration range: approximately 2 to 1600 IU/mL. See graph.) By linear regression:

$$(IML\ 2000) = 0.96 (IML) + 10.0\ IU/mL$$
$$r = 0.988$$

Means:

225 IU/mL (IMMULITE 2000)

223 IU/mL (IMMULITE)

References

1) Barbee RA, et al. Distribution of IgE in a community population sample: correlation with age, sex and allergen skin test reactivity. J Allergy Clin Immunol 1981;68:106–11. 2) Geha RS. Human IgE. J Allergy Clin Immunol 1984;74:109–20. 3) Gordon RR, et al. Immunoglobulin E and the eczema-asthma syndrome in early childhood. Lancet 1982;1: 72–74. 4) Halpern GM. Markers of human allergic disease. J Clin Immunoassay 1983 Jun;6(2):131–39. 5) Hamilton RG, Adkinson NF. Clinical laboratory methods in allergic disease. Lab Management 1983 Dec;21(12):37–50. 6) Hamilton RG, Adkinson NF. Quantitation of allergen-specific IgE in serum using the radioallergosorbent test. J Clin Immunoassay 1983 Jun;6(12):147–54. 7) Homberger HA, Yunginger JW. Laboratory testing in the diagnosis and management of allergic diseases. Clin Lab Annual 1983;2:351–88. 8) Kjellman N-IM, Johansson SGO, Roth A. Serum IgE levels in healthy children. Clin Allergy 1976;6:51–59. 9) Mandy FF, Perelmutter L. Laboratory measurement of total human serum IgE. J Clin Immunoassay 1983 Jun;6(2):140–46.

10) Norman PS. The clinical significance of IgE. Hosp Pract 1975 Aug;10(8):41–9. 11) Yunginger JW. Clinical significance of IgE. In: Middleton E, et al, editors. Allergy principles and practice. Vol 2. 3rd ed. St. Louis: C.V. Mosby, 1988:849–60. 12) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 13) Bazara M, Hamburger RN. Standardization and stability of immunoglobulin E (IgE). J Allerg Clin Immunol 1972;49:189–91. 14) Nye L, Merrett TG, Landon J, White RJ. A detailed investigation of circulating IgE levels in a normal population. Clin Allergy 1975;1:13–24. 15) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:358. 16) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377–82, 387–8. 17) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 18) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

In the United States, contact
Siemens Healthcare Diagnostics
Technical Services department.
Tel: 877.229.3711. Outside the United
States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare
Diagnostics Products Ltd. is certified to ISO
13485.

Tables and Graphs


Precision (IU/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	50	1.65	3.3%	2.54	5.1%
2	107	4.78	4.5%	7.20	6.7%
3	223	9.85	4.4%	15.0	6.7%
4	600	24.4	4.1%	38.3	6.4%

Linearity (IU/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	14	—	—
	4 in 8	6.4	7.0	91%
	2 in 8	3.8	3.5	109%
	1 in 8	1.8	1.8	100%
2	8 in 8	21	—	—
	4 in 8	11	11	100%
	2 in 8	5.4	5.3	102%
	1 in 8	3.1	2.6	119%
3	8 in 8	195	—	—
	4 in 8	106	98	108%
	2 in 8	47	49	96%
	1 in 8	26	24	108%
4	8 in 8	267	—	—
	4 in 8	133	134	99%
	2 in 8	68	67	102%
	1 in 8	34	33	103%
5	8 in 8	411	—	—
	4 in 8	207	206	101%
	2 in 8	97	103	94%
	1 in 8	46	51	90%
6	8 in 8	823	—	—
	4 in 8	415	412	101%
	2 in 8	225	206	109%
	1 in 8	121	103	118%
7	8 in 8	838	—	—
	4 in 8	381	419	91%
	2 in 8	201	210	96%
	1 in 8	107	105	102%
8	8 in 8	1050	—	—
	4 in 8	536	525	102%
	2 in 8	255	263	97%
	1 in 8	134	131	102%

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
All rights reserved.
CONTROL AB