# IMMULITE®

## Anti-TPO Ab

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**SIEMENS**

# IMMULITE®/IMMULITE 1000 Anti-TPO Ab

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of antibodies against thyroid peroxidase (TPO) in serum and EDTA plasma, as an aid in the clinical diagnosis of thyroid diseases.

Catalog Number: **LKTO1** (100 tests)

Test Code: **ATA** Color: **Aqua**

### Summary and Explanation

Antithyroid peroxidase (TPO) antibodies are autoantibodies directed against the thyroid peroxidase enzyme. This enzyme catalyzes the iodination of tyrosine in thyroglobulin during the biosynthesis of T3 and T4.<sup>2</sup> Historically, these antibodies were referred to as antimicrosomal antibodies (AMA) because the antibodies bind to the microsomal part of the thyroid cells. Recent research has identified thyroid peroxidase as the primary antigenic component of microsomes.<sup>2,3,13</sup>

Thyroid autoimmune disease is the major factor underlying hypothyroidism and hyperthyroidism and tends to occur in a genetically predisposed population.

The major thyroid autoimmune diseases are: Hashimoto's thyroiditis and Graves' disease.<sup>13</sup> In virtually all cases of Hashimoto's disease and in the majority of Graves' disease cases, TPO autoantibodies are elevated.<sup>13</sup> High levels of TPO autoantibodies, in the context of the clinical presentation of hypothyroidism, confirms the diagnosis of Hashimoto's disease.

### Principle of the Procedure

IMMULITE/IMMULITE 1000 Anti-TPO Ab is a solid-phase, enzyme-labeled, chemiluminescent sequential immunometric assay.

**Incubation Cycles:** 2 × 30 minutes.

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Anti-TPO Ab has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 50 µL of the prediluted patient sample (1-in-101 dilution). (Sample cup must contain at least 100 µL more than the total volume required.)

**Storage:** 2 days at 2–8°C. For longer storage, aliquot and freeze at –20°C.

**Patient Sample Dilution:** Samples must be prediluted 1-in-101 in Thyroid Autoantibody Sample Diluent (LAAZ4), e.g. by adding 10 µL sample to 1,000 µL Diluent.

Note that printed results have been corrected for the 1-in-101 dilution. Samples expected to have values outside the calibration range should be further diluted. Correction for further dilution must be entered manually.

## Warnings and Precautions

For *in vitro* diagnostic use.

### CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.<sup>16-18</sup>

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**Danger!** Toxic in contact with skin. Harmful if swallowed. Toxic to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.

H311, H302, H411

P280, P273, P301 + P312, P302 + P312, P501

**Contains:** sodium azide, 2-methyl-2H-isothiazol-3-one; Anti-TPO Ab Adjustors

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### Anti-TPO Ab Test Units (LTO1)

Each barcode-labeled unit contains one bead coated with highly purified human TPO. Stable at 2–8°C until expiration date.

**LKTO1:** 100 units

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

### Anti-TPO Ab Reagent Wedges (LTOA, LTOB)

With barcodes. **LTOA:** one wedge (7.5 mL) of buffer matrix, with preservatives. **LTOB:** one wedge (7.5 mL) alkaline phosphatase (bovine calf intestine) conjugated to a monoclonal murine anti-human IgG antibody in buffer, with preservatives. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

**LKTO1:** 1 set

### Anti-TPO Ab Adjustors (LTOL, LTOH)

Two vials (Low and High), of lyophilized TPO autoantibodies in a human serum/buffer matrix, with preservative. At least 30 minutes before use, reconstitute each vial by adding **4.0 mL** distilled or deionized water. Mix by *gentle* swirling or inversion until the lyophilized material is fully dissolved. (No further dilution is required.) Aliquot and freeze in tightly capped plastic containers. Avoid repeated freeze/thaw cycles. Stable at 2–8°C for 7 days after reconstitution, or for 6 months (aliquotted) at –20°C.

**LKTO1:** 1 set

## Thyroid Autoantibody Sample Diluent (LAAZ4)

For the manual dilution of patient samples. 100 mL anti-TG/anti-TPO antibody-free buffer matrix, with preservative. Stable at 2–8°C for 6 months after opening.

**LKTO1:** 1 vial

The specific Thyroid Autoantibody Sample Diluent (LAAZ4) lot should always be used in combination with the reagents from the kit lot with which it was provided.

## Kit Components Supplied Separately

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash Module

**LKPM:** Probe Cleaning Kit

**LCHx-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)

**LSCC:** Sample Cup Caps (optional)

**LAACM:** IMMULITE Thyroid Autoantibody Control Module

Also Required

Sample transfer pipets, pipets for predilution of patient samples, distilled or deionized water, controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, adjustment, assay and quality control procedures.

Note that both Reagent Wedges A and B must be loaded on the carousel to run this assay.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

**Recommended Adjustment Interval:**

2 weeks.

**Quality Control Samples:** Use controls or sample pools with at least two levels (low and high) of anti-TPO Ab.

## Expected Values

Thyroid peroxidase autoantibodies are unique for each patient; therefore, each patient sample will show its own individual

binding characteristics. Normal levels for most patients have been defined as less than 35 IU/mL for the IMMULITE Anti-TPO Ab assay. A numerical result greater than or equal to 35 IU/mL indicates an *elevated* anti-TPO level in serum, whereas a result less than 35 IU/mL indicates a *normal* anti-TPO level in serum. However, due to individual variability it is also possible that a normal patient with no history of thyroid disease may have reproducible elevated values greater than 35 IU/mL.

Two studies are presented below that demonstrate the variability between two random populations.

Using the IMMULITE Anti-TPO Ab assay, the levels of TPO autoantibodies were determined for 75 apparently healthy subjects (26 males and 49 females), with ages ranging from 13 to 87 years. All had normal TSH levels, and none had thyroid autoimmune disorders. Using the recommended 35 IU/mL cutoff, 79% (59/75) of the subjects had normal anti-TPO values. The "Expected Values 1 ( $n = 75$ )" histogram (see Tables and Graphs) shows a distribution of subjects for different levels of TPO autoantibodies.

In another study, the IMMULITE Anti-TPO Ab procedure was used to determine the TPO autoantibody levels of 130 apparently healthy subjects. Using the recommended 35 IU/mL cutoff, 92% (118/130) of the subjects had normal anti-TPO values. The "Expected Values 2 ( $n = 130$ )" histogram (see Tables and Graphs) shows a distribution of subjects for different levels of TPO autoantibodies.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

## Limitations

For diagnostic purposes, anti-TPO results should be used in conjunction with other test results, the overall clinical presentation to the physician, and all other appropriate information. Autoantibodies may be found in less than 10% of the normal population at low levels and in patients with non-thyroidal illnesses, such as the inflammatory rheumatic diseases.<sup>14</sup>

Patient TPO autoantibodies and the antibodies used for calibrating the assay may differ in their binding properties. This

may lead to deviations from strict dilutional parallelism for certain patient samples.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in IU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Calibration Range:** Up to 1,000 IU/mL (WHO 1st IRP 66/387).

**Analytical Sensitivity:** 7 IU/mL

### High-dose Hook Effect:

None up to 18,500 IU/mL.

### Intraassay Precision (Within-Run):

Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

### Interassay Precision (Run-to-Run):

Statistics were calculated for samples assayed in 20 different runs. (See "Interassay Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.) Due to the heterogeneity of TPO autoantibodies, certain samples may not maintain dilutional parallelism in the IMMULITE Anti-TPO Ab procedure.

**Recovery:** Samples spiked 1 to 19 with three anti-TPO Ab solutions (2,320, 4,740 and 11,480 IU/mL) were assayed. (See "Recovery" table for representative data.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Biotin:** Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

**Hemolysis:** Presence of packed red blood cells in concentrations up to 30  $\mu$ L/mL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** Blood was collected from 10 laboratory volunteers into plain and EDTA vacutainer tubes. All samples were spiked with various amounts of TPO autoantibodies, with a spiking ratio not exceeding 1 to 10, and then assayed by the IMMULITE Anti-TPO Ab procedure, after the recommended 1-in-101 dilution.

(EDTA) = 0.97 (Serum) + 12 IU/mL  
 $r = 0.997$

Means:  
325 IU/mL (Serum)  
327 IU/mL (EDTA)

In another study, blood was collected from 18 laboratory volunteers into plain, heparinized and Becton Dickinson SST<sup>®</sup> vacutainer tubes. All samples were spiked with various amounts of TPO autoantibodies, with a spiking ratio not exceeding 1 to 10, and then assayed by the IMMULITE Anti-TPO Ab procedure, after the recommended 1-in-101 dilution.

(Heparin) = 0.94 (Serum) + 8.3 IU/mL  
 $r = 0.984$

(SST) = 1.0 (Plain Tubes) + 2.1 U/mL  
 $r = 0.977$

Means:  
281 IU/mL (Serum)  
273 IU/mL (Heparin)  
286 U/mL (SST)

**Method Comparison:** The assay was compared to an enzyme immunometric assay on 94 samples. (Concentration range: from 0 to > 2,000 IU/mL, as measured by Kit A.) (See "Method Comparison" table.) The clinical classification of 89 out of 91 samples agreed by both methods.

**References**

1) Becker W, Reimers C, Brner W. Immunological criteria for the differential diagnosis of thyroid disease. *Therapiewoche* 1985;35:1167-76. 2) Czarnocka B, Ruf J, Ferrand M, et al. Purification of the human thyroid peroxidase and its identification as the microsomal antigen involved in autoimmune thyroid diseases. *FEBS lett.* 1985;190:147-52. 3) Czarnocka B, Ruf J, Ferrand M, et al. Interaction of highly purified thyroid peroxidase with anti-microsomal antibodies in autoimmune thyroid diseases. *J Endocrinol Invest* 1986;9:135-8. 4) Ruf J, Czarnocka B, Ferrand M, et al. Thyroid peroxidase is the organ-specific 'microsomal' antigen involved in thyroid autoimmunity. *Acta Endocrinol (Copenh)* 1987;Suppl. 281: 49-55. 5) Ruf J, Czarnocka B, Ferrand M, et al. Relationship between immunological structure and biochemical properties of human thyroid peroxidase. *Endocrinol* 1988;125 (3):1211-8. 6) Horster FA. The importance of MAK, TAK, TRAK and thyroglobulin in the diagnosis of thyroid disease. *Internist* 1988;29:538-40. 7) Schatz H, Lbig H. The diagnostic and prognostic importance of antibodies against thyroid microsomes and thyroglobulin. *Akt Endokr Stoffw* 1989;10:146-53. 8) Hashimoto H. Understanding the lymphatic changes of the thyroid. *Arch Klin Chir* 1912;97:219-48. 9) Libert F, Ruel J, Ludgate M, et al. Thyroperoxidase, an auto-antigen with a mosaic structure made of nuclear and mitochondrial gene modules. *EMBO J* 1987;6:4193-6. 10) McKenzie JM, Zakarija M, Sato A. Humoral immunity in Graves' disease. *Clin Endocrinol Metab* 1978;7:31-45. 11) Scherbaum WA, Berg PA. The importance of autoantibodies in the diagnosis of endocrine disease. *Dtsch Med Wschr* 1981;106:308-13. 12) Scherbaum WA, Berg PA. Thyroid autoimmune disease – new aspects in development and diagnosis. *Dtsch Med Wschr* 1984;109:1574-81. 13) Feldt-Rasmussen U. Analytical and clinical performance goals for testing autoantibodies to thyroperoxidase, thyroglobulin, and thyrotropin receptor. *Clin Chem* 1996;42:160-3. 14) Wild D, editor. *The immunoassay handbook*. Great Britain: Stockton Press, 1994: 342-3. 15) National Committee for Clinical Laboratory Standards. *Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard*. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 16) Centers for Disease Control. Update:

Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377–82, 387–8. 17) Clinical and Laboratory Standards Institute (formerly NCCLS). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 18) Federal Occupational Safety and Health Administration, *Bloodborne Pathogens Standard*, 29 CFR 1910.1030.

**Technical Assistance**

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

---

**Tables and Graphs**

---

**Intraassay Precision (IU/mL)**

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	54	2.3	4.3%
2	185	6.5	3.5%
3	825	46	5.6%

**Interassay Precision (IU/mL)**

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	57	6.0	10.5%
2	183	16	8.7%
3	345	27	7.8%
4	619	65	10.5%


### Linearity (IU/mL)

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	1:100	119	—	—
	1:200	62	60	103%
	1:400	32	30	107%
	1:800	14	15	93%
2	1:100	307	—	—
	1:200	174	154	113%
	1:400	83	77	108%
	1:800	33	38	87%
3	1:100	716	—	—
	1:200	396	358	111%
	1:400	191	179	107%
	1:800	94	90	104%
	1:1,600	46	45	102%


### Recovery (IU/mL)

	Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	23	—	—
	A	138	138	100%
	B	232	259	90%
	C	497	596	83%
2	—	34	—	—
	A	140	148	95%
	B	264	269	98%
	C	519	606	86%
3	—	71	—	—
	A	164	184	89%
	B	276	305	90%
	C	554	642	86%
4	—	146	—	—
	A	264	255	104%
	B	375	376	100%
	C	662	713	93%

### Expected Values 1 (n = 75)

### Expected Values 2 (n = 130)

### Method Comparison (IU/mL)

Kit A <sup>1</sup>	IMMULITE Anti-TPO Ab			
	Normal <35 <sup>2</sup>	Elevated ≥35 <sup>2</sup>	Relative Sensitivity <sup>4</sup>	Relative Specificity <sup>5</sup>
Normal <sup>6</sup> < 100	43	2		
Borderline <sup>7</sup> 100–150	1	2	100%	95.6%
Elevated <sup>8</sup> > 150	0	46		

Agreement: 97.8%<sup>9</sup>

95% Confidence Limits for Relative Sensitivity and Specificity, respectively: 92.3% – 100% and 84.9% – 99.5%. (Borderline results omitted in analysis.)

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2008 Siemens Healthcare Diagnostics. All rights reserved.
**CONTROL AB**