IMMULITE®

Progesterone

For use on the IMMULITE®
and IMMULITE 1000 systems

SIEMENS

IMMULITE®/IMMULITE 1000 Progesterone

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® and IMMULITE 1000 Analyzers — for the quantitative measurement of progesterone in serum as an aid in the diagnosis and treatment of disorders of the ovaries or placenta.

Catalog Number: **LKPW1** (100 tests)

Test Code: **PRG** Color: **Light green**

Summary and Explanation

Progesterone is a steroid hormone which plays an important role in the preparation for and maintenance of pregnancy. It is synthesized from cholesterol via pregnenolone — then rapidly metabolized to pregnenediol, for the most part, in the liver.^{1,8,12} The ovary and placenta are the major production sites; but a small amount is also synthesized by the adrenal cortex in both men and women.

Circulating progesterone levels, which are characteristically low during the follicular phase, increase sharply during the luteal phase of menstrual cycles, reaching a maximum some 5 to 10 days after the midcycle LH peak.¹¹ Unless pregnancy occurs, a steep decline to follicular levels sets in about 4 days before the next menstrual period. This pattern constitutes the rationale behind the well-established use of serum progesterone measurements as a simple and reliable method for ovulation detection.^{2,3,15}

There is growing literature on luteal phase defects.^{4,7,10,16,17} Daily progesterone levels are considered the most accurate means for documenting a defective luteal phase.^{4,7} However, some investigators have found that three samples or even a single sample^{13,17} (if well-timed) can provide valuable information on the adequacy of the luteal phase.

Measurements of serum progesterone have also been used to check the effectiveness of ovulation induction,¹⁵ to monitor progesterone replacement therapy¹⁵ and to detect and evaluate patients at risk for abortion during the early weeks of pregnancy.^{5,6,14}

On the other hand, although progesterone levels increase throughout pregnancy, they are not considered a suitable means for monitoring fetal well-being during the third trimester.⁸

Principle of the Procedure

IMMULITE/IMMULITE 1000 Progesterone is a solid phase, competitive immunoassay using enzyme-labeled chemiluminescent technology. The solid phase (bead) is coated with polyclonal rabbit anti-progesterone antibody. The liquid phase consists of alkaline phosphatase (bovine calf intestine) conjugated to progesterone.

The patient sample and the reagent are incubated together with the coated bead for 30 minutes. During this time, progesterone in the sample competes with enzyme-conjugated progesterone in the reagent for a limited number of antibody binding sites on the bead. Unbound sample and enzyme conjugate are then removed by centrifugal washes. Finally, chemiluminescent substrate is added to the test unit containing the bead and the signal is generated in proportion to the bound enzyme.

Incubation Cycles: 1 × 30 minutes

Time to First Result: 42 minutes

Specimen Collection

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on tube materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Progesterone

has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

EDTA: Because EDTA would have a significant effect on results, it should not be used as an anticoagulant.

Gel Barrier Tubes: Time-dependent decreases in progesterone levels have been reported when serum samples are collected and stored in gel barrier tubes.¹⁸⁻²⁰

Lipemia: Gross lipemia has been shown to cause a decrease in the apparent progesterone concentration. The use of an ultracentrifuge is recommended to clear lipemic samples.

Volume Required: 25 µL serum. (Sample cup must contain at least 100 µL more than the total volume required.)

Storage: 7 days at 2–8°C or 6 months at –20°C.

Dilution of High Samples: All samples expected to have levels above the assay's calibration range should be manually diluted before assay. Correction for dilution must be entered manually.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.²²⁻²⁴

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Caution: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Safety data sheets (MSDS) available on www.siemens.com/diagnostics.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The bar-code labels are needed for the assay.

Progesterone Test Units (LPW1)

Each bar-code-labeled unit contains one bead coated with polyclonal rabbit anti-progesterone. Stable at 2–8°C until expiration date.

LKPW1: 100 units

Allow the temperature of the Test Unit bags to reach room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

Progesterone Reagent Wedge (LPW2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to progesterone, in buffer. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKPW1: 2 wedges

Progesterone Adjustors (LPWL, LPWH)

Two 3 mL vials (Low and High) of progesterone in processed human serum with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKPW1: 1 set

Kit Components Supplied Separately

Progesterone Sample Diluent (LPWZ)

For the manual dilution of patient samples. 25 mL of processed progesterone-free human serum with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (bar-coded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

Tri-level, multi-constituent Lyphochek® Control available from Bio-Rad.

Also Required

Sample transfer pipets, distilled or deionized water, controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:
2 weeks

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of progesterone.

Follow government regulations or accreditation requirements for quality control frequency.

Expected Values

A study involving women in apparent good health (age: 21–51 years), who volunteered

to have blood samples drawn on a daily basis throughout one complete ovulatory cycle yielded the following results:

Progesterone, ng/mL			
Ovulatory Cycles	n*	Median	Absolute Range
Follicular Phase	31 (174)	0.32	ND–13.4
Midfollicular Days 5 to 11	31 (80)	0.29	ND–1.33
Midcycle	18 (18)	1.06	0.48–2.01
Luteal Phase	31 (122)	7.70	ND–26.6
Midluteal, Days 7 to 8 of Luteal Phase	17 (18)	11.3	2.83–26.1

*Number of subjects (total number of results)

Progesterone, nmol/L			
Ovulatory Cycles	n*	Median	Absolute Range
Follicular Phase	31 (174)	1.02	ND–42.6
Midfollicular Days 5 to 11	31 (80)	0.92	ND–4.23
Midcycle	18 (18)	3.37	1.53–6.39
Luteal Phase	31 (122)	24.5	ND–84.6
Midluteal, Days 7 to 8 of Luteal Phase	17 (18)	35.9	9.00–83.0

*Number of subjects (total number of results)

See Menstrual Cycle Graph (see Tables and Graphs).

Another study yielded the following results:

Progesterone, ng/mL			
Group	n	Median	Absolute Range
Adult Males	100	0.44	ND–1.20
Adult Females:			
Follicular Phase	31	0.32	ND–12.8
Luteal Phase	31	7.70	ND–25.3
Postmenopausal	55	ND	ND–11.0
Oral Contraceptives	50	0.48	ND–8.11
Pregnant Females:			
First Trimester	29	24.4	12.5–32.7
Second Trimester	50	35.0	15.4–139
Third Trimester	51	60.1	10.1–255

ND: not detectable

Group	n	Progesterone, nmol/L	
		Median	Absolute Range
Adult Males	100	1.40	ND–3.82
Adult Females:			
Follicular Phase	31	1.02	ND–40.7
Luteal Phase	31	24.5	ND–80.5
Postmenopausal	55	ND	ND–35.0
Oral Contraceptives	50	1.53	ND–25.8
Pregnant Females:			
First Trimester	29	77.6	39.8–104
Second Trimester	50	111.3	49.0–442
Third Trimester	51	191.1	32.1–811

ND: not detectable

In pregnancy, the general tendency is for values to increase. There is considerable interpersonal variation in progesterone values, particularly in groups associated with elevated levels. (Note that measurement of progesterone levels is generally considered unsuitable for monitoring fetal well-being in the later weeks of pregnancy.⁸)

Progesterone concentrations were measured in pediatric serum samples using the IMMULITE/IMMULITE 1000 assay. The observed values and ranges for each population represented are shown below:

Group	Age (yr)	n	Progesterone, ng/mL	
			Median	Absolute Range
Females				
	0.1–0.4	11	1.49	ND–76.4
	0.5–1.0	21	0.76	ND–0.90
	1.1–9.0	45	ND	ND–2.42
Males				
	0.1–0.4	18	0.83	0.24–34.1
	0.5–1.0	24	ND	ND–1.04
	1.1–9.0	45	ND	ND–1.36
Combined				
	0.1–0.4	29	0.84	ND–76.4
	0.5–1.0	45	ND	ND–1.04
	1.1–9.0	90	ND	ND–2.42

ND: not detectable

Group	Age (yr)	n	Progesterone, nmol/L	
			Median	Absolute Range
Females				
	0.1–0.4	11	4.74	ND–243
	0.5–1.0	21	2.42	ND–2.86
	1.1–9.0	45	ND	ND–7.70
Males				
	0.1–0.4	18	2.64	0.76–108
	0.5–1.0	24	ND	ND–3.31
	1.1–9.0	45	ND	ND–4.32
Combined				
	0.1–0.4	29	2.67	ND–243
	0.5–1.0	45	ND	ND–3.31
	1.1–9.0	90	ND	ND–7.70

ND: not detectable

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

The supplement Dehydroepiandrosterone (DHEA) may cause falsely elevated progesterone results in immunoassays. At an initial concentration of Progesterone of 1.20 ng/mL, a 91.7% change in concentration was observed at the supraphysiologic level of DHEAS (DHEA metabolite) of 20,000 ng/mL. For patients being treated with DHEA, an alternate method that is not expected to show cross reactivity to DHEAS (DHEA metabolite), such as Liquid Chromatography-Mass Spectrometry (LC-MS), should be used.

With the advent of new steroid based medications (analogues) with similar chemical structures to progesterone, there is the possibility of cross-reactivity and falsely elevated results. For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings. If the progesterone results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.

Performance Data

See Tables and Graphs for data representative of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

ng/mL $\times$ 3.18 $\rightarrow$ nmol/L

Reportable Range: 0.20–40 ng/mL
(0.64–127 nmol/L)

Limit of Detection: 0.20 ng/mL
(0.64 nmol/L)

Limit of Blank: 0.05 ng/mL
(0.16 nmol/L)

Functional Sensitivity: 0.46 ng/mL
(1.46 nmol/L)

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See *Precision* table.)

Linearity: Samples were assayed under various dilutions. (See *Linearity* table for representative data.)

Recovery: Samples spiked 1 to 19 with three progesterone solutions (40, 80 and 160 ng/mL) were assayed. (See *Recovery* table for representative data.)

Specificity: The antibody is highly specific for progesterone. (See *Specificity* table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Gross lipemia has been shown to cause a decrease in the apparent

progesterone concentration.
(See *Lipemia* table.)

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 40 volunteers into plain serum, heparinized, EDTA and Becton Dickinson SST vacutainer glass tubes. Equal volumes of the matched samples were spiked with various concentrations of Progesterone, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE/IMMULITE 1000 Progesterone procedure.

(Lithium Heparin) = 1.02 (Serum) – 0.03 ng/mL
 $r = 0.994$

(EDTA) = 1.21 (Serum) + 0.24 ng/mL
 $r = 0.999$

(SST) = 1.06 (Plain Tubes) – 0.04 ng/mL
 $r = 0.999$

Means:

4.16 ng/mL (Serum)
4.21 ng/mL (Heparin)
5.28 ng/mL (EDTA)
4.38 ng/mL (SST)

Method Comparison: The IMMULITE/IMMULITE 1000 Progesterone procedure was compared the IMMULITE 2000 Progesterone on 172 patient samples (Concentration range: approximately 0.20 to 40 ng/mL. See *Method Comparison* graph.)
By linear regression:

(IML) = 0.99 (IML 2000) + 0.31 ng/mL
 $r = 0.991$

Means:

7.81 ng/mL (IMMULITE)
7.56 ng/mL (IMMULITE 2000)

References

- 1) Aufrere MB, Benson H. Progesterone: an overview and recent advances. *J Pharm Sci* 1976; 65:783-800.
- 2) Bauman J. Basal body temperature: unreliable method of ovulation detection. *Fertil Steril* 1981; 36:729-33.
- 3) Brown JB. Timing of ovulation. *Med J Austral* 1977; ii:780-83.
- 4) Gautray JP, et al. Clinical investigation of the menstrual cycle: clinical, endometrial and endocrine aspects of luteal defect. *Fertil Steril* 1981; 35:296-303.
- 5) Hensleigh PA, Fainstat T. Corpus luteum dysfunction: serum progesterone levels in diagnosis and assessment of therapy for recurrent and threatened abortion. *Fertil Steril* 1979; 32:396-9.
- 6) Hernandez Horta JL, et al. Direct evidence of luteal insufficiency in women with habitual abortion. *Obstet Gynecol* 1977; 49:705-8.
- 7) Jones G. Luteal phase defects. In: Behrman SJ, Kistner RW, editors. *Progress in infertility*.

Boston: Little & Brown, 2nd Edition, 1975: 299-324. 8) Klopfer A, Fuchs F. Progestagens. In: Fuchs F, Klopfer A, editors. Endocrinology of pregnancy. Hagerstown: Harper & Row, 1977: 99-122. 9) Lehmann F, Bettendorf G. The endocrine shift from a normal cycle to anovulation. In: Insler V, Bettendorf G, editors. Advances in diagnosis and treatment of infertility. Amsterdam: Elsevier/North Holland, 1981: 105-13. 10) March CM. Luteal phase defects. In: Mishell DR, Davajan V, editors. Reproductive endocrinology, infertility and contraception. Philadelphia: F.A. Davis, 1979: 469-76. 11) March CM, Goebelsmann U, Nakamura RM and Mishell DR. Roles of estradiol and progesterone in eliciting the midcycle luteinizing hormone and follicle-stimulating hormone surges. *J Clin Endocrinol Metab* 1979; 49:507-13. 12) Progesterone (Rochester: Bioeducational Publications, 1981). A BIO-ED slide/seminar educational program. 13) Radwanska E, et al. Plasma progesterone and oestradiol estimations in the diagnosis and treatment of luteal insufficiency in menstruating infertile women. *Acta Eur Fertil* 1976; 7:39-47. 14) Radwanska E, et al. Plasma progesterone levels in normal and abnormal early human pregnancy. *Fertil Steril* 1978; 30:398-402. 15) Radwanska E, et al. Single midluteal progesterone assay in the management of ovulatory infertility. *J Reprod Med* 1981; 26:85-89. 16) Sheehan KL, Casper RF, Yen SSC. Luteal phase defects induced by an agonist of luteinizing hormone-releasing factor: a model for fertility control. *Science* 1982; 215:170-72. 17) Wentz A. Pathophysiology of luteal phase inadequacy. In: Tozzini RI, Reeves G and Pineda RL, editors. Endocrine physiopathology of the ovary. Amsterdam: Elsevier/North Holland, 1980, 257-74. 18) Hilborn S, Krahn J. Effect of time of exposure to gel-barrier tubes on results for progesterone and some other endocrine tests. *Clin Chem* 1987; 33:204. 19) Reimers TJ, et al. Effect of storage times and temperature on T3, T4, LH, prolactin, insulin, cortisol, and progesterone concentrations in blood samples from cows. *J Anim Sci* 1983; 57:683-691. 20) Smith RL. Effect of serum-separating gels on progesterone assays. *Clin Chem* 1985; 31:1239. 21) Clinical and Laboratory Standards Institute (formerly NCCLS) *Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. Approved Guideline Fourth ed.* NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 22) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377-82, 387-8. 23) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 24) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR

1910.1030.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	%CV ⁵	SD	%CV
1	1.04	0.13	12.5%	0.14	13.2%
2	1.73	0.13	7.8%	0.15	8.5%
3	4.28	0.29	6.9%	0.33	7.6%
4	7.51	0.41	5.5%	0.48	6.4%
5	18.9	1.11	5.9%	1.19	6.3%
6	31.4	1.74	5.5%	2.03	6.5%

Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	9.64	—	—
	7 in 8	8.49	8.43	101
	6 in 8	7.36	7.23	102
	5 in 8	6.33	6.02	105
	4 in 8	5.01	4.82	104
	3 in 8	3.62	3.61	100
	2 in 8	2.49	2.41	103
	1 in 8	1.16	1.20	96
2	8 in 8	7.71	—	—
	7 in 8	6.74	6.75	100
	6 in 8	6.23	5.78	108
	5 in 8	5.10	4.82	106
	4 in 8	4.32	3.86	112
	3 in 8	3.49	2.89	121
	2 in 8	2.32	1.93	121
	1 in 8	1.07	0.96	111
3	8 in 8	18.89	—	—
	7 in 8	17.54	16.79	104

	6 in 8	15.33	14.63	105
	5 in 8	12.80	12.06	106
	4 in 8	10.52	9.68	109
	3 in 8	8.12	7.34	111
	2 in 8	5.31	4.87	109
	1 in 8	2.58	2.42	107
4	8 in 8	18.16	—	—
	7 in 8	15.96	15.89	100
	6 in 8	14.34	13.62	105
	5 in 8	12.40	11.35	109
	4 in 8	10.15	9.08	112
	3 in 8	7.61	6.81	112
	2 in 8	5.07	4.54	112
	1 in 8	2.37	2.27	104
5	8 in 8	35.08	—	—
	7 in 8	32.07	30.69	105
	6 in 8	26.52	26.30	101
	5 in 8	22.13	21.92	101
	4 in 8	18.09	17.54	103
	3 in 8	13.75	13.15	105
	2 in 8	8.98	8.77	102
	1 in 8	4.35	4.38	99

Dilution ¹		Observed ²	Expected ³	%O/E ⁴
6	8 in 8	34.40	—	—
	7 in 8	32.27	30.10	107
	6 in 8	26.88	25.80	104
	5 in 8	22.22	21.50	103
	4 in 8	17.97	17.20	104
	3 in 8	13.22	12.90	102
	2 in 8	9.13	8.60	106
	1 in 8	4.53	4.30	105

Specificity

Compound ¹	ng/mL Added ²	Apparent ng/mL ³	%Cross-reactivity ⁴
Androstenedione	100	ND	ND
Androstenediol	100	ND	ND
Corticosterone	100	1.93	1.90
Cortisol	1000	ND	ND
Cortisone	500	ND	ND
Clomiphene citrate	100	ND	ND
Danazol	1000	ND	ND
DHEA sulfate*	5000	0.44	0.01
11-Deoxy-corticosterone	100	0.83	0.80
11-Deoxycortisol	500	ND	ND
Estrone	10	ND	ND
Estradiol	20	ND	ND
Estriol	100	ND	ND
Testosterone	500	0.47	0.10
Medroxy-progesterone	1000	0.28	0.00
Pregnenolone	200	1.15	0.60
17a-Hydroxy-progesterone	100	0.64	0.60
20a-Dihydro-progesterone	100	1.00	1.00
17a, 20a-Dihydroxy-progesterone	200	ND	ND
17a, 20b-Dihydroxy-progesterone	200	ND	ND
5a-Pregnan-3b, 20b-diol-11-one	200	ND	ND
5b-Pregnane-3a, 20a-diol-3-glucuronide	200	ND	ND



Compound ¹	ng/mL Added ²	Apparent ng/mL ³	%Cross-reactivity ⁴
5b-Dihydro-progesterone	20	1.31	6.60
Prednisolone	200	ND	ND

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IMMULITE is a trademark of Siemens Healthcare Diagnostics. Becton Dickinson SST is a trademark of Becton, Dickinson and Company.
© 2011 Siemens Healthcare Diagnostics. All rights reserved.
ADJUSTOR AB