# IMMULITE®

## Free T3

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**SIEMENS**

# IMMULITE®/IMMULITE 1000 Free T3

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE® and IMMULITE 1000 Analyzers — for the quantitative measurement of free T3 levels in serum, as an aid in the clinical assessment of thyroid status.

Catalog Number: **LKF31** (100 tests)

Test Code: **T3F** Color: **Aqua**

### Summary and Explanation

Under normal physiological conditions triiodothyronine (T3) represents approximately 5 percent of the thyroid hormone in plasma. Although present in lower concentration than thyroxine (T4), T3 has greater metabolic activity, faster turnover and a larger volume of distribution. It is produced largely through extrathyroidal conversion of T4. Like T4, it circulates almost entirely bound to the carrier proteins TBG, pre-albumin and albumin. Free T3 constitutes only about 0.25% of the total T3 in circulation.

Measurement of total T3 by immunoassay has a number of well-established uses. In the presence of elevated free or total T4, total T3 measurements help confirm a hyperthyroidism diagnosis. Abnormal elevation of total T3 may also occur when the total T4 concentration is normal — a condition known as "T3 toxicosis".

For the most part, free T3 levels correlate closely with total T3 levels. Total T3, however, depends not only on thyroid status and the peripheral conversion of T4 to T3, but also on the concentration of thyroid hormone-binding proteins. Free T3, on the other hand, is largely unaffected by variations in these carrier proteins. Thus, the TBG elevations typical of pregnancy, oral contraceptive use, and estrogen therapy effects an increase in the total T3 level while leaving the free T3 concentration basically unchanged.

The free T3 concentration typically reflects a patient's actual thyroid status more reliably than the total T3 concentration.

### Principle of the Procedure

IMMULITE/IMMULITE 1000 Free T3 is a competitive, analog-based immunoassay.<sup>25,26</sup>

**Incubation Cycles:** 2 × 30 minutes

The IMMULITE Free T3 procedure is a *direct* or *single test* assay,<sup>2</sup> in the sense that its results are not calculated as a function of total T3, but interpolated from a (stored) standard curve calibrated in terms of free T3 concentrations. In this respect it differs from classic equilibrium dialysis methods and from so-called free T3 index determinations as well. It requires neither a pre-incubation step nor preliminary isolation of the free fraction by dialysis or column chromatography.

The assay employs several features to preserve the equilibrium between free and protein-bound T3 and to measure accurately the unbound fraction. First, optimal concentrations of blocking agents prevent ligand-labeled T3 analog from binding to endogenous proteins (including albumin) while maintaining native T3 binding characteristics. The blocking agents also minimize artifacts arising from abnormal levels of albumin or free fatty acids. Second, the T3 analog has a nondetectable binding affinity for TBG. Third, the assay's specific antibody binds T3 about as tightly as albumin binds T3 and thereby avoids stripping the hormone from thyroid binding proteins. Finally, the procedure operates at physiological temperature, pH and ionic strength.

### Specimen Collection

Because EDTA would affect results, it should not be used as an anticoagulant.

Lipemic or grossly contaminated samples may give erroneous results. The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin,

ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Free T3 has not been tested with all possible variations of tube types.

**Volume Required:** 100  $\mu$ L serum (Sample cup must contain at least 250  $\mu$ L more than the total volume required.)

**Storage:** 2 days at 2–8°C or 2 months at –20°C.

## Warnings and Precautions

For *in vitro* diagnostic use.

### CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.<sup>28-30</sup>

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

|                   |                                                                                                                                                                                           |
|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>H412</b>       | Harmful to aquatic life with long lasting effects.                                                                                                                                        |
| <b>P273, P501</b> | Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations.<br><b>Contents:</b> sodium azide; Free T3 Adjustors |

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### Free T3 Test Units (LF31)

Each barcode-labeled unit contains one bead coated with monoclonal murine anti-T3 antibody. Stable at 2–8°C until expiration date.

**LKF31:** 100 units

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

### Free T3 Reagent Wedges (LF3A, LF3B)

With barcodes. **LF3A:** Two wedges.

7.5 mL ligand-labeled T3 analog in buffer, with preservative. **LF3B:** One wedge.

7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to anti-ligand in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

**LKF31:** 1 set

### Free T3 Adjustors (LF3L, LF3H)

Two vials (Low and High) of lyophilized free T3 in processed human serum, with preservative. Reconstitute each vial with **4.0 mL** distilled or deionized water. Mix by

gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution. For longer storage, aliquot and freeze at –20°C.

**LKF31:** 1 set

## Kit Components Supplied Separately

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash Module

**LKPM:** Probe Cleaning Kit

**LCHx-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)

**LSCC:** Sample Cup Caps (optional)

Also Required

Sample transfer pipets, distilled or deionized water, controls

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for preparation, setup, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Note that both Reagent Wedges A and B must be loaded on the carousel to run this assay.

The sample cup must contain at least 250 µL more than the total volume required.

**Recommended Adjustment Interval:**  
2 weeks

**Quality Control Samples:** Follow government regulations or accreditation requirements for quality control frequency.

Use controls or sample pools with at least two levels (low and high) of free T3.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined

by an appropriate internal laboratory quality control scheme.

## Expected Values

A reference range study for IMMULITE Free T3 was performed on 212 apparently healthy adults yielding a central 95% range of 1.5–4.1 pg/mL for euthyroid adults.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

## Limitations

The free T3 assay serves primarily as a confirmatory test for hyperthyroidism and is of very limited value in the context of hypothyroidism. There may be some overlap between the euthyroid and hyperthyroid ranges.

The interpretation of free T3 results is complicated by a variety of drugs which can affect the binding of T3 to the thyroid hormone carrier proteins.

In severe nonthyroidal illness, the assessment of thyroid status becomes especially difficult. Since some patients in this category may suffer from concomitant primary hypothyroidism or from a compensatory secondary hypothyroidism, TSH immunoassay has been recommended as a confirmatory test in this context.

In rare conditions associated with extreme variations in the albumin-binding capacity — such as familial dysalbuminemic hyperthyroxinemia (FDH) — direct free thyroid hormone assays may yield misleading results.

Circulating autoantibodies to T3, and hormone-binding inhibitors, may interfere with the assay.<sup>2,3,12</sup>

Heparin has been reported to have both *in vivo* and *in vitro* effects on free thyroid hormones.<sup>8</sup> Hence, samples should not be collected during or soon after the administration of this anticoagulant.

Since dilution shifts the equilibrium between free and protein-bound T3, the assay cannot be expected to maintain linearity under dilution. Accordingly, no attempt should be made to dilute samples with high free T3 results. In particular, any sample with a result greater than the

upper limit of the assay's calibration range should simply be reported as such.

Because EDTA would have a significant effect on results, it should not be used as an anticoagulant.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data representative of the assay's performance. Results are expressed in pg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

### Conversion Factor:

pg/mL  $\times$  1.536  $\rightarrow$  pmol/L

**Calibration Range:** 1–40 pg/mL  
(1.5–61 pmol/L)

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

**Analytical Sensitivity:** 1.0 pg/mL  
(1.5 pmol/L)

**Precision:** Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Specificity:** The antibody is highly specific for free T3. (See "Specificity" table.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Biotin:** Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

**Hemolysis:** Presence of hemoglobin in concentrations up to 381 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Highly lipemic samples may interfere with the assay. (See "Effect of Lipemia" table.)

**Effect of Albumin:** Samples spiked with charcoal-absorbed human serum albumin (10, 20 and 50 mg/mL) were assayed. The assay reports increased free T3 values in the presence of albumin at 50 mg/mL. (See "Effect of Albumin" table.)

**Effect of Nonesterified Fatty Acids (NEFA):** Samples spiked with oleic acid (2.5 and 5.0 mmol/L) were assayed. Results show no significant effect at 2.5 mmol/L, but decreased values at 5.0 mmol/L. Increased levels of NEFA may be encountered in late pregnancy and in heparin therapy. (See "Effect of Nonesterified Fatty Acids" table.)

**Other Compounds:** Results for samples assayed after spiking with phenylbutazone (10, 100 and 1000  $\mu$ g/dL), or phenytoin (5, 10, 25, and 50  $\mu$ g/mL), or salicylate (10, 25, 50 and 100 mg/dL) showed no clinically significant effect from these compounds. (See "Effect of Phenylbutazone", "Effect of Phenytoin" and "Effect of Salicylate" tables.)

**Effect of TBG:** The assay's zero calibrator was spiked with T3-stripped TBG resin (20, 40 and 80  $\mu$ g/mL) and assayed. (See "Effect of TBG" table.) The results indicate that TBG does not interfere with the assay.

**Method Comparison:** The assay was compared to a commercially available chemiluminescent immunoassay for free T3 (Kit A) on 237 samples. (Concentration range: approximately 1.3 to 19 pg/mL. See graph.) By linear regression:

(IML) = 0.99 (Kit A) – 0.11 pg/mL  
 $r = 0.930$

Means:

3.13 pg/mL (IMMULITE)

3.28 pg/mL (Kit A)

## References

- 1) Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982.
- 2) Beck-Peccoz P, Romelli PB, et al. Free T4 and free T3 measurement in patients with anti-iodothyronine autoantibodies. In: Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982: 231-8.
- 3) Bhagat CI, et al. Interference in radioimmunoassay of total serum thyroxin and free thyroxin due to thyroxin-binding autoantibodies. Clin Chem 1983;29:1324-5.
- 4) Borst GC, Eil C, Burman KD. Euthyroid hyperthyroxinemia; a review. Ann Intern Med 1983;98:366-78.
- 5) Burger AG, Lambert MJ, Cullen M. Pharmacological influence on free thyroid hormones. In: Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982: 303-14.
- 6) Caroyan P, et al. Assessment and clinical significance of free thyroid hormone radioimmunoassays. In: Ekins RP, et al, editors. Free thyroid hormones. Amsterdam: Excerpta Medica, 1979: 181-93.
- 7) Carayon P, et al. Free thyroid hormones in euthyroidism: evaluation of normal ranges in various physiological conditions. In: Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982: 179-85.
- 8) Diagnostic Products Corporation. Coat-A-Count Free T4, package insert (PITKF4).
- 9) Ekins RP. Methods for the measurement of free thyroid hormones. In: Ekins RP, editor. Free thyroid hormones. Amsterdam: Excerpta Medica, 1979:72-106.
- 10) Ekins RP. The radioimmunoassay of free thyroid hormones in blood. In: Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982: 73-90.
- 11) Ellis SM, Ekins RP. The radioimmunoassay of serum free triiodothyronine and thyroxine. In: Pasternak CA, editor. Radioimmunoassay in clinical biochemistry. London: Heyden, 1975: 187-94.
- 12) John R. Autoantibodies to thyroxine and interference with free thyroxine assay. Clin Chem 1983;29:581-2.
- 13) Larsen PR. The use of serum T3 measurements by radioimmunoassay in the diagnosis of thyroid disease. In: Abraham GE, editor. Radioassay systems in clinical endocrinology. New York: Marcel Dekker, 1981: 117-29.
- 14) Parslow ME, Oddie TH, Fisher DA. Evaluation of serum triiodothyronine and adjusted triiodothyronine (free triiodothyronine index) in pregnancy. Clin Chem 1977;23:490-2.
- 15) Lin T, Nankin HR. Direct measurement of free thyroxine in patients with T3-toxicosis. Horm Metab Res 1980;12:611-4.
- 16) Mardell R, Gerson M. A method of assessing serum triiodothyronine concentrations that is independent of subject's age and variations in concentrations of binding proteins in serum. Clin Chem 1978;24:1792-6.
- 17) Pennisi F, Romelli PB, Vancheri L. Measurement of free thyroid hormones in serum by column adsorption chromatography and radioimmunoassay. In: Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982: 107-12.
- 18) Roux F, et al. Clinical merits in the radioimmunoassay of free serum tetraiodothyronine (FT4) and free serum triiodothyronine (FT3). Int J Nuclear Med Biol 1980;7:386-90.
- 19) Salmon D, et al. Chemical hyperthyroidism: serum triiodothyronine levels in clinically euthyroid individuals treated with levothyroxine. Arch Intern Med 1982;142:571-3.
- 20) Sawin CT, Chopra D, et al. The free triiodothyronine (T3) index. Ann Intern Med 1978;88:474-7.
- 21) Siegel L, McDonald LJ, Robin NI. Estimation of free triiodothyronine in serum: a new method and its clinical relevance. Clin Chem 1978;24:1891-4.
- 22) Smals AGH, Ross AH, Kloppenborg PWC. Dichotomy between serum free triiodothyronine and free thyroxine concentrations in familial thyroxine-binding globulin deficiency. J Clin Endocrinol Metab 1981;53:917-22.
- 23) Walfish PG. Thyroid physiology and pathology. In: Collu R, Ducharme JR, Guyda H, editors. Pediatric Endocrinology. New York: Raven Press, 1981: 357-431.
- 24) Weeke J, et al. A longitudinal study of serum TSH, and total and free iodothyronines during normal pregnancy. Acta Endocrinol 1982;101:531-7.
- 25) Witherspoon LR, El Shami AS, et al. Chemically blocked analog assays for free thyronines. Clin Chem 1988;34:9-16 and 17-23.
- 26) Diagnostic Products Corporation. Albumin-Related Artifacts in Free Thyroid Hormone Assays. Los Angeles, 1987 (Document number: ZJ026).
- 27) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4. Wayne, PA: NCCLS, 1998.
- 28) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8.
- 29) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
- 30) NCCLS Document M29-A3. 30) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

## Technical Assistance

In the United States, contact  
Siemens Healthcare Diagnostics  
Technical Services department.  
Tel: 877.229.3711. Outside the United  
States, contact your National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of Siemens Healthcare  
Diagnostics Products Ltd. is certified to  
ISO 13485.

## Tables and Graphs

### Precision (pg/mL)

|   | Mean <sup>3</sup> | Within-Run <sup>1</sup> |                 |  | Total <sup>2</sup> |      |
|---|-------------------|-------------------------|-----------------|--|--------------------|------|
|   |                   | SD <sup>4</sup>         | CV <sup>5</sup> |  | SD                 | CV   |
| 1 | 1.7               | 0.24                    | 14.1            |  | 0.29               | 17.1 |
| 2 | 2.9               | 0.32                    | 11.0            |  | 0.41               | 14.1 |
| 3 | 3.7               | 0.31                    | 8.4             |  | 0.39               | 10.5 |
| 4 | 4.6               | 0.34                    | 7.4             |  | 0.42               | 9.1  |
| 5 | 5.6               | 0.34                    | 6.1             |  | 0.42               | 7.5  |
| 6 | 10.2              | 0.59                    | 5.8             |  | 0.78               | 7.6  |

### Specificity

| Compound <sup>1</sup>                | pg/mL<br>Added <sup>2</sup> | Apparent<br>Conc.<br>(pg/mL) <sup>3</sup> | %<br>Cross-<br>reactivity <sup>4</sup> |
|--------------------------------------|-----------------------------|-------------------------------------------|----------------------------------------|
| L-Thyroxine (T4)                     | 50,000                      | 0.71                                      | 0.0014%                                |
| D-Thyroxine                          | 50,000                      | 0.70                                      | 0.0014%                                |
| Triiodothyroacetic<br>acid           | 5000                        | 2.62                                      | 0.052%                                 |
| Tetraiodothyroacetic<br>acid         | 200,000                     | 0.89                                      | 0.0004%                                |
| Diiodo-L-tyrosine                    | 10,000,000                  | 2.31                                      | 0.0002%                                |
| Methimazole                          | 10,000,000                  | 2.4                                       | 0.0002%                                |
| 5,5'-<br>Diphenylhydantoin           | 10,000,000                  | ND                                        | ND                                     |
| Phenylbutazone                       | 10,000,000                  | 1.04                                      | 0.00001%                               |
| 6- <i>n</i> -Propyl-2-<br>thiouracil | 10,000,000                  | 0.47                                      | 0.000005%                              |
| Triiodo-L-tyrosine                   | 10,000,000                  | 1.59                                      | 0.00002%                               |

ND: not detectable<sup>5</sup>

### Effect of Albumin (pg/mL)

|    | Unspiked <sup>2</sup> | Albumin (mg/mL) <sup>1</sup> |      |      |
|----|-----------------------|------------------------------|------|------|
|    |                       | 10                           | 20   | 50   |
| 1  | 2.7                   | 3.0                          | 3.3  | 4.0  |
| 2  | 3.4                   | 3.8                          | 3.8  | 4.6  |
| 3  | 3.7                   | 3.8                          | 3.7  | 4.3  |
| 4  | 3.8                   | 3.6                          | 4.1  | 4.7  |
| 5  | 3.9                   | 3.6                          | 3.9  | 4.9  |
| 6  | 8.0                   | 10.3                         | 9.6  | 10.1 |
| 7  | 8.7                   | 7.6                          | 8.8  | 10.2 |
| 8  | 10.5                  | 10.1                         | 11.6 | 12.5 |
| 9  | 11.3                  | 13.1                         | 12.6 | 15.0 |
| 10 | 13.6                  | 13.9                         | 14.1 | 16.8 |
| 11 | 14.9                  | 13.6                         | 15.9 | 18.2 |

### Effect of Nonesterified Fatty Acids (pg/mL)


|   | Unspiked <sup>2</sup> | Oleic acid (mmol/L) <sup>1</sup> |     |
|---|-----------------------|----------------------------------|-----|
|   |                       | 2.5                              | 5.0 |
| 1 | 2.2                   | 2.2                              | 1.4 |
| 2 | 3.0                   | 3.1                              | 2.5 |
| 3 | 3.4                   | 3.1                              | 1.9 |
| 4 | 3.5                   | 3.5                              | 2.5 |
| 5 | 4.1                   | 4.1                              | 2.9 |
| 6 | 6.0                   | 6.7                              | 4.7 |
| 7 | 7.8                   | 8.0                              | 6.0 |
| 8 | 8.3                   | 8.2                              | 6.5 |
| 9 | 10.0                  | 11.1                             | 7.0 |

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IMMULITE and Coat-A-Count are trademarks of Siemens Healthcare Diagnostics.
All rights reserved.