Anti-TPO Ab

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Anti-TPO Ab

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of antithyroid peroxidase (TPO) antibodies in serum and EDTA plasma, as an aid in the clinical diagnosis of thyroid diseases.

Catalog Number: **L2KTO2** (200 tests)
L2KTO6 (600 tests)

Test Code: **ATA** Color: **Aqua**

Summary and Explanation

Antithyroid peroxidase (TPO) antibodies are autoantibodies directed against the thyroid peroxidase enzyme. This enzyme catalyzes the iodination of tyrosine in thyroglobulin during the biosynthesis of T3 and T4.²

Historically, these antibodies were referred to as antimicrosomal antibodies (AMA) because the antibodies bind to the microsomal part of the thyroid cells. Recent research has identified thyroid peroxidase as the primary antigenic component of microsomes.^{2,3,13}

Thyroid autoimmune disease is the major factor underlying hypothyroidism and hyperthyroidism and tends to occur in a genetically predisposed population.

The major thyroid autoimmune diseases are: Hashimoto's thyroiditis and Graves' disease.¹³ In virtually all cases of Hashimoto's disease and in the majority of Graves' disease cases, TPO autoantibodies are elevated.¹³ High levels of TPO autoantibodies, in the context of the clinical presentation of hypothyroidism, confirms the diagnosis of Hashimoto's disease.

Principle of the Procedure

IMMULITE 2000 Anti-TPO Ab is a solid-phase, enzyme-labeled, chemiluminescent sequential immunometric assay.

Incubation Cycles: 2 × 30 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Anti-TPO Ab has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 5 µL sample

Storage: 2 days at 2–8°C, or for longer periods frozen at –20°C.

Automatic Predilution Factor: 100

Dispense required volume of diluent (L2AAZ) to a suitable test tube with barcode label applied.

Warnings and Precautions

For *in vitro* diagnostic use.

**CAUTION! POTENTIAL BIOHAZARD**

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.¹⁵⁻¹⁷

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H311, H302, H411

P280, P273, P301 + P312, P302 + P312, P501

Danger! Toxic in contact with skin. Harmful if swallowed. Toxic to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection.

Avoid release to the environment. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 2-methyl-2H-isothiazol-3-one, sodium azide; Anti-TPO Ab Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Anti-TPO Ab Bead Pack (L2TO12)

With barcode. 200 beads, coated with highly purified human TPO. Stable at 2–8°C until expiration date.

L2KTO2: 1 pack **L2KTO6:** 3 packs

Anti-TPO Ab Reagent Wedge (L2TOA2)

With barcode. 11.5 mL buffer matrix, with preservative. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-human IgG in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KTO2: 1 wedge **L2KTO6:** 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Anti-TPO Ab Adjustors (LTOL, LTOH)

Two vials (Low and High), of lyophilized TPO autoantibodies in a human serum/buffer matrix, with preservative. At least 30 minutes before use, reconstitute each vial by adding **4.0 mL** distilled or deionized water. Mix by *gentle* swirling or inversion until the lyophilized material is fully dissolved. (No further dilution is required.) Aliquot and freeze in tightly capped plastic containers. Avoid repeated freeze/thaw cycles. Stable at 2–8°C for 7 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KTO2: 1 set **L2KTO6:** 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Thyroid Autoantibody Sample Diluent (L2AAZ)

Test code: **TAD**. For the on-board dilution of patient samples. 50 mL concentrated (ready-to-use) anti-TG/anti-TPO antibody-free buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or 6 months (aliquotted) at –20°C.

L2KTO2: 1 vial **L2KTO6**: 3 vials

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2KTO2: 3 labels **L2KTO6**: 5 labels

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LAACM: Bi-level control module

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks

Quality Control Samples: Use controls or serum pools with at least two levels (low and high) of TG and TPO autoantibodies.


Expected Values

Based on its relationship to IMMULITE Anti-TPO Ab (see Method Comparison), the assay can be expected to have essentially the same reference ranges.


Thyroid peroxidase autoantibodies are unique for each patient; therefore, each patient sample will show its own individual binding characteristics. Normal levels for most patients have been defined as less than 35 IU/mL for the IMMULITE Anti-TPO Ab assay. A numerical result greater than or equal to 35 IU/mL indicates an elevated anti-TPO level in serum, whereas a result less than 35 IU/mL indicates a *normal* anti-TPO level in serum. However, due to individual variability it is also possible that a normal patient with no history of thyroid disease may have reproducible elevated values greater than 35 IU/mL.

Two studies are presented below that demonstrate the variability between two random populations.

Using the IMMULITE Anti-TPO Ab assay, the levels of TPO autoantibodies were determined for 75 apparently healthy subjects (26 males and 49 females), with ages ranging from 13 to 87 years. All had normal TSH levels, and none had thyroid autoimmune disorders. Using the recommended 35 IU/mL cutoff, 79% (59/75) of the subjects had normal anti-TPO values. The following histogram shows the distribution of subjects for different levels of TPO autoantibodies.

In another study, IMMULITE Anti-TPO Ab procedure was used to determine the TPO autoantibody levels of 130 apparently healthy subjects. Using the recommended 35 IU/mL cutoff, 92% (118/130) of the subjects had normal anti-TPO values. The distribution of subjects for different levels of TPO autoantibodies is presented in the following histogram.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

For diagnostic purposes, anti-TPO results should be used in conjunction with other test results, the overall clinical presentation to the physician, and all other appropriate information. Autoantibodies may be found in less than 10% of the normal population at low levels and in patients with non-thyroidal illnesses, such as the inflammatory rheumatic diseases.¹⁴

Patient TPO autoantibodies and the antibodies used for calibrating the assay may differ in their binding properties. This may lead to deviations from strict dilutional parallelism for certain patient samples.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in IU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: up to 1,000 IU/mL (WHO 1st IRP 66/387)

Analytical Sensitivity: 5.0 IU/mL

Hook-Effect: None up to 16,624 IU/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three anti-TPO Ab solutions (63, 207 and 520 IU/mL) were assayed. (See "Recovery" table for representative data.)

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 20 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of TPO autoantibodies, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 Anti-TPO Ab procedure.

(Heparin) = 0.84 (Serum) + 17 IU/mL
r = 0.984

(EDTA) = 0.89 (Serum) + 15 IU/mL
r = 0.981

(SST) = 0.88 (Plain Tubes) + 17 IU/mL
r = 0.984

Means:

299 IU/mL (Serum)
268 IU/mL (Heparin)
282 IU/mL (EDTA)
279 IU/mL (SST)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 381 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison: The assay was compared to IMMULITE Anti-TPO Ab on 246 samples. (Concentration range: approximately 10 to 1,000 IU/mL. See graph.) By linear regression:

$$(IML\ 2000) = 0.97 (IML) - 2.6\ IU/mL$$

$$r = 0.987$$

Means:

188 IU/mL (IMMULITE 2000)

197 IU/mL (IMMULITE)

References

- 1) Becker W, Reimers C, Börner W. Immunological criteria for the differential diagnosis of thyroid disease. *Therapiewoche* 1985;35:1167-76.
- 2) Czarnocka B, Ruf J, Ferrand M, et al. Purification of the human thyroid peroxidase and its identification as the microsomal antigen involved in autoimmune thyroid diseases. *FEBS Lett.* 1985;190:147-52.
- 3) Czarnocka B, Ruf J, Ferrand M, et al. Interaction of highly purified thyroid peroxidase with anti-microsomal antibodies in autoimmune thyroid diseases. *J Endocrinol Invest* 1986;9:135-8.
- 4) Ruf J, Czarnocka B, Ferrand M, et al. Thyroid peroxidase is the organ-specific 'microsomal' antigen involved in thyroid autoimmunity. *Acta Endocrinol (Copenh)* 1987;Suppl. 281: 49-55.
- 5) Ruf J, Czarnocka B, Ferrand M, et al. Relationship between immunological structure and biochemical properties of human thyroid peroxidase. *Endocrinol* 1988;125 (3):1211-8.
- 6) Horster FA. The importance of MAK, TAK, TRAK and thyroglobulin in the diagnosis of thyroid disease. *Internist* 1988;29:538-40.
- 7) Schatz H, Löbig H. The diagnostic and prognostic importance of antibodies against thyroid microsomes and thyroglobulin. *Akt Endokr Stoffw* 1989;10: 146-53.
- 8) Hashimoto H. Understanding the lymphatic changes of the thyroid. *Arch Klin Chir* 1912;97:219-48.
- 9) Libert F, Ruel J, Ludgate M, et al. Thyroperoxidase, an auto-antigen with a mosaic structure made of nuclear and mitochondrial gene modules. *EMBO J* 1987;6:4193-6.
- 10) McKenzie JM, Zakarija M, Sato A. Humoral immunity in Graves' disease. *Clin Endocrinol Metab* 1978;7:31-45.
- 11) Scherbaum WA, Berg PA. The importance of autoantibodies in the diagnosis of endocrine disease. *Dtsch Med Wschr* 1981;106:308-13.
- 12) Scherbaum WA, Berg PA. Thyroid autoimmune disease – new aspects in development and diagnosis. *Dtsch Med Wschr* 1984;109:1574-81.
- 13) Feldt-Rasmussen U. Analytical and clinical performance goals for testing autoantibodies to thyroperoxidase, thyroglobulin, and thyrotropin receptor. *Clin*

Chem 1996;42:160-3.

- 14) Wild D, editor. The immunoassay handbook. Great Britain: Stockton Press, 1994: 342-3.
- 15) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377-82, 387-8.
- 16) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.
- 17) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics


The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

Tables and Graphs

Precision (IU/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	55.54	2.90	5.21%	3.41	6.13%
2	91.15	4.52	4.96%	6.00	6.58%
3	122.96	7.72	6.27%	8.83	7.18%
4	544.14	30.31	5.57%	35.20	6.47%
5	636.29	33.68	5.29%	51.85	8.15%

Method Comparison

(IML 2000) = 0.97 (IML) – 2.6 IU/mL
r = 0.987

Linearity (IU/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	219	—	—
	4 in 8	109	110	99%
	2 in 8	57	55	104%
	1 in 8	27	27	100%
2	8 in 8	394	—	—
	4 in 8	218	197	111%
	2 in 8	103	99	104%
	1 in 8	52	49	106%
3	8 in 8	474	—	—
	4 in 8	253	237	107%
	2 in 8	131	119	110%
	1 in 8	59	59	100%
4	8 in 8	547	—	—
	4 in 8	255	274	93%
	2 in 8	129	137	94%
	1 in 8	62	68	91%
5	8 in 8	651	—	—
	4 in 8	336	326	103%
	2 in 8	166	163	102%
	1 in 8	88	81	109%
6	8 in 8	932	—	—
	4 in 8	478	466	103%
	2 in 8	223	233	96%
	1 in 8	110	117	94%

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.

© 2009 Siemens Healthcare Diagnostics.
All rights reserved.

Made in: UK

Siemens Healthcare
Diagnostics Products Ltd.
Glyn Rhonwy, Llanberis,
Gwynedd LL55 4EL
United Kingdom

2015-04-20

PIL2KTO – 23

Changes in this Edition:

cc#EU22513: 1) Under Warnings and Precautions, added biohazard symbol and potential biohazard statement, caution statement, required GHS pictograms/signal words/hazard and precautionary codes and statements, offending chemicals and sources. 2) Under References, added citations from the potential biohazard statement. 3) Changed "Origin: UK" to "Made in: UK". 4) Updated the manufacturer address to include Glyn Rhonwy. 5) In Understanding the Symbols, revised the symbol for catalog number, added a GHS pictogram for health hazard, and updated existing symbols with GHS pictograms and terminology for exclamation mark, corrosion, skull and crossbones, and environment.

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