Growth Hormone (hGH)

(Recombinant 98/574)

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000

Growth Hormone (hGH)

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of human growth hormone in serum, as an aid in patient management.

Catalog Number: **L2KGRH2** (200 tests)

Test Code: **GRH** Color: **Red**

Summary and Explanation

Human growth hormone (hGH, somatotropin) is a polypeptide originating in the anterior pituitary. It is 191 amino acids in length and has a molecular mass of approximately 22,000 daltons. Its metabolic effects are primarily anabolic. It promotes protein conservation and engages a wide range of mechanisms for protein synthesis. It also enhances glucose transport and facilitates the buildup of glycogen stores. Another family of peptide hormones, the somatomedins, mediates its cascade of growth-promoting actions.¹

Measurement of hGH is primarily of interest in the diagnosis and treatment of various forms of inappropriate growth hormone secretion. Clinical disorders of hyposecretion include dwarfism and unattained growth potential. Hypersecretion is associated with gigantism and acromegaly.

Caution must be exercised in the clinical interpretation of growth hormone levels. These vary throughout the day, making it difficult to define a reference range or to judge an individual's status based on single determinations. Many factors are known to influence the rate of growth hormone secretion, including periods of sleep and wakefulness, exercise, stress, hypoglycemia, estrogens, corticosteroids, L-dopa, and others.¹ Because of its similarity to prolactin and placental lactogen, earlier growth hormone immunoassays were often plagued with

falsely high values in pregnant and lactating women.

Because not all acromegalic individuals have elevated baseline levels, suppression tests based on glucose loading are of value in this context. In spite of the induced hyperglycemia, there is rarely a decrease from baseline levels in acromegaly.²

Growth hormone-deficient individuals have fasting/resting levels similar to those found in healthy individuals. Various challenge tests have therefore been devised to differentiate these groups. Thus, with the onset of deep sleep or after 15 to 20 minutes of vigorous exercise, growth hormone levels normally show a rise. Other tests of growth hormone responsiveness are based on the administration of L-dopa, arginine and insulin. Propanolol or estrogen are sometimes given in conjunction with the primary stimulus to accentuate the response.^{4,5}

A small number of cases of dwarfism have been documented in which both the basal level and the response to challenge testing were normal. Such cases may involve tissue insensitivity to either growth hormone or the somatomedins, or the presence of antibodies or immunoreactive but biologically inactive growth hormone.⁴

Principle of the Procedure

IMMULITE 2000 Growth Hormone is a solid-phase, two-site chemiluminescent immunometric assay.

The solid phase (bead) is coated with murine monoclonal anti-hGH antibody. The reagent contains alkaline phosphatase (bovine calf intestine) conjugated to a rabbit anti-hGH polyclonal antibody. The reagent and hGH in the sample are incubated together with a bead coated with a murine anti-hGH monoclonal antibody to form an antibody-sandwich complex. Unbound enzyme conjugate is then removed by a centrifugal wash. Finally, chemiluminescent substrate is added to

the reaction tube and signal is generated in proportion to the bound enzyme.

Incubation Cycles: 1 × 30 minutes

Time to First Result: 35 minutes

Specimen Collection

Patient must be fasting and at complete rest 30 minutes before blood collection.^{6,7}

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Growth Hormone has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 25 µL serum

Storage: Refrigerate immediately. Stable at 2–8°C for 8 hours,^{6,7} or for 2 months (aliquotted) at –20°C. Avoid repeated thawing and freezing.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412	Harmful to aquatic life with long lasting effects.
P273, P501	Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations.
	Contents: sodium azide; GH Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

GH Bead Pack (L2GRH12)

With barcode. 200 beads, coated with murine monoclonal anti-hGH antibody. Stable at 2–8°C until expiration date.

L2KGRH2: 1 pack

GH Reagent Wedge (L2GRHA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to rabbit polyclonal anti-hGH antibody in buffer. Stable at 2–8°C until expiration date.

L2KGRH2: 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

GH Adjustors (LGHL, LGHH)

Two vials (Low and High) containing lyophilized hGH in a nonhuman serum, with preservative. Reconstitute each vial with **3.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KGRH2: 1 set

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Multi-Diluent 2 (L2M2Z, L2M2Z4)

For the on-board dilution of high samples. One vial of concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M2Z: 3 labels **L2M2Z4:** 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

Also Required

Distilled or deionized water; test tubes; controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:

2 weeks

Quality Control Samples: Follow government regulations or accreditation requirements for quality control frequency.

Use controls or serum pools with at least two levels (low and high) of human growth hormone.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

Expected Values

A study performed on 67 healthy adult males and 149 healthy adult females yielded the following reference intervals:

Males: Up to 3 ng/mL

Females: Up to 8 ng/mL

Secretion of growth hormone by the pituitary gland is episodic and pulsatile, and transient levels up to 40 ng/mL have been observed in healthy subjects. Because serum levels are rather low between pulses in healthy individuals, immunoassays with sufficient sensitivity are necessary to distinguish patients with abnormally low values from healthy subjects who have values that happen to fall in the low-normal range. Thus, basal or random levels of hGH provide little diagnostic information. In some individuals, spontaneous growth hormone secretion is better monitored by using a continuous withdrawal pump, or by drawing specimens for hGH assay every 20 to 30 minutes over a 12-to-24-hour period. A number of provocative tests designed to stimulate or suppress release are also used to assess growth hormone secretion; the definition of abnormal responses varies with the laboratory.¹

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Patients on hGH therapy may develop human antibodies to hGH that will interfere in the assay and cause falsely low values.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

ng/mL \times 3.0 \rightarrow mIU/L
WHO NIBSC 2nd IS 98/574

Reportable Range: 0.05 to 40 ng/mL
(0.15–120 mIU/L) WHO NIBSC 2nd IS 98/574

Analytical Sensitivity: 0.01 ng/mL

High Dose Hook Effect: None up to 17,000 ng/mL

Precision: Five samples were processed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three hGH solutions (120, 240 and 480 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for human growth hormone. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 20 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of hGH, to obtain values throughout the reportable range of the assay, and then assayed by the IMMULITE 2000 Growth Hormone procedure.

(Heparin) = 0.99 (Serum) + 0.07 ng/mL
 $r = 0.997$

(EDTA) = 1.0 (Serum) – 0.01 ng/mL
 $r = 0.997$

(SST) = 0.98 (Plain Tubes) + 0.02 ng/mL
 $r = 0.999$

Means:

4.15 ng/mL (Heparin)

4.12 ng/mL (EDTA)

4.14 ng/mL (Serum)

4.09 ng/mL (SST)

Method Comparison: The assay was compared to IMMULITE Growth Hormone (LKGRH, WHO NIBSC 2nd IS 98/574) on 72 samples. (Concentration range: approximately 0.1 to 28 ng/mL. See graph.) By linear regression:

(IML 2000) = 1.03 (IML) + 0.06 ng/mL
 $r = 0.997$

Means:

6.75 ng/mL (IMMULITE 2000)

6.48 ng/mL (IMMULITE)

References

1) Whitley RJ, Meikle AW, Watts NB. Endocrinology. Part 2: Protein hormones. In: Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: Saunders, 1994: 1665-70. 2) Daughaday WM, Cryer PE. Growth hormone hypersecretion and acromegaly. Hospital Practice. no. 8 August 1978;13:75-80. 3) Raite S. The standards for human growth hormone assays. In: Laron Z, Butenandt O, editors. Evaluation of Growth Hormone Secretion. Basel: Karger, 1983: 162-9. 4) Van Wyk JJ. Diagnosis of growth hormone deficiency and other causes of growth failure. The Endocrine Society: Syllabus for the 35th Annual Postgraduate Assembly October 1983: 90-106. 5) Van Wyk JJ, Underwood L. Growth hormone, somatomedins and growth failure. Hospital Practice no. 8 August 1978;13:57-67. 6) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: W.B. Saunders, 1995: 300. 7) Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: W.B. Saunders, 1994: 1670.

Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	2.6	0.09	3.5%	0.17	6.5%
2	3.7	0.17	4.6%	0.21	5.7%
3	5.3	0.18	3.4%	0.29	5.5%
4	7.9	0.23	2.9%	0.33	4.2%
5	17	0.72	4.2%	1.13	6.6%

Specificity

Compound ¹	ng/mL Added ²	% Cross reactivity ³
FSH	10,000	0.001%
Free Beta HCG	280	ND
HPL	100	0.03%
LH	3663	0.5%
Prolactin	40,000	ND
TSH	11,125	0.01%

ND: nondetectable⁴


Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	0.95	—	—
	4 in 8	0.48	0.47	102%
	2 in 8	0.25	0.24	104%
	1 in 8	0.12	0.12	100%
2	8 in 8	2.09	—	—
	4 in 8	1.06	1.04	102%
	2 in 8	0.55	0.52	106%
	1 in 8	0.27	0.26	104%
3	8 in 8	3.53	—	—
	4 in 8	1.83	1.76	104%
	2 in 8	0.86	0.88	98%
	1 in 8	0.43	0.44	98%
4	8 in 8	6.07	—	—
	4 in 8	2.81	3.04	92%
	2 in 8	1.39	1.52	91%
	1 in 8	0.72	0.76	95%
5	8 in 8	30.1	—	—
	4 in 8	14.6	15.1	97%
	2 in 8	7.17	7.53	95%
	1 in 8	3.78	3.76	101%
6	8 in 8	36.9	—	—
	4 in 8	16.3	18.4	89%
	2 in 8	8.02	9.22	87%
	1 in 8	4.15	4.61	90%

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IMMULITE is a trademark of
All rights reserved.