IMMULITE®

TBG

For use on the IMMULITE®
and IMMULITE® 1000 systems

SIEMENS

IMMULITE®/IMMULITE 1000 TBG

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® and IMMULITE 1000 Analyzers — for the quantitative measurement of thyroxine-binding globulin (TBG) in serum, as an aid in the diagnosis and treatment of thyroid diseases.

Catalog Number: **LKTB1** (100 tests),

Test Code: **TBG** Color: **Dark Blue**

Summary and Explanation

Thyroxine-binding globulin (TBG) is an acidic glycoprotein consisting of a single polypeptide chain with a molecular weight of 54,000 daltons.^{11,12,13} It is one of three principal carrier proteins for both thyroxine (T4) and 3,5,3'-triiodothyronine (T3); the other thyroid hormone carrier proteins are thyroxine-binding prealbumin (TBPA) and albumin. TBG, though present in significantly lesser amounts than TBPA and albumin, has a higher affinity for the thyroid hormones and is therefore the chief carrier protein.^{1,2} In healthy subjects, under 0.05 percent of the total serum T4 is present as unbound (free) hormone. The bound fraction is distributed among the carrier proteins as follows: TBG, 70–75%; TBPA, 15–20%; and albumin, 5–10%.²

Both bound and unbound T4 are present in an equilibrium that tends to reassert itself in the face of altered levels of the carrier proteins by inducing a corresponding alteration in the *total* level of T4 in circulation, while leaving the free T4 level relatively unchanged.² Hence the free T4 concentration may be expected to correlate more closely than the total T4 concentration with clinical thyroid status,³ for an abnormal total T4 result may signify *either* an abnormality in thyroid function *or* simply a variation (physiological or pathological) in the carrier proteins.¹

Thus, for example, the TBG elevations typical of pregnancy,⁴ oral contraceptives and estrogen therapy will cause the total T4 level to increase, often beyond the limits of normal, without inducing a corresponding elevation in the free T4 level.^{5,6} Again, alterations in the TBG level

sometimes mask the effects of abnormal thyroid function by elevating the total T4 level of a hypothyroid patient — or lowering that of a hyperthyroid patient⁴ — into the euthyroid reference range. Here, too, the free T4 concentration will typically reflect the patient's actual thyroid status more reliably than the total T4 concentration.

Historically, estimation of free T4, often via the free thyroxine index (FT4I), has been the most frequently used test for thyroid dysfunction.⁷ Calculation of the FT4I typically involves multiplying a total T4 result by a T3-uptake (or T4-uptake) result. The T3 uptake test yields a relative measure of the unsaturated binding sites on TBG, rather than a direct, quantitative estimate of the concentration of this carrier protein. The TBG saturation index (TBG-SI) has been used, under various names, e.g. thyroid hormone binding ratio (THBR), as a substitute for the FT4I.¹⁴ The TBG-SI is essentially just the ratio of total T4 to TBG, multiplied by some factor (depending on the units) so as to achieve a standard reporting scale. When total T4 and TBG are both expressed in molar units as nmol/L and the factor is taken as 100, the TBG-SI represents the fractional occupation of TBG binding sites by T4, expressed as a percent.¹⁴

Principle of the Procedure

IMMULITE/IMMULITE 1000 TBG is a solid-phase, competitive chemiluminescent enzyme immunoassay.

Incubation Cycles: 1 × 30 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of

samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants.

IMMULITE/IMMULITE 1000 TBG has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 20 μ L serum. (Sample cup must contain at least 100 μ L more than the total volume required.)

Storage: 2 days at 2–8°C or 1 month at –20°C.⁸

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.¹⁶⁻¹⁸

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

TBG Test Units (LTB1)

Each barcode-labeled unit contains one bead coated with monoclonal murine anti-TBG. Stable at 2–8°C until expiration date.

LKTB1: 100 units

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

TBG Reagent Wedge (LTB2)

With barcode. 7.5 mL containing alkaline phosphatase (bovine calf intestine) conjugated to human TBG in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKTB1: 1 wedge

TBG Adjustors (LTBL, LTBH)

Two vials (Low and High), 4 mL each, of human TBG in buffer, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKTB1: 1 set

Kit Components Supplied Separately

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

LTBCM: A bi-level control module, containing human TBG

Also Required

Sample transfer pipets, distilled or deionized water, controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:
1 week

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of TBG.

Calculation of TBG Saturation Index

The TBG Saturation Index (TBG-SI) represents the fractional occupation of TBG binding sites by T4, expressed as a percent.¹⁴ This is a natural measure since TBG has a single binding site for thyroid hormones, with a higher affinity for T4 than for T3; and T4 circulates at a much higher concentration than T3.

The TBG-SI is defined as 100 times the *molar* ratio of total T4 to TBG:

$$\text{TBG-SI} = \frac{\text{Total T4 (nmol/L)}}{\text{TBG (nmol/L)}} \times 100$$

The definition makes the TBG-SI independent of the units used to *report* total T4 and TBG results, since it requires putting both results into the same concentration units (e.g. nmol/L) before computing their ratio. For example, given total T4 and TBG concentrations expressed in mass units of µg/dL and µg/mL, respectively, multiply the T4 result by 12.87 to convert from µg/dL to nmol/L, multiply the TBG result by 18.5 to convert µg/mL to nmol/L, then multiply by 100 as in the formula above. Alternatively, multiply the T4/TBG ratio in the mass units indicated by $(12.87 / 18.5) \times 100 = 69.5$:

$$\text{TBG-SI} = \frac{\text{Total T4 (}\mu\text{g/dL)}}{\text{TBG (}\mu\text{g/mL)}} \times 69.5$$

Expected Values

A reference range study was performed with IMMULITE TBG kit on 75 serum samples. These were obtained from adult men and nonpregnant women in good health; all samples had normal total T4 and free T4 index (FT4I) levels. The TBG results showed an absolute range of 13 to 39 µg/mL, a central 95% range of 14 to 31 µg/mL, and a median of 19 µg/mL. This is in good agreement with the literature, which indicates a reference range of 15 to 34 µg/mL for TBG as determined by RIA.⁸

TBG Saturation Index (TBG-SI) values were calculated for these 75 samples, using results obtained with the IMMULITE TBG and IMMULITE Total T4 kits. The TBG-SI values ranged from 17 to 39% (central 95%), with a median of 27%. This is likewise in agreement with the literature, which indicates a reference range of 2 to 5 for the T4/TBG ratio calculated as T4 (in µg/dL) times 10 divided by TBG (in µg/mL);⁸ this corresponds to a range of 14 to 35% for the TBG-SI.

Twenty-four samples obtained from women in various stages of pregnancy, when assayed by the IMMULITE TBG kit, yielded TBG levels ranging from 27 to 66 µg/mL, with a median of 43 µg/mL. The literature suggests 47 to 59 µg/mL as the range of TBG values expected for pregnant women in the third trimester.⁸

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Many physiological, pharmacological, pathological, and genetic factors affect the interpretation of TBG results. (Refer to the Summary and Explanation section.)

Increased TBG levels may occur with estrogen therapy (including oral contraceptives), pregnancy, some cases of hepatitis, and as an inherited condition. Decreased TBG levels may be due to the presence of androgens or anabolic steroids; nephrotic syndrome; severe liver

disease; increased protein catabolism with fever, wasting diseases, surgery, stress, and as an inherited condition.^{10,11}

The TBG saturation index (TBG-SI), like other total T4/TBG ratios, does not necessarily convey the same information as a free thyroxine index (FT4I).¹⁵

This assay has not been validated for use with neonates and children.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in µg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

µg/mL × 18.5 → nmol/L
(WHO 1st IS 88/638)

Calibration Range: 3.5–80 µg/mL
(65–1480 nmol/L)

Analytical Sensitivity: 1.1 µg/mL
(20 nmol/L)

Intraassay Precision (Within-Run): See "Intraassay Precision" table for statistics based on at least 20 degrees of freedom per sample.

Interassay Precision (Run-to-Run):

Statistics were calculated for samples assayed in at least 20 different runs. (See "Interassay Precision" table.)

Specificity: The antibody is highly specific for TBG. (See "Specificity" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three TBG solutions (308, 529 and 1068 µg/mL) were assayed. (See "Recovery" table for representative data.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of packed red blood cells in concentrations up to 30 µL/mL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 17 volunteers into plain, heparinized, EDTA and Becton Dickinson SST® vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of TBG, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE TBG procedure.

(Heparin) = 0.96 (Serum) + 0.35 µg/mL
r = 0.95

(EDTA) = 1.0 (Serum) + 1.1 µg/mL
r = 0.96

(SST) = 1.0 (Plain Tubes) – 1.3 µg/mL
r = 0.96

Means:

26 µg/mL (Serum)
26 µg/mL (Heparin)
26 µg/mL (SST)
27 µg/mL (EDTA)

Method Comparison: The assay was compared to an immunoradiometric assay for TBG (Kit A) on a total of 99 samples. Seventy-five of these samples were obtained from adult men and nonpregnant

women in good health; the other 24 were obtained from women in various stages of pregnancy. (Concentration range: approximately 14 to 58 µg/mL. See graph.) By linear regression:

$$(IML) = 1.1 (\text{Kit A}) - 2.1 \mu\text{g/mL} \\ r = 0.94$$

Means:

26 µg/mL (IMMULITE)

25 µg/mL (Kit A)

The assay was also compared to Coat-A-Count® T3 Uptake on the same 99 samples: a correlation coefficient (r) of 0.94 was observed.

References

1) Gershengorn MC, Glinoe D, Robbins J. Transport and metabolism of thyroid hormones. In: De Visscher M, editor. The thyroid gland. New York: Raven Press, 1980: 81-121. 2) Refetoff S. Thyroid hormone transport. In: DeGroot LJ, et al., editors. Endocrinology, vol. 1. San Francisco: Grune & Stratton, 1979: 347-56. 3) Carayon P, Castanas E, Guibout M, Codaccioni J-L. Assessment and clinical significance of free thyroid hormone radioimmunoassays. In: Ekins R, et al., editors. Free Thyroid Hormones—Proceedings of the International Symposium held in Venice, Dec 1978. Princeton: Excerpta Medica, 1979: 181-93. 4) Caplan RH, Pagliara AS, Wickus G. Laboratory diagnosis of hyperthyroidism: a reappraisal. Postgraduate Medicine 1979;66: 75-90. 5) Wenzel, KW. Pharmacological interference with *in vitro* tests of thyroid function. Metabolism 1981;7:17-32. 6) Robbins J, Rall JE. The iodine-containing hormones. In: Gray CH, James VHT, editors. Hormones in blood, 3rd edition. San Francisco: Academic Press, 1979: 576-688. 7) Witherspoon LR, El Shami AS, Shuler SE, Neely H, Sonnemaker R, Gilbert SS, Alyea K. Chemically blocked analog assays for free thyronines. I. The effect of chemical blockers on T4 analog and T4 binding by albumin and by thyroxine-binding globulin. Clin Chem 1988;34:9-23. 8) Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: Saunders, 1994: 2208. 9) National Committee for Clinical Laboratory Standards, Procedures for the collection of diagnostic blood specimens by venipuncture, 3rd ed, NCCLS Document H3-A3, 1991. 10) Davidsohn I, Henry JB, editors. Todd—Sanford clinical diagnosis by laboratory methods. 14th ed. Philadelphia: Saunders, 1969: 632. 11) Bartalena L. Studies on thyroxine-binding globulin. J Endocrinol Invest 1993;16:353-71. 12) Bristow AF, Gaines-Das RE, et al. The international standard for thyroxine binding globulin. Clin Endocrinol 1993;38:361-6. 13) Gershengorn MC, et al. Characterization of human thyroxine-binding

globulin; evidence for a single polypeptide chain. J Biol Chem 1977;252:8713-8. 14) Glinoe D, et al. Regulation of maternal thyroid during pregnancy. J Clin Endocrinol Metab 1990;71:276-87. 15) Wilke TJ. Free thyroid hormone index, thyroid hormone / thyroxine-binding globulin ratio, triiodothyronine uptake, and thyroxine-binding globulin compared for diagnostic value regarding thyroid function. Clin Chem 1983;29:74-79. 16) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8. 17) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 18) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Intraassay Precision (µg/mL)

	Mean ¹	SD ²	CV ³
1	6.7	0.8	12%
2	13	1.2	9.2%
3	21	1.6	7.6%
4	27	2.0	7.4%
5	42	3.0	7.1%
6	52	3.8	7.3%

Interassay Precision (µg/mL)

	Mean ¹	SD ²	CV ³
1	6.7	1.0	15%
2	13	1.4	11%
3	21	2.0	9.5%
4	34	2.9	8.5%
5	54	4.4	8.1%

Recovery (µg/mL)


	Spiking Solution ¹	O Observed ²	E Expected ³	%O/E ⁴
1	—	16	—	—
	A	30	31	97%
	B	43	42	102%
	C	71	69	103%
2	—	17	—	—
	A	29	31	94%
	B	45	42	107%
	C	73	69	106%
3	—	20	—	—
	A	33	34	97%
	B	47	45	104%
	C	73	72	101%
4	—	23	—	—
	A	36	37	97%
	B	52	48	108%
	C	82	75	109%

Specificity

Compound ¹	Amount Added ²	Apparent µg/mL ³	% Cross reactivity ⁴
Thyroglobulin	300 ng/mL	ND	ND
Transferrin	4 mg/mL	ND	ND
SHBG	80 nmol/L	ND	ND
AFP	58,000 ng/mL	ND	ND
HSA	5 mg/dL	ND	ND

ND: not detectable⁵

Method Comparison

(IML) = 1.1 (Kit A) – 2.1 µg/mL
r = 0.94

Linearity (µg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	42	—	—
	4 in 8	24	21	114%
	2 in 8	11	11	100%
	1 in 8	4.6	6.0	77%
2	8 in 8	46	—	—
	4 in 8	23	23	100%
	2 in 8	11	11	100%
	1 in 8	5.8	5.7	102%
3	8 in 8	70	—	—
	4 in 8	32	35	91%
	2 in 8	17	17	100%
	1 in 8	8.1	8.7	93%
4	8 in 8	80	—	—
	4 in 8	40	40	100%
	2 in 8	19	20	95%
	1 in 8	9.5	10.0	95%

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IMMULITE and Coat-A-Count are trademarks of Siemens Healthcare Diagnostics.
All rights reserved.