Free T4

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Free T4

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of non-protein-bound thyroxine (free T4) in serum and heparinized plasma. Measurements of free T4 are used in the diagnosis and treatment of thyroid disease.

Catalog Numbers:

L2KFT42 (200 tests)

L2KFT46 (600 tests)

Test Code: **F4** Color: **Dark Green**

Summary and Explanation

The principal thyroid hormone, thyroxine (T4), circulates almost entirely bound to carrier proteins, chief of which is thyroxine-binding globulin (TBG). Altered carrier protein concentrations induce changes in total T4 levels, and free T4 concentrations tend to stay within a tight range. For this reason, total T4 measurements do not always reflect thyroid status. TBG levels may vary under different physiological conditions, such as during pregnancy, oral contraceptive use, and estrogen therapy.^{1,2} Total T4 levels may increase above the normal range while free T4 remains normal. Alternatively, patients with a dysfunctional thyroid gland and altered TBG levels can have normal total T4 levels, masking the illness. Since abnormal T4 levels may signify either abnormal thyroid function or carrier protein variation (physiological or pathological), free T4 measurements more highly correlate with thyroid status than total T4 measurements.^{3,4}

Principle of the Procedure

IMMULITE 2000 Free T4 is a solid-phase, enzyme-labeled chemiluminescent competitive immunoassay. The solid phase (bead) is coated with monoclonal murine anti-T4 antibody. The liquid phase consists of alkaline phosphatase (bovine calf intestine) conjugated to T4.

The patient sample and the reagent are incubated together with the coated bead for 30 minutes. During this time, free T4 in the sample competes with enzyme conjugated T4 in the reagent for a limited number of antibody binding sites on the bead. Unbound patient sample and enzyme conjugate are then removed by centrifugal washes. Finally, chemiluminescent substrate is added to the reaction tube containing the bead and the signal is generated in proportion to the bound enzyme.

The IMMULITE 2000 Free T4 procedure is a *direct* or *single test* assay, in the sense that its results are not calculated as a function of total T4, but interpolated from a (stored) standard curve calibrated in terms of free T4 concentrations.⁶ In this respect it differs from so-called free T4 index determinations. Unlike the classic equilibrium dialysis methods, it requires neither a pre-incubation step nor preliminary isolation of the free fraction by dialysis or column chromatography.

Incubation Cycles: 1 × 30 minutes

Time to First Result: 35 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Free T4 has not been tested with all possible variations of tube types.

Volume required: 10 µL serum or heparinized plasma.

Storage. 2 days at 2–8°C, or 2 months frozen at –20°C.¹

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.⁹⁻¹¹

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**H302 + H312,
H412**

**P280, P273,
P301 + P312,
P302 + P312,
P501**

Warning! Harmful if swallowed or in contact with skin. Harmful to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: sodium azide; Free T4 Adjustors

For professional use.

For prescription use only.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: avoid contamination and exposure to direct sunlight (see insert).

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Free T4 Bead Pack (L2FT412)

With barcode. 200 beads, coated with monoclonal murine anti-T4 antibody. Stable at 2–8°C until expiration date.

L2KFT42: 1 pack

L2KFT46: 3 packs

Free T4 Reagent Wedge (L2FT4A2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to T4 in buffer. Stable at 2–8°C until expiration date.

L2KFT42: 1 wedge

L2KFT46: 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Free T4 Adjustors (LFT4L, LFT4H)

Two vials (Low and High) of lyophilized free T4 in processed human serum, with preservative. Reconstitute each vial with **2.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KFT42: 1 set

L2KFT46: 2 sets

Before running the adjustors, place the appropriate Aliquot Labels (supplied with the kit) on test tubes, so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

Also Required

Sample transfer pipets; distilled or deionized water; controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks

Quality Control: Follow government regulations or accreditation requirements for quality control frequency.

To monitor system performance and chart trends, as a minimum requirement, quality control materials with at least two levels (low and high) of Free T4 should be assayed on each day that samples are analyzed. Quality control samples should also be assayed when performing adjustment. Treat all quality control samples the same as patient samples.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory

level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

If the quality control results do not fall within the Expected Values or within the laboratory's established range, do not report patient results. Take the following actions:

- Verify that the materials are not expired.
- Verify that required instrument maintenance was performed.
- Verify that the assay was performed according to the instructions for use.
- Re-run the assay with fresh quality control samples.
- If necessary, contact technical support provider for assistance.

Expected Values

Adult: Based on its relationship to ADVIA Centaur FrT₄ assay (see Method Comparison), IMMULITE 2000 Free T4 can be expected to have essentially the same reference ranges:

	FT4 Range (ng/dL)	FT4 Range (pmol/L)
Euthyroid	0.89–1.76	11.5–22.7
Hypothyroid	< 0.89	< 11.5
Hyperthyroid	> 1.76	> 22.7

Pediatric: Reference intervals for the pediatric population (children and adolescents) were established for the IMMULITE Free T4 assay in accordance with CLSI guideline EP28-A3C.¹²

For analysis of data, the population was divided into three age subgroups:

- Infants: subjects aged 1–23 months
- Children: subjects aged 2–12 years
- Adolescents: subjects aged 13–20 years

A non-parametric approach was used to establish the reference intervals for children and adolescents where the 2.5 and 97.5 percentiles of the distribution of values were calculated. For the infant population, a robust measure of location and spread, as developed by Horn and

Pesce, was used to estimate the 2.5 and the 97.5 percentile reference intervals, accommodating the smaller sample size.^{12, 13, 14}

The reference intervals detailed by age group and number of samples are presented in the Reference Intervals table.

**IMMULITE 2000 Free T4
Pediatric Reference Intervals**

Age Group	n	Conventional (ng/dL)	SI (pmol/L)
Infants (1–23 Months)	81	0.80–1.27	10.3–16.3
Children (2–12 Years)	197	0.74–1.28	9.5–16.5
Adolescents (13–20 Years)	148	0.75–1.27	9.7–16.3

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

The interpretation of free T4 results is complicated by a variety of drugs which can affect the binding of T4 to the thyroid hormone-carrier proteins or interfere with its metabolism to T3.⁴

In severe nonthyroidal illness, the assessment of thyroid status becomes especially difficult. Since some patients in this category may suffer from concomitant primary hypothyroidism or from a compensatory secondary hypothyroidism, TSH immunoassay has been recommended as a confirmatory test in this context.³

In rare conditions associated with extreme variations in the albumin-binding capacity for T4 — such as familial dysalbuminemic hyperthyroxinemia (FDH) — direct free thyroid hormone assays may yield misleading results.

Circulating autoantibodies to T4, and hormone-binding inhibitors may interfere.

Heparin has been reported to have *in vivo* and *in vitro* effects on free T4 assays. Hence, samples should not be collected during or soon after the administration of this anticoagulant.

Because EDTA would have a significant effect on results, it should not be used as an anticoagulant.

Heterophilic antibodies in human serum or plasma can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data representative of the assay's performance. Results are expressed in ng/dL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

ng/dL × 12.87 → pmol/L

Reportable Range: 0.3–6 ng/dL

(3.9–77.2 pmol/L)

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity:

Limit of Blank (highest value expected for a sample with no analyte; determined in accordance with CLSI EP17-A⁷): 0.11 ng/dL (1.42 pmol/L)

Limit of Detection (lowest detectable concentration; determined in accordance with CLSI EP17-A⁷): 0.22 ng/dL (2.83 pmol/L)

Functional Sensitivity: (concentration with 20% coefficient of variation (CV) determined in accordance with CLSI EP5-A2⁸): 0.31 ng/dL (3.99 pmol/L)

Precision: Samples were processed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Specificity: The antibodies are highly specific for free T4. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 634 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 1000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 49 volunteers into plain, heparinized and Becton Dickinson SST® plastic vacutainer tubes. Ten of the matched samples were spiked with various concentrations of T4, to obtain values throughout the reportable range of the assay, and then assayed by the IMMULITE 2000 Free T4 procedure.

(Heparin) = 1.05 (Serum) – 0.018 ng/dL
 $r = 0.989$

(SST) = 1.04 (Serum) – 0.044 ng/dL
 $r = 0.997$

Means:

1.53 ng/dL (Serum)
1.59 ng/dL (Heparin)
1.55 ng/dL (SST)

Method Comparison: The IMMULITE 2000 Free T4 assay was compared to ADVIA Centaur FrT₄ on 282 patient samples. (Concentration range: approximately 0.3 to 5.2 ng/dL. See graph.) By linear regression:

(IML 2000) = 1.06 (Centaur) – 0.001 ng/dL
 $r = 0.981$

Means:

1.43 ng/dL (IMMULITE 2000)
1.35 ng/dL (Centaur)

References

- 1) Hay ID, Bayer MF, et al. American Thyroid Association assessment of current free thyroid hormone and thyrotropin measurements and guidelines for future clinical assays. Clin Chem 1991;37:2002-8.
- 2) Lindstedt G, et al. Clinical use of laboratory thyroid tests and investigations. JIFCC 1994;6:136-41.
- 3) Mandel SJ, Brent GA, Larsen PR. Levothyroxine therapy in patients with thyroid disease. Ann Intern Med 1993;119:492-502.

- 4) Singer PA, Cooper DS, et al. Treatment guidelines for patients with hyperthyroidism and hypothyroidism. JAMA 1995;273:808-12.
- 5) Witherspoon LR, El Shami AS, et al. Chemically blocked analog assays for free thyronines. Clin Chem 1988;34:9-16 and 17-23.
- 6) Wosilait WD. A theoretical analysis of the distribution of thyroxine among sites on thyroid binding globulin, thyroid binding prealbumin and serum albumin. Res Commun Chem Pathol Pharmacol 1977;16:541-8.
- 7) CLSI. Protocols for the Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. CLSI document EP17-A Vol 24 (No 34). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2004.
- 8) Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition. CLSI document EP5-A2. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2004.
- 9) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8.
- 10) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
- 11) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.
- 12) Clinical and Laboratory Standards Institute. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Guideline EP28-A3C.
- 13) Horn PS, Pesce AJ, Reference Intervals. A User's Guide, Washington, DC: AACC Press; 2005.
- 14) Reed AH, Henry RJ, Mason WB. Influence of statistical method used on the resulting estimate of normal range. 1971; 17:275-284.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (ng/dL)


	Within Run ¹			Total ²	
	Mean ³	SD ⁴	CV ⁵	SD	CV
1	0.51	0.040	7.8%	0.052	10.2%
2	0.85	0.038	4.5%	0.060	7.1%
3	1.13	0.067	5.9%	0.072	6.4%
4	1.49	0.072	4.8%	0.090	6.0%
5	2.91	0.103	3.5%	0.104	3.6%
6	4.82	0.144	3.0%	0.172	3.6%

Specificity

ng/dL Added ¹	Apparent ng/dL ²	%Cross-reactivity ³
Albumin		
50	ND	ND
D-T4		
1000	ND	ND
3,5-Diiodo-L-thyronine		
10,000	ND	ND
3,5-Diiodo-DL-thyrosine		
100,000	ND	ND
L-T3		
100,000	1.15	ND
5,5-Diphenylhydantoin		
4,000,000	0.70	ND
3-Monoiodo-L-tyrosine		
100,000	ND	ND
Salicylic Acid		
50	0.63	1.26%
3,3',5,5'-Tetraiodothyroacetic Acid		
2000	ND	ND

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
All rights reserved.