IMMULITE®

IGF-I

For use on the IMMULITE®
and IMMULITE® 1000 systems

SIEMENS

IMMULITE®/IMMULITE 1000 IGF-I

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® and IMMULITE 1000 Analyzers — for the quantitative measurement of insulin-like growth factor I (IGF-I) in serum or heparinized plasma, as an aid in the evaluation of growth disorders.

Catalog Number: **LKIGF1** (100 tests)

Test Code: **IGF** Color: **Dark Blue**

Summary and Explanation

Insulin-like growth factor I is a single polypeptide chain with three intra-molecule disulfide bonds. It consists of 70 amino acid residues with a molecular mass of 7649 daltons.¹ It is structurally homologous to IGF-II and insulin. IGF-I circulates primarily in a high molecular weight tertiary complex with IGF-binding protein-3 (IGFBP-3) and acid-labile subunit.^{2,3} *In vivo*, IGF-I synthesis is stimulated by growth hormone (GH) and nutritional intake.⁴

In humans, plasma IGF-I levels are barely detectable at birth, rise gradually during childhood, peak during mid-puberty until approximately 40 years of age, then decline gradually.⁵ Maternal plasma levels increase during pregnancy. In the diagnosis of growth disorders, measurements of IGF-I are a useful indicator of growth hormone (GH) secretion. A normal plasma or serum IGF-I concentration is strong evidence against GH deficiency. A low IGF-I value implies GH deficiency and requires additional testing to determine whether GH secretion is subnormal. Measurement of IGF-I is also useful in assessing change of nutritional status.

Measurement of IGF-I in serum is complicated by the presence of acid-labile components and binding proteins.¹ Acid treatment is necessary to release IGF-I to ensure accurate quantitation.^{1,6}

Principle of the Procedure

IMMULITE/IMMULITE 1000 IGF-I is a solid-phase, enzyme-labeled chemiluminescent immunometric assay.

Incubation Cycles: 1 × 60 minutes.

Specimen Collection

Because EDTA would affect results, it should not be used as an anticoagulant.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on tube materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 IGF-I has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 25 µL of the pretreated serum or heparinized plasma sample (1-in-10 dilution with IGF-I Pretreatment Solution). (Sample cup must contain at least 100 µL more than the total volume required.)

Storage: Samples stored at 2–8°C or at room temperature (22°C) are stable for up to 24 hours. Samples stored at –25°C are stable for 12 months.⁸

Patient Sample Pretreatment: Samples must be pretreated manually by diluting 1-in-10 with IGF-I Pretreatment Solution (LGFA), e.g. by adding 25 µL sample to 225 µL solution.


The pretreated samples, both serum and heparinized plasma, are stable refrigerated at 2–8°C for 24 hours.⁹

Note that printed results have been corrected for the 1-in-10 dilution. Samples expected to have values outside the calibration range should be further diluted. Correction for further dilution must be entered manually.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

IGF-I Test Units (LIGF1)

Each barcode-labeled unit contains one bead coated with monoclonal murine anti-IGF-I. Stable at 2–8°C until expiration date.

LKIGF1: 100 units

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

IGF-I Reagent Wedge (LIGF2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal rabbit anti-IGF-I in buffer. Store capped and refrigerated: stable at 2–8°C until expiration date.

Recommended usage is within 30 days after opening when stored as indicated.

LKIGF1: 1 wedge

IGF-I Adjustors (LGFL, LGFH)

Two vials (Low and High), 4 mL each, of IGF-I in a protein/buffer matrix. (No predilution is required.) Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKIGF1: 1 set

IGF-I Pretreatment Solution (LGFA)

For the manual pretreatment of patient samples. 30 mL IGF-I-free buffer matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKIGF1: 1 vial

Kit Components Supplied Separately

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

LGCOCM: Bi-level IGF control module in a protein/buffer matrix

Also Required

Sample transfer pipets, pipets for pretreatment of patient samples, distilled or deionized water, controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:

2 weeks

Quality Control Samples: Use controls or serum pools with at least two levels (low and high) of IGF-I.

Expected Values



A reference range study was performed with the IMMULITE IGF-I kit on 1805 pediatric and adult samples from an apparently healthy population, yielded the following results, presented in the tables and graphs below.

IGF-I Pediatric Reference Ranges

IGF-I, ng/mL			
Age, y	Gender	Median	Central 95% Range
0–3	Male	41	< 15–129
4–6	Male	82	20–196
7–9	Male	125	39–251
10–11	Male	159	65–300
12–13	Male	286	133–471
14–15	Male	304	163–464
16–18	Male	246	151–359
0–3	Female	62	< 15–166
4–6	Female	92	25–211
7–9	Female	129	47–265
10–11	Female	230	106–431
12–13	Female	279	149–473
14–15	Female	288	167–431
16–18	Female	256	168–367

IGF-I Adult Reference Ranges

IGF-I, ng/mL		
Age, y	Median	Central 95% Range
19–21	208	113–332
22–24	171	96–299
25–29	162	82–269
30–34	142	69–243
35–39	125	61–229
40–44	121	56–224
45–49	115	51–221
50–54	108	46–216
55–59	99	42–214
60–64	108	40–222
65–69	101	37–225
70–79	87	32–214
80–90	83	28–205

IGF-I Pediatric Data**IGF-I Adult Data**

IGF-I Tanner Stage Reference Ranges

Female, $n = 234$

IGF-I, ng/mL		
Tanner Stage	Median	Central 95% Range
1	184	70–386
2	285	120–499
3	309	155–512
4	288	161–441
5	247	155–365

Male, $n = 250$

IGF-I, ng/mL		
Tanner Stage	Median	Central 95% Range
1	142	57–255
2	224	105–386
3	269	146–456
4	248	145–393
5	223	141–336

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. Unless otherwise noted, all were generated on serum samples.

Conversion Factor: $\text{ng/mL} \times 1 \rightarrow \mu\text{g/L}$

$\text{ng/mL} \times 0.13 \rightarrow \text{nmol/L}$

Calibration Range: 15–1000 ng/mL

The IMMULITE 1000 IGF-I assay is traceable to an internal standard.

Detection Limits:

Limit of Detection (LoD) = 14.4 ng/mL

Limit of Quantitation (LoQ) = 19.7 ng/mL

LoD and LoQ were determined in accordance with Clinical and Laboratory Standards Institute (CLSI) EP17-A2.¹⁰

High-dose Hook Effect:

None up to 60,000 ng/mL

Precision: Samples were processed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates per CLSI guideline EP05-A3.¹¹ (See "Precision" table.)

Linearity: Linearity was assessed in accordance with CLSI guideline EP06-A.¹² The assay was demonstrated to be linear from 15–1000 ng/mL.

Recovery: Mixing recovery studies were carried out using native patient samples. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for IGF-I, as assessed per CLSI guideline EP07-A2.¹³ (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay. Assessed per CLSI guideline EP07-A2.¹³

Hemolysis: Presence of hemoglobin may cause a depression in apparent IGF-I concentrations. (See "Hemolysis" table.) Assessed per CLSI guideline EP07-A2.¹³

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay. Assessed per CLSI guideline EP07-A2.¹³

Alternate Sample Types: To assess the effect of alternate sample types, blood was collected from 43 volunteers into plastic serum tubes, plastic heparinized tubes and plastic gel barrier tubes (SST®). All tubes were from Becton Dickinson. The samples were assayed by the IMMULITE 1000 IGF-I procedure. Results were analyzed by Passing and Bablok

linear regression. Because EDTA would affect results, it should not be used as an anticoagulant.

(Heparin Plastic) = 0.99 (Serum Plastic) – 0.7 ng/mL
(SST Plastic) = 0.98 (Serum Plastic) – 0.3 ng/mL

Method Comparison 1:

The IMMULITE 1000 IGF-I assay was compared to the IMMULITE 2000 IGF-I assay (WHO NIBSC 1st IS 02/254) on 163 patient samples. (Concentration range: approximately 26–862 ng/mL. See graph.) Results were analyzed by Passing and Bablok linear regression.

(IML) = 1.00 (IML 2000) – 9.9 ng/mL

Method Comparison 2:

The IMMULITE 1000 IGF-I assay was compared to a commercially available IGF-I assay (IDS iSYS) on 164 patient samples. (Concentration range: approximately 23–900 ng/mL. See graph.) Results were analyzed by Passing and Bablok linear regression.

(IML) = 0.83 (IDS iSYS IGF-I) – 5.2 ng/mL

References

1) Daughaday WH, Rotwein P. Insulin-like growth factors I and II. Peptide, messenger ribonucleic acid and gene structures, serum, and tissue concentrations. *Endocr Rev* 1989;10: 68-91. 2) Baxter RC, Martin JL, Beniac VA. High molecular weight insulin-like growth factor binding protein complex. Purification and properties of the acid-labile subunit from human serum. *J Biol Chem* 1989;264:11843-8. 3) Rechler MM. Insulin-like growth factor binding proteins. *Vitam Horm* 1993;47:1-114. 4) Lee PD, Rosenfeld RG. Clinical utility of insulin-like growth factor assays. *Pediatrics* 1987;14: 154-61. 5) Hall K, Enberg G, Ritzen M, et al. Somatomedin A levels in serum from healthy children and from children with growth hormone deficiency or delayed puberty. *Acta Endocrinol (Copenh)* 1980;94:155-65. 6) Underwood LE, Murphy MG. Radioimmunoassay of the somatomedins. In: Patrono C, Peskar BA, editors. *Radioimmunoassay in basic and clinical pharmacology. Handbook of experimental pharmacology. vol. 82.* New York: Springer-Verlag; 1987: 561-74. 7) National Committee for Clinical Laboratory Standards. *Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed.* NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 8) Elmlinger MW, Zwirner, M, Kühnel, W. Stability of Insulin-Like Growth Factor (IGF)-I and IgF Binding Protein (IGFBP)-3 Measured by the Immulite® Automated Chemi-Luminescence Assay System in Different Blood Specimens. *Clin Lab* 2005; 51: 145-152. 9) Elmlinger MW,

Kühnel W, Weber MM, Ranke MB. Reference ranges for two automated chemiluminescent assays for serum insulin-like growth factor I (IGF-I) and IGF-binding protein 3 (IGFBP-3). *Clin Chem Lab Med* 2004;42:654-64. 10) EP17-A2: Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition, CLSI, Wayne, PA, 2012. 11) EP05-A3: Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition, CLSI, Wayne, PA, 2014. 12) EP06-A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline, CLSI, Wayne, PA, 2003. 13) EP07-A2: Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition, CLSI, Wayne, PA, 2005.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

siemens.com/healthcare

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Repeatability ¹		Within-Lab ²	
		SD ⁴	CV ⁵	SD	CV
1	50	1.7	3.4%	2.1	4.2%
2	99	4.8	4.8%	5.7	5.7%
3	223	5.8	2.6%	6.3	2.8%
4	501	10.2	2.0%	16.9	3.4%
5	718	20.7	2.9%	22.5	3.1%
6	381	7.7	2.0%	11.9	3.1%

Hemolysis

	Unspiked ²	Hemoglobin ¹		
		150 mg/dL	300 mg/dL	500 mg/dL
1	183	171	150	140
2	359	353	304	285

Recovery (ng/mL)


Sample Mixture ¹	Observed ²	Expected ³	% Bias ⁴
1	429	470	-9%
2	465	480	-3%
3	460	480	-4%
4	404	408	-1%
5	376	389	-3%

Specificity


Compound ¹	Added ²	% Cross-reactivity at 0 ng/mL ³	% Cross-reactivity at 60 ng/mL ⁴
Proinsulin (ng/mL)	155,000	ND	0.00%
Insulin (μIU/mL) μIU/mL × 0.042 → ng/mL	18,000	ND	-0.29%
LH (mIU/mL)	70,000	ND	0.07%
TSH (mIU/mL)	310	ND	0.00%
IGF-II (ng/mL)	6000	0.25%	0.50%

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
All rights reserved.