IMMULITE®

***H. pylori* IgG**

**For the Detection of IgG
Antibodies to *H. pylori*
in Human Serum**

**For use on the IMMULITE®
and IMMULITE® 1000 systems**

SIEMENS

IMMULITE®/IMMULITE® 1000 *H. pylori* IgG

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the qualitative detection of IgG antibodies to *Helicobacter pylori* in human serum from symptomatic adults, as an aid in the diagnosis of *Helicobacter pylori* infection.

Catalog Number: **LKHPG1** (100 tests)

Test Code: **HPG** Color: **Dark Green**

Summary and Explanation

In 1983, Warren and Marshall described a curved bacillus associated with gastric mucosa in cases of chronic gastritis, tentatively drawing comparisons between that organism and the genus *Campylobacter*. More recently, investigators have shown a correlation between colonization with this organism (now termed *Helicobacter pylori*) and gastric and duodenal ulcers and chronic gastritis. *H. pylori* colonization is chronic in nature and appears to cause histological inflammation of the gastric mucosa. When *H. pylori* is eliminated from the gastric mucosa, the inflammation tends to decrease. If the organism recolonizes, the inflammation increases in severity and coincides with the appearance of gastrointestinal symptoms. A lack of gastrointestinal symptoms in the presence of *H. pylori* indicates colonization rather than infection. The most recent clinical trials have confirmed *H. pylori* to be the causative agent for most cases of chronic gastritis and ulcers. Evidence has been presented that *H. pylori* is also associated with gastric carcinoma.

Current procedures for detecting *H. pylori* rely on isolation of the organism from tissue obtained from endoscopic biopsy. The tissue is also tested by histology, gram staining, and/or testing for the enzyme urease, which is produced by the organism in abundance. Each of the methods, including histology and gram staining, require biopsy material taken from multiple sites. In addition to the necessity for performing an invasive gastroscopy, each of these methods

suffers from other inadequacies.⁶ The presence of *H. pylori* has also been detected with a urea breath test (using radioactive or nonradioactive isotopes) and by serological methods.¹⁴⁻¹⁶

A positive serological response to *H. pylori* antigens has been determined in individuals with duodenitis, chronic gastritis, and gastric or duodenal ulcer. Further, many people without clinical symptoms are seropositive for *H. pylori* antibodies, with the prevalence increasing with age. Therefore, while *H. pylori* serology is a sensitive method to determine colonization, differentiation between colonization and active disease is not possible.

Principle of the Procedure

IMMULITE/IMMULITE 1000 *H. pylori* IgG is a solid-phase, chemiluminescent two-step immunometric assay.

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 *H. pylori* IgG has not been tested with all possible variations of tube types.

Volume Required: 10 μ L (prediluted) sample. (Sample cup must contain at least 100 μ L more than the total volume required.)

Storage: 3 days at 2–8°C or 6 months at –20°C.¹⁸

Patient Sample Dilution: Samples must be prediluted 1-in-21 in *H. pylori* IgG Sample Diluent (LHPGZ1). This dilution must be performed manually, e.g. by adding 20 μ L sample to 400 μ L Sample Diluent, (LHPGZ1, supplied in the kit).

Alternatively, the IMMULITE 1000 Windows® system performs the dilution on board, automatically, using LIDW (provided in the ID2 Sample Diluent Module, L1KIDW1).

Printed results have been corrected for the 1-in-21 dilution.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.²¹⁻²³

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

The bead is coated with *inactivated* *H. pylori* antigen. However, caution is advised because of the possible presence of residual organism when working with, or disposing of, the materials supplied.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

H. pylori IgG Test Units (LHPG1)

Each barcode-labeled unit contains one bead coated with inactivated, partially purified *H. pylori* antigen. Stable at 2–8°C until expiration date.

LKHPG1: 100 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

H. pylori IgG Reagent Wedges (LHPGA, LHPGB)

With barcodes. **LHPGA:** 7.5 mL of a buffer solution. **LHPGB:** 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-human IgG antibody in buffer. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKHPG1: 1 set.

H. pylori IgG Adjustors (LHPGL, LHPGH)

Two vials (Low and High), 4 mL each, of human serum with IgG reactive to *H. pylori*, in buffer, with preservative. (No predilution is required.) Stable at 2–8°C

for 14 days after opening, or for 6 months (aliquotted) at –20°C.

LKHPG1: 1 set.

***H. pylori* IgG Controls (LHPGC1, LHPGC2, LHPGC3)**

Three vials (Negative, Low Positive and Positive) containing 2 mL each. **LHPGC1 (Negative Control):** One vial containing human serum with IgG nonreactive to *H. pylori*, with preservative. **LHPGC2, LHPGC3 (Low Positive Control, Positive Control):** Two vials containing human serum with IgG reactive to *H. pylori*, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

LKHPG1: 1 set.

IMMULITE: *H. pylori* IgG controls require the same **manual** 1-in-21 predilution as patient samples, with *H. pylori* IgG Sample Diluent (LHPGZ1).

IMMULITE 1000: The IMMULITE 1000 software performs automatic on-board dilution of patient and control samples. Enter controls as controls.

Refer to the control insert for concentration levels.

***H. pylori* IgG Sample Diluent (LHPGZ1)**

For the manual dilution of patient samples and controls. Buffer solution, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKHPG1: 2 × 25 mL

For on-board dilution of patient samples using IMMULITE 1000 Windows®, refer to the Kit Components Supplied Separately section.

Kit Components Supplied Separately

ID2 Sample Diluent Module (L1KIDW1)

Required for the on-board dilution of patient samples and controls using IMMULITE 1000 Windows®. Sample Diluent wedges (LIDW) containing concentrates (2.5X) buffer solution with preservative. Stable at 2–8°C until expiration date.

L1KIDW1: 2 Wedges.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

DCHS: Dilution Cup Holders (barcoded) for IMMULITE 1000 Windows® on-board dilutions.

LHPGCM: Tri-level *H. pylori* IgG Control Module

Also Required

Sample transfer pipets, distilled or deionized water.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Note that both Reagent Wedges A and B must be loaded on the carousel to run this assay. For IMMULITE 1000 Windows®, the Sample Diluent Wedge must also be loaded on the carousel.

Recommended Adjustment Interval: 1 week.

Quality Control Samples: Use the control(s) supplied with the kit.

Interpretation of Results

Positive: A result greater than or equal to 1.1 U/mL is considered to be "positive", and indicates that *H. pylori* IgG antibodies were detected in the sample.

Indeterminate: A result greater than or equal to 0.9 U/mL and less than 1.1 U/mL is considered to be "Indeterminate". Indeterminate samples should be retested. Samples still yielding "indeterminate" values should be examined by an alternate method, or a second sample should be collected, if possible, within a reasonable time period (e.g., one week).

Negative: A result less than 0.9 U/mL is considered to be "negative", and indicates that *H. pylori* IgG antibodies were not detected in the sample. Negative results by this test do not preclude recent primary infection.

The presence of IgG antibodies to *H. pylori* is an indication of previous exposure to the organism. A single specimen can only be used to determine the immune status of the individual.

The anti- *H. pylori* IgG results determined for a given specimen with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. Therefore, the results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE/IMMULITE 1000 *H. pylori* IgG EIA. Results obtained from other manufacturers' assay methods may not be used interchangeably."

Expected Values

The majority of individuals exposed to *H. pylori* possess IgG antibodies to the organism. In addition, the presence of *H. pylori* antibodies is a function of age, race, geography and clinical condition.¹⁶ Age-specific rates for *H. pylori* infection are similar for males and females.

The prevalence of *H. pylori* infection is between 30 and 40% in the US and Canada, about 20% in Australia, 70% in Europe, and between 70 and 90% in Africa, South America and Asia.¹⁸

H. pylori has also been identified as a risk factor for gastric cancer. Nearly 100% of patients with duodenal ulcer, 85% of those with gastric ulcer,¹⁹ and 50 – 100% of patients with gastric cancer are infected with *H. pylori*.²⁰

A relatively large proportion of patients who have positive levels of antibody are asymptomatic, even though they are colonized with the *H. pylori* organism. Therefore, antibody levels do not necessarily correlate with the severity of clinical symptoms.¹⁶

Limitations

For the determination of seroconversion from *nonreactive* to *reactive*, two serum samples should be drawn three to four weeks apart, during the acute and

convalescent stages of the infection. The acute phase sample should be stored and tested in parallel with the convalescent sample.

Individuals with an acute *H. pylori* infection may not exhibit any detectable IgG antibodies at the early stage of infection.

The results in HIV patients, in patients undergoing immunosuppressive therapy, or in patients with other disorders leading to immunosuppression, should be interpreted with caution.

The performance characteristics of this assay have not been established for use with specimens from cord blood, neonates, pediatric patients, or pretransplant patients.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in U/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: 0.4 – 8.0 U/mL

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity: 0.4 U/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Specificity: The assay is highly specific for *H. pylori* IgG antibodies, with no observed crossreactivity to *Campylobacter jejuni*, *Campylobacter fetus* and *Campylobacter coli* microorganisms.

Effect of Gel Barrier Tubes: Gel Barrier tubes have no clinically significant effect on the measurement of *H. pylori* IgG antibodies.

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 504 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison: The assay was compared to clinical diagnosis by biopsy on 155 randomly selected retrospective specimens from patients suspected of *H. pylori* infection. A specimen was considered clinically positive if the culture test, or both CLO and histology tests, were positive. A specimen was considered clinically negative if none of the biopsy tests were positive.

IMMULITE

Clinical	Pos	Ind	Neg	Sensitivity	Specificity
Diagnosis					
Pos	64	3	2	97.0%	98.8%
Neg	1	1	84		

Agreement: 98.0%

Central 95% Confidence Limits for Sensitivity and Specificity, respectively: 89.5% – 99.6% and 93.6% – 100%.

The assay was also compared to a commercially available enzyme immunoassay for *H. pylori* IgG (Kit A) on 70 patient samples.

IMMULITE

Kit A	Pos	Ind	Neg
Positive	50	3	2
Indeterminate	1	0	0
Negative	0	0	14

Agreement: 91.4%

References

- 1) Antonescu CG, Marshall BJ. *Helicobacter pylori*: A potentially curable form of peptic ulcer disease. Gastroenterol J Club 1990;1:3-12.
- 2) Barthel JS, Everett ED. Diagnosis of *campylobacter pylori* infections: the "Gold Standard" and the alternatives. Rev Infect Dis 1990;12 (suppl.1):107-14.
- 3) DeCross AJ, Peura DA. Role of *H. pylori* in peptic ulcer disease. Contemp Gastroenterol 1992;5:18-28.
- 4) George LL, Borody TJ, Andrews P, Devine M, Moore-Jones D et al. Cure of duodenal ulcer after eradication of *H. pylori*. Med J Australia 1990;153.
- 5) Goodwin CS, Armstrong JA, Marshall BJ. *Campylobacter pylori*, gastritis, and peptic ulceration. J Clin Pathol 1986;39:353-65.
- 6) Graham DY, Klein PD. What you should know about the methods, problems, interpretations, and uses of urea breath tests. Amer J Gastroenterol 1991;86 (9):1118-22.
- 7) Kosunen TU, Sepälä K, Sarna S, Sipponen P. Diagnostic value of decreasing IgG, IgA, and IgM antibody titers after eradication of *Helicobacter pylori*. Lancet 1992;339:893-5.
- 8) Parsonnet J, Friedman GD, Vandersteen DP, Chang Y, Vogelman JH et al. *Helicobacter pylori* infection and the risk of gastric carcinoma. New Engl J Med 1991;325:1127-31.
- 9) Rauws EA, Tytgat GNJ. Cure of duodenal ulcer associated with eradication of *Helicobacter pylori*. Lancet 1990;335:1233-5.
- 10) Talley NJ, Newell DG, Ormand JE, Carpenter HA, Wilson WR et al. Serodiagnosis of *Helicobacter pylori*: Comparison of enzyme-linked immunosorbent assays. J Clin Microbiol 1991;29:1635-9.
- 11) Jones DM et al. Antibody to the gastric *Campylobacter*-like organism (*Campylobacter pylori*) — clinical correlations and distribution in the normal population. J Med Microbiol 1986;22:57-62.
- 12) Thomas JE, Whatmore AM, Barer MR, Eastham EJ, Kehoe MA. Serodiagnosis of *Helicobacter pylori* infection in childhood. J Clin Microbiol 1990;28:2641-6.
- 13) Warren JR, Marshall BJ. Unidentified curved bacillus on gastric epithelium in active chronic gastritis. Lancet 1983;1:1273-5.
- 14) Fennerty MB. A review of tests for the diagnosis of *Helicobacter pylori* infection. Lab Med 1998;29:561-6.
- 15) Lichtenstein DR. Peptic ulcer disease: clinical options for detecting and eradicating *H. pylori* infection. Mod Med 1994;62.
- 16) Dunn BE, et al. *Helicobacter pylori*. Clin Microbiol Rev 1997;720-41.
- 17) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:358.
- 18) Bonney RC. The developing market for *Helicobacter pylori* testing. PJB Publications Ltd, 1997:43-4.
- 19) Marshall BJ. *Helicobacter pylori*: the etiological agent for peptic ulcer. Journal of the American Medical Association 1995; 274:1064-1067.
- 20) Vanagunas A. The link between *Helicobacter pylori* and gastric cancer. Infect Med 1998;15:644-640,656.
- 21) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human

immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377–82, 387–8. 22) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 23) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (U/mL)

	Mean ³	Within-Run ¹			Total ²	
		SD ⁴	CV ⁵	SD	CV	
1	<0.4*					
2	0.76	0.021	2.8%	0.033	4.4%	
3	1.25	0.027	2.2%	0.062	4.9%	
4	2.74	0.058	2.1%	0.106	3.9%	
5	5.14	0.157	3.0%	0.251	4.9%	
6	7.23	0.213	2.9%	0.403	5.6%	

* Results were consistently below 0.4 U/mL.⁶

Dansk. ¹I samme kørsel, ²Totalt, ³Middelværdi, ⁴Standardafvigelse, ⁵Varianskoefficient, ⁶Resultaterne var konstant under 0,4 U/mL. **Eesti.** ¹Mõõtmisseeria sisene, ²Seeriaveheline, ³Keskmine väärtused, ⁴Standardhälve, ⁵Variatsioonikoeffitsient, ⁶Tulemused olid püsivalt alla 0,4 U/mL. **Latviski.** ¹Sērijas robežās, ²Kopā, ³Vidējā vērtība, ⁴SD, ⁵Variāciju koeficients, ⁶Rezultāti bija zem 0,4 U/mL. **Lietuviškai.** ¹Vieno tyrimo metu, ²Bendras, ³Vidurkiai, ⁴SD, ⁵CV, ⁶Rezultatai nuolat buvo žemesni už 0,4 U/ml. **Svenska.** ¹Inom körning, ²Totalt, ³Medel, ⁴SD, ⁵CV, ⁶Resultaten var konsekvent under 0,4 U/mL.

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IMMULITE® is a trademark of Siemens Healthcare Diagnostics.
© 2009 Siemens Healthcare Diagnostics. All rights reserved.