IL2R

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 IL2R

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of soluble interleukin 2 receptor (IL2R) in serum, as an aid in the study of inflammatory diseases.

Catalog Number: **L2KIP2** (200 tests)

Test Code: **I2R** Color: **Light Blue**

Summary and Explanation

The receptor of the cytokine interleukin 2 (IL-2) plays a crucial role in the regulation of the immune response. Binding of IL-2 to its receptor (IL2R) on the surface of T-lymphocytes triggers a series of intracellular signaling events that result in the activation and proliferation of resting T cells and, ultimately, in the generation of helper, suppressor and cytotoxic T cells, which mediate immune reactions.

The IL-2 receptor is unique among growth factor receptors in that it is made up of at least three distinct membrane components: the α chain (IL2R α) — the so-called Tac antigen (MW 55kDa), a β -chain (IL2R β ; MW 70-75 kDa), and the γ -chain (IL2R γ ; MW 64 kDa). Different combinations of these three components give rise to the generation of various forms of the IL-2 receptor, each of which manifests different binding affinities to IL-2.

Most resting T cells, B cells, large granular lymphocytes (LGLs) and monocytes do not express significant numbers of this receptor on their surfaces. Upon activation, receptor molecules are expressed on the surface of the cells, and a soluble form (sIL2R) is released, which is about 10 kDa smaller than the membrane-bound protein.

It has been found that sIL2R is present at low levels in serum of healthy individuals and at significantly elevated levels in a broad range of disorders such as neoplastic diseases, autoimmune diseases, organ allograft rejection and different infections. Thus, it appears that

sIL2R can serve as a marker for diagnosis, therapeutic evaluation and management of cancer, as well as an indicator of a wide spectrum of disorders involving immune activation.

Principle of the Procedure

IMMULITE 2000 IL2R is a solid-phase, two-site chemiluminescent immunometric assay.

The solid phase (a bead coated with murine monoclonal anti-IL-2R antibody) is mixed with reagent (containing alkaline phosphatase conjugated to a rabbit anti-IL-2R polyclonal antibody) and sample. Solid phase, reagent, and sample are incubated together. Soluble IL-2R in the sample binds to the solid phase antibody and reagent antibody forming a sandwich complex on the bead. Unbound enzyme conjugate is then removed by a centrifugal wash. Finally, chemiluminescent substrate is added to the bead and signal is generated in proportion to the bound enzyme.

Incubation Cycles: 1 × 30 minutes

Specimen Collection

EDTA plasma should not be used as a sample type.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants.

IMMULITE 2000 IL2R has not been tested with all possible variations of tube types.

Volume Required: 15 μ L serum

Storage: 2 days at 2–8°C or for longer at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**H302 + H312,
H412**

**P280, P273,
P301 + P312,
P302 + P312,
P501**

Warning! Harmful if swallowed or in contact with skin. Harmful to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: sodium azide; IL2R Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

IL2R Bead Pack (L2IP12)

With barcode. 200 beads, coated with a murine monoclonal anti-IL2R antibody. Stable at 2–8°C until expiration date.

L2KIP2: 1 pack

IL2R Reagent Wedge (L2IPA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to rabbit polyclonal anti-IL2R antibody in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KIP2: 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

IL2R Adjustors (LIPL, LIPH)

Two vials (Low and High) of lyophilized IL2R in a protein-based buffer, with preservative. Reconstitute each vial with **3.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 3 months after reconstitution. Do not freeze. Following reconstitution, the Adjustors should be aliquoted, to avoid contamination from repeated usage.

L2KIP2: 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Multi-Diluent 2 (L2M2Z, L2M2Z4)

For the on-board dilution of patient samples. One vial of concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcode can be read by the on-board reader.

L2M2Z: 3 labels **L2M2Z4:** 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LILCM: IMMULITE Cytokine Control Module (Bi-level, human serum-based)

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
4 weeks

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of IL2R.

Expected Values

A study performed with IMMULITE 2000 IL2R on 100 apparently healthy adults yielded a median of 333 U/mL and a 95% reference range of 158 to 623 U/mL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitation

Dilution of samples to concentrations below 200 U/mL is not recommended.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all

immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in U/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Reportable Range: 50 to 7,500 U/mL

Analytical Sensitivity: 5 U/mL

High-dose Hook Effect: None up to 1,030,000 U/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three IL2R solutions (5,000, 25,000 and 50,000 U/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for IL2R. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison: The assay was compared to IMMULITE IL2R on 131 spiked and unspiked patient samples. (Concentration range: approximately 130 to 6,990 U/mL. See graph.) By linear regression:

$$(IML\ 2000) = 1.03 (IML) + 5.5\ U/mL$$

$$r = 0.995$$

Means:

767 U/mL (IMMULITE 2000)

736 U/mL (IMMULITE)

References

1) Rubin LA, et al. Identification and characterization of a released form of the Interleukin-2 Receptor. Leukocytes and host defense 1986;95:102. 2) Rubin LA, et al. Soluble interleukin-2 receptors are released from activated human lymphoid cells *in vitro*. J Immunol 1985;135:3172-7. 3) Rubin LA, et al. The released interleukin-2 receptor binds interleukin-2 efficiently. J Immunol 1986;137:3841-4. 4) Pizzolo G. The soluble Interleukin-2 receptor as a new biological marker in diseases. Immunol Clin 1988;VII:1. 5) Pui Ch-H. Serum Interleukin-2 Receptor: Clinical and biological Implications. Leukemia 1989;3:323-7. 6) Rubin LA. The soluble Interleukin-2 Receptor: biology function, and clinical application. Annals Internal Med 1990;113:619-27. 7) Zerler B. The soluble Interleukin-2 Receptor as a marker for Human Neoplasia and Immune status. Cancer Cells 1991;3:471-9. 8) Taniguchi T, Minami Y. The IL-2/IL-2 Receptor System: A current overview. Cell 1993;73:5-8. 9) Waldmann ThA. The IL-2/IL-2 receptor system: a target for rational immune intervention. Immunol Today 1993;14: 264-70. 10) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

Available outside the United States only.
For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (U/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	207	5.7	2.8 %	8.0	3.9%
2	482	16.3	3.4%	21.2	4.4%
3	957	30.7	3.2%	34.1	3.6%
4	1,635	40.9	2.5%	59.9	3.7%
5	3,602	91	2.5%	133	3.7%
6	6,077	174	2.9%	250	4.1%

Specificity

Compound ¹	U/mL Added ²	Apparent U/mL ³	% Cross-reactivity ⁴
IL-1 α	10,000	ND	ND
IL-2	10,000	ND	ND
IL-4	10,000	ND	ND
IL-6	10,000	ND	ND
IL-8	10,000	ND	ND
IL-1 β	10,000	ND	ND
Interferon γ	10,000	ND	ND
TNF α	10,000	ND	ND

ND: not detectable⁵


Linearity (U/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	568	—	—
	4 in 8	271	284	95%
	2 in 8	126	142	89%
	1 in 8	56	71	78%
2	8 in 8	619	—	—
	4 in 8	300	310	97%
	2 in 8	135	155	87%
	1 in 8	62	77	80%
3	8 in 8	632	—	—
	4 in 8	295	316	93%
	2 in 8	133	158	84%
	1 in 8	57	79	72%
4	8 in 8	676	—	—
	4 in 8	316	338	93%
	2 in 8	148	169	88%
	1 in 8	63	85	75%
5	8 in 8	764	—	—
	4 in 8	360	382	94%
	2 in 8	170	191	89%
	1 in 8	74	96	78%
6	8 in 8	807	—	—
	4 in 8	384	404	95%
	2 in 8	171	202	85%
	1 in 8	79	101	78%

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2008 Siemens Healthcare Diagnostics. All rights reserved.