IMMULITE®

Homocysteine

For use on the IMMULITE®
and IMMULITE® 1000 systems

SIEMENS

IMMULITE®/IMMULITE 1000 Homocysteine

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® and IMMULITE 1000 Analyzers — for the quantitative determination of L-homocysteine in human plasma or serum. This device can assist in the diagnosis and treatment of patients suspected of having hyperhomocysteinemia or homocystinuria.

Catalog Number: **LKHO1** (100 tests)

Test Code: **HCY** Color: **Dark Gray**

Warning: Specimens from patients who are on drug therapy involving S-adenosyl-methionine may show falsely elevated levels of homocysteine.

Results on specimens obtained from patients taking methotrexate, carbamazepine, phenytoin, nitrous oxide, anti-convulsants, or 6-azauridine triacetate should be interpreted with caution as these substances interfere with homocysteine determination.

Summary and Explanation

Total homocysteine (tHcy) has emerged as an important risk factor in the assessment of cardiovascular disease.^{3-7,12} Hcy, a thiol-containing amino acid, is produced by the intracellular demethylation of methionine. Hcy thereby serves as a pool that can be later scavenged for use in the remanufacture of either methionine through the action of the folate-dependent enzyme methionine synthase or cysteine using the B6 dependent transsulfuration pathway.^{1,2} Hcy in plasma is found primarily in a protein bound form but free, oxidized and disulfide forms are also present.

Highly elevated levels of tHcy are found in patients with homocystinuria, a rare genetic disorder of the enzymes involved in Hcy metabolism.^{1,2,3} Homocystinuria patients exhibit arterial thromboembolism, mental retardation and early

arteriosclerosis.¹ Less severe genetic defects are also associated with moderate levels of Hcy.^{3,4,5}

Homocysteine has been identified as an indicator of cardiovascular disease. A meta-analysis of 27 epidemiological studies has suggested that a 5 µmol/L increase in tHcy could be associated with an odds ratio for coronary artery disease (CAD) of 1.6 for men and 1.8 for women, the same increase in risk as a 0.5 mmol/L increase in cholesterol.⁷ In addition, patients with chronic renal disease complicated by arteriosclerotic cardiovascular disease show elevated tHcy due to the inability of the kidney to remove Hcy from the blood.^{1,2,3}

Principle of the Procedure

Competitive Immunoassay.

IMMULITE/IMMULITE 1000

Homocysteine involves a preliminary manual sample pretreatment step. Homocysteine in the patient plasma or serum sample is released from its binding proteins and converted to S-adenosyl-homocysteine (SAH) by an off-line 30-minute incubation at 37°C in the presence of S-adenosyl-L-homocysteine hydrolase and dithiothreitol (DTT).

The treated sample and alkaline phosphatase-labeled anti-SAH antibody are simultaneously introduced into a test unit containing an SAH-coated polystyrene bead. During a 30-minute incubation, the converted SAH from the patient sample competes with the immobilized SAH for binding the alkaline phosphatase labeled-anti-SAH antibody conjugate. Unbound enzyme conjugate is removed by centrifugal wash. Substrate is added, and the procedure continues as described for typical immunoassays in the Operator's Manual.

Incubation Cycle: 1 × 30 minutes

Time to first result: 42 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Heparinized and EDTA plasma are the samples of choice, but serum is also suitable for use. It is important to separate the plasma or serum from the cells as soon as possible after collection, as synthesis of HCY will take place in red blood cells after sampling. **Samples should be stored on ice between the time of sampling and centrifugation.** Note that storage on ice makes the use of serum samples particularly difficult.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Homocysteine has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required – untreated sample:
15 µL untreated plasma or serum

Volume Required – treated sample:
75 µL treated sample. (Sample cup must contain at least 100 µL more than the total volume required.)

Storage: 14 days at 2–8°C, or 6 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412	Harmful to aquatic life with long lasting effects.
P273, P501	Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations. Contents: zinc chloride, 2-methyl-2H-isothiazol-3-one; Homocysteine Reagent Wedge

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight (see insert).

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

Homocysteine Test Units (LHO1)

Each barcode-labeled unit contains one bead coated with S-adenosyl-L-homocysteine (SAH). Stable at 2–8°C until expiration date.

LKHO1: 100 units

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

Homocysteine Reagent Wedge (LHO2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to murine monoclonal anti-SA-H in buffer. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is

within 30 days after opening when stored as indicated.

LKHO1: 1 wedge

Pretreatment A Solution (LHOA)

20 mL bovine S-adenosyl-L-homocysteine hydrolase in buffer, with preservative.

Stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKHO1: 1 vial

Pretreatment B Solution (LHOB)

5.5 mL dithiothreitol (DTT) in buffer. Stable at 2–8°C until expiration date.

Recommended usage is within 30 days after opening when stored as indicated.

LKHO1: 1 vial

Pretreatment B Solution must be diluted 1-in-10 with distilled or deionized water before combining it with Pretreatment A Solution. See section on Preparation of Sample Pretreatment Working Solution.

Homocysteine Adjustors (LHOL, LHOH)

Two amber vials (Low and High), 2 mL each, of synthetically-derived S-adenosyl-L-homocysteine (SAH) in a protein/buffer matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKHO1: 1 set

Kit Components Supplied Separately

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM : Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

Also Required

Sample transfer pipets, distilled or deionized water, controls,

variable micropipets: 10–500 µL

Preparation of Sample Pretreatment Working Solution

- 1 Dilute sufficient volume of Pretreatment B Solution 1-in-10 with distilled or deionized water before combining it with Pretreatment A Solution in Step 2.
- 2 Combine equal volumes of Pretreatment A Solution and diluted

Pretreatment B Solution in sufficient quantities for the number of samples to be processed.

For example, see chart below:

No. of Samples	Wkg Sol req. mL	H ₂ O mL	Pre-B µL	Pre-A mL	Total mL
5	1.5	0.90	100	1	2
10	3	1.8	200	2	4
15	4.5	2.25	250	2.5	5
20	6.0	3.15	350	3.5	7
25	7.5	3.6	400	4	8

Use Sample Pretreatment Working Solution within two hours of preparation.

Sample Pretreatment

Adjustors, Controls and patient samples must go through the Sample Pretreatment Procedure prior to analysis.

For Controls and Patient Samples

- 1 Label one test tube for each sample to be processed.
- 2 Add 300 µL of Sample Pretreatment Working Solution into the tubes prepared.
- 3 Add 15 µL of patient untreated plasma or serum to each tube.
- 4 Cover the tubes and mix well.
- 5 Incubate for 30 minutes at 37°C in a water bath or oven.
- 6 Transfer all of the treated sample into an IMMULITE Sample Cup.

For Adjustors

Because Adjustors are processed in replicates of 4, additional volume is required.

At Step 2, use 400 µL of Sample Pretreatment Working Solution.

At Step 3, use 20 µL of each of the Adjustor levels. Proceed with remaining pretreatment steps.

Treated samples, both serum and plasma, are stable at room temperature (15–28°C) or refrigerated at 2–8°C for 1 hour prior to assay.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:
4 weeks

Quality Control Samples: Follow government regulations or accreditation requirements for quality control frequency.

Use controls or sample pools with at least two levels (low and high) of homocysteine.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

Expected Values

Homocysteine levels can vary with age, gender, geographical area and genetic factors, therefore it is important for laboratories to establish their own reference ranges based on their local populations. Literature suggests a reference range of 5–15 $\mu\text{mol/L}$ for adult males and females,¹¹ but also notes that men tend to have higher levels than women, and that postmenopausal women tend to have higher levels than premenopausal.¹²

Based on its relationship to IMMULITE 2000 Homocysteine procedure (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

One hundred and twenty samples from apparently healthy adult male and female volunteers age 22–66 were analyzed using the IMMULITE 2000 Homocysteine procedure. The samples were drawn in

heparinized plasma tubes and kept on ice prior to separation of the plasma from the cells. The median value was 7.7 $\mu\text{mol/L}$, with a central 95% range of 5.0–12 $\mu\text{mol/L}$.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Specimens from patients who are on drug therapy involving S-adenosyl-methionine may show falsely elevated levels of homocysteine.

Although analysis of the parent compounds of carbamazepine, phenytoin, 6-azauridine, and anthopterin indicate no cross-reactivity, specimens obtained from patients treated with these drugs, as well as with methotrexate, nitrous oxide, and other anticonvulsants, should be interpreted with caution as these substances have been shown to interfere in some homocysteine assays.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in $\mu\text{mol/L}$. (Unless otherwise noted, all were generated on heparinized plasma samples.)

Reportable Range: 2–50 $\mu\text{mol/L}$

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity: 1.2 µmol/L

Precision: Samples were assayed in quadruplicate over the course of several days, for a total of 20 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three homocysteine solutions (105, 207 and 411 µmol/L) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for homocysteine. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To determine whether serum can be used in the IMMULITE Homocysteine procedure, blood was collected on ice from 20 laboratory volunteers into plain, heparinized and EDTA vacutainer tubes. Samples were separated from the cells, and matched samples were spiked with homocysteine and then assayed by the IMMULITE Homocysteine procedure, with the following results.

(EDTA) = 1.08 (Heparin) – 1.10 µmol/L
 $r = 0.983$

(Serum) = 1.11 (Heparin) – 1.79 µmol/L
 $r = 0.992$

Means:
9.18 µmol/L (Heparin)
8.82 µmol/L (EDTA)
8.38 µmol/L (Serum)

To assess the effect of storage temperature, heparinized, EDTA and plain vacutainer tubes were collected from five

volunteers for each tube type. Some tubes were stored at room temperature, while other tubes were kept on ice for various time periods *prior to separation*. The graphs below show the effect of storage time and temperature for heparin, EDTA and serum. (See graphs 1–3.)

Method Comparison: The assay was compared to IMMULITE 2000 Homocysteine on 257 plasma samples. (Concentration range: approximately 4 to 50 µmol/L, as measured by IMMULITE 2000. See graph.) By linear regression:

(IML) = 0.97 (IML 2000) + 0.86 µmol/L
 $r = 0.975$

Means:
13.5 µmol/L (IMMULITE)
13.1 µmol/L (IMMULITE 2000)

References

- 1) Selhub J, Miller JW. The pathogenesis of homocysteinemia: interruption of the coordinate regulation by S-adenosylmethionine of the remethylation and transsulfuration of homocysteine. *Am J Clin Nutr* 1992;55:131-8.
- 2) Jacques PF, Bostom AG, Williams RR, et al. Relation between folate status, a common mutation in methylenetetrahydrofolate reductase and plasma homocysteine concentrations. *Circulation* 1996;93:7-9.
- 3) Malinow MR. Plasma homocysteine and arterial occlusive diseases: A mini review. *Clin Chem* 1995;40:173-76.
- 4) Clarke R, Daly L, Robinson K, Naughten E, et al. Hyperhomocysteinemia: An independent risk factor for vascular disease. *New Eng J Med* 1991;324:1149-55.
- 5) Deloughery TG, Evans A, Sadeghi A, et al. Common mutation in methylenetetrahydrofolate reductase: Correlation with homocysteine metabolism and late-onset vascular disease. *Circulation* 1996;94:3074-78.
- 6) Bostom AG, Selhub J. Homocysteine and arteriosclerosis: Subclinical and clinical disease associations. *Circulation* 1999;99:2361-63.
- 7) Boushey CJ, Beresford SAA, et al. A quantitative assessment of plasma homocysteine as a risk factor for vascular disease. *JAMA* 1995;274(13):1049-57.
- 8) Guttormsen AB, Svarstad E, et al. Elimination of homocysteine in subjects with end-stage renal failure. *Irish J Med Sci* 1995;164:8.
- 9) Bostom AG, Lathrop L. Hyperhomocysteinemia in end-stage renal disease (ESDR): Prevalence, etiology and potential relationship to arteriosclerotic outcomes. *Kidney Int* 1997;52:10-20.
- 10) Refsum H, Ueland PM. Clinical significance of pharmacological modulation of homocysteine metabolism. *TIPS* 1990;11:411-16.
- 11) Ueland PM, Refsum H, Stabler SP, et al. Total homocysteine in plasma or serum: Methods and clinical applications. *Clin Chem* 1993;39:1764-79.
- 12) Nehler MR, Taylor LM, Porter JM.

Homocysteine as a risk factor for atherosclerosis: A Review. Cardiovasc Pathol 1997;6:1-9.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (μmol/L)

	Mean	Within-Run			Total	
		SD	CV	SD	CV	
1	11.05	1.01	9.1%	1.09	9.9%	
2	19.06	1.37	7.2%	1.82	9.5%	
3	25.33	1.79	7.1%	1.93	7.6%	
4	27.49	1.85	6.7%	3.00	10.9%	
5	35.99	2.25	6.3%	3.48	9.7%	

Specificity

Compound ¹	μmol/L Added ²	% Cross reactivity ³
Adenosine	5,000	ND
S-adenosyl-L-methionine	500	0.6%
Cystathionine	500	5.7%
L-Cysteine	100,000	ND
Glutathione	100,000	ND

ND: Not detectable.⁴

Linearity (μmol/L)


	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8	7.28	—	—
	4 in 8	3.92	3.64	108%
	2 in 8	1.63	1.82	90%
	1 in 8	0.73	0.91	80%

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
2	8 in 8	8.68	—	—
	4 in 8	4.13	4.34	95%
	2 in 8	1.87	2.17	86%
	1 in 8	0.94	1.09	86%
3	8 in 8	28.33	—	—
	4 in 8	12.68	14.17	89%
	2 in 8	7.03	7.08	99%
	1 in 8	2.81	3.54	79%
4	8 in 8	32.39	—	—
	4 in 8	16.88	16.20	104%
	2 in 8	8.17	8.10	101%
	1 in 8	3.94	4.05	97%
5	8 in 8	36.91	—	—
	4 in 8	17.50	18.46	95%
	2 in 8	8.62	9.23	93%
	1 in 8	4.18	4.61	91%

Recovery (μmol/L)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	6.09	—	—
	A	12.65	11.04	115%
	B	17.39	16.14	108%
	C	30.02	26.34	114%
2	—	9.38	—	—
	A	15.75	14.16	111%
	B	20.88	19.26	108%
	C	32.87	29.46	112%
3	—	20.88	—	—
	A	26.89	25.09	107%
	B	30.09	30.19	100%
	C	40.77	40.39	101%
4	—	30.25	—	—
	A	31.51	33.99	93%
	B	40.02	39.09	102%
	C	49.77	44.29	101%
5	—	34.60	—	—
	A	38.32	38.12	101%
	B	42.65	43.22	99%
	C	57.27	53.42	107%

Method Comparison


$$(IML) = 0.97 (IML\ 2000) + 0.86\ \mu\text{mol/L}$$


$$r = 0.975$$

Effect of Storage Temperature on Different Sample Types¹

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
Homocysteine: Homoisteina. **Effect of Storage Temperature on Different Sample Types:**

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2018 Siemens Healthcare Diagnostics. All rights reserved.