# IMMULITE®

## Intact PTH

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**SIEMENS**

# IMMULITE®/IMMULITE 1000 Intact PTH

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE® and IMMULITE 1000 Analyzers — for the quantitative measurement of intact parathyroid hormone (parathyrin, PTH) in EDTA plasma or serum, as an aid in the differential diagnosis of hypercalcemia and hypocalcemia.

Catalog Number: **LKPP1** (100 tests)

Test Code: **IPT** Color: **Dark Green**

## Summary and Explanation

Parathyroid hormone (parathyrin, PTH), a single-chain polypeptide containing 84 amino acids, exerts significant influence in the maintenance of optimal calcium ion concentrations. PTH raises serum ionized calcium levels through direct action on bone and the kidneys: it increases the rate of calcium ion flow from bone to the extracellular fluid, and increases both the renal tubular reabsorption of ionized calcium and the renal excretion of phosphate. Long-term regulation of total body calcium by PTH occurs through its stimulation of vitamin D metabolism, which results in enhanced intestinal absorption of ionized calcium.<sup>1</sup>

In healthy individuals, PTH is secreted in response to circulating calcium ion levels. Any dip below an individual's normal level triggers a pronounced increase in PTH secretion. Calcium levels returning to normal exert a negative feedback effect, thus inhibiting PTH secretion by the parathyroid glands.<sup>1</sup>

PTH undergoes proteolysis to a lesser extent in the parathyroid glands but mostly peripherally — especially in the liver but also in the kidneys and bone — to yield N-terminal fragments and longer lived C-terminal and midregion fragments. The N-terminal fragment contains the region that confers bioactivity. C-terminal and N-terminal fragments are initially generated in equivalent amounts, but the N-terminal fragments disappear rapidly. The C-terminal fragment has a half-life of several hours. In renal failure, C-terminal

fragment clearance by glomerular filtration is impaired, so that high levels are found. C-terminal assays (as well as midregion assays) are consequently likely to be especially unreliable in chronic renal failure, where increased PTH is typically just a reflection of impaired renal clearance.<sup>1,2</sup>

For the intact hormone, the *in vivo* half-life is 2 to 5 minutes.<sup>3</sup> Intact PTH clearance is accomplished by both peritubular uptake and glomerular filtration followed by reabsorption. In normal renal function, intact PTH is the greatest part of circulating PTH-like bioactivity<sup>4</sup> and is present in the circulation at concentrations of  $10^{-11}$  to  $10^{-12}$  mol/L.<sup>2</sup>

In hypercalcemia due to primary hyperparathyroidism or to ectopic PTH production (pseudohyperparathyroidism), the majority of patients have elevated PTH levels. By contrast, in hypercalcemia due to malignancy or other causes, the concentration of PTH in circulation is typically low or within normal reference range limits. PTH levels are also characteristically high in secondary hyperparathyroidism — usually associated with renal failure — as a result of constant stimulation of the parathyroid gland by low calcium levels. Hypocalcemia accompanied by a low PTH level, on the other hand, is to be expected in hypoparathyroidism, either postsurgical or idiopathic.<sup>2,5,6</sup>

Immunoassays specific for various PTH fragments have been developed. Most rely on antisera specific for a discrete region: the C-terminal, N-terminal, or midmolecule. The antisera employed in such assays recognize not only the specific region, but similar fragments as well.<sup>1,4</sup>

Recent assays for intact PTH have the necessary sensitivity for detecting circulating intact PTH in normals and for discriminating between normals and those with primary hyperparathyroidism. These assays also seem to discriminate better between primary hyperparathyroidism and hypercalcemia of malignancy compared with previous assays, and do so virtually without any significant overlap between these groups.<sup>4</sup>

Much improved clinical sensitivity is reported for PTH assays when dynamic reference intervals (based on a range of serum PTH values obtained by acute modification of serum calcium concentrations in healthy subjects) are used, rather than a gaussian reference range (based on PTH values seen in normocalcemic individuals). Using an immunoradiometric assay for intact PTH and applying a dynamic reference range, Lepage, *et al.* obtained average clinical sensitivity of up to 100 percent with primary hyper- and hypoparathyroid samples. Moreover, only the intact PTH assay allowed complete separation between primary hyperparathyroid and nonparathyroidal hypercalcemic patients.<sup>7</sup>

## Principle of the Procedure

IMMULITE/IMMULITE 1000 Intact PTH is a solid-phase, two-site chemiluminescent enzyme-labeled immunometric assay.

**Incubation Cycles:** 1 × 60 minutes

## Specimen Collection

Because of the nocturnal rise in intact PTH levels, samples should be collected in the morning, after 7 a.m., preferably after an overnight fast. (A study by Logue, *et al.* suggests that sampling after 10 a.m. may optimize discrimination between normal subjects and patients with mild primary hyperparathyroidism.<sup>8</sup>)

Collect blood by venipuncture into plastic EDTA plasma tubes or plain plastic serum tubes (with or without gel barrier), avoiding hemolysis. Both plasma and serum should be separated from the cells as soon as possible.

**For EDTA plasma**, keep specimens cold (2–8°C) throughout the collection and separation process. Separate the plasma from the cells, using a refrigerated centrifuge, if possible. Note that EDTA collection tubes must be filled to their capacity. Failure to completely fill the tube will result in excess concentration of EDTA which will interfere with the assay, causing a false depression of values.

**For serum**, allow the specimens to clot at room temperature. Separate the serum from the cells, using a refrigerated centrifuge, if possible.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Lipemic, hemolyzed, icteric or grossly contaminated samples may give erroneous results.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Intact PTH has been tested with plastic BD Vacutainer™ tubes (plain serum, SST, and EDTA). It has not been tested with all possible variations of tube types.

**Volume Required:** 50 µL serum or plasma. (Sample cup must contain at least 100 µL more than the total volume required.)

**Storage:** Samples can be stored at 2–8°C for up to 8 hours after collection. For longer storage, aliquot and freeze up to 2 months at –20°C.

## Warnings and Precautions

*For in vitro* diagnostic use.

### CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.<sup>21-23</sup>

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**H311, H302,  
H411**

**P280, P273,  
P301 + P312,  
P302 + P312,  
P501**

**Danger!** Toxic in contact with skin. Harmful if swallowed. Toxic to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.

**Contains:** sodium azide;  
Intact PTH Adjustors

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws. Bags of Test Units should be allowed to come to room temperature (15–28°C) before opening to avoid condensation.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### Intact PTH Test Units (LPP1)

Each barcode-labeled unit contains one bead coated with affinity-purified murine monoclonal anti-PTH (44–84) antibody.

Stable at 2–8°C until expiration date.

**LKPP1:** 100 units

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

### Intact PTH Reagent Wedge (LPP2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to affinity-purified goat polyclonal anti-PTH (1–34) antibody in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

**LKPP1:** 1 wedge

### Intact PTH Adjustors (LPHL, LPHH)

Two vials (Low and High) of lyophilized, synthetic human intact PTH in a buffered matrix. Reconstitute each vial with **2.0 mL** distilled or deionized water, then place Adjustors on ice immediately. Mix by gentle, intermittent swirling. Use *only freshly reconstituted* Adjustors for each assay. *Do not freeze.*

**LKPP1:** 4 sets

## Kit Components Supplied Separately

### Intact PTH Sample Diluent (LPHZ)

For the manual dilution of samples. One vial containing 25 mL of a PTH-free buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash Module

**LKPM:** Probe Cleaning Kit

**LCHx-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)

**LSCC:** Sample Cup Caps (optional)

**LPHCM:** Bi-level Intact PTH control, in a buffered matrix

Also Required:

Sample transfer pipets; distilled or deionized water; controls; 2.0 mL volumetric pipet ( $\pm 0.02$  mL) for reconstitution of Adjustors; ice bath for keeping samples and Adjustors cold before assay; controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for preparation, setup, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

### Recommended Adjustment Interval:

4 weeks

**Quality Control Samples:** Use controls or patient sample pools with at least two levels (low and high) of intact PTH.

## Expected Values

A reference range study was conducted on matched EDTA plasma and serum samples from 88 apparently healthy volunteers, collected in BD plastic Vacutainer™ tubes. The matched samples were analyzed with the IMMULITE/IMMULITE 1000 Intact PTH and with the IMMULITE 2000 Intact PTH assays. Analysis of the data indicated no statistically significant difference between platforms, although there was a clear difference in reference ranges between EDTA plasma and serum. The reference ranges suggested by this study for IMMULITE/IMMULITE 1000 Intact PTH and IMMULITE 2000 Intact PTH, for both EDTA plasma and serum, are shown in the following table.

	Median	Central 95% Range (pg/mL)
EDTA Plasma	38	16 – 87
Serum	30	11 – 67

The reference range for serum is in good agreement with a previous multi-site reference range study, conducted with IMMULITE/IMMULITE 1000 Intact PTH on 194 serum samples from apparently healthy subjects. This population yielded a median of 31 pg/mL. The reference range suggested by this study is:

10 – 69 pg/mL (1.1 – 7.3 pmol/L).

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

## Limitations

The IMMULITE/IMMULITE 1000 Intact PTH kit is intended strictly as an aid in the differential diagnosis of hypercalcemia and hypocalcemia. It is *not* intended for the diagnosis or management of malignancy.

Because of the physiological relationship between circulating calcium and PTH, it is always important to interpret PTH results in the light of total or ionized calcium levels.<sup>9,10</sup> The finding of a persistently high-normal calcium accompanied by a high-normal PTH (alternatively, a low-normal calcium accompanied by a low-normal PTH) warrants further investigation; for the PTH level, though itself within normal limits, may still be inappropriately high (or inappropriately low) relative to the circulating calcium level.<sup>11</sup>

Indices of renal function — e.g. creatinine levels; measurement of albumin, as an adjunct to measurement of *total* calcium levels;<sup>12–15</sup> and determinations of phosphorus,<sup>14,16</sup> chloride,<sup>12</sup> nephrogenous cyclic AMP<sup>14,17,18</sup> and possibly calcitonin<sup>14</sup> — may also (in certain circumstances) aid in the interpretation of PTH and calcium results. It should also be remembered that hypercalcemia and hypocalcemia may be secondary to disordered vitamin D metabolism.<sup>19</sup>

Lipemic, hemolyzed, icteric or grossly contaminated samples may give erroneous results.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this

assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL.

Unless otherwise specified, all performance data were generated with EDTA plasma samples.

### Conversion Factor:

pg/mL × 0.1053 → pmol/L

**Calibration Range:** Up to 2,500 pg/mL (263 pmol/L)

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

**Analytical Sensitivity:** 3 pg/mL (0.3 pmol/L)

**High-dose Hook Effect:** None up to 500,000 pg/mL (50,000 pmol/L)

**Intraassay Precision (Within-Run):** Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

**Interassay Precision (Run-to-Run):** Statistics were calculated for samples assayed in 10 different runs. (See "Interassay Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three intact PTH solutions (900, 3,600 and 7,200 pg/mL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The IMMULITE Intact PTH assay is highly specific for intact PTH, with particularly low crossreactivity to most PTH fragments, as well as to other naturally occurring compounds which may be present in patient samples.

The antibodies used in the IMMULITE Intact PTH kit were purified by affinity chromatography to achieve specificity for well-defined regions of the intact PTH molecule. The antibodies immobilized to the solid phase (coated bead) are specific for the C-terminal region (44-84) and have no detectable crossreactivity to the N-terminal region

(1-34). Conversely, the enzyme-labeled antibody recognizes only the N-terminal region, and has no detectable crossreactivity to C-terminal or midmolecule regions. Accordingly, the assay, which requires binding by both enzyme-labeled and solid-phase antibodies, is able to recognize only intact PTH and very large PTH fragments that are nearly as long as intact PTH itself. One such fragment, PTH 7–84, exhibits significant crossreactivity (48.3%)<sup>20</sup> in the IMMULITE PTH assay.

**Bilirubin:** Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L may cause a depression of values. (See "Bilirubin" table.)

**Biotin:** Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

**Hemolysis:** Presence of hemoglobin in concentrations up to 513 mg/dL may cause a depression of values. (See "Hemolysis" table.)

**Lipemia:** Presence of triglycerides in concentrations up to 3,000 mg/dL may cause a depression of values. (See "Lipemia" table.)

**Method Comparison:** The assay was compared to IMMULITE 2000 Intact PTH on 140 serum samples. (Concentration range: approximately 5 to 1,000 pg/mL. See graph.) By linear regression:

(IML) = 1.08 (IML 2000) – 1.21 pg/mL  
r = 0.997

Means:  
145 pg/mL (IMMULITE 2000)  
156 pg/mL (IMMULITE)

## References

- 1) Armitage EK. Parathyrin (parathyroid hormone): metabolism and methods for assay. Clin Chem 1986;32:418-24. 2) Blind E, Schmidt-Gayk H, Scharla S, Flentje D, Fischer S, Goehring U, Hitzler W. Two-site assay of intact parathyroid hormone in the investigation of primary hyperparathyroidism and other disorders of calcium metabolism compared with a midregion assay. J Clin Endocrinol Metab 1988;67:353-60. 3) Kao PC, Grant CS, Klee GG, Khosla S. Clinical performance of parathyroid hormone immunometric assays. Mayo Clin Proc 1992;67:637-45. 4) Measuring the PTH level. Lancet 1988 (Jan 16);94-5. 5) Christensen MS.

Radioimmunoassay of human parathyroid hormone. *Dan Med Bull* 1979;26:157-73.

6) Hawker CD, Di Bella FP. Radioimmunoassays for intact and carboxyl-terminal parathyroid hormone: clinical interpretation and diagnostic significance. *Ann Clin Lab Sci* 1980;10:76-87.

7) Lepage R, Whittom S, Bertrand S, Bahsali G, D'Amour P. Superiority of dynamic over static reference intervals for intact, midmolecule, and C-terminal parathyrin in evaluating calcemic disorders. *Clin Chem* 1992;38:2129-35.

8) Logue FC, Fraser WD, O'Reilly DSJ, Beasall GH. The circadian rhythm of intact parathyroid hormone (1-84) and nephrogenous cyclic adenosine monophosphate in normal men. *J Endocrinol* 1989;121:R1-R3.

9) Bowers GN, Brassard C, Sena SF. Measurement of ionized calcium in serum with ion-selective electrodes. *Clin Chem* 1986;32:1437-47.

10) Halvorsen JF, Gautvik VT, Gautvik KM. Improved diagnosis of primary hyperparathyroidism by defining the inverse relationship between serum immunoreactive parathyroid hormone and calcium. *Scand J Clin Lab Invest* 1986;46:435-42.

11) Hawker CD, Clark SW, Martin KJ, Siatopolsky E, Di Bella FP. Radioimmunoassay of parathyroid hormone: clinical utility and interpretation. In: Thompson NW, Vinik AI, editors. *Endocrine Surgery Update*. New York: Grune & Stratton, 1983: 321-40.

12) Boyd JC, Ladenson JH. Value of laboratory tests in the differential diagnosis of hypercalcemia. *Am J Med* 1984;77:863-72.

13) Groth TL, Ljunghall S, De Verdier C-H. Optimal screening for patients with hyperparathyroidism with use of serum calcium observation; a decision-theoretical analysis. *Scand J Clin Lab Invest* 1983;43:699-707.

14) Kleeman K, Kleeman CR. Parathyroid hormone and calcitonin. *Contemp Endocrinol* 1985;2:247-344.

15) Orwoll ES, Meier D. Alterations in calcium, vitamin D, and parathyroid hormone physiology in normal men with aging: relationship to the development of senile osteopenia. *J Clin Endocrinol Metab* 1986;63:1262-9.

16) Kao PC. Parathyroid hormone assay. *Mayo Clin Proc* 1982;57:596-7.

17) Heath H. Tests of parathyroid function: utility and limitations. *American Association for Clinical Chemistry: Continuing Education in Endocrinology and Metabolism*. *Endo* (Feb) 1984;2(8):10 pages.

18) Johnson DR, Howarth AT. Differential diagnosis of hypercalcemia. *CRC Crit Rev Clin Lab Sci* 1984;21:51-97.

19) Audran M, Kumar R. The physiology and pathophysiology of vitamin D. *Mayo Clin Proc* 1985;60:851-66.

20) UK National External Quality Assurance Scheme (NEQAS) for PTH, ACTH, and Calcitonin. *Distribution* 65, 2001.

21) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377-82, 387-8.

22) Clinical and Laboratory Standards Institute (formerly NCCLS). *Protection of Laboratory Workers From*

Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 23) Federal Occupational Safety and Health Administration, *Bloodborne Pathogens Standard*, 29 CFR 1910.1030.

## Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

## Tables and Graphs

### Intraassay Precision (pg/mL)

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	57	3.6	6.3%
2	236	13	5.5%
3	662	40	6.3%

### Interassay Precision (pg/mL)

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	59	5.1	8.6%
2	433	34	7.9%

### Specificity

Compound <sup>1</sup>	pg/mL Added <sup>2</sup>	% Cross reactivity <sup>3</sup>
PTH (1-34)	300	ND
PTH (1-44)	100,000	ND
PTH (44-68)	100,000	ND
PTH (53-84)	100,000	ND
Calcitonin	10,000	ND

ND: Not detectable<sup>4</sup>

## Linearity (pg/mL)


	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8 <sup>5</sup>	175	—	—
	4 in 8	85	88	97%
	2 in 8	43	44	98%
	1 in 8	22	22	100%
2	8 in 8	208	—	—
	4 in 8	98	104	94%
	2 in 8	53	52	102%
	1 in 8	26	26	100%
3	8 in 8	811	—	—
	4 in 8	410	406	101%
	2 in 8	225	203	111%
	1 in 8	111	101	110%

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2008 Siemens Healthcare Diagnostics. All rights reserved.
#### **CONTROL AB**