# IMMULITE®

## Total T3

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**SIEMENS**

# IMMULITE®/IMMULITE 1000 Total T3

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE® and IMMULITE 1000 Analyzers — for the quantitative measurement of total circulating triiodothyronin (T3) in serum as an aid in the clinical assessment of thyroid status.

Catalog Number: **LKT31** (100 tests)

Test Code: **T3** Color: **Violet**

### Summary and Explanation

Under normal physiological conditions, T3 represents approximately 5 percent of the thyroid hormones in serum. Although present in lower concentration, T3 has a greater intrinsic metabolic activity, faster turnover and larger volume of distribution than circulating T4.<sup>1</sup> Reports that thyrotoxicosis may be caused by abnormally high concentrations of T3 — rather than T4 — have reinforced the importance of T3 measurements.<sup>2-5</sup> In addition, T3 determination is an important tool for monitoring hypothyroid patients receiving sodium liothyronine therapy. Unlike "T3 Uptake" tests, which estimate the saturation of thyroid hormone binding proteins, total T3 analysis measures circulating levels of triiodothyronine. Most reports indicate that T3 levels distinguish clearly between euthyroid and hyperthyroid subjects, but provide a less clear-cut separation between hypothyroid and euthyroid subjects.<sup>1</sup>

Numerous conditions unrelated to thyroid disease may cause abnormal T3 values.<sup>1,6-9</sup> Consequently, total T3 values should not be used on their own in establishing the thyroid status of an individual. The levels of serum T4, TBG (thyroid binding globulin), TSH (thyroid stimulating hormone) and other clinical findings must be considered as well.

### Principle of the Procedure

IMMULITE/IMMULITE 1000 Total T3 is a solid-phase, competitive chemiluminescent enzyme immunoassay.

**Incubation Cycles:** 1 × 30 minutes

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed or grossly contaminated samples may give erroneous results.

EDTA plasma is not recommended for use.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Total T3 has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 25 µL serum (Sample cup must contain at least 100 µL more than the total volume required.)

**Storage:** 7 days at 2–8°C or 2 months at –20°C.

### Warnings and Precautions

For *in vitro* diagnostic use.

#### **CAUTION! POTENTIAL BIOHAZARD**

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should

---

---

be handled using good laboratory practices and universal precautions.<sup>11-13</sup>

---

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### Total T3 Test Units (LT31)

Each barcode-labeled unit contains one bead coated with monoclonal murine anti-T3. Stable at 2–8°C until expiration date.

**LKT31:** 100 units

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

### Total T3 Reagent Wedge (LT32)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to T3, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

**LKT31:** 1 wedge

### Total T3 Adjustors (LT3L, LT3H)

Two vials (Low and High), 2.0 mL each, of T3 in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**LKT31:** 1 set

## Kit Components Supplied Separately

### Total T3 Sample Diluent (LT3Z)

For the manual dilution of samples. One vial containing 25 mL of processed, T3-free human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash Module

**LKPM:** Probe Cleaning Kit

**LCHx-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)

**LSCC:** Sample Cup Caps (optional)

Also Required

Sample transfer pipets, distilled or deionized water, controls

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

**Recommended Adjustment Interval:**  
2 weeks

**Quality Control Samples:** Follow government regulations or accreditation requirements for quality control frequency.

Use controls or sample pools with at least two levels (low and high) of T3.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system,

or within an established range determined by an appropriate internal laboratory quality control scheme.

## Expected Values

139 serum samples from apparently healthy subjects were processed by the IMMULITE Total T3 assay.

Non-parametric analysis of the results yielded a median value of 126 ng/dL and a central 95% range of

81–178 ng/dL (1.2–2.7 nmol/L)

A cross-sectional study of pediatric values performed with IMMULITE Total T3 at a "wellness" clinic in the southwestern United States yielded the following results.

Age (years)	Median (ng/dL)	Central 95% Range
1.0	139	33 – 256
2.0	141	36 – 254
3.0	141	42 – 248
4.0	137	44 – 236
5.0	134	44 – 231
6.0	131	44 – 225
7.0	128	43 – 220
8.0	126	42 – 216
9.0	124	42 – 212
10.0	122	42 – 209
11.0	120	41 – 205
12.0	118	41 – 202

Reference ranges were generated using serum samples collected in Becton Dickinson plastic vacutainer tubes without gel barriers.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

## Limitations

Many physiological, pharmacological, pathological, and genetic factors affect the interpretation of total T3 results. (Refer to the Summary and Explanation section.)

More than 20% of patients with non-thyroidal disease who are critically ill have low circulating total T3 levels.

EDTA plasma has an effect on the measurement of T3 in the procedure.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/dL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

### Conversion Factor:

ng/dL  $\times$  0.01536  $\rightarrow$  nmol/L

**Calibration Range:** 40–600 ng/dL (0.61–9.2 nmol/L)

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

**Analytical Sensitivity:** 35 ng/dL (0.54 nmol/L)

**Precision:** Run-to-run (interassay) statistics were determined for a representative spectrum of samples assayed in at least 50 runs. Also tabulated are the within-run statistics expected at these concentrations for samples assayed in replicates of 10 or 20 in a single run. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three T3 solutions (1395, 2576 and 3318 ng/dL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for T3. (See "Specificity" table.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Biotin:** Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

**Hemolysis:** Gross hemolysis has some effect on the IMMULITE Total T3 assay, causing approximately 10–20% depression of values. (See "Hemolysis" table.)

**Lipemia:** Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** To assess the effect of alternate sample types, blood was collected from 20 volunteers into plain, heparinized, and Becton Dickinson SST<sup>®</sup> vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of T3, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE Total T3 procedure.

(Heparin) = 1.05 (Serum) – 5.1 ng/dL  
r = 0.998

(SST) = 0.99 (Plain Tubes) + 1.4 ng/dL  
r = 0.994

Means:

189 ng/dL (Heparin)  
183 ng/dL (SST)  
184 ng/dL (Serum)

In another experiment, blood was collected from 65 laboratory volunteers into plain and EDTA vacutainer tubes.

(EDTA) = 1.15 (Serum) + 21.7 ng/dL  
r = 0.837

Means:

168 ng/dL (EDTA)  
127 ng/dL (Serum)

**Method Comparison:** The assay was compared to Coat-A-Count<sup>®</sup> Total T3 on 118 samples. (Concentration range: approximately 56 to 580 ng/dL. See graph.) By linear regression:

(IML) = 0.96 (CAC) – 0.12 ng/dL  
r = 0.979

Means:

154 ng/dL (IMMULITE)  
162 ng/dL (Coat-A-Count)

## References

- 1) Hollander CS, Shenkman L. Radioimmunoassays for triiodothyronine and thyroxine. In: Rothfeld B, editor. Nuclear medicine *in vitro*. Philadelphia: Lippincott, 1974, 136-49.
- 2) Ahmed M, Doe RP, Nuttall FQ. Triiodothyronine thyrotoxicosis following iodide ingestion: a case report. *J Clin Endocrinol Metab* 1974;38:574-6.
- 3) Hollander CS, Nihei N, Burday SZ, Mitsuma T, Shenkman L, Blum M. Clinical and laboratory observations in cases of triiodothyronine toxicosis confirmed by radioimmunoassay. *Lancet* 1972;i:609-11.
- 4) Mitsuma T, Owens R, Shenkman L, Reiter E, Hollander CS. T3 toxicosis in childhood: hyperthyroidism due to isolated hypersecretion of triiodothyronine. *J Pediatr* 1972;81:982-4.
- 5) Sterling K, Refetoff S, Selenkow HA. Thyrotoxicosis due to elevated serum triiodothyronine levels. *JAMA* 1970;213:571-5.
- 6) Bates HM. *Clin Lab Prod* 1974;3:16.
- 7) Utiger RD. Serum triiodothyronine in man. *Annu Rev Med* 1974;25:289-302.
- 8) Larson PR. Triiodothyronine: review of recent studies of its physiology and pathophysiology in man. *Metabolism* 1972;21:1073-92.
- 9) Oppenheimer JH. Role of plasma proteins in the binding, distribution and metabolism of the thyroid hormones. *N Engl J Med* 1968;278:1153-62.
- 10) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4. Wayne, PA: NCCLS, 1998.
- 11) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377–82, 387–8.
- 12) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.
- 13) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

## Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor. [siemens.com/healthcare](http://siemens.com/healthcare)

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

# Tables and Graphs

## Precision (ng/dL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Run-to-Run <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	62	8.2	13.2%	9.7	15.6%
2	90	9.3	10.3%	12	13.3%
3	159	12	7.5%	16	10.1%
4	186	13	7.0%	18	9.7%
5	296	17	5.7%	21	7.1%
6	392	21	5.4%	30	7.7%

## Linearity (ng/dL)

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8	313	—	—
	4 in 8	156	157	99%
	2 in 8	76	78	97%
	1 in 8	40	39	103%
2	8 in 8	324	—	—
	4 in 8	153	162	94%
	2 in 8	90	81	111%
	1 in 8	44	41	107%
	3	8 in 8	500	—
4 in 8		257	250	103%
2 in 8		123	125	98%
1 in 8		65	63	103%

## Recovery (ng/dL)

	Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	142	—	—
	A	204	205	100%
	B	260	264	98%
	C	321	301	107%
2	—	245	—	—
	A	314	303	104%
	B	330	362	91%
	C	378	399	95%
3	—	391	—	—
	A	443	441	100%
	B	486	500	97%
	C	554	537	103%

## Specificity


Compound <sup>1</sup>	ng/dL Added <sup>2</sup>	% Cross reactivity <sup>3</sup>
Triiodo-L-thyronine (T3)	—	100%
Reverse T3	500,000	0.007%
Triiodo-D-thyronine	134	100%
L-Thyroxine (T4)	24,000	0.24%
D-Thyroxine	5000	1.0%
Triiodothyroacetic acid	10,000	0.7%
	1000	0.6%
Tetraiodothyroacetic acid	20,000	1.3%
Monoiodotyrosine	1,000,000	ND
Diodo-L-tyrosine	1,000,000	0.005%
Methimazole	1,000,000	0.006%
5,5'-Diphenylhydantoin	1,000,000	ND
Phenylbutazone	1,000,000	0.005%
6- <i>n</i> -Propyl-2-thiouracil	1,000,000	0.006%

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IMMULITE and Coat-A-Count are trademarks of Siemens Healthcare Diagnostics.
All rights reserved.
#### **CONTROL AB**