Digitoxin

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE[®] 2000 Digitoxin

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Systems Analyzers — for the quantitative measurement of digitoxin in serum, as an aid in the therapeutic administration of this cardioglycoside.

Catalog Numbers: **L2KDG2** (200 tests)

Test Code: **DGT** Color: **Red**

Summary and Explanation

Digitoxin is a long-acting cardiac glycoside used in the treatment of congestive heart failure and certain arrhythmias. It is structurally identical to digoxin except for the absence of a single hydroxyl group. Being less polar, digitoxin is absorbed more completely, and circulates to a much greater extent protein-bound. In the liver, it is extensively metabolized to both active and inactive forms. Digitoxin has a comparatively long half-life of 5 to 7 days and, except in the nephrotic syndrome, is largely unaffected by renal dysfunction.

A transient peak is associated with the distribution phase; thereafter the circulating level exhibits a steady and very gradual decline. Accordingly, blood samples for digitoxin analysis should be collected just before the next dose — or in any case not less than 6 hours after the previous dose — to allow time for plasma and tissue levels to reach equilibrium.

Serum determinations can assist in verifying compliance and in adjusting the dose to achieve a desired plasma concentration. Monitoring circulating digitoxin levels is especially valuable when there are enterohepatic or renal complications, and during concomitant therapy with drugs like quinidine, which can affect the pharmacokinetics of digitoxin.

It must be remembered, however, that the relation between circulating level and clinical response varies greatly from one individual to another, due to the many factors which are at work. This is reflected in the considerable overlap between

therapeutic and toxic ranges summarized in the literature.

Thus, circulating levels of approximately 10 to 25 ng/mL are observed in the majority of patients exhibiting a satisfactory therapeutic response to digitoxin, while levels greater than 45 ng/mL are typical of patients with symptoms of digitoxin toxicity.^{2,4,16} But concentrations higher than the conventional therapeutic range may be required to control certain arrhythmias. Conversely, circulating levels well below 45 ng/mL may give rise to symptoms of toxicity in patients rendered more sensitive to the drug's action by, for example, a lowered potassium level. Even so, the higher the circulating level, the more probable it is that clinical symptoms of toxicity can be attributed to the drug — and the lower the level, the less likely.

Serum determinations should supplement clinical judgment, not substitute for it. An isolated value should not be used as the sole criterion for the assessment of efficacy or toxicity; rather it should always be interpreted in the total clinical context.

Principle of the Procedure

IMMULITE 2000 Digitoxin is a solid-phase, competitive chemiluminescent enzyme immunoassay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

Collect blood into plain (without anticoagulant) or heparinized tubes, and separate the serum from the cells. The time of collection should be noted.

To avoid misleadingly high digitoxin results, patient samples should be drawn just before the next dose — or in any case no less than six hours after the administration of this cardiac glycoside.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Digitoxin has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 10 µL serum.

Storage: 7 days at 2–8°C or 2 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.¹⁷⁻¹⁹

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were

tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Digitoxin Bead Pack (L2DG12)

With barcode. 200 beads, coated with polyclonal rabbit anti-digitoxin antibody. Stable at 2–8°C until expiration date. **L2KDG2:** 1 pack.

Digitoxin Reagent Wedge (L2DGA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to digitoxin in buffer. Stable at 2–8°C until expiration date. **L2KDG2:** 1 wedge.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Digitoxin Adjustors (LDGL, LDGH)

Two vials (Low and High), 4 mL each, of digitoxin in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C. **L2KDG2:** 1 set.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate
L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit
LRXT: Reaction Tubes (disposable)

SDCM: Bi-level, multi-constituent control

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of digitoxin.

Therapeutic Range

The following results are reported in the literature for therapeutic and toxic circulating levels of digitoxin.¹⁶

	ng/mL	nmol/L
Therapeutic Range	10–25	13–33
Toxic Level	>45	>59

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Many physiological, pharmacological, pathological and genetic factors may affect the interpretation of results (see the Summary and Explanation section).

In cases of overdosage, specimens should be obtained *before* administering digoxin immune Fab antidote preparations (e.g. Digibind®) since these can be expected to interfere with any digitoxin immunoassay procedure. Once such a preparation has been administered, accurate determinations of digitoxin can be made only by specialized techniques (not routinely available in most clinical settings) until the Fab fragments have been

eliminated from the body: this takes several days for patients with intact renal function, and a week or longer for patients with renal impairment. In addition to *in vitro* effects, these preparations, administered intravenously, cause a rapid rise in pharmacologically inactive, Fab-bound, total circulating digitoxin that may be detected immunologically.¹⁵

The presence of Uzara® in the serum of patients can result in falsely elevated digitoxin results.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

ng/mL × 1.307 → nmol/L

Calibration Range: 5 – 60 ng/mL
(6.5 – 78 nmol/L).

Analytical Sensitivity: 1 ng/mL
(1.3 nmol/L)

Precision: Samples were assayed in duplicate over the course of multiple days, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with 3 digitoxin solutions (193, 385 and 770 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for digitoxin. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of hemoglobin in concentrations up to 600 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 20 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer plastic tubes. Equal volumes of the matched samples were spiked with various concentrations of digitoxin, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 Digitoxin procedure.

(Heparin) = 0.90 (Serum) + 0.26 ng/mL
 $r = 0.987$

(EDTA) = 0.89 (Serum) + 0.39 ng/mL
 $r = 0.973$

(SST) = 1.07 (Plain Tubes) – 0.07 ng/mL
 $r = 0.975$

Means:

7.96 ng/mL (Serum)
7.46 ng/mL (Heparin)
7.50 ng/mL (EDTA)
8.48 ng/mL (SST)

Method Comparison 1: The assay was compared to IMMULITE Digitoxin on 86 patient samples. (Concentration range: approximately 7 to 57 ng/mL. See graph.) By linear regression:

(IML 2000) = 0.97 (IML) – 0.4 ng/mL
 $r = 0.98$

Means:

24.0 ng/mL (IMMULITE 2000)
25.1 ng/mL (IMMULITE)

Method Comparison 2: The assay was also compared to Coat-A-Count Digitoxin radioimmunoassay on 89 patient samples. (Concentration range: approximately 5 to 55 ng/mL. See graph.) By linear regression:

(IML 2000) = 1.04 (Coat-A-Count) – 1.1 ng/mL
 $r = 0.98$

Means:

24.2 ng/mL (IMMULITE 2000)
24.1 ng/mL (Coat-A-Count)

References

- 1) Aronson JK. Clinical pharmacokinetics of cardiac glycosides in patients with renal dysfunction. *Clin Pharmacokinet* 1983;8:155-78.
- 2) Baer DM. Key points for interpretation of drug concentrations: digoxin and digitoxin. In: Baer DM, Dito WR, editors. Interpretations in therapeutic drug monitoring. Chicago: American Society of Clinical Pathologists, 1982: 351-4.
- 3) Bresnahan JF, Vlietstra RE. Digitalis glycosides. *Mayo Clin Proc* 1979;54:675-84.
- 4) Dito WR. Cardiac glycosides and therapeutic drug monitoring. In: Baer DM, Dito WR, editors. Interpretations in therapeutic drug monitoring. Chicago: American Society of Clinical Pathologists, 1981: 61-73.
- 5) Fenster PE, et al. Digitoxin-quinidine interaction: pharmacokinetic evaluation. *Ann Int Med* 1980;93:698-701.
- 6) Garty M. Digitoxin elimination reduced during quinidine therapy. *Ann Int Med* 1981;94:35-7.
- 7) Kuemmel HC, et al. Studies to determine individual toxic doses of digitoxin. *Munchener Med Wochenschr* 1982;124:545-9.
- 8) Kuhlmann J. Sources of digitalis intoxication. *Med Klinik* 1982;77:39-44.
- 9) Larsen A, Storstein L. Digitoxin kinetics and renal excretion in children. *Clin Pharmacol Ther* 1983;33:717-27.
- 10) Maier C, Kolenda KD. Digitoxin in the presence of hepatorenal failure. *Deutsche Med Wochenschr* 1983;108:1475-80.
- 11) Rietbrock N, Alken RG. The therapy of heart failure with digitalis. *Deutsche Med Wochenschr* 1980;105:1622-8.
- 12) Smith TW. Digitalis in the management of heart failure. *Hosp Pract* 1984 Mar;19(3):67-92.
- 13) Smith TW, et al. Digitalis glycosides: mechanisms and manifestations of toxicity. *Prog Cardiovasc Dis* 1984;26:413-58 and 495-540.
- 14) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.
- 15) Digoxin immune Fab. In: McEvoy GK, editor. *AHFS Drug Information* 92. American Society of Hospital Pharmacists, 1992: 1977-81.
- 16) National Academy of Clinical Biochemistry. Guidelines for therapeutic drug monitoring services. Printed by Durik Advertising, Inc., 1999: 82-83.
- 17) Centers

for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377–82, 387–8.

18) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.

19) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	8.9	0.73	8.2%	1.04	11.7%
2	16.4	1.11	6.7%	1.28	7.8%
3	24.1	1.54	6.4%	2.11	8.8%
4	31.9	1.73	5.4%	3.18	10.0%
5	44.1	2.68	6.1%	4.14	9.4%

Linearity (ng/mL)



	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	2 in 2	18.0	—	—
	1 in 2	10.3	9.0	114%
2	2 in 2	26.3	—	—
	1 in 2	12.1	13.2	92%
3	4 in 4	33.1	—	—
	2 in 4	16.7	16.6	101%
	1 in 4	8.78	8.28	106%
4	4 in 4	39.3	—	—
	2 in 4	19.0	19.7	96%
	1 in 4	10.4	9.83	106%
5	8 in 8	49.7	—	—
	4 in 8	27.0	24.9	108%
	2 in 8	13.2	12.4	106%
	1 in 8	5.79	6.21	93%

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IMMULITE® and Coat-A-Count® are trademarks of Siemens Healthcare Diagnostics.
© 2008 Siemens Healthcare Diagnostics. All rights reserved.
CONTROL AB