 IMMULITE[®] 2000

Thyroid Stimulating Immunoglobulins (TSI)

For use on IMMULITE[®] 2000 systems

SIEMENS

IMMULITE® 2000 TSI

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative determination of thyroid stimulating autoantibodies specific to thyroid stimulating hormone receptor (TSHR) in human serum or EDTA and lithium heparin plasma. The determination of thyroid stimulating autoantibodies, in conjunction with other clinical and laboratory findings, is used as an aid in the assessment of patients with suspected Graves' disease.

Catalog Number: **L2KSI2** (200 tests)

SMN: **10713448**

Test Code: **TSI** Color: **Dark Blue**

Summary and Explanation

Toxic diffuse goiter, also known as Graves' disease (GD), is an autoimmune disorder and the most common cause of hyperthyroidism. In GD, thyroid stimulating immunoglobulins (TSI) bind to the TSH receptor (TSHR) and mimic TSH stimulation of the thyroid gland.¹ Because TSI induced thyroid hormone secretion is not controlled by negative feedback, such stimulation leads to GD hyperthyroidism.² The TSH receptor contains a large extracellular domain that presents epitopes for a variety of autoantibodies, including TSI and Thyroid Blocking Immunoglobulins (TBI).³ In contrast to TSI, TBI bind to the TSH receptor and inhibit TSH stimulation of thyroid cells, leading to hypothyroidism.^{4,5} TSHR autoantibody (TRAb) assays do not distinguish between TSI and TBI.⁶⁻⁸ The IMMULITE 2000 TSI assay utilizes recombinant human TSH receptors (hTSHR) for the specific detection of thyroid stimulating autoantibodies.

The clinical utility of TSI measurement includes a determination of the autoimmune etiology of thyrotoxicosis,⁹ monitoring Graves' patient therapy,¹⁰ prediction of remission or relapse,¹¹ confirmation of Graves' ophthalmopathy,¹² and prediction of hyperthyroidism in neonates.¹³

Principle of the Procedure

The IMMULITE 2000 TSI assay is an automated, two-cycle, chemiluminescent immunoassay. It employs a pair of recombinant hTSHR constructs in a bridging immunoassay format. The capture receptor is immobilized on the solid phase (polystyrene bead). The signal receptor is an alkaline phosphatase labeled recombinant hTSHR in a buffer solution.

In the first cycle, sample is incubated with the solid phase for 30 minutes, allowing the TSI in the sample to bind through one arm to the capture receptor. Next, centrifugal washes remove residual sample.

In the second cycle, the signal receptor is added to the reaction tube and incubated for 30 minutes. The complexed TSI bind the signal receptor through the second arm, forming a bridge. Unbound signal receptor is then removed by centrifugal washes.

Finally, chemiluminescent substrate is added to the reaction tube and a signal is generated in direct relation to the amount of TSI in the sample.

Incubation Cycles: 2 × 30 minutes

Time to first result: 65 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 TSI has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Short drawn EDTA plasma samples may result in under-recovery of IMMULITE 2000 TSI results.

Volume Required: 50 µL serum or plasma.

Storage: 24 hours at 20-25°C; 7 days at 2–8°C; or at –20°C for 12 months.

The handling and storage information provided is based on data or references maintained by the manufacturer. It is the responsibility of the individual laboratory to use all available references and/or its own studies when establishing alternate stability criteria to meet specific needs.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.^{14–16}

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Caution: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested by FDA approved methods or equivalent and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

TSI Bead Pack (L2S112)

With barcode. 200 beads, coated with recombinant human TSH receptor. Stable at 2–8°C until expiration date.

L2KS12: 1 pack

TSI Reagent Wedge (L2SIA2)

With barcode. 11.5 mL buffer, with preservative. 11.5 mL alkaline phosphatase labeled recombinant human TSH receptor in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KS12: 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

TSI Adjustors (L2SIL, L2SIH)

Two vials (low and high) containing lyophilized human monoclonal thyroid stimulating antibody in a bovine serum matrix, with preservative. Reconstitute each vial with **2.0 mL** distilled or deionized water. Mix by gentle swirling or inversion on a rotator for at least 30 minutes until the lyophilized material is fully dissolved. Stable at 2–8°C for 90 days after reconstitution, or for 4 months (aliquotted) at –20°C. Avoid repeated freeze/thaw cycles.

L2KS12: 1 set

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

TSI Controls (L2SIC1, L2SIC2)

Two vials each of 2 levels (low, and high) containing lyophilized human monoclonal thyroid stimulating antibody in a human serum matrix, with preservative.

Reconstitute each vial with **4.0 mL** distilled or deionized water. Mix by gentle swirling or inversion on a rotator for at least 30 minutes until the lyophilized material is fully dissolved. (No further dilution is required.) Stable at 2–8°C for 90 days after reconstitution, or for 6 months (aliquotted) at –20°C. Avoid repeated freeze/thaw cycles.

L2KSI2: 1 set

For the current Control ranges in IU/L, please refer to the Control insert.

Aliquot Labels with barcodes are supplied with the kit, for use with the Adjustors and Controls. Before use, place the appropriate Aliquot Labels on test tubes, so the barcodes can be read by the barcode reader on the IMMULITE 2000.

Kit Components Supplied Separately

Multi-Diluent 2 (L2M2Z, L2M2Z4)

For the on-board dilution of patient samples exceeding the analytical measuring range of the assay. One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M2Z: 3 labels **L2M2Z4:** 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

Also Required

Distilled or deionized water; test tubes.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:

4 weeks

Quality Control Samples: Follow government regulations or accreditation requirements for quality control frequency.

Interpretation of Results

Standardization: This assay is traceable to World Health Organization (WHO) 2nd International Standard for Thyroid Stimulating Antibody, National Institute for Biological Standards and Control (NIBSC) Code: 08/204.

Assay Cut-Off: The sera from 164 treated Graves' disease patients and 270 individuals with other thyroid or autoimmune diseases were evaluated to establish a cut-off for hyperthyroid GD patients. Based on a receiver operating characteristic (ROC) curve, an optimal cut-off of 0.55 IU/L was selected.

Expected Values

A total of 842 serum samples from apparently healthy males (n = 151), non-pregnant females (n = 155), first trimester (n = 169), second trimester (n = 191), and third trimester (n = 176) pregnant donors were analyzed using the IMMULITE 2000 TSI assay. The results from this study suggest a nonparametric 97.5th percentile of < 0.10 IU/L.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all

immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in IU/L. (Unless otherwise noted, all were generated on serum samples.)

Detection Limits:

Limit of Blank (LoB) = 0.03 IU/L

Limit of Detection (LoD) = 0.06 IU/L

Limit of Quantitation (LoQ) = 0.10 IU/L

LoB, LoD, and LoQ were determined in accordance with Clinical and Laboratory Standards Institute (CLSI) EP17-A2.¹⁷

Reportable Range: 0.10–40 IU/L

Samples with results > 40 IU/L may be diluted on-board with Multi-Diluent 2 (L2M2Z, L2M2Z4) and reassayed. Because of the heterogeneity of TSHR autoantibodies, certain patient samples may exhibit a non-linear dilution.

Samples with results < 0.10 IU/L will be reported as "< 0.10 IU/L" by the instrument.

High-dose Hook Effect: No effect up to 800 IU/L

Precision: Precision was evaluated according to CLSI guideline EP5-A2.¹⁸ Serum samples were assayed in duplicate over the course of 20 days, 2 runs per day, for a total of 40 runs and 80 replicates (see Precision table).

Cross Reactivity: The cross reactivity of the IMMULITE 2000 TSI assay with FSH, LH, TSH, hCG, anti-TG, and anti-TPO was determined by adding these hormones to TSI-free human serum (see Cross Reactivity table for representative data).

Clinical Sensitivity and Specificity:

Serum samples from 236 treated and untreated Graves' disease patients, 138 individuals with other thyroid or autoimmune diseases (10 Addison's Disease, 20 Sjögren's disease, 9 chronic thyroiditis, 2 multinodular goiter, 12 thyroid cancer, 10 Type 1 diabetes, 10 Systemic Lupus Erythematosus, 10 Rheumatoid Arthritis, 20 Celiac disease, 15 Crohn's disease, 10 Ulcerative Colitis, 10 Pernicious Anemia) and 200 apparently healthy individuals were evaluated. The TSI values for the patients with other thyroid or autoimmune diseases had an upper limit of 0.39 IU/L. At 0.55 IU/L cut-off, the clinical sensitivity and specificity were 98.3% (232/236) and 99.7% (337/338) respectively (see Sensitivity and Specificity table).

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 40 mg/dL has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin at a concentration of 200 mg/dL may affect recovery.

Lipemia: Presence of triglycerides in concentrations up to 1000 mg/dL has no effect on results, within the precision of the assay.

Rheumatoid Factor: Samples containing up to 200 IU/mL of Rheumatoid Factors cause no interference.

Alternate Sample Type: To assess the effect of alternate sample types, 40 matched sets of samples from patients with Graves' disease were drawn into lithium heparin, EDTA, and BD SST Vacutainer® tubes. All samples were analyzed by the IMMULITE 2000 TSI assay, with the following results:

(Li Heparin) = 0.98 (Serum) – 0.00 IU/L

r = 0.995

(EDTA) = 1.04 (Serum) – 0.01 IU/L

r = 0.995

(SST) = 0.99 (Serum) – 0.01 IU/L

r = 0.991

This study showed 100% positive and negative agreements between all samples collected in serum, lithium heparin, EDTA and SST tubes.

Method Comparison: The IMMULITE 2000 TSI assay was compared to the Thyretain™ TSI Reporter BioAssay on 244 serum samples from patients with Graves' and other thyroid or autoantibody diseases (see Method Comparison table).

References

1) B. R. Smith, R. Hall, Thyroid-stimulating immunoglobulins in Graves' disease, *Lancet*, 304:427-430, 1974. 2) B. Rapoport, G.D. Chazenbalk, J. C. Jaume, S. M. McLachlan, The thyrotropin receptor: interaction with thyrotropin and autoantibodies, *Endocr. Rev.*, 19:673-716, 1998. 3) M. Evans, J. Sanders, T. Tagami, P. Sanders, S. Young, E. Roberts, J. Wilmot, X. Hu, K. Kabelis, J. Clark, S. Holl, T. Richards, A. Collyer, J. Furmaniak, B. R. Smith, Monoclonal autoantibodies to the TSH receptor, one with stimulating activity and one with blocking activity, obtained from the same blood sample, *Clin. Endocrin.*, 73:404-412, 2010. 4) J. Orgiazzi, D. E. Williams, I. J. Chopra, D. H. Solomon, Human thyroid adenyl cyclase-stimulating activity in immunoglobulin G of patients with Graves' disease, *J. Clin. Endocrin. Metab.*, 42: 341-354, 1976. 5) K. Endo, K. Kasagi, J. Konishi, K. Ikekubo, O. Tatsuyo, Y. Takeda, T. Mori, K. Torizuka, Detection and properties of TSH binding inhibitor immunoglobulins in patients with Graves' disease and Hashimoto's thyroiditis. *J. Clin. Endocrin. Metab.*, 46:734-739, 1978. 6) C. Massart, R. Sapin, J. Gibassier, A. Agin, M. d'Herbomez, Intermethod variability in TSH-receptor antibody measurement: Implication for the diagnosis of Graves disease and for the follow-up of Graves ophthalmopathy, *Clin. Chem.*, 55:183-186, 2009. 7) K. Kamijo, H. Murayama, T. Uzu, K. Togashi, P. D. Olivo, G. J. Kahaly, Similar clinical performance of a novel chimeric thyroid-stimulating hormone receptor bioassay and an automated thyroid-stimulating hormone receptor binding assay in Graves' disease, *Thyroid*, 21:1-5, 2011. 8) N. G. Morgenthaler, S. C. Ho, W. B. Minich, Stimulating and blocking TSH-receptor autoantibodies from patients with Graves' disease and autoimmune hypothyroidism have very similar concentration, TSH-receptor affinity and binding sites, *J. Clin. Endocrin. Metab.*, 92:1058-1065, 2007. 9) J. Ginsberg, Diagnosis and management of Graves' disease, *Canad. Med. Assoc. J.*, 168(5):575-585, 2003. 10) P. Laurberg, G. Wallin, L. Tallstedt, M. Abraham-Nordling, G. Lundell, O. Tørring, TSH-receptor autoimmunity in Graves' disease after therapy with anti-thyroid drugs, surgery, or radioiodine: a 5-year prospective randomized study, *Euro. J. Endocrinol.*, 158:69-75, 2008. 11) C. Giuliani, D. Cerrone, N. Harii, M. Thornton, L.D. Kohn, N.M. Dagia, I. Bucci, M. Carpentieri, B. Di Nanno, A. Di Blasio, P. Vitti, F. Monaco, G. Napolitano, A TSHR-LH/CGR chimera that measures functional thyroid-stimulating autoantibodies (TSAb) can predict remission or recurrence in Graves' patients

undergoing antithyroid drug (ATD) treatment, *J. Clin. Endocrinol. Metab.*, 97(7):E1080-E1087, 2012. 12) A.K. Eckstein, M. Plicht, H. Lax, M. Neuhauser, K. Mann, S. Lederbogen, C. Heckmann, J. Esser, N. G. Morgenthaler, Thyrotropin receptor autoantibodies are independent risk factors for Graves' ophthalmopathy and help to predict severity and outcome of the disease, *J. Clin. Endocrinol. Metab.*, 91(9):3464-3470, 2006. 13) M.R. Bjørgaas, H. Farstad, S.C. Christiansen, H-G.K. Blaas, Impact of thyrotropin receptor antibody levels on fetal development in two successive pregnancies in a woman with Graves' disease, *Horm. Res. Paediatr.*, 79:39-43, 2012. 14) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377-82, 387-8. 15) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 16) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030. 17) EP17-A2: Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition, CLSI, Wayne, PA, 2012. 18) EP5-A2: Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition, CLSI, Wayne, PA, 2004.

This product or portions thereof is manufactured under license from RSR Ltd., Cardiff, UK, under European Patent Number 1021721 and US Patent Number 6844162 and foreign equivalents of these patent rights.

Technical Assistance

Outside the United States, contact your National Distributor.

siemens.com/healthcare

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (IU/L)

Sample ⁶	Mean ³ (IU/L)	Repeatability ¹		Within-Lab ²	
		SD ⁴	%CV ⁵	SD	%CV
1	0.34	0.02	7.0	0.03	8.3
2	0.69	0.03	4.1	0.03	5.0
3	1.57	0.07	4.4	0.08	5.3
4	4.43	0.18	4.0	0.26	5.9
5	7.80	0.27	3.5	0.42	5.4
6	29.09	1.91	6.6	2.11	7.3

Cross Reactivity

Cross Reactivity Compound ¹	IU/L Added ²	Apparent Dose ³ (IU/L)
Follicle Stimulating Hormone (FSH)	750	0.01
Luteinizing Hormone (LH)	500	ND
Thyroid-Stimulating Hormone (TSH)	0.14	ND
Human Chorionic gonadotropin (hCG)	100,000	0.01
Anti-Thyroglobulin (aTG)	2,000,000	ND
Anti-Thyroid Peroxidase (aTPO)	2,000,000	ND

ND: Not Detected⁴

Sensitivity and Specificity

IMMULITE 2000 TSI	Diagnostic Accuracy ¹	
	Positive ²	Negative ³
	Positive	232
Negative	4	337

Overall:⁴ 99.1%
Sensitivity:⁵ 98.3%
Specificity:⁶ 99.7%

Method Comparison

IMMULITE 2000 TSI	Thyretain™ TSI Reporter BioAssay	
	Positive ¹	Negative ²
	Positive	128
Negative	1	106

Overall Agreement:³ 95.90%
Positive % Agreement:⁴ 99.22%
Negative % Agreement:⁵ 92.17%

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2018 Siemens Healthcare Diagnostics. All rights reserved.
CONTROL AB