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Vasc Biol 1997;17: 989-95.

7. Malinow MR, Bostom AG, Krauss RM. Homocyst(e)line, diet, and cardiovascular diseases a statement for healthcare professionals from the Nutrition Committee, American Heart Association. *Circulation* 1999; 99:178-82.

8. Nygard O, Nordrehaug JE, Refsum H, Ueland PM, Farstad M, Vollset, SE. Plasma homocysteine levels and mortality in patients with coronary artery disease. *N Engl J Med* 1997;337:230-6.

9. Graham IM, Daly LE, Refsum HM, et al. Plasma homocysteine as a risk factor for vascular disease. *JAMA*. 1997;277(22):1775-1781.

10. Konecny N, Malinow MR, Tunick PA, et al. Correlation between plasma homocyst(e)line and aortic atherosclerosis. *Am Heart J*. 1997;133(5):534-540.

11. Wang H, Yoshizumi M, Lai K, et al. Inhibition of growth and p21ras methylation in vascular endothelial cells by homocysteine but not cysteine. *J Biol Chem*. 1997;272(40):25380-25385.

12. Hogeveen EK, Kostense PJ, Jakobs C, et al. Hyperhomocysteinemia increases risk of death, especially in type-2 diabetes: 5-year follow-up of the Hoorn Study. *Circulation*. 2000;101(13):1506-1511.

13. Still RA, McDowell IF. ACP Broadsheet No. 152: March 1998. Clinical implications of plasma homocysteine measurement in cardiovascular disease. *J Clin Pathol*. 1998;51(3): 183-188.

14. Ubbink JB, Vermaak WJ, van der Merwe A, Becker PJ. The effect of blood sample aging and food consumption on plasma total homocysteine levels. *Clin Chim Acta*. 1992;207(1-2):119-128.

15. Clinical and Laboratory Standards Institute. Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.

16. Clinical and Laboratory Standards Institute. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.

17. Clinical and Laboratory Standards Institute. Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.

18. Clinical and Laboratory Standards Institute. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.

19. Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem*. 1999;45(7):942-956.

20. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem*. 1992;38(9):1737-1742.

21. Clinical and Laboratory Standards Institute. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document

EP28-A3c.

22. Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Proposed Guideline. Wayne, PA: Clinical and Laboratory Standards Institute; 1986. NCCLS Document EP7-P.

23. Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

24. Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.

25. Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.

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Siemens Healthcare Diagnostics Inc.

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