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1) Beischer W. Proinsulin and C-peptide in humans. In: Fotherby

K, Pal S, editors. H Hormones in Normal and Abnormal Human Tissues. Vol 3. Berlin: Walter DeGruyter, 1983: 1-43.

2) Beyer J, Krause U, Cordes U. C-peptide: its biogenesis, structure, determination and clinical significance, *Giornale Ital Chem Clin* 1979; 4 Suppl 1:9-22.

3) Bonser A, Garcia-Webb P. C-Peptide measurement: methods and clinical utility. *CRC Crit Rev Clin Lab Sci* 1984; 19:297-352.

4) Gonen G, Rubenstein AH, Horwitz DL, Blix PM. Clinical significance of C-peptide. In: Baba S, Kaneko T, Yanaihara N, editors, *Proinsulin, insulin, C-Peptide*. Amsterdam: Excerpta Medica, 1979: 246-53.

5) Horwitz D, et al. Proinsulin, insulin and C-Peptide concentrations in human portal and peripheral blood. *J Clin Invest* 1975; 55: 1278-83.

6) Horwitz D, Kuzuya H, Rubenstein AH. Circulating serum C-peptide. *N Engl J Med* 1976; 295: 207-9.

7) Polonsky KS, Rubenstein AH. C-peptide as a measure of the secretion and hepatic extraction of insulin; pitfalls and limitations. *Diabetes* 1984; 33: 486-93.

8) Rendell M. Expanding clinical use of C-peptide radioimmunoassay. *Acta Diabetol Lat* 1983; 20: 105-13.

9) Rubenstein AH, Kuzuya H, Horwitz DL. Clinical significance of circulating C-peptide in diabetes mellitus and hypoglycemic disorders. *Arch Unrestricted Med* 1977; 137: 625-32.

10) Turkington RW, Estkowski A, Link M. Secretion of insulin or connecting peptide; a predictor of insulin dependence of obese diabetics, *Arch Unrestricted Med* 1982; 142: 1102-5.

11) Data on file.

12) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostics blood specimens by venipuncture; approved standard, 4th ed, NCCLS Document H3-A4, Wayne, PA; NCCLS, 1998.

13) CLSI. Protocols for the Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. CLSI document EP17-A Vol. 24 (No. 34).CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2004.

14) Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5-A2. CLSI, 940 W EST Valley Road, Suite 1400 Wayne, Pennsylvania 19087-1898, USA 2004.

15) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377-82, 387-8.

16) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.

17) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

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Siemens Healthcare Diagnostics Products
Limited

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