


IT Information Guide

Introduction

This Technical Guide provides an overview of the Information Technology aspects of the **epoc® Blood Analysis System**. This guide includes information about following topics:

- epoc Enterprise Data Manager (EDM) server requirements;
- VPN remote access for implementation, technical support, and troubleshooting;
- epoc Host Wi-Fi connectivity requirements;
- HL7 Interface between the EDM and the Laboratory or Hospital Information System;
- Installation and continued IT support for the epoc Blood Analysis System.

The following diagram shows the basic components of the epoc Blood Analysis System:



IT Information Guide

epoc Enterprise Data Manager (EDM)

The EDM Software is a web-based software application used to collect test results from the epoc Host devices and to help manage the entire epoc System. The EDM Software is comprised of the following components:

- User interface accessed using Microsoft® Internet Explorer®, Firefox®, or Chrome™;
- Microsoft SQL Server Database Engine.

Note: This software resides on a dedicated server or virtual server built to the minimum specification described in Appendix A.

epoc Blood Analysis System

The epoc Blood Analysis System consist of three components: epoc Reader, epoc Host, and epoc Test Card.

epoc Reader

The epoc Reader is a portable device that reads information from the epoc BGEM Test Card and, using a Bluetooth wireless connection, sends raw measurement data to the epoc Host. The epoc Host uses this data to calculate test results and present them to the caregiver.

epoc Host

The epoc Host is a Microsoft Windows® Mobile-based handheld computer that wirelessly connects to the epoc Reader using Bluetooth® technology to provide a caregiver with the ability to perform patient blood testing. The epoc Host provides the following:

- A user interface from which users can enter patient data, view test results, and perform QA testing;
- The ability to view previously obtained test results stored on the epoc Host and print them to a portable printer;
- A barcode scanner to scan patient and operator IDs;
- A Bluetooth connection to the epoc Reader to perform testing;
- Transmission of test results to the EDM by connecting to an existing Wireless LAN (WLAN) using industry standard wireless security protocols.


Note: See Appendix B for detailed hardware information about the epoc Host.

IT Information Guide

epoc HL7 Interface

The epoc HL7 Interface is a customizable interface used to send patient test results from the EDM to Laboratory or Hospital Information Systems. The following interface capabilities are supported:

- Sending test results to the LIS/HIS with standard acknowledgements;
- Optionally receiving Enhanced Acknowledgements initiated by the LIS/HIS indicating success or failure to process a test result;
- Optionally receiving ADT information containing patient demographic information.



IT Information Guide

Appendix A: EDM Server Requirements

System Requirements				
32/64-bit Operating Systems:	Microsoft Windows Server 2016 R2		Microsoft Windows 7	
	Microsoft Windows Server 2012 R2		Microsoft Windows 8.1	
	Microsoft Windows Server 2008 R2		Microsoft Windows 10	
Supported Languages:	Czech	Estonian	Hungarian	Portuguese (Portugal)
	Danish (Denmark)	Finnish (Finland)	Japanese (Japan)	Romanian (Romania)
	Dutch (Netherlands)	French (France)	Italian (Italy)	Russian (Russian)
	English (US)	German (Germany)	Norwegian (Bokmål)	Simplified Chinese
	Estonian	Greek (Greece)	Polish	Spanish (Spain)
				Swedish (Sweden)
System Type:	Dedicated hardware or running in a virtual environment.			
Processor Type:	Intel® Core™ i5, i7, Intel Xeon™, or equivalent			
Number of CPUs:	1 (4 cores per processor)			
Memory:	4 GB or more			
Display Resolution:	1 024x768 minimum			
Network Interface	Ethernet 10BASE-T, 100BASE-T, or 1000BASE-T			
Estimated Storage Requirements				
Primary Drive:	60GB or more			
Secondary Drive:	40GB or more (optional but recommended for larger installations)			
Database Size:	5MB initial size growing by approximately 10MB per 1000 Test Records			
Software Requirements				
Database Engine:	Microsoft SQL Server 2008 R2, 2012, 2014, 2016, 2017			
Required Software:	Microsoft Internet Explorer 11, Firefox, Chrome			
	Microsoft .NET Framework 4.6.1 or newer.			
	Microsoft Internet Information Services (IIS) (version dependent on OS)			
	Java™ Runtime 8.0 (including the latest service packs)			
	Microsoft SQL Server Management Studio			
	Adobe® Reader® (latest version)			
	Remote Desktop Enabled			
Data Protection				
Virus Protection:	Customer is responsible for installing anti-virus and anti-malware protection on the EDM Server.			
Backup and Recovery:	Customer is responsible for EDM database backup.			
Microsoft updates:	Customer is responsible to installing Microsoft updates and security patches according to their policies and procedures.			
VPN and Remote Access				
VPN:	VPN access is required for Installation and continued system support.			
Remote Control:	Microsoft Remote Desktop is used to remotely access the EDM Server.			
Administration:	A login account with local administrator privileges is required for installation, interface development, and system maintenance.			
User Access				
EDM User Access:	Users can access the EDM Web application using the following:			
	• Microsoft Internet Explorer 11, Firefox, Chrome			
	• Minimum display resolution of 1024 x 768			

IT Information Guide

Appendix B: epoc Host² Specifications

epoc Host ² Hardware	
Mobile Computer:	Zebra MC55A0-HC, mobile computer
Operating system:	Microsoft Windows Embedded Handheld 6.5 (previously Windows Mobile 6.5 Classic)
Processor:	Marvell™ PXA 320 @ 806 MHz processor
Memory:	256 MB RAM; 1 GB Flash
Display:	PenTile™ 3.5 inch Color VGA, super bright 650+ NITS
Battery:	Rechargeable Lithium Ion 3.7V 3600mAh Smart Battery
Wi-Fi Interface	
Purpose:	Transmit test results to a data manager using existing Wireless infrastructure
Network Interface:	Tri-mode IEEE™ 802.11a (5 GHz) 802.11b/g (2.4GHz)
Output Power:	14.5 dBm (typical)
Range:	Up to 30m
Wi-Fi Security:	WEP (40/128), WPA (TKIP), WPA2 (AES), TLS, TTLS (MS-CHAP), TTLS (MS-CHAP v2), TTLS (CHAP), TTLS-MD5, TTLS-PAP, PEAP-TLS, PEAP (MS-CHAP v2), LEAP, FAST
Wi-Fi Certification:	CCXv4 certified FIPS 140-2 certified (NIST FIPS 140-2 Certificate #1210)
Bluetooth® Interface	
Purpose:	Communicate with epoc Readers and Bluetooth Printers
Bluetooth Version:	Class II, v 2.1 Enhanced Data Rate (EDR)
Bluetooth Security:	128-bit encryption with PIN authentication.
Output Power:	Maximum output power: 2.5mW (4dBm) Nominal output power: 1mW (0dBm) Minimum output power: 0.25 mW (-6dBm)
Range:	Up to 10m (Data Rate: 3Mbit/s)
Barcode Scanner	
Purpose:	To scan Operator ID badges and Patient bar-coded wristbands
Barcode Scanner:	2D Imager (SE4500-DL)
Barcode Formats:	1D: UPC-A, UPC-E, UPC-E1, EAN-8, EAN-13, Code-128, UCC EAN-128, Code-39 Full ASCII Conversion, Code-39, ISBT-1281, Code-93, Interleaved 2 of 5, Discrete 2 of 5, Codabar, MSI 2D: PDF417, MICROPDF, MACROPDF, MaxiCode, DataMatrix, QRCode, MACROMICROPDF, Composite AB, Composite, TLC39, AZTEC, MICROQR

IT Information Guide

Appendix B: epoc Host² Specifications (continued)

Certification and Compliance	
Electrical Safety:	Certified to UL/cUL 60950-1, IE C/EN60950-1
Environmental:	RoHS-compliant; RoHS2 Directive 2011/65/EU and REACH
Wi-Fi and Bluetooth:	USA: FCC Part 15.247, 15.407 Canada: RSS-210, RSS 310 EU: EN 300 328, EN 300 440-2, EN 301 893 Japan: ARI B STD-T33, ARI B STD-T66, ARI B, STD-T70 & T71 Australia: AS/NZS 4268
RF Exposure:	USA: FCC Part 2, FCC OET Bulletin 65, Supplement C Canada: RSS-102 EU: EN 50392 Australia: Radio communications Standard 2003
EMI/RFI	North America: FCC Part 15, Class B Canada: ICES 003 Class B EU: EN55022 Class B, EN 301 489-1, EN 301 489-17, EN 301 489-19, EN 60601-1-2
Laser Safety	IE C/Class 2/FDA II in accordance with IE C60825-1/EN 60825-1
Other	
Operating Temperature	-10°C -50°C [14°F -122°F]
Charging Temperature	0° C - 40° C [32° F - 104° F]
Shipping/Storage Temperature	-40° C - 70° C [-40° F - 158° F]
Operating/Shipping/Storage Humidity	Up to 85% relative humidity, non-condensing
Water Ingress	IP64

Questions

If you have any questions, please contact your Siemens Customer Care Center or your local Siemens technical support representative.

Manufactured by:
Epocal Inc.
 2060 Walkley Rd.
 Ottawa, ON K1G 3P5 Canada
 Tel.: +1-613-738-6192
 Fax: +1-613-738-6195
siemens.com/epoc

Authorized Representative:
Emergo Europe
 Prinsessegracht 20, 2514 AP
 The Hague, Netherlands
service@emergogroup.com

© 2019 Siemens Healthcare Diagnostics Inc. All rights reserved.
 epoc and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates.

51002506 Rev.: 13 2019-05