Dimension® clinical chemistry system

Flex® reagent cartridge

DGNA

See shaded sections: Updated information from 2018-09 version.


Issue Date 2019-04-01

Digoxin

Intended Use: The DGNA method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended to measure digoxin, a cardiovascular drug, in human serum and plasma. DGNA test results are used in the diagnosis and treatment of digoxin overdose and in monitoring levels of digoxin to ensure appropriate therapy.

Summary: Digoxin, a cardiac glycoside, is used as an antiarrhythmic agent, both alone and in conjunction with other drugs. Absorption from the gastrointestinal tract is variable: 60 – 80% of the administered dose is absorbed. Digoxin is excreted by the kidney almost entirely unchanged. Therefore the patient's renal function is an important consideration in determining dosage. In persons with normal kidney function the half-life is about 1.5 days. The most serious complications of digoxin toxicity are ventricular arrhythmias: ventricular tachycardia and ventricular fibrillation.^{1,2}

Principles of Procedure: The digoxin method uses an immunoassay technique in which free and digoxin-bound antibody-enzyme species are separated using magnetic particles. The DGNA chemistry is optimized for measurement of β -galactosidase activity. Magnesium acetate is included to activate the enzyme and N-2-Hydroxyethylpiperazine-N'-1-ethanesulfonic acid (HEPES) buffer to provide optimum pH. The methodology for DGNA involves Antibody Conjugate reagent mixing with patient's serum or plasma. The Antibody Conjugate reagent utilizes the F(ab')₂ fragment of the antibody to eliminate interference from rheumatoid factor. Digoxin in the sample is bound by the F(ab')₂- β -galactosidase in the Antibody Conjugate reagent. Magnetic particles coated with the digoxin analog ouabain are added to bind free (unbound) antibody-enzyme conjugate. The reaction mixture is then separated magnetically. Following separation, the supernatant containing the digoxin-antibody-enzyme complex is transferred and mixed with a substrate. The β -galactosidase (β -gal) portion of the Digoxin-F(ab')₂- β -galactosidase complex catalyzes the hydrolysis of chlorophenol- β -D-galactopyranoside (CPRG) to chlorophenol red (CPR). The change in absorbance at 577 nm due to the formation of CPR is directly proportional to β -galactosidase activity. Since β -galactosidase is not present in serum, its activity is directly proportional to digoxin in the patient's sample and is measured using a bichromatic (577, 700 nm) rate technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1, 2	Liquid	Antibody Conjugate		Rabbit
3, 4	Tablets	Reagent and stabilizers	c	
		Ouabain		
		Magnetic Particles	0.3%	
5, 6	Tablets	CPRG	7mM	
7	Liquid	Substrate Diluent, Buffer	100mM	

- a. Wells are numbered consecutively from the wide end of the cartridge.
b. Nominal value per test at manufacture.
c. Antibody titre and conjugate activity vary from lot to lot.

Risk and Safety:

H319, H317
P280, P305 + P351 + P338, P501

Warning!

Causes serious eye irritation. May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: Chromium oxide; 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: Hydrating, diluting and mixing are automatically performed by the instrument.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6

10 days for well 7

Specimen Collection and Handling: Serum and heparinized plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.³

Follow the instructions provided with your specimen collection device for use and processing.⁴

Each laboratory must determine its appropriate sample timing/collection protocol.

Complete clot formation should take place before centrifugation.⁵ Specimens should be free of particulate matter.

Separated specimens are stable for 8 hours at 20 – 25 °C, 7 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C or colder for up to 6 months.⁶

The purpose of specimen storage information is to provide guidance to users; however, users may validate their own procedures for storing patient samples.

Procedure

Materials Provided

DGNA Flex® reagent cartridge, Cat. No. DF35A

Materials Required But Not Provided

Reaction vessels (required for HM only) RXV1A

Drug Calibrator, Cat. No. DC22B

Quality Control Materials

Test Steps

Sampling,⁴ reagent delivery, mixing, separation, processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

- d. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

Test Conditions

Test conditions for separation module or heterogeneous immunoassay module (HM):

Reaction Vessel

Sample Size	30 μ L
Antibody Conjugate Reagent Volume	100 μ L
Ouabain Particle Volume	75 μ L
Diluent Volume	115 μ L
Incubation Temperature	42 °C*

*Dimension® EXL™ with LOCI® Module and Non-HM Dimension® systems: 37 °C

Cuvette

Transfer Volume	60 μ L
Substrate Reagent Volume	175 μ L
Diluent Volume	165 μ L
Temperature	37 °C
Wavelength	577 and 700 nm
Type of Measurement	Bichromatic rate

Calibration

Assay Range	0.20 – 5.00 ng/mL [0.26 – 6.41 nmol/L] ^a
Calibration Material	Drug Calibrator, Cat. No. DC22B
Calibration Scheme	5 levels, n = 2
Units	ng/mL [nmol/L] (ng/mL x 1.28) = [nmol/L]
Typical Calibration Levels	0.00, 0.60, 1.20, 2.50, 5.00 ng/mL [0.00, 0.77, 1.54, 3.20, 6.41 nmol/L]
Calibration Frequency	Every month for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C_0 1.0 C_1 108.0 C_2 -2.0 C_3 1.77 C_4 0.5

e. Système International d'Unités [SI Units] are in brackets.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of a Quality Control (QC) material with known digoxin concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of digoxin in ng/mL [nmol/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.20 – 5.00 ng/mL [0.26 – 6.41 nmol/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 5.00 ng/mL [6.41 nmol/L] should be repeated on dilution.

Manual Dilution: Dilute the sample with a digoxin-free serum or Level 1 (0 ng/mL) of Drug Calibrator to produce a sample concentration within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): Autodilution is not recommended for this method.

DGNA results less than 0.20 ng/mL [0.26 nmol/L] should be reported as "less than 0.20 ng/mL [0.26 nmol/L]" instead of the numerical value.

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
0.60 ng/mL [0.77 nmol/L]	>0.06 ng/mL [0.08 nmol/L]
2.40 ng/mL [3.07 nmol/L]	>0.09 ng/mL [0.12 nmol/L]

Lower concentrations of 0.5 – 1.2 ng/mL [0.64 – 1.54 nmol/L] have been found to be more appropriate in certain populations such as chronic heart failure patients.^{7,8} Digoxin toxicity is commonly associated with serum levels > 2.0 ng/mL [2.6 nmol/L] but may occur with lower digoxin levels. Significant overlap of toxic and nontoxic values has been reported. Consequently, analysis of serum concentrations alone is not sufficient for optimization of digoxin therapy. Additional factors such as age, thyroid condition, electrolyte balance, hepatic and renal functions, and other clinical symptoms must be considered.⁹

Interfering Substances

Endogenous, digoxin-like immunoreactive factors (DLIF) have been detected in the serum and plasma of neonates, pregnant women, and patients in renal and hepatic failure. Several studies have established that these factors can cause falsely elevated digoxin measurements when assayed by commercially available immunoassays.¹⁰

Digoxin bound to Fab fragments of antidigoxin antibodies (ex. Digibind®), as found in serum and plasma of individuals being treated for digoxin toxicity, may result in misleading digoxin values.¹¹

Lipemia (Intralipid®) at 3000 mg/dL [33.9 mmol/L] tripped a test report message; therefore the magnitude of the interference could not be determined.

Digibind® is a registered trademark of Glaxo SmithKline, Inc., Research Triangle Park, NC
Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

The assay has also been designed to minimize interference from antibodies to β-galactosidase. These antibodies can be encountered in samples as a consequence of bacterial infection and may produce falsely elevated results which may not be consistent with clinical evaluation. In very rare instances, immunoassays may produce falsely elevated or decreased results due to other patient-specific interferents.¹²

Therapeutic Range: Therapeutic digoxin concentrations vary significantly, depending on the individual. A range of 0.90 – 2.00 ng/mL [1.15 – 2.56 nmol/L] includes effective serum concentrations for many patients; however, some individuals are best treated at concentrations outside this range. Concentrations greater than 2.00 ng/mL [2.56 nmol/L] are often associated with toxic symptoms.^{13,14}

Specific Performance Characteristics¹

Material	Precision ^a		Standard Deviation (% CV)	
	Mean ng/mL [nmol/L]	Within-run	Total	
Dade® TDM PLUS XL				
Level 1	0.56 [0.72]	0.02 [0.03] (4.1)	0.03 [0.04] (5.2)	
Level 2	1.58 [2.02]	0.02 [0.03] (1.1)	0.03 [0.04] (2.0)	
Level 3	3.08 [3.94]	0.03 [0.04] (1.0)	0.05 [0.06] (1.8)	

- f. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).
- g. Specimens at each level were analyzed in duplicate, two runs/day for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

Dade® is the registered trademark of Siemens Healthcare Diagnostics, Deerfield, IL 60015.

Method Comparison				
Regression Statistics ^b				
Comparative Method	Slope	Intercept ng/mL [nmol/L]	Correlation Coefficient	n
DGN method				
on aca® analyzer	1.03	0.16 [0.20]	0.976	197 ⁱ
Abbott Tdx®	1.10	0.02 [0.03]	0.985	206 ^j
Stratus®	1.03	0.12 [0.15]	0.978	148 ^k

h. Model equation for regression statistics is: [results of Dimension® system] = slope x [comparative method results] + intercept.

i. Range of samples: 0.21 – 3.61 ng/mL [0.3 – 4.6 nmol/L]

j. Range of samples: 0.26 – 4.47 ng/mL [0.3 – 5.7 nmol/L]

k. Range of samples: 0.31 – 3.93 ng/mL [0.4 – 5.0 nmol/L]

Tdx® is the registered trademark of Abbott Diagnostics, Abbott Park, IL 60064.

Specificity

Hemolysis, Icterus, Lipemia (HIL) Interference

The DGNA method was evaluated for interference according to CLSI/NCCLS EP7-P. Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Test Concentration SI Units	DGNA Concentration ng/mL [nmol/L]	Bias % ^l
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L] (monomer)	2.16 [2.76]	<10
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	2.16 [2.76]	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L] 3000 mg/dL [33.9 mmol/L]	2.33 [2.98] 2.33 [2.98]	<10 m

l. Analyte results should not be corrected based on this bias.

m. Lipemia (Intralipid®) of 3000 mg/dL [33.9 mmol/L] and above tripped an error flag on this method, so the magnitude of the interference could not be determined.

Non-Interfering Substances

The following substances have no measurable effect on the DGNA method at the concentrations indicated:

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1323 µmol/L
Acetylsalicylic Acid	50 mg/dL	1.67 mmol/L
Aldosterone	1 µg/mL	8.32 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Amobarbital	10 mg/dL	442 µmol/L
Ampicillin	4 mg/dL	114 µmol/L
Ascorbic Acid	3 mg/dL	170 µmol/L
Beta-Carotene	0.6 mg/dL	11.2 µmol/L
Caffeine	10 mg/dL	515 µmol/L
Carbamazepine	12 mg/dL	508 µmol/L
Chloramphenicol	25 mg/dL	774 µmol/L
Chlordiazepoxide	2 mg/dL	67 µmol/L
Chlorpromazine	5 mg/dL	157 µmol/L
Cimetidine	10 mg/dL	396 µmol/L
Codeine	10 mg/dL	334 nmol/L
Cortisol	3 µg/mL	8.28 nmol/L
Cortisone	3 µg/mL	8.28 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 75	2500 mg/mL	33 µmol/L
Diazepam	2 mg/dL	70 µmol/L
Erythromycin	20 mg/dL	273 µmol/L
17-Estradiol	3 µg/mL	11 µmol/L
17-β-Estradiol	3 µg/mL	11 µmol/L
Estriol	3 µg/mL	10.4 µmol/L
Estrone	3 µg/mL	11.1 µmol/L
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 µmol/L
Furosemide	2 mg/dL	60 µmol/L
Gentamicin	12 mg/dL	221 µmol/L
Ibuprofen	40 mg/dL	1939 µmol/L
Lidocaine	6 mg/dL	256 µmol/L
Lithium	3.5 mg/dL	5.07 mmol/L
Nicotine	2 mg/dL	123 µmol/L
Nortriptyline	1000 ng/mL	3797 nmol/L
Penicillin	80 mg/mL	2247 µmol/L
Pentobarbital	10 mg/dL	442 µmol/L
Phenobarbital	15 mg/dL	646 µmol/L
Phenytoin	10 mg/dL	396 µmol/L
Prednisone	3 µg/mL	10.8 µmol/L
Primidone	10 mg/dL	458 µmol/L
Progesterone	3 µg/mL	9.45 µmol/L
Propoxyphene	0.4 mg/dL	12 µmol/L
Protein	4 g/dL	40 g/L
Protein	12 g/dL	120 g/L
Rheumatoid Factor	707 IU/L	707 IU/L
Salicylic Acid	50 mg/dL	3.62 mmol/L
Secobarbital	5 mg/dL	210 µmol/L
Sodium Fluoride	1 mg/dL	0.01 g/L
Testosterone	2 µg/mL	6.93 µmol/L
Theophylline	25 mg/dL	1388 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3472 µmol/L

Cross-reactivity

The following substances have been shown to cross-react with the DGNA method:

Substance	Level Tested	% Cross Reactivity ⁿ
Digoxigenin	5 ng/mL [12.80 nmol/L]	142.6
Gitoxin	25 ng/mL [32.02 nmol/L]	7.6
Digitoxin	25 ng/mL [32.75 nmol/L]	11.2
Acetyldigitoxin	25 ng/mL [30.75 nmol/L]	10.7
Dihydrodigoxin	25 ng/mL [30.50 nmol/L]	9.4

n. % Cross Reactivity = apparent digoxin concentration divided by the concentration of the cross-reactant multiplied by 100.

Analytical Sensitivity: 0.20 ng/mL [0.26 nmol/L]

The analytical sensitivity represents the lowest concentration of digoxin that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations at the 0.00 ng/mL [nmol/L] level 1 Drug Calibrator.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, Flex®, EXL™ and LOGI® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics

All rights reserved.

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Finn AL, Taylor WJ. Individualizing Drug Therapy, Practical Applications of Drug Monitoring. New York: Gross, Townsend, Frank, Inc., 1981: 2-6.
2. Burtis CA, Ashwood ER, Bruns DE. Tietz Textbook of Clinical Chemistry and Molecular Biology, Fourth Edition, Elsevier Saunders, St. Louis, MO; 2006, p. 1256.
3. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
4. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
6. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, 3rd ed. Washington, DC: AACC Press, 2007, p. 343.
7. ESC Guidelines for diagnosis and treatment of patients with heart failure 2008. Eur. J Heart Failure 10: 933-989, 2008.
8. 2009 Focused Update Incorporated into the ACC/AHA 2005 Guidelines for the Diagnosis and Management of Heart Failure in Adults: A Report of the American College of Cardiology Foundation/American Heart Association Task Force on Practice Guidelines: Developed in Collaboration With the International Society for Heart and Lung Transplantation, Circulation, 2009; 119: e391 – e479.
9. Dec, William G., Digoxin remains useful in the management of chronic heart failure, Med Clin N Am 87: 317-337, 2003.
10. Stone JA, and Soldin SJ. Update On Digoxin, Clin Chem, 1989, Vol 35, pp 1326-1331.
11. Reiney. Effects of Digoxin Immune Fab, American Journal of Clin Path, 1990, Vol 92, pp 779-786.
12. Ismail AA, Walker, PL, Cawood ML and Arth JH. Interference in immunoassay is an underestimated problem. Ann. Clin. Biochem 2002; 39:1737-1742.
13. Evans WE, Schentag JJ, Jusko WJ. Applied Pharmacokinetics, Applied Therapeutics Inc., San Francisco, CA, 1980, p 330.
14. Huffman DH, Crow JW, Pentikainen P, Azarnoff DL. Association between clinical cardiac status, laboratory parameters and digoxin use, Amer Heart Journal 1976; 91:28-34.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFIS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare