Dimension® clinical chemistry system

Flex® reagent cartridge

PTN

See shaded sections: Updated information from 2017-12 version.

Issue Date 2019-04-01

Phenytoin


Intended Use: The PTN method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended to measure phenytoin (dilantin, diphenylhydantoin), an anti-epileptic drug, in human serum and plasma. PTN test results are used in the diagnosis and treatment of phenytoin overdose and in monitoring levels of phenytoin to ensure appropriate therapy.

Phenytoin is widely used and effective for all types of seizure disorders except absence seizures.¹ It is of value in the treatment of elementary partial (focal) or complex partial epilepsy (psychomotor, temporal lobe) seizures, but ineffective in petit mal epilepsy.² Occasionally, it is used in the treatment of cardiac arrhythmias.¹

Because of considerable inter-individual variation and the limited capacity for the liver to metabolize phenytoin, blood concentrations should be monitored to obtain maximal anti-seizure effect. Once metabolism is saturated, small dosage changes may result in disproportionate changes in blood concentration and can cause wide variations in dosage requirements among patients.¹

Summary: The phenytoin method is based on a particle enhanced turbidimetric inhibition immunoassay (PETINIA) technique.

Principles of Procedure: The phenytoin method is a homogeneous particle enhanced turbidimetric inhibition immunoassay (PETINIA) technique which uses a synthetic particle-phenytoin reagent (PR) and phenytoin-specific monoclonal antibody (AB). Phenytoin present in the sample competes with the particles for the antibody, thereby decreasing the rate of aggregation. Hence, the rate of aggregation is inversely proportional to the concentration of phenytoin in the sample. The rate of aggregation is measured using bichromatic turbidimetric readings at 340 nm and 700 nm.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1, 2, 3	Liquid	Antibody	40 µL ^c	Mouse, monoclonal
4, 5, 6	Liquid	Particle Reagent	40 µL ^c	
7, 8	Liquid	Buffer	290 µL	

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value in final reaction mixture.

c. The antibody and particle reagent are a matched pair. The antibody titer and particle reagent concentration may vary from lot to lot.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501

Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone.

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6
7 days for wells 7 – 8

Specimen Collection and Handling: Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.³

Separated samples drawn in non-gel collection tubes are stable for 24 hours at room temperature, 2 days at 2 – 8 °C and 5 months stored at -20 °C.⁴

Some studies have shown the stability of phenytoin to be related to storage conditions and can vary with time, temperature and blood collection tube.^{5,6,7} Studies have shown variable rates of absorption of phenytoin by the gel barrier in several gel formulations available.

The amount of absorption and the subsequent clinical significance of this phenomenon are directly related to the length of time the serum is stored upon the gel, the volume of sample above the gel, the surface area of the gel and the temperature of the sample storage.^{6,7,8}

Follow instructions with your specimen collection device for use and processing.⁸

Complete clot formation should take place before centrifugation.⁹ Specimens should be free of particulate matter.

Heparin, citrate and EDTA anticoagulants do not interfere with the PTN method at concentrations normally found in blood collection tubes.

Each laboratory must determine its appropriate sample timing/collection protocol.

Procedure

Materials Provided

PTN Flex® reagent cartridge, Cat. No. DF64

Materials Required But Not Provided

Drug Calibrator, Cat. No. DC22B

Quality Control Materials

Test Steps

Sampling,⁴ reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

d. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

Test Conditions

Sample Size	4 µL
Buffer Volume	290 µL
Particle Reagent Volume	40 µL
Reagent 3 Volume	40 µL
Diluent Volume	111 µL
Temperature	37 °C
Wavelength	340 and 700 nm
Type of Measurement	Turbidimetric rate

Calibration

Assay Range	0.5 – 40.0 µg/mL [2.0 – 158.4 µmol/L] ^a
Calibration Material	Drug Calibrator, Cat. No. DC22B
Calibration Scheme	5 levels, n=2
Units	µg/mL [µmol/L] (µg/mL x 3.96) = [µmol/L]
Typical Calibration Levels	0.0, 5.0, 10.0, 20.0, 40.0 µg/mL [0.0, 19.8, 39.6, 79.2, 158.4 µmol/L]
Calibration Frequency	Every 30 days for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ 600.0 C ₁ -682.0 C ₂ -1.90 C ₃ 17.1 C ₄ 0.500

e. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known phenytoin concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of phenytoin in µg/mL [µmol/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.5 – 40.0 µg/mL [2.0 – 158.4 µmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

• Samples with results in excess of 40.0 µg/mL [158.4 µmol/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilutions with level 1 (0 µg/mL) of Drug Calibrator or drug-free serum to obtain result within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): Refer to your Dimension® Operator's Guide.

Results less than 0.5 mg/mL [2.0 µmol/L] should be reported as "less than 0.5 mg/mL [2.0 µmol/L]" instead of the numerical value.

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide for troubleshooting.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
10.0 µg/mL [39.6 µmol/L]	> 1.0 µg/mL [3.96 µmol/L]
20.0 µg/mL [79.2 µmol/L]	> 1.7 µg/mL [6.7 µmol/L]

Interfering Substances

Lipemia (Intralipid®) at 1000 mg/dL [11.3 mmol/L] and above tripped a test report message; therefore the magnitude of the interference could not be determined.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Therapeutic Range: Therapeutic phenytoin concentrations vary significantly depending on the individual patient. A range of 10.0 – 20.0 µg/mL [39.6 – 79.2 µmol/L] for peak drug levels indicates effective plasma or serum levels for many patients; however, some individuals are best treated at concentrations outside this range. Concentrations greater than 30.0 µg/mL [118.9 µmol/L] are often associated with toxic symptoms.¹⁰

Specific Performance Characteristics¹

Material	Precision ^a		
	Mean	Standard Deviation (% CV)	
	µg/mL [µmol/L]	Within-run	Total
Serum pool	4.6 [18.2]	0.14 [0.55] (3.1)	0.32 [1.26] (7.1)
Serum pool	32.4 [128.4]	0.60 [2.37] (1.9)	1.17 [4.63] (3.6)
Liquichek® TDM I	6.3 [24.9]	0.15 [0.59] (2.4)	0.42 [1.66] (6.7)
Liquichek® TDM II	13.3 [52.7]	0.37 [1.46] (2.8)	0.68 [2.69] (5.2)

- f. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).
- g. Specimens at each level were analyzed in duplicate, two runs/day for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

Liquichek® is a registered trademark of Bio-Rad, ECS Division, Anaheim, CA 92806.

Comparative Method	Regression Statistics ^b			n
	Slope	Intercept	Correlation Coefficient	
		µg/mL [µmol/L]		
aca® discrete clinical analyzer Serum/Plasma	1.03	-0.58 [-2.3]	0.991	137
TDx® Serum/Plasma	0.98	-2.14 [-8.5]	0.989	161

- h. Model equation for regression statistics is: results of Dimension® system = [slope x (comparative method results)] + intercept.

TDx® is a registered trademark of Abbott Diagnostics, Abbott Park, IL 60064.

Specificity

HIL Interference

The PTN method was evaluated for interference according to CLSI/NCCLS EP7-P. Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Test Concentration SI Units	Phenytoin Concentration µg/mL [µmol/L]	Bias ⁱ %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L] (monomer)	28.30 [112.1]	<10
Billirubin (unconjugated)	80 mg/dL [1368 µmol/L]	30.13 [119.3]	<10
Lipemia (Intralipid®)	600 mg/dL [6.78 mmol/L]	30.35 [120.2]	<10

- i. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances do not interfere with the PTN method when present in serum or plasma in the amounts indicated. Systematic inaccuracies (bias) due to these substances are less than 2 µg/mL [7.9 µmol/L] at a phenytoin concentration of 20.0 µg/mL [79.2 µmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	200 µg/mL	1324 µmol/L
Acetylsalicylic Acid	500 µg/mL	1.67 mmol/L
Amikacin	150 µg/mL	256 µmol/L
Amobarbital	100 µg/mL	442 µmol/L
Ampicillin	50 µg/mL	143 µmol/L
Ascorbic Acid	30 µg/mL	170 µmol/L
Butabarbital	100 µg/mL	472 µmol/L
Caffeine	100 µg/mL	515 µmol/L
Carbamazepine	120 µg/mL	508 µmol/L
Chloramphenicol	250 µg/mL	774 µmol/L
Chlordiazepoxide	20 µg/mL	67 µmol/L
Chlorpromazine	50 µg/mL	157 µmol/L
Cholesterol	500 mg/dL	13.0 mmol/L
Cimetidine	100 µg/mL	397 µmol/L
Diazepam	20 µg/dL	70 µmol/L
Digoxin	3.23 ng/mL	4.13 nmol/L
Ethanol	750 mg/dL	163 mmol/L
Ethosuximide	300 µg/mL	2.13 mmol/L
Erythromycin	200 µg/mL	272 µmol/L
5-ethyl-5-phenylhydantoin	30 µg/mL	147 µmol/L
Fluoride	10 µg/mL	26 µmol/L
Gentamicin	120 µg/mL	221 µmol/L
Hydroxyphenyl-phenylhydantoin (HPPH)	30 µg/mL	112 µmol/L
Ibuprofen	400 µg/mL	1939 µmol/L
Mephenytoin	250 µg/mL	1147 µmol/L
Mephobarbital	150 µg/mL	610 µmol/L
Methsuximide	75 µg/mL	369 µmol/L
Nicotine	20 µg/mL	123 µmol/L
Nortriptyline	1000 ng/mL	3.8 nmol/L
2-phenyl-2-ethylmalonamide (PEMA)	100 µg/mL	485 µmol/L
Penicillin	800 µg/mL	2247 µmol/L
Pentobarbital	100 µg/mL	442 µmol/L
Phenobarbital	150 µg/mL	647 µmol/L
Primidone	100 µg/mL	459 µmol/L
Salicylic Acid	50 mg/dL	3620 µmol/L
Secobarbital	50 µg/mL	210 µmol/L
Theophylline	250 µg/mL	1389 µmol/L
Valproic Acid	500 µg/mL	3467 µmol/L

Analytical Sensitivity: 0.4 µg/mL [2.0 µmol/L]

The sensitivity represents the lowest concentration of phenytoin that can be distinguished from zero. **The sensitivity is defined as the mean value (n = 20) plus two standard deviations above the level 1 (0.0 µg/mL [0 µmol/L]) Drug Calibrator.**

Symbols Key: See adjacent panel.

Dimension®, aca® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografia

1. Baer DM and Dito WR. Interpretations in Therapeutic Drug Monitoring, Educational Products Division, ASCP, Chicago, IL, 1981, pp 49, 141 (monitoring of anti-epileptic drugs) and pp 362-363 (key points for interpretation of drug concentrations).
2. Finn AL, and Taylor WJ. *Individualizing Drug Therapy*, Practical Applications of Drug Monitoring, Gross, Townsend, Frank, Inc., New York, NY, 1981, pp 63-85.
3. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
4. Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 3rd ed. Washington, DC: AACC Press, 2007: p 714.
5. Parish RC and Alexander T. Stability of Phenytoin in Blood Collected in Vacuum Blood Collection Tubes, *Ther Drug Monit* 1990; 12: 85-89.
6. Dasgupta A, Dean R, Saldana S, Kinnaman G, and McLawhon R. Absorption of Therapeutic Drugs by Barrier Gels in Serum Separator Blood Collection Tubes, *AJCP* 1994; 101:4, 456-461.
7. Bergqvist Y, Eckerborn S, and Funding L. Effect of use of Gel-Barrier Sampling Tubes on the Determination of Some Anti-epileptic Drugs in Serum, *Clin Chem* 1984; 30: 465-466.
8. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
9. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
10. Rall, TW and Schleifer, LS, Goodman and Gilman's *The Pharmacological Basis of Therapeutics*, MacMillan Publishing Co., New York, NY 1980, pp.448-455.

Symbols Key
Symbol Schlüssel
Explication des Symboles
Interpretazione simboli
Clave de los Símbolos

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabricante / Fabricante
	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos
	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
	Contents / Inhalt / Contenu / Contenido / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
	Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare