Dimension® clinical chemistry system

Flex® reagent cartridge

T4

See shaded sections: Updated information from 2017-08 version.


Issue Date 2019-04-01

Thyroxine

Intended Use: The T4 method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of total thyroxine in human serum and plasma.

Summary: The thyroxine method is an adaptation of the EMIT® homogeneous enzyme immunoassay technology.

Principles of Procedure: The thyroxine method is a homogeneous enzyme immunoassay for the measurement of total thyroxine (T4) (protein bound and free). Sodium hydroxide and sodium salicylate are used to release T4 from binding proteins. The released T4 binds to anti-T4 antibody, thereby reducing the amount of T4-glucose-6-phosphate dehydrogenase conjugate (T4-G6PDH) that can be bound to the antibody. The T4-G6PDH which is not bound to the antibody catalyzes the oxidation of glucose-6-phosphate with the simultaneous reduction of NAD⁺ to NADH which absorbs at 340 nm. In the absence of T4, T4-G6PDH binds with the antibody and enzyme activity is reduced. The increase in absorbance at 340 nm due to the formation of NADH over a 60 second measure period is proportional to the activity of the T4-G6PDH. The T4 concentration is measured using a bichromatic (340, 383 nm) rate technique.

where: TBP = Thyroxine binding proteins
Ab = Antithyroxine antibody
T4 = Thyroxine
T4-G6PDH = Thyroxine-glucose-6-phosphate dehydrogenase conjugate

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1, 2	Liquid	T4-G6PDH Buffer Microbial Inhibitor	c	Sheep
3, 4	Tablet ^d	NAD G6P	9.6 mM 14.4 mM	
5, 6	Liquid	Antibody Buffer Microbial Inhibitor	c	
7	Liquid	Buffer		
8	Liquid	Sodium hydroxide Sodium salicylate	19.2 mM 1.2 mM	

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per test at manufacture.

c. The antibody and conjugate are a matched pair. The antibody titer and the conjugate activity may vary from lot to lot.

d. Tablet contains excipients.

Risk and Safety:

H290, H319, H315
P280, P305 + P351 + P338, P501

Warning!

May be corrosive to metals. Causes serious eye irritation. Causes skin irritation.

Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: Sodium hydroxide

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Contains sodium azide (<0.1%) as a preservative. Sodium azide can react with copper or lead pipes in drain lines to form explosive compounds. Dispose of properly in accordance with local regulations.

Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Use Limitations: Performance of this assay has not been established with neonatal specimens.

Reagent Preparation: Hydrating, diluting and mixing are automatically performed by the Dimension® system.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6
15 days for wells 7 – 8

Specimen Collection and Handling: Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.¹

Follow the instructions provided with your specimen collection device for use and processing.²

Complete clot formation should take place before centrifugation.³

Specimens should be free of particulate matter.

Separated samples may be stored for 7 days at 20 – 25 °C or 2 – 8 °C, and for 30 days at -20 °C or colder.⁴

The use of EDTA as an anticoagulant is not recommended.

Lithium heparin and sodium fluoride in blood collection tubes do not affect the T4 result.

Corvac® and SST® collection tubes do not affect the T4 result.

Corvac® is a registered trademark of Monoject, Division of Sherwood Medical, St. Louis, MO. SST® is a registered trademark of Becton-Dickinson, Rutherford, NJ.

Procedure

Materials Provided

T4 Flex® reagent cartridge, Cat. No. DF65

Materials Required But Not Provided

T4 Calibrator, Cat. No. DC13

Quality Control Materials

Test Steps

Sampling,⁵ reagent delivery, mixing, processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

e. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus the dead volume. Precise container filling is not required.

Test Conditions

Sample Size	16 µL
Reagent 1 Volume	20 µL
Reagent 2 Volume	140 µL
Reagent 3 Volume	50 µL
Reagent 4 Volume	50 µL
Reagent 5 Volume	50 µL
Diluent Volume	174 µL
Temperature	37 °C
Wavelength	340 and 383 nm
Type of Measurement	Bichromatic rate

Calibration

Assay Range	0.5 – 24 µg/dL [6.4 – 309 nmol/L] ^f
Calibration Material	T4 Calibrator, Cat. No. DC13
Calibration Scheme	5 levels, n = 2
Units	µg/dL [nmol/L] (µg/dL) x 12.87 = [nmol/L]
Typical Calibration Levels	1, 4.5, 8, 14, 20 µg/dL [13, 58, 103, 180, 257 nmol/L]
Calibration Frequency	Every 2 months for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ 190.0 C ₁ 170.0 C ₂ -3.2 C ₃ 23.0 C ₄ 0.5

f. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known total thyroxine concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of thyroxine in µg/dL [nmol/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.5 – 24 µg/dL
[6.4 – 309 nmol/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 24 µg/dL [309 nmol/L] should be repeated on dilution.
Manual Dilution: Make appropriate dilution with Reagent grade water to obtain result within assay range.
Enter dilution factor. Reassay. Resulting readout is corrected for dilution.
Autodilution: Refer to your Dimension® Operator's Guide.
Results less than 0.5 µg/dL [6.4 nmol/L] should be reported as "less than 0.5 µg/dL [6.4 nmol/L]".

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide for troubleshooting.
Venipuncture should occur prior to sulfasalazine administration due to the potential for falsely elevated results.
A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
5.1 µg/dL [66 nmol/L]	>0.8 µg/dL [10 nmol/L]
13.9 µg/dL [179 nmol/L]	>1.0 µg/dL [13 nmol/L]

United States Federal law restricts this device to sale by or on the order of a physician or other practitioner licensed by the laws of the State in which he practices, to use or order the use of the device.

Interfering Substances
Lipemia (Intralipid®) of 600 mg/dL [6.78 mmol/L] and above at a T4 concentration of 4.4 µg/mL [59 nmol/L] tripped a test report message; therefore the magnitude of the interference could not be determined.
Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Expected Values
Adult Population: 4.7 – 13.3 µg/dL [60 – 171 nmol/L]

Males: 4.5 – 12.1 µg/dL [58 – 156 nmol/L]
Females: 4.8 – 13.9 µg/dL [62 – 179 nmol/L]

This hospitalized reference population consisted of the following:
71 males
100 females

The reference interval was calculated non-parametrically and represents the central 95% of the population.

Pediatric Population:

Pediatric Age Group	Number of Samples	Reference Intervals	
		[µg/dL]	[nmol/L]
Infants (01 – 23 months)	75	6.6 – 13.4	85 – 173
Children (02 – 12 years)	186	5.8 – 11.8	75 – 152
Adolescents (13 – 20 years)	147	5.4 – 10.6	70 – 136

Reference intervals for the pediatric population (infants, children, and adolescents) were established for the Dimension® T4 method in accordance with the CLSI guideline EP28-A3c.⁵ Samples were collected prospectively from apparently healthy (euthyroid) pediatric subjects, using predefined inclusion criteria.

The reference interval for infants for the T4 method was calculated by a robust measure of location and spread as developed by Horn and Pesce.⁶ A non-parametric approach based on the CLSI guideline was used to establish the reference intervals for children and adolescents. The 2.5th and 97.5th percentiles of the distribution of values were calculated for each age group.

Each laboratory should establish its own reference interval for thyroxine as performed on the Dimension® system.

Specific Performance Characteristics^{a,h}

Material	Mean µg/dL [nmol/L]	Precision ⁱ	
		Standard Deviation (% CV)	
		Within-run	Between-day
QCS™ Control			
Serum Assayed			
Normal	7.0 [90]	0.25 [3.2] (3.5)	0.51 [6.6] (7.4)
Abnormal	18.7 [241]	0.48 [6.2] (2.6)	1.24 [16] (6.6)

g. The sensitivity or lowest detectable limit (lower assay range) is 0.5 µg/dL [6.4 nmol/L] as determined by dilution in linearity studies.
h. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).
i. Specimens at each level were analyzed in triplicate for 20 days. The within-run and between-day standard deviations were calculated by the analysis of variance method.
QCS™ is a trademark of Gifford, Irvine, CA 92714.

Method Comparison				
Regression Statistics ^j				
Comparative Method	Slope	Intercept µg/dL [nmol/L]	Correlation Coefficient	n
T4 method on the aca® discrete clinical analyzer	1.08	-0.04 [-0.5]	0.939	186

j. Model equation for regression statistics is: [results of Dimension® system] = slope x [comparative method results] + intercept.

Specificity

HIL Interference			
The T4 method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference."			
Substance Tested	Test Concentration SI Units	T4 Concentration µg/dL [nmol/L]	Bias ^k %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L] (monomer)	4.4 [56]	<10
Bilirubin (unconjugated)	5 mg/dL [85.5 µmol/L]	4.5 [58]	<10
	20 mg/dL [342 µmol/L]	4.5 [58]	11
Lipemia (Intralipid®)	200 mg/dL [2.26 mmol/L]	4.4 [56]	<10

k. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances do not interfere with the T4 method when present in serum or plasma in the amounts indicated. Systematic inaccuracies (bias) due to these substances are less than 10% at thyroxine concentrations of 7.5 µg/dL [96 nmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	1000 µg/mL	6620 µmol/L
Ampicillin	20 µg/mL	57 µmol/L
Diazepam	20 µg/dL	70 µmol/L
Digoxin	20 ng/mL	25.6 nmol/L
Ethanol	800 mg/dL	174 mmol/L
Gentamicin	16 µg/mL	29.4 µmol/L
Methimazole	10 µg/dL	876 nmol/L
Nortriptyline	1000 ng/mL	3797 nmol/L
Phenobarbital	80 µg/mL	344 µmol/L
Phenytoin	30 µg/mL	119 µmol/L
Propylthiouracil	1000 µg/dL	58744 nmol/L
Salicylate	100 mg/dL	7.24 mmol/L
Theophylline	100 µg/mL	555 µmol/L
Potassium oxalate	500 mg/dL	5 g/L

Recovery

The effect of various thyroxine congeners on the recovery of the T4 method is shown below:

Compound	Concentration	% Recovery ⁱ
D-Thyroxine	10 µg/dL [129 nmol/L]	210
L-Triiodothyronine	2.5 µg/dL [38 nmol/L]	97
Diiodotyrosine	10 µg/dL [231 nmol/L]	95
Monoiodotyrosine	10 µg/dL [326 nmol/L]	99
Triiodothyroacetic Acid (TRIAC)	5 µg/dL [80 nmol/L]	109
Tetraiodothyroacetic Acid (TETRAC)	5 µg/dL [67 nmol/L]	123

i. At a thyroxine concentration of approximately 8.5 µg/dL [109 nmol/L].

Analytical Sensitivity: 0.5 µg/dL [6.4 nmol/L]

The analytical sensitivity represents the lowest concentration of thyroxine that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations of the 0.00 µg/dL [nmol/L] level.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, Flex® and aca® are trademarks of Siemens Healthcare Diagnostics.
©2008 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
2. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
3. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
4. Tietz NW. *Clinical Guide to Laboratory Tests*. 4th ed., 2006, Saunders WB Co, Philadelphia, PA., pp 1050.
5. Clinical and Laboratory Standards Institute (formerly NCCLS). *Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition*. CLSI Document EP28-A3c. CLSI 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2010.
6. Horn PS, Pesce AJ, Reference Intervals. A User's Guide, Washington, DC: AACC Press; 2005.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

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Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare