Dimension® clinical chemistry system

Flex® reagent cartridge

EC02

See shaded sections: Updated information from 2016-03 version.


Issue Date 2019-04-01

Enzymatic Carbonate

Intended Use: The EC02 method for the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended to quantitatively measure total carbon dioxide in human serum or heparinized plasma.

Summary: The enzymatic carbonate (EC02) method for the Dimension® system employs a phosphoenolpyruvate carboxylase-malate dehydrogenase coupled enzymatic reaction and a stable analog of the cofactor NADH.

Principles of Procedure: The bicarbonate anion reacts with phosphoenolpyruvate in the presence of phosphoenolpyruvate carboxylase (PEPC) and Mg⁺⁺ to form oxaloacetate and inorganic phosphate (P_i). The oxaloacetate is reduced to malate by malate dehydrogenase (MDH) with simultaneous oxidation of the reduced form of an analog (aNADH)* of the cofactor, NADH.

The reduction in absorbance of aNADH is proportional to the total CO₂ concentration in the sample and is measured bichromatically at wavelengths of 405 nm (primary) and 700 nm (secondary).

*U.S. patent #5,801,006

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1 – 6	Liquid	Potassium Phosphoenolpyruvate	15.6 mM	
		Magnesium Chloride	19.0 mM	
		Sodium Oxamate	1.8 mM	
		3-Acetylpyridine Adenine Dinucleotide	1 mM	
		Phosphoenolpyruvate carboxylase	494 U/L	Bacterial
		Malate Dehydrogenase	5062 U/L	Porcine Heart
		Buffers and stabilizers		

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Represents the nominal value in the final reaction mixture.

Risk and Safety

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Contains sodium azide (<0.1%) as a preservative. Sodium azide can react with copper or lead pipes in drain lines to form explosive compounds. Dispose of properly in accordance with local regulations.

Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 2 days for wells 1 – 6

Specimen Collection and Handling: Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.¹

Follow the instructions provided with your specimen collection device for use and processing.²

Complete clot formation should take place before centrifugation.³

Samples should be analyzed as promptly as possible after collection and centrifugation of the blood in the unopened tube.⁴

Unopened, separated samples may be stored for 8 hours at room temperature, 2 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C or colder for 6 months.⁵

Total carbon dioxide concentration may be lowered by as much as 6 mmol/L when uncapped specimens are exposed to the air for one hour.⁴

Underfilling of vacutainers may account for low total carbon dioxide results of up to 3 mmol/L.⁶

Procedure

Materials Provided

EC02 Flex® reagent cartridge, DF137

Materials Required But Not Provided

EC02 Calibrator, DC137 or CHEM III Calibrator, Cat. No. DC130

Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® clinical chemistry system. For details of this processing, refer to your Dimension® Operator's Guide.

Test Conditions

Sample Size	5 µL
Reagent 1 Volume	100 µL
Temperature	37° C
Wavelength	405 and 700 nm
Type of Measurement	Bichromatic rate

Calibration

Assay range	5 – 45 mmol/L
Calibration Material	EC02 Calibrator, Cat. No. DC137 or CHEM III Calibrator, Cat. No. DC130
Calibration Scheme	3 levels, n = 3
Units	mmol/L
Typical Calibration Levels	0, 25, 50 mmol/L
Calibration Frequency	Every 90 days for any one lot
A new calibration is required	<ul style="list-style-type: none"> For each new lot of Flex® reagent cartridges After major maintenance or service, if indicated by quality control results As indicated in laboratory quality control procedures When required by government regulations

Assigned Coefficients

C₀ 0.000

C₁ 1.000

Quality Control

At least once each day of use, analyze two levels of a Quality Control material (QC) with known carbon dioxide concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of carbon dioxide in mmol/L using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 5 – 45 mmol/L

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 45 mmol/L should be repeated on dilution.

Manual Dilution: Make appropriate dilution with de-ionized water to obtain result within the assay range. Enter dilution factor. Reassay within 15 minutes. Resulting readout is corrected for dilution.

Results less than 5 mmol/L should be reported as "less than 5.0 mmol/L" instead of the numerical value.

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions.

Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
25 mmol/L	> 1.2 mmol/L
50 mmol/L	> 2.0 mmol/L

Interfering Substances

In rooms with poor ventilation, an open Flex® reagent cartridge well can absorb CO₂ which may cause results to be elevated by up to 30%.

Hemoglobin (hemolysate) of 1000 mg/dL [0.62 mmol/L]^c (monomer) decreases an EC02 result of 13 mmol/L by 21%.

Lipemia (Intralipid®) of 3000 mg/dL [33.9 mmol/L] decreases an EC02 result of 13 mmol/L by 16%.

c. Système International d'Unités [SI Units] are in brackets.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Expected Values: 21 – 32 mmol/L^d

Each laboratory should establish its own reference interval for total carbon dioxide as performed on the Dimension® system.

d. This reference interval applies to serum and plasma samples from 120 healthy males and females ages 18 and higher, analyzed using the Dimension® EC02 method.

Specific Performance Characteristics*

Material	Precision ^a		
	Mean mmol/L	Standard Deviation (% CV)	
		Within-run	Total
Dade® Moni-trol® TOTAL Control			
Level 1	13.1	0.5 (4.0)	0.7 (5.0)
Level 2	30.5	0.6 (2.0)	0.9 (3.0)
Plasma Pool	24.5	0.7 (3.0)	1.1 (4.5)
Serum Pool	24.4	0.7 (3.0)	1.1 (4.5)

e. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).

f. Reproducibility testing was done in accordance with the CLSI/NCCLS Approved Guideline for Evaluation of Precision Performance of Clinical Chemistry Devices (EP5-A, Feb. 1999).

g. Specimens at each level were analyzed in duplicate, once a day, for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

Dade® is a registered trademark of Siemens Healthcare Diagnostics, Deerfield, IL 60015.

Method Comparison

Regression Statistics^h

Comparative Method	Slope	Intercept mmol/L	Correlation Coefficient	n
Dimension® Rxl TC02	0.99	1.6	0.993	267 ⁱ

h. Model equation for regression statistics is: Result of Dimension® system = [slope x (comparative method result)] + Intercept.

i. The range of total carbon dioxide values in the correlation study was 5.8 – 43.4 mmol/L.

Specificity

HIL Interference

The ECO2 method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	ECO2 Concentration mmol/L	Bias ⁱ %
Hemoglobin (hemolysate)	500 mg/dL [0.31 mmol/L] (monomer)	13	<10
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	14	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	13	<10

j. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances do not interfere with the ECO₂ method when present in serum and plasma at the concentrations indicated. Systemic inaccuracies (bias) due to these substances are less than 10% at nominal total carbon dioxide concentrations of 26 to 36 mmol/L.

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1323 µmol/L
Albumin	6 g/dL	60 g/L
Amikacin	15 mg/dL	256 µmol/L
Ascorbic Acid	3 mg/dL	170 µmol/L
Caffeine	10 mg/dL	515 µmol/L
Carbamazepine	12 mg/dL	508 µmol/L
Chloramphenicol	25 mg/dL	774 µmol/L
Chloridiazepoxide	2 mg/dL	67 µmol/L
Chlorpromazine	5 mg/dL	157 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	10 mg/dL	396 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 75	2500 mg/dL	333 µmol/L
Diazepam	2 mg/dL	70 µmol/L
Digoxin	5 ng/mL	6.4 nmol/L
Erythromycin	20 mg/dL	273 µmol/L
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 µmol/L
Furosemide	2 mg/dL	61 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Heparin	8 U/mL	8000 U/L
Ibuprofen	40 mg/dL	1939 µmol/L
Immunoglobulin G	6 g/dL	60 g/L
Lidocaine	6 mg/dL	256 µmol/L
Lithium	3.5 mg/dL	5.04 mmol/L
Nicotine	2 mg/dL	123 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	442 µmol/L
Phenobarbital	15 mg/dL	646 µmol/L
Phenytoin	10 mg/dL	396 µmol/L
Primidone	10 mg/dL	458 µmol/L
Propoxyphene	0.4 mg/dL	12 µmol/L
Protein (low)	8 g/dL	80 g/L
Protein (high)	13.9 g/dL	139 g/L
Rheumatoid factor	550 IU/mL	550 IU/mL
Salicylic Acid	67 mg/dL	3.62 mmol/L
Theophylline	25 mg/dL	1388 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic acid	50 mg/dL	3467 µmol/L

Analytical Sensitivity: 5.0 mmol/L

The analytical sensitivity represents the lowest concentration of carbon dioxide that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations of the level 1 (0.0 mmol/L) ECO₂ Calibrator.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
2. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN-1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
3. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN-1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
4. Tietz NW. Textbook of Clinical Chemistry, W. B. Saunders Co., Philadelphia, PA 1999, pp 1084 (reference intervals) and pp1371 (specimen handling).
5. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests.3rd ed. Washington, DC: AACC Press, 2007: p171.
6. Herr RD, Swanson T. Serum Bicarbonate Declines with Sample Size in Vacutainer Tubes. Am J Clin Pathol. 1992; 97: 213-216.

**Symbols Key
Symbolschlüssel
Explication des Symboles
Interpretazione simboli
Clave de los Símbolos**

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
LOT	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
REF	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante
EC REP	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenido suficiente per "n" saggi / Contenido suficiente para <n> ensayos
IVD	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
CE	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
CONTENTS	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
LEVEL	Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFIS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare