Dimension® clinical chemistry system Heterogeneous Immunoassay Module

Flex® reagent cartridge

TSH

See shaded sections: Updated information from 2017-08 version

Issue Date 2019-04-30


Thyroid Stimulating Hormone

Intended Use: The TSH method for the Dimension® clinical chemistry system with the heterogeneous immunoassay module is an *in vitro* diagnostic test intended to quantitatively measure human thyroid stimulating hormone (TSH, thyrotropin) in human serum and plasma.

Summary: The TSH method is a sensitive, colorimetric immunoassay that measures thyroid stimulating hormone, a glycoprotein secreted by the anterior lobe of the pituitary gland. TSH stimulates the normal thyroid gland to synthesize and secrete thyroxine (T4) and triiodothyronine (T3).

Principles of Procedure: The TSH method is a one step enzyme immunoassay based on the "sandwich" principle. Sample is incubated with chromium dioxide particles coated with monoclonal antibodies specific for the intact TSH molecule, and a conjugate reagent [alkaline phosphatase (ALP) labeled monoclonal antibodies specific for the TSH beta subunit] to form a particle/TSH/conjugate sandwich. Unbound conjugate is removed by magnetic separation and washing. After separation and washing, the particle/TSH/conjugate sandwich is transferred to the cuvette where the sandwich bound alkaline phosphatase (ALP) triggers an amplification cascade. * ALP dephosphorylates synthetic flavin adenine dinucleotide phosphate (FADP) to give FAD. FAD binds to apo d-amino acid oxidase and converts it to active holo d-amino acid oxidase. Each molecule of holo d-amino acid oxidase then produces multiple molecules of hydrogen peroxide (H₂O₂) which, in the presence of horseradish peroxidase (HRP), convert 3, 5-dichloro-2-hydroxybenzenesulfonic acid (DCHBS) and 4-aminoantipyrine (4-AAP) to a colored product that absorbs at 510 nm. The color change measured is directly proportional to the concentration of TSH present in the patient sample.

* Technology licensed from London Biotechnology, Ltd., London, U.K.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1, 2	Tablets ^c	FADP	0.055 mg/mL	
		D-Proline	69 mg/mL	
3, 4	Tablets ^c	Apo d-amino acid oxidase	0.23 mg/mL	Porcine Kidney or Recombinant, E. coli
		Horseradish peroxidase	0.19 mg/mL	Horseradish
		DCHBS	3.9 mg/mL	
		4-AAP	0.94 mg/mL	
6	Tablets ^c	Antibody-CrO ₂	1.3 mg/mL ^d	Mouse
7	Liquid	Buffer and Stabilizers		
8	Liquid	TSH Ab-Alkaline Phosphatase	0.02 mg/mL ^d	Mouse/bovine
		Buffer and Stabilizers		

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value in hydrated cartridge.

c. Tablets contain excipients, buffers, and stabilizers.

d. Antibody titer and conjugate activity vary from lot to lot.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501
Warning!
May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.
Contains: 2-chloroacetamide

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Contains sodium azide (< 0.1%) as a preservative. Sodium azide can react with copper or lead pipes in drain lines to form explosive compounds. Dispose of properly in accordance with local regulations.

Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

Reagent Preparation: Hydrating, diluting, and mixing are automatically performed by the Dimension® system.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 4
7 days for wells 6 – 8

Specimen Collection and Handling: Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.¹

Follow the instructions provided with your specimen collection device for use and processing.²

Specimens should be free of particulate matter. To prevent the appearance of fibrin in serum samples, complete clot formation should take place before centrifugation.³ If clotting time is increased due to thrombolytic or anticoagulant therapy, the use of plasma specimens will allow for faster sample processing and reduce the risk of particulate matter.

Blood collected in the presence of oxalate potentially can cause clumping of the chrome particles and should not be used.

Separated samples are stable for 1 day at ambient temperature and 7 days at 2 – 8 °C. For longer storage, samples should be frozen at -20 °C for 1 month.^{4,5} Avoid repeated freezing and thawing.

Frozen plasma samples with insufficient anticoagulants potentially can cause clumping of chrome particles and should not be used.

Procedure

Materials Provided

TSH Flex® reagent cartridge, Cat. No. RF412

Materials Required But Not Provided

Reaction Vessels, Cat. No. RXV1A
Chemistry Wash, Cat. No. RD701
Probe Cleaner, Cat. No. RD702
Sample Probe Cleaner, Cat. No. RD703
Sample Diluent, Cat. No. 791092901
Thyroid Calibrator, Cat. No. RC410
Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, separation, processing and printing of results are automatically performed by the Dimension® system with the heterogeneous immunoassay module. For details of this processing, refer to your Dimension® Operator's Guide.

Test Conditions

Reaction Vessel		
Sample Size	60 µL	
Antibody-CrO ₂	20 µL	
Antibody-ALP	50 µL	
Incubation Temperature	42 °C	
Incubation Period	8.5 minutes	
Cuvette		Blanking
Transfer Volume	65 µL	0 µL
FADP Reagent Volume	30 µL	30 µL
APD Reagent Volume	30 µL	30 µL
Diluent Volume	325 µL	390 µL
Temperature	37.0 °C	
Wavelength	510 and 700 nm	
Type of Measurement	Bichromatic rate	

Calibration

Assay Range	0.01 – 50.00 µIU/mL [mIU/L] ^e
Calibration Material	Thyroid Calibrator, Cat. No. RC410
Calibration Scheme	Level 1, n = 4 Level 2 & 3, n = 3 Level 4 & 5, n = 2
Units	µIU/mL [mIU/L]
Typical Calibration Levels	0.00, 1.00, 4.00, 20.00 and 55.00 µIU/mL [mIU/L]
Calibration Frequency	Every 60 days for any lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ -3300 C ₁ 16000 C ₂ 2.0 C ₃ 53 C ₄ 0.5

e. Système International d' Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a quality control (QC) material with known TSH concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of TSH in $\mu\text{IU/mL}$ [mIU/L].

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.01 – 50.00 $\mu\text{IU/mL}$ [mIU/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 50.00 $\mu\text{IU/mL}$ [mIU/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilution with Sample Diluent to obtain result within assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution: Refer to your Dimension® Operator's Guide.

Results less than 0.01 $\mu\text{IU/mL}$ [mIU/L] should be reported as "less than 0.01 $\mu\text{IU/mL}$ [mIU/L]".

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

Performance of this assay has not been established with neonatal specimens.

Patient samples may contain heterophilic antibodies that could react in immunoassays to give falsely elevated or depressed results. This assay has been designed to minimize interference from heterophilic antibodies.⁸ Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed.⁷

As with any immuno-recognition measurement of a peptide, extremely rare genetic variants may exhibit varying degrees of detection.

A system malfunction may exist if the following 5-test precision is observed:

TSH Concentration	SD
0.97 $\mu\text{IU/mL}$ [mIU/L]	> 0.05 $\mu\text{IU/mL}$ [mIU/L]
3.98 $\mu\text{IU/mL}$ [mIU/L]	> 0.16 $\mu\text{IU/mL}$ [mIU/L]

Expected Values: 0.34 – 4.82 $\mu\text{IU/mL}$ [mIU/L]

The TSH reference interval was determined from populations of apparently healthy individuals. A total of 293 samples were tested. This group was comprised of 124 males, ages 17 – 81, and 169 females, ages 19 – 81. The reference interval for TSH was calculated nonparametrically and represents the central 95% of this population.

Each laboratory should establish its own reference interval for TSH as performed on the Dimension® system.

Specific Performance Characteristics¹

Material	Precision ^{a,b}		Standard Deviation (% CV)
	Mean $\mu\text{IU/mL}$ [mIU/L]	Within-run	Total
Dade® Immunoassay Control			
Level 1	0.77	0.03 (4.0)	0.05 (6.2)
Level 2	9.41	0.21 (2.2)	0.45 (4.8)
Level 3	40.02	0.91 (2.3)	2.49 (6.2)
Serum Pool			
Level 1	0.08	0.01 (6.9)	0.01 (12.4)
Level 2	1.59	0.05 (3.2)	0.08 (5.2)

- f. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).
- g. Reproducibility testing was done in accordance with the CLSI/NCCLS Tentative Guideline for User Evaluation of Precision Performance of Clinical Chemistry Devices (EP5-T2, 1992).
- h. Specimens at each level were analyzed in duplicate, once a day, for 20 days. The within-run and total standard deviations were calculated by analysis of variance method.

Method Comparison

Regression Statistics¹

Comparative Method	Slope	Intercept $\mu\text{IU/mL}$ [mIU/L]	Correlation Coefficient	n ¹
AxSYM® Ultrasensitive TSH	0.96	-0.23	0.99	99

- i. Model equation for regression statistics is: Result of Dimension® analyzer = [Slope x (comparative method result)] + Intercept.
- j. The range of TSH values in the correlation study was 0.02 – 44.70 $\mu\text{IU/mL}$ [mIU/L].

AxSYM® is a registered trademark of Abbott Laboratories, Abbott Park, IL 60064.

Specificity

Non-Interfering Substances

The following substances do not interfere with the TSH method when present in serum and plasma at the concentrations indicated. Systemic inaccuracies (bias) due to these substances are less than 5% at a TSH level of 4.65 $\mu\text{IU/mL}$ [mIU/L]:

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1323 $\mu\text{mol/L}$
Amikacin	15 mg/dL	256 $\mu\text{mol/L}$
Ampicillin	5 mg/dL	143 $\mu\text{mol/L}$
Ascorbic acid	3 mg/dL	170 $\mu\text{mol/L}$
Billirubin	20 mg/dL	342 $\mu\text{mol/L}$
Caffeine	10 mg/dL	515 $\mu\text{mol/L}$
Carbamazepine	12 mg/dL	508 $\mu\text{mol/L}$

Chloramphenicol	25 mg/dL	774 $\mu\text{mol/L}$
Chlordiazepoxide	2 mg/dL	67 $\mu\text{mol/L}$
Chlorpromazine	5 mg/dL	157 $\mu\text{mol/L}$
Cimetidine	10 mg/dL	396 $\mu\text{mol/L}$
Creatinine	30 mg/dL	2652 $\mu\text{mol/L}$
Dextran 75	2500 mg/dL	333 $\mu\text{mol/L}$
Diazepam	2 mg/dL	70 $\mu\text{mol/L}$
Digoxin	5 ng/mL	6.4 nmol/L
Erythromycin	20 mg/dL	273 $\mu\text{mol/L}$
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 $\mu\text{mol/L}$
Furosemide	2 mg/dL	61 $\mu\text{mol/L}$
Gentamicin	12 mg/dL	251 $\mu\text{mol/L}$
Hemoglobin	1000 mg/dL	0.62 mmol/L (monomer)
Heparin (sodium)	8 U/mL	8000 U/L
Ibuprofen	40 mg/dL	1939 $\mu\text{mol/L}$
Lidocaine	6 mg/dL	256 $\mu\text{mol/L}$
Lipemia	3000 mg/dL	33.9 mmol/L (triglyceride)
Lithium chloride	3.5 mEq/L	5.07 mmol/L
Nicotine	2 mg/dL	123 $\mu\text{mol/L}$
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	442 $\mu\text{mol/L}$
Phenobarbital	15 mg/dL	646 $\mu\text{mol/L}$
Phenytoin	10 mg/dL	396 $\mu\text{mol/L}$
Primidone	10 mg/dL	458 $\mu\text{mol/L}$
Propoxyphene	0.4 mg/dL	12 $\mu\text{mol/L}$
Protein, human albumin	6 g/dL	60 g/L
Protein, human IgG	6 g/dL	60 g/L
Salicylic acid	50 mg/dL	3.62 mmol/L
Theophylline	25 mg/dL	1388 $\mu\text{mol/L}$
Urea	500 mg/dL	83.3 mmol/L
Uric acid	20 mg/dL	1.2 mmol/L
Valproic acid	50 mg/dL	3467 $\mu\text{mol/L}$

The following substances do NOT interfere with the TSH assay when tested at TSH level of 5.16 $\mu\text{IU/mL}$ [mIU/L] at the concentrations indicated:

Substance	Concentration
FSH	1000 mIU/mL [IU/L]
LH	1000 mIU/mL [IU/L]
hCG	200,000 mIU/mL [IU/L]

Hook Effect

One step sandwich immunometric assays are susceptible to a high-dose "hook effect", where an excess of antigen prevents simultaneous binding of the capture and detection antibodies to a single analyte molecule.⁸ Such samples must be diluted and reassayed prior to reporting results (see dilution information). The TSH method shows no hook effect up to 300 $\mu\text{IU/mL}$ [mIU/L].

Recovery

Known amounts of TSH were added to human serum and plasma samples with baseline TSH values of 0.7 – 2.76 $\mu\text{IU/mL}$ [mIU/L]. The sample concentrations were measured and the calculated percent recovery ranged from 90% to 110% with a mean recovery of 100%.

% Recovery: $\frac{\text{Value obtained} - \text{Baseline}}{\text{Amount Added}} \times 100$

Ten serum samples with TSH values from 14.9 – 53.7 $\mu\text{IU/mL}$ [mIU/L] were serially diluted with aca® plus/Dimension® Sample Diluent and assayed in triplicate by the Dimension® TSH assay.

The observed values were plotted against the expected values for each sample and slopes of 0.97 to 1.02 were obtained. The linear relationship also appeared acceptable by visual evaluation.

Analytical Sensitivity: 0.01 $\mu\text{IU/mL}$ [mIU/L]

The analytical sensitivity represents the lowest concentration of TSH that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations of the level 1 (0 $\mu\text{IU/mL}$ [mIU/L]) Thyroid Calibrator.

Functional Sensitivity: 0.04 $\mu\text{IU/mL}$ [mIU/L]

The Dimension® TSH functional sensitivity (at 20% CV) was determined using the study protocol recommended by the National Academy of Clinical Biochemistry.⁹

Conducted over 8 weeks, the study involved testing two lots of reagents, each calibrated with a different lot of calibrator. The first reagent lot was used during the first four weeks of testing and the second lot of reagent was used during the last four weeks of testing.

Ten patient pools whose values span the TSH assay range were assayed using one replicate per sample, twice a week for 8 weeks. All tests were run in random order. The inter-assay coefficient of variation was plotted by nonlinear regression versus concentration. Based on this study, the Dimension® TSH functional sensitivity was determined to be 0.04 $\mu\text{IU/mL}$ [mIU/L].

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics

©2009 Siemens Healthcare Diagnostics
All rights reserved.

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Dimension® clinical chemistry system Heterogeneous Immunoassay Module

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Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía:

1. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
2. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
3. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
4. Tietz NW. *Clinical Guide to Laboratory Tests*. 4th ed., 2006, Saunders WB Co, Philadelphia, PA., pp 1040.
5. Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 3rd ed. Washington, DC: AACC Press, 2007: p 867.
6. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG, *Clin Chem* 1992; 38:1737-1742.
7. Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. *Clin Chem* 1999; 45:7,942-956.
8. Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the “hook effect” in two-site immunoradiometric assays, *Anal Biochem* 1982;127:308-15
9. Standards of Laboratory Practice, Laboratory support for the diagnosis & monitoring of thyroid disease, National Academy of Clinical Biochemistry, 2002.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFIS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare