Dimension® clinical chemistry system

Flex® reagent cartridge

UCFP

See shaded sections: Updated information from 2017-12 version.

Issue Date 2019-04-01

Urinary/Cerebrospinal Fluid Protein


Intended Use: The UCFP method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of total protein in human urine and cerebrospinal fluid.

Summary: The UCFP method allows direct quantitation of proteins in urine and cerebrospinal fluid specimens. The UCFP method is an adaptation of pyrogallol red-molybdate method by Y. Fujita, I. Mori and S. Kitano.¹

Measurement of the protein content in urine is used in diagnosis and treatment of kidney diseases.

Measurement of the protein content in cerebrospinal fluid is used in the diagnosis and treatment of central nervous system diseases.

Principles of Procedure: In the reaction sequence, pyrogallol red combines with sodium molybdate to form a red complex with maximum absorbance at 470 nm. The protein in the sample reacts with this complex in acid solution to form a bluish-purple colored complex, which absorbs at 600 nm. The absorbance at 600 nm is directly proportional to the concentration of protein in the sample. The analyte concentration is determined by calculation using a logit curve fit on a previously stored calibration curve.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1 – 5	Liquid	Pyrogallol red in methanol Sodium molybdate Stabilizers and Surfactants	0.20 mM 0.35 mM
7 – 8 ^c	Liquid	NaOH, 0.5N	

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value in final reaction mixture.

c. Sodium hydroxide is used as a probe cleaning solution and is not used in the reaction.

Risk and Safety:

H290, H302, H319, H315, H371
P280, P301 + P312, P305 + P351 + P338, P501

Warning!

May be corrosive to metals. Harmful if swallowed. Causes serious eye irritation. Causes skin irritation. May cause damage to organs.

Wear protective gloves/protective clothing/eye protection/face protection.
IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: Methanol; Sodium Hydroxide

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 5 days for wells 1 – 5, 7, 8

Specimen Collection and Handling: Normal procedures for collecting urine and cerebrospinal fluid may be used for samples to be analyzed by this method.

Specimens stored at 4 °C with no additives are stable for at least three days.²

Hemolyzed samples should be avoided since hemolysis increases UCFP results (at 25 mg/dL hemoglobin).

Blood collection tubes containing sodium citrate should not be used.

Random urine specimens may be used but timed 24-hr specimens are preferred. No preservative is required during 24-hr collection, but thereafter urine aliquots should be stored at 2 – 4 °C for < 72 hours or frozen at -20 °C for up to 1 year. Urine specimens must be free of any particulate matter before analysis.^{3,4}

Cerebrospinal fluid specimens should be collected with care to avoid contamination with plasma proteins. Blood present in the cerebrospinal fluid invalidates the protein values since it reflects contamination with plasma proteins. Analyze fresh or store at 4 °C for < 72 hours or frozen for up to 6 months.⁴

Procedure

Materials Provided

UCFP Flex® reagent cartridge, Cat. No. DF26

Materials Required But Not Provided

UCFP Calibrator, Cat. No. DC45

Quality Control Materials

Test Steps

Sampling,⁴ reagent delivery, mixing, processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

d. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus the dead volume. Precise container filling is not required.

Test Conditions

Sample Size	10 µL
Reagent Volume	350 µL
Diluent Volume	50 µL
Temperature	37 °C
Wavelength	600 and 700 nm
Type of Measurement	Bichromatic endpoint

Calibration

Assay Range	6 – 250 mg/dL [60 – 2500 mg/L] ^a
Calibration Material	UCFP Calibrator, Cat. No. DC45
Calibration Scheme	5 levels, n = 2
Units	mg/dL [mg/L] (mg/dL x 10) = [mg/L]
Typical Calibration Levels	6.0, 30.0, 60.0, 135.0, 270.0 mg/dL [60, 300, 600, 1350, 2700 mg/L]
Calibration Frequency	Every 2 months for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C_0 -193.89 C_1 1337.5 C_2 -2.595 C_3 462.58 C_4 0.5

e. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known UCFP concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of UCFP in mg/dL [mg/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 6 – 250 mg/dL [60 – 2500 mg/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 250 mg/dL [2500 mg/L] should be repeated on dilution.

Manual Dilution:

Make appropriate dilution with Reagent grade water to obtain result within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Automated Urine (AD): For urine and CSF specimens: Refer to your Dimension® Operator's Guide.

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
32 mg/dL	1.2 mg/dL
142 mg/dL	2.0 mg/dL

Interfering Substances

Samples containing amikacin, gentamicin, kanamycin, and tobramycin should be avoided since these substances falsely increase UCFP results.

Neomycin sulfate at 15 µg/mL increases UCFP results by 11% and at 7.5 µg/mL the interference is less than 5%.

Expected Values

Urine < 11.9 mg/dL, < 149.1 mg/day [< 119 mg/L]

Spinal Fluid³ 15 – 45 mg/dL [150 – 450 mg/L]

This urine reference population consisted of 115 adults (50 females, 65 males). The reference interval was calculated non-parametrically and represents 95% of the population.

Each laboratory should establish its own reference interval for urine and cerebrospinal fluid protein as performed on the Dimension® system.

Specific Performance Characteristics¹

Material	Precision ^a		Standard Deviation (% CV)
	Mean mg/dL		Within-run Total
Bio-Rad Urine Control			
Level 1	20.4	0.60 (2.9)	1.51 (7.4)
Level 2	69.8	0.62 (0.9)	1.19 (1.7)
Ciba-Corning			
CSF/LCR Control	26.9	0.58 (2.2)	1.34 (5.0)
UP/CFP Calibrator	142.3	0.88 (0.6)	1.24 (0.9)

- f. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).
- g. Control product samples were analyzed in duplicate twice a day for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

Method Comparison					
Comparative Method	Regression Statistics ^b		Correlation Coefficient	n	Range mg/dL
	Slope	Intercept			
UP method on the aca® analyzer Urine	1.07	5.0	0.993	130	6 – 240
CFP method on the aca® analyzer Cerebrospinal fluid	0.97	6.6	0.986	142	10 – 200

- h. Model equation for regression statistics is: results of Dimension® system = [slope x (comparative method results)] + intercept.

Specificity

Non-Interfering Substances

The following substances have shown less than 10% interference on the UCFP results at the concentrations indicated.

Substance	Test Concentration
Acetaminophen	20 mg/dL
Acetazolamide	135 µg/mL
Acetylsalicylic acid	50 mg/dL
Ammonia	150 mg/dL
Amobarbital	10 mg/dL
Ampicillin	5 mg/dL
Amphotericin B	4 µg/dL
Ascorbic acid	200 mg/dL
Bilirubin	20 mg/dL
Boric acid	300 mg/dL
Caffeine	10 mg/dL
Calcium	100 mg/dL
Carbamazepine	12 mg/dL
Citric acid	25 mg/dL
Chloride	100 mmol/L
Chloramphenicol	25 mg/dL
Chlordiazepoxide	2 mg/dL
Chlorpromazine	4 mg/dL
Cimetidine	10 mg/dL
Codeine	10 mg/dL
Copper	5 mg/dL
Creatinine	30 mg/dL
Cyclosporine A	1500 ng/mL
Dextran 75	2500 mg/dL
Diazepam	2 mg/dL
Digoxin	3.0 ng/mL
EDTA	1 g/L
Erythromycin	20 mg/dL
Ethanol	350 mg/dL
Ethosuximide	30 mg/dL
Furosemide	2 mg/dL

Glucose	5000 mg/dL
Heparin	1000 U/L
Ibuprofen	40 mg/dL
Imipramine	25 µg/mL
Iron	10 mg/dL
Lactic acid	25 mmol/L
Lidocaine	6 mg/dL
Lithium	3.5 mg/dL
Magnesium	25 mg/dL
Mannitol	19.4 µg/mL
Methicillin	180 µg/dL
Metrizamide	20 mg/dL
Methotrexate	9 mg/dL
Morphine	800 ng/mL
Nicotine	2 mg/dL
Penicillin V	80 mg/dL
Pentobarbital	10 mg/dL
Phenacetin	250 µg/mL
Phenobarbital	15 mg/dL
Phenytoin	100 µg/mL
Phosphorus	200 mg/dL
Primidone	10 mg/dL
Potassium	17 mmol/L
Procaine	12 µg/mL
Propoxyphene	0.4 mg/dL
Salicylic acid	50 mg/dL
Secobarbital	5 mg/dL
Sodium	600 mg/dL
Sodium fluoride	1 mg/dL
Sulfadiazine	200 mg/dL
Sulfisoxazole	200 mg/dL
Sulindac	200 mg/dL
Streptomycin	6 mg/dL
Tartaric acid	25 mg/dL
Theophylline	250 µg/mL
Thymol	150 mg/dL
Tolbutamide	1200 µg/mL
Urea	500 mg/dL
Uric acid	7.5 mg/dL
Valproic acid	50 mg/dL
Vancomycin	100 µg/mL

Recovery

Recovery of various proteins by the UCFP method is shown below.

Protein	% Recovery
Albumin	89
Gamma Globulin	57
Transferin	81
Orosomucoid	44
β-Lipoprotein	63
Lysosome	5

Analytical Sensitivity: 6 mg/dL [60 mg/L]

The analytical sensitivity represents the lowest concentration of urine or cerebrospinal fluid protein that can be distinguished from zero. This sensitivity is defined as the low assay range.

Symbols Key: See Adjacent Panel

Bibliography: See Adjacent Panel

Dimension®, aca® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2009 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía:

1. Fujita Y, Mori I, Kitano S. Color reaction between Pyrogallol Red – Molybdenum (VI) complex and protein. *Bunseki Kagaku* 1983; 32:379-386.
2. Gadsden RH. et al. An evaluation of the urinary protein method for the aca® discrete analyzer, DuPont, Wilmington, DE (October 1982).
3. Burtis CA, Ashwood ER. *Tietz Textbook of Clinical Chemistry*, WB Saunders Co., Philadelphia, PA, 1994: pp 58-69 (specimen collection), pp 75-76 (separation and storage of specimens), pp 717-730 (proteins in urine and cerebrospinal fluid).
4. Tietz NW. *Clinical Guide to Laboratory Tests*. 4th ed., 2006, Saunders WB Co, Philadelphia, PA., pp 216-218.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

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Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare