Dimension® clinical chemistry system

Flex® reagent cartridge

CRE2

See shaded sections: Updated information from 2016-12 version.

Issue Date 2019-04-16

Creatinine


Intended Use: The CRE2 method is an *in vitro* diagnostic test for the quantitative measurement of creatinine in human serum, plasma, and urine on the Dimension® clinical chemistry system. Creatinine measurements are used in the diagnosis and treatment of certain renal disease, in monitoring renal dialysis, and as a calculation basis for measuring other urine analytes.

Summary: The creatinine method employs a modification of the kinetic Jaffe reaction. This method has been reported to be less susceptible than conventional methods to interference from non-creatinine, Jaffe-positive compounds.¹ Creatinine is generally regarded as the most useful endogenous substance to measure for the assessment of kidney function.¹

In most chronic kidney diseases (CKD), excretion, endocrine and metabolic functions decline altogether. The glomerular filtration rate (GFR) is used to measure the filtering capacity of the kidneys and is widely accepted as the best index of kidney function. There are equations to estimate GFR (eGFR_{creat}) from serum/plasma creatinine that are applicable for most clinical circumstances for diagnosis, staging and tracking the progression of CKD. For the initial assessment of CKD the recommendation is to report an eGFR_{creat} in addition to a serum/plasma creatinine result.²

Additionally, eGFR_{creat} may be utilized for estimating kidney function and to establish effective therapeutic dosing levels of pharmaceuticals. Calculated creatinine clearance (CRCL) values may also be used for drug-dosing levels when limitations for use of eGFR_{creat} occur.³

Principles of Procedure: The CRE2 method uses a modified kinetic Jaffe technique. In the presence of a strong base such as NaOH, picric acid reacts with creatinine to form a red chromophore. The rate of increasing absorbance at 510 nm due to the formation of this chromophore is directly proportional to the creatinine concentration in the sample and is measured using a bichromatic (510, 600 nm) rate technique. Bilirubin is oxidized by potassium ferricyanide⁴ to prevent interference.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1 – 3 (Reagent 1)	Liquid	Lithium Picric acid	125 mM
4 – 6 (Reagent 2)	Liquid	NaOH K ₃ Fe(CN) ₆	2000 mM 2.7 mM

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

Risk and Safety:

H290, H314
P280, P301 + P310 + P331, P303 + P361 + P353 + P310, P305 + P351 + P338,
P390, P501

Danger!

May be corrosive to metals. Causes severe skin burns and eye damage.

Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. Do NOT induce vomiting. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Immediately call a POISON CENTER or doctor/physician. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Absorb spillage to prevent material damage. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: Sodium hydroxide

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use.

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8°C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6

Specimen Collection and Handling

Recommended specimen types: serum, plasma (lithium heparin) and urine.

Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁵

Follow the instructions provided with your specimen collection device for use and processing.⁶

For serum, complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.⁷

Separated serum and plasma specimens are stable for 24 hours at room temperature, 7 days at 2 – 8°C. For longer storage, specimens may be frozen for up to 3 months at -20°C or colder.⁸

Urine Creatinine measurements require no special patient preparation for specimen collection. Urine (random or 24 hour) specimens should be free of particulate matter.

Collect random and timed urine specimens using recommended procedures for collection, transportation, and preservation of urine specimens.⁹ Random or 24 hour urine collections require no preservatives. Specimens previously preserved with 6N HCl or Boric Acid are acceptable.

Urine (random or 24 hour collections) should be stored at 2 – 8°C and analyzed within 4 days. Freeze for longer storage.¹⁰

Thawed frozen specimens or urine specimens which are turbid must be clarified by centrifugation prior to testing.

The purpose of specimen handling and storage information is to provide guidance to users; however, users may validate their own procedures for handling and storing patient samples.

Procedure

Materials Provided

CRE2 Flex® reagent cartridge, Cat. No. DF33B

Materials Required But Not Provided

CHEM I Calibrator, Cat. No. DC18C

Quality Control Materials

Test Steps

Sampling¹¹, reagent delivery, mixing, and processing are automatically performed by the Dimension® clinical chemistry system. For details of this processing, refer to your Dimension® Operator's Guide.

c. The sample container must contain sufficient quantity to accommodate the sample volume plus dead volume.

Precise container filling is not required.

Test Conditions

Sample Volume (delivered to the cuvette)	20 µL
Reagent 1 Volume	74 µL
Reagent 2 Volume	18 µL
Diluent Volume	258 µL
Temperature	37°C
Reaction Time (reagent 1 to final read)	2.0 minutes
Wavelength	510 nm and 600 nm
Type of Measurement	Bichromatic rate

Calibration

Assay Range for serum and plasma	0.15 – 20.00 mg/dL [13 – 1768 µmol/L] ¹⁴
Assay Range for urine	13.00 – 400.00 mg/dL [1149 – 35360 µmol/L]
Calibration Material	Dimension® CHEM I CAL, Cat. No. DC18C
Calibration Scheme	Three levels (n = 3)
Units	mg/dL [µmol/L]
Conversion Factor	(mg/dL x 88.4) = [µmol/L]
Decimal places	2 for mg/dL and 0 for [µmol/L]
Typical Calibration Levels	Level 1: 0.00 mg/dL [0 µmol/L] Level 2: 11.00 mg/dL [972 µmol/L] Level 3: 22.00 mg/dL [1945 µmol/L]
Calibration Frequency	Every 90 days for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ -0.3866 C ₁ 0.0823

d. Système International d'Unités [SI Units] are in brackets.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of a Quality Control (QC) material with known creatinine concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the concentration of creatinine in mg/dL [µmol/L] using the calculation scheme described in your Dimension® Operator's Guide. To correct for non-specific reaction due to the presence of proteins, serum and plasma measurements are corrected by -0.05 mg/dL [-4.42 µmol/L].

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Serum and Plasma

Analytical Measurement Range (AMR): 0.15 – 20.00 mg/dL [13 – 1768 µmol/L]

Urine

Analytical Measurement Range (AMR): 13.00 – 400.00 mg/dL [1149 – 35360 µmol/L]

- Automated Urine Dilution (AUD): The automated urine dilution sample volume is 15 µL. Urine samples are automatically diluted 1:20 with the Dimension® System water (1 part urine sample and 19 parts water) resulting in an analytical measurement range of 13.00 – 400.00 mg/dL [1149 – 35360 µmol/L].

The serum, plasma, and urine AMRs are the ranges of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Serum and plasma samples with results in excess of 20.00 mg/dL [1768 µmol/L] are reported as "Above Assay Range" and should be repeated on dilution.

- Autodilution (AD): The autodilute sample volume is 10 µL and recommended dilution factor is 1:2 for serum and plasma. This extends the serum and plasma reportable range to 40 mg/dL [3536 µmol/L]. Refer to your Dimension® Operator's Guide.

Urine samples with results in excess of 400.00 mg/dL [35360 µmol/L] are reported as "Above Assay Range" and should be repeated on dilution.

- Autodilution (AD): Autodilution is not available for urine samples.
- Manual Dilution: Urine samples can be manually diluted with Reagent grade water to obtain results within the analytical measurement range. The recommended dilution factor is 1:2. This extends the urine reportable range to 800.0 mg/dL [70720 µmol/L]. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.

Serum and plasma samples with results less than 0.15 mg/dL [13 µmol/L] will produce an instrument flag "below assay range" and should be reported as "less than 0.15 mg/dL [13 µmol/L]".

Urine samples with results less than 13.00 mg/dL [1149 µmol/L] will produce an instrument flag "below assay range" and should be reported as "less than 13.00 mg/dL [1149 µmol/L]".

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in CRE2 results. Refer to your Dimension® Operator's Guide for the meaning of report flags and comments. Any creatinine result containing flags and/or comments should be addressed according to your laboratory's procedure manual.

A system malfunction may exist if the following 5 test precision is observed:

Concentration	SD
1.0 mg/dL [88 µmol/L]	> 0.10 mg/dL [9 µmol/L]
18.6 mg/dL [1644 µmol/L]	> 0.50 mg/dL [44 µmol/L]

Interfering substances

The CRE2 method was evaluated for interference according to CLSI EP07-A2.¹¹ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	Creatinine mg/dL [µmol/L]	Bias* %
Serum/plasma	Acetone	9.375 mg/dL [0.016 mmol/L]	1.50 [133]
		18.75 mg/dL [0.032 mmol/L]	1.50 [133]
		37.5 mg/dL [0.065 mmol/L]	5.00 [442]
		75 mg/dL [0.13 mmol/L]	5.00 [442]
		20 mg/dL [342 µmol/L]	1.50 [133]
Bilirubin (unconjugated)	Cefoxitin	3.75 mg/dL [87.8 µmol/L]	1.50 [133]
		5.00 mg/dL [117 µmol/L]	1.50 [133]
		5.00 mg/dL [117 µmol/L]	5.00 [442]
		10 mg/dL [250 µmol/L]	1.50 [133]
		12.5 mg/dL [312 µmol/L]	1.50 [133]
Cephalothin		15 mg/dL [375 µmol/L]	5.00 [442]
		400 mg/dL [22.0 mmol/L]	1.50 [133]
		500 mg/dL [27.6 mmol/L]	1.50 [133]
		600 mg/dL [33.1 mmol/L]	5.00 [442]
		1500 mg/dL [17.0 mmol/L]	1.50 [133]
Lipemia (Intralipid®)	Pyruvate	1.32 mg/dL [0.150 mM]	1.50 [133]
		5.26 mg/dL [0.6 mM]	5.00 [442]
		10.5 mg/dL [1.2 mM]	5.00 [442]
		2500 mg/dL [28.3 mmol/L]	1.50 [133]
		3000 mg/dL [34.0 mmol/L]	1.50 [133]
Triglycerides		3000 mg/dL [34.0 mmol/L]	5.00 [442]
			< 10
Urine	Ascorbic Acid	0.1 g/dL [3.0 µmol/L]	40.00 [3536]
		0.15 g/dL [4.5 µmol/L]	40.00 [3536]
		0.2 g/dL [6.0 µmol/L]	175.00 [15470]
		0.3 g/dL [9.0 µmol/L]	175.00 [15470]
			+12

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

* Analyte results should not be corrected based on this bias.

Expected Values:

Serum and Plasma

Males: 0.70 – 1.30 mg/dL [62 – 115 µmol/L]¹²

Females: 0.55 – 1.02 mg/dL [49 – 90 µmol/L]¹³

Urine

Males: 0.95 – 2.49 g/24 hr [8.4 – 22.0 mmol/24 hr]¹⁴

Females: 0.60 – 1.80 g/24 hr [5.3 – 15.9 mmol/24 hr]¹⁵

The reference intervals were transferred from literature references in accordance with CLSI EP28-A3c¹⁶.

Each laboratory should establish its own expected values for CRE2 as performed on the Dimension® clinical chemistry system.

Specific Performance Characteristics

The following data represent typical performance for the Dimension® clinical chemistry system and were collected on the Dimension® EXL™ 200 Systems.

Precision^{17, *}

Material	Mean mg/dL [µmol/L]	Standard Deviation (%CV)	
		Repeatability	Within-Lab
Serum			
	Serum pool 1	1.32 [116.7]	0.04 [3.5] (3.0)
	Serum pool 2	15.79 [1395.8]	0.19 [16.8] (1.2)
BioRad® Multiquel®			
	Level 1	0.71 [62.8]	0.03 [2.7] (4.7)
	Level 2	1.79 [158.2]	0.04 [3.5] (2.1)
Urine			
	Level 3	7.04 [622.3]	0.07 [6.2] (1.0)
			0.09 [8.0] (1.3)
Urine			
	Urine pool 1	39.31 [3475.0]	1.52 [134.4] (3.9)
	Urine pool 2	339.56 [30017.1]	4.28 [378.4] (1.3)
			1.53 [135.3] (3.9)
BioRad® Liquichek™			
	Level 1	62.28 [5505.6]	0.61 [53.9] (1.0)
	Level 2	142.80 [12623.5]	1.56 [137.9] (1.1)
			3.39 [299.7] (2.4)

e. CLSI EP05-A2¹⁷ was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

BioRad® is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

Multiquel® is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

Liquichek™ is a trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

Method Comparison ¹⁸ Regression Statistics ⁱ				
Comparative Method	Slope	Intercept mg/dL [µmol/L]	Correlation Coefficient	n
Serum				
	Dimension® CREA	1.00	-0.08 [-7.1]	0.999
	Dimension® EZCR	1.02	+0.02 [+1.8]	0.999
	IDMS* primary reference method	1.04	+0.02 [+1.8]	0.999
Urine				
	Dimension® CREA	1.04	-3.58 [-316.5]	0.996

* IDMS: Isotope dilution mass spectrometry

f. CLSI EP09-A2-IR¹⁸ was used. The method used to fit the linear regression line was ordinary least-squares.

g. The range of creatinine values in the correlation study was 0.4 – 19.8 [35 – 1750].

h. The range of creatinine values in the correlation study was 0.53 – 14.11 [47 – 1247].

i. The range of creatinine values in the correlation study was 0.18 – 6.32 [16 – 559].

j. The range of creatinine values in the correlation study was 13.5 – 372.7 [1193 – 32947].

Serum and Plasma Equivalency				
No clinically significant difference was observed between serum (x) and plasma (y) as shown in the following ordinary least-squares regression.				
Sample Comparison	Slope	Intercept mg/dL [µmol/L]	Correlation Coefficient	n
Lithium heparin plasma vs.				
Serum	1.05	-0.02 [-1.8]	0.998	56 ^k

k. The range of creatinine values in the correlation study was 0.50 – 17.35 [44 – 1529.3].

Specificity

Hemolysis, Icterus, Lipemia (HIL) Interference

The CRE2 method was evaluated for interference according to CLSI EP07-A2.¹¹ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	Creatinine mg/dL [µmol/L]	Bias* %
Hemoglobin (hemolysate)	500 mg/dL [0.31 mmol/L]	1.50 [133]	< 10
	1000 mg/dL [0.62 mmol/L]	1.50 [133]	-11.1
	1000 mg/dL [0.62 mmol/L]	5.00 [442]	< 10
Bilirubin (conjugated)	20 mg/dL [342 µmol/L]	1.50 [133]	< 10
	40 mg/dL [684 µmol/L]	1.50 [133]	-17.2
	40 mg/dL [684 µmol/L]	5.00 [442]	< 10
Bilirubin (unconjugated)	10 mg/dL [171 µmol/L]	1.50 [133]	< 10
	20 mg/dL [342 µmol/L]	1.50 [133]	-20.2
	40 mg/dL [684 µmol/L]	5.00 [442]	< 10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	1.50 [133]	< 10
	1500 mg/dL [17.0 mmol/L]	1.50 [133]	+11.3
	2000 mg/dL [22.6 mmol/L]	5.00 [442]	< 10

* Analyte results should not be corrected based on this bias.

Non-Interfering Substances – Serum and Plasma

The following substances do not interfere with the CRE2 method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than or equal to 10% at creatinine concentrations of 1.50 mg/dL [133 µmol/L] and 5.00 mg/dL [442 µmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1324 µmol/L
Acetoacetate	20 mg/dL	2.0 mmol/L
Amikacin	8 mg/dL	137 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	6 mg/dL	342 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Cephalexin	25 mg/dL	720 µmol/L
Cephapirin	25 mg/dL	562 µmol/L
Cephadrine	25 mg/dL	769 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	503 mg/dL	13.3 mmol/L
Cimetidine	2 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.51 mg/dL	18 µmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
EDTA	200 mg/dL	2 g/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	1 mg/dL	21 µmol/L
Heparin (196 Units/mg)	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	24.2 mmol/L
Immunoglobulin G (IgG)	5000 mg/dL	50 g/L
Isopropanol	1.0 g/dL	100 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Nortriptyline	1000 ng/mL	3797 nmol/L
Penicillin G (1654 Units/mg)	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Potassium oxalate	500 mg/dL	5 g/L

Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Protein: Albumin	6000 mg/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Sodium fluoride	400 mg/dL	4 g/L
Theophylline	4 mg/dL	222 µmol/L
Urea	500 mg/dL	83 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L
Vancomycin	10 mg/dL	69 µmol/L

Non-Interfering Substances – Urine

The following substances do not interfere with the CRE2 method when present in urine at the concentrations indicated. Inaccuracies (biases) due to these substances are less than or equal to 10% at creatinine concentrations of 40.00 mg/dL [3536 µmol/L] and 175.00 mg/dL [15470 µmol/L].

Substance	Test Concentration	SI Units
50% Acetic Acid	25 mL/24 hr collection	25 mL/24 hr collection
6N HCl	0.6%	0.6%
6N HNO3	0.6%	0.6%
Acetone	100 mg/dL	0.172 mol/L
Bilirubin, conjugated	2.0 mg/dL	34.2 µmol/L
Boric acid	1% w/v	162 mmol/L
Ethanol	1.0 g/dL	217 mmol/L
Gamma globulin	0.5 g/dL	5.0 g/L
Glucose	2.0 g/dL	0.11 mol/L
Hemoglobin	115 mg/dL	0.07 mmol/L
Human serum albumin	0.5 g/dL	5.0 g/L
Oxalic acid	0.1 g/dL	110 mmol/L
Sodium carbonate	5g/24 hr collection	5g/24 hr collection
Sodium fluoride	1% w/v	238 mmol/L

Recovery

Certified standard serum reference materials from the National Institute of Standards and Technology (NIST) SRM 967a with known amounts of creatinine at concentrations of 0.85 and 3.88 mg/dL [75 and 343 µmol/L] were run as unknowns. The recovered creatinine concentrations on these samples were then measured and the following calculated percent bias recovery was observed.

NIST SRM 967a Level	IDMS Target Value mg/dL [µmol/L]	CRE2 Value mg/dL [µmol/L]	Bias Recovery
Level 1	0.85 [75]	0.89 [79]	+5%
Level 2	3.88 [343]	3.86 [341]	-1%

% Recovery: $\frac{[\text{Value obtained}] - [\text{IDMS Target Value}]}{[\text{IDMS Target Value}]}$ x 100

Limit of Detection and Limit of Blank^{18,1}

The Limit of Detection (LoD) for CRE2 for serum and plasma is 0.10 mg/dL [9 µmol/L] and for urine is 2.00 mg/dL [177 µmol/L], determined consistent with CLSI EP17-A2¹⁹ and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 240 determinations, with 120 blank and 120 low level samples. The Limit of Blank (LoB) for serum and plasma is 0.05 mg/dL [4 µmol/L] and for urine is 1.00 mg/dL [88 µmol/L].

I. LoD is the lowest concentration of creatinine that can be detected reliably. LoB is the highest concentration of creatinine that is likely to be observed for a blank sample.

Limit of Quantitation^{19,m}

The Limit of Quantitation (LoQ) for CRE2 for serum and plasma is 0.15 mg/dL [13 µmol/L] and based on allowable total error of 0.15 mg/dL [13 µmol/L], determined consistent with CLSI EP17-A2¹⁹.

The Limit of Quantitation (LoQ) for CRE2 for urine is 13.00 mg/dL [1149 µmol/L] and based on allowable total error of 3.00 mg/dL [265 µmol/L], determined consistent with CLSI EP17-A2¹⁹.

m. LOQ is the lowest amount of creatinine that can be quantitatively determined within a defined total error.

Standardization

The Dimension® CRE2 method is standardized to the Isotope dilution mass spectrometry (IDMS) reference method via correlation of patient samples and is verified using the IDMS traceable National Institute of Standards and Technology (NIST) standard reference material SRM 967.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

© 2014 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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© 2014 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Burtis CA, Ashwood ER. Tietz Fundamentals of Clinical Chemistry, 5th ed. Philadelphia: W.B. Saunders, 2001:23-25, 419-420.
2. Kidney Disease: Improving Global Outcomes (KDIGO) CKD Work Group. KDIGO 2012 Clinical Practice Guideline for the Evaluation and Management of Chronic Kidney Disease. Kidney inter., Suppl. 2013; 3: 1-150.
3. National Kidney Disease Education Program NKDEP: Improving the understanding, detection, and management of kidney disease. CKD and Drug Dosing: Information for Providers Retrieved June 2013, available: <http://www.nkdep.nih.gov/resources/CKD-drug-dosing.shtml>.
4. Knapp ML, and Mayne PD. Development of an automated kinetic Jaffe method designed to minimize bilirubin interference in plasma creatinine assays. Clin Chim Acta 1987; 168: 239-246.
5. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. CLSI document H3-A6 [ISBN 1-56238-650-6]. Wayne, PA: CLSI; 2007.
6. Clinical and Laboratory Standards Institute (CLSI). *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. CLSI document H01-A6 [ISBN 1-56238-740-5]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2010.
7. Clinical and Laboratory Standards Institute (CLSI). *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI document H18-A4 [ISBN 1-56238-724-3]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2010.
8. Young DS. Effects of Preamerical Variables on Clinical Laboratory Tests, 3rd ed. Washington, D.C. AACC Press, 2007: p. 316.
9. Clinical and Laboratory Standards Institute (CLSI). *Urinalysis; Approved Guideline—Third Edition*. CLSI document GP16-A3 [ISBN 1-56238-687-5] Wayne, PA: CLSI; 2009.
10. Wu A, ed. Tietz Clinical Guide to Laboratory Tests, 4th ed. Philadelphia: W.B. Saunders Company, 2006: p. 316-320.
11. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI document EP07-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 2500, Wayne, PA 19087 USA, 2005.
12. Burtis CA, Ashwood ER. Tietz Fundamentals of Clinical Chemistry, 5th ed. Philadelphia: W.B. Saunders, 2001:23-25, 422.
13. Ceriotti F, Boyd JC., *Reference Intervals for Serum Creatinine concentrations: Assessment of Available Data for Global Application*. Clinical Chemistry 54:3 (2008) 559-566.
14. Junge W., Wilke B., *Determination of reference intervals for serum creatinine, creatinine excretion and creatinine clearance with an enzymatic and a modified Jaffe method*. Clinica Chimica Acta 344 (2004) 137-148.
15. Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia: W.B. Saunders, 1999:1809.
16. Clinical and Laboratory Standards Institute (CLSI). *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. CLSI document EP28-A3c [ISBN 1-56238-682-4]. CLSI, 940 West Valley Road, Suite 2500, Wayne, PA 19087 USA, 2008.
17. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI document EP05-A2 [ISBN 1-56238-542-9]. Wayne, PA: CLSI; 2004.
18. Clinical and Laboratory Standards Institute (CLSI). *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition (Interim Revision)*. CLSI document EP09-A2-IR [ISBN 1-56238-731-6]. Wayne, PA: CLSI; 2010.
19. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI document EP17-A2 [ISBN 1-56238-795-2]. Wayne, PA: CLSI; 2012.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_ERGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare