Dimension® clinical chemistry system

Flex® reagent cartridge

GENT

See shaded sections: Updated information from 2017-04 version.


Issue Date 2019-04-01

Gentamicin

Intended Use: The GENT method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended to measure gentamicin, an aminoglycoside antibiotic in human serum or plasma. GENT test results may be used in the diagnosis and treatment of gentamicin overdose and in monitoring levels of gentamicin to ensure appropriate therapy.

Summary: Gentamicin is an antibiotic effective against gram negative aerobic bacteria. It has a wide spectrum of antibiotic activity and relatively low toxicity. Gentamicin is a naturally occurring antibiotic produced by the organism *Micromonospora purpurea*. Gentamicin is administered either intramuscularly or intravenously. Peak concentrations are reached 60 minutes after intramuscular injection and after completion of intravenous injection.^{1,2}

Principles of Procedure: The methodology for GENT involves a homogeneous particle enhanced turbidimetric inhibition immunoassay (PETINIA) technique which uses a synthetic particle-gentamicin conjugate (PR) and gentamicin-specific monoclonal antibody (Ab). Gentamicin present in the sample competes with the particles for the antibody, thereby decreasing the rate of aggregation. Hence, the rate of aggregation is inversely proportional to the concentration of gentamicin in the sample. The rate of aggregation is measured using bichromatic turbidimetric readings at 340 nm and 700 nm.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1, 2	Liquid	Particle Reagent	3.3 g/L ^{c,d}	Mouse, monoclonal
3, 4	Liquid	Buffer	25.7 mM	
5, 6	Liquid	Antibody	0.027 g/L ^{c,d}	
8	Liquid	0.5 N NaOH	^e	

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

c. Antibody titer and conjugate activity may vary from lot to lot.

d. Contains buffers, stabilizers and preservatives.

e. Sodium hydroxide is used as a probe cleaning solution and is not used in the reaction.

Risk and Safety:

H290, H319, H315, H317
P280, P305 + P351 + P338, P501

Warning!

May be corrosive to metals. Causes serious eye irritation. Causes skin irritation. May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: Sodium hydroxide; 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6
30 days for well 8

Specimen Collection and Handling: Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.³

Follow the instructions provided with your specimen collection device for use and processing.⁴ Complete clot formation should take place before centrifugation.⁵

Separated specimens are stable for 8 hours at room temperature, 2 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C.⁵

Heparin, oxalate, and EDTA anticoagulants do not interfere with the GENT method at concentrations normally found in blood collection tubes.

Procedure

Materials Provided

GENT Flex® reagent cartridge, Cat. No. DC49D

Materials Required But Not Provided

Drug Calibrator II, Cat. No. DC49D

Quality Control Materials

Test Steps

Sampling¹, reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

f. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus dead volume: Precise container filling is not required.

Test Conditions

Sample Size	3 µL
Buffer Volume	115 µL
Particle Reagent Volume	64 µL
Antibody Volume	64 µL
Diluent Volume	216 µL
Temperature	37 °C
Wavelength	340 and 700 nm
Type of Measurement	Turbidimetric rate

Calibration

Assay Range	0.0 – 12.0 µg/mL [0.0 – 25.92 µmol/L] ^g
Calibration Material	Drug Calibrator II, Cat. No. DC49D
Calibration Scheme	5 levels, n = 2
Units	µg/mL [µmol/L]

Typical Calibration Levels

Calibration Frequency

A new calibration is required

Assigned Coefficients

C ₀	471.674
C ₁	-468.119
C ₂	-2.524
C ₃	4.885
C ₄	0.500

g. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known gentamicin concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of GENT in µg/mL [µmol/L] using the calculation scheme outlined in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.0 – 12.0 µg/mL [0.0 – 25.92 µmol/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 12.0 µg/mL [25.92 µmol/L] should be repeated on dilution.

Manual Dilution:

Make appropriate dilution with Level 1 (0 µg/mL [µmol/L]) Drug Calibrator II or drug-free serum to obtain results within assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution:

Refer to your Dimension® Operator's Guide.

GENT results less than 0.20 µg/mL [0.4 µmol/L] should be reported as "less than 0.20 µg/mL [0.4 µmol/L]" instead of the numerical value.

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
1.5 µg/mL [3.24 µmol/L]	> 0.17 µg/mL [0.37 µmol/L]
6.0 µg/mL [12.96 µmol/L]	> 0.21 µg/mL [0.45 µmol/L]

Interfering Substances

Lipemia (Intralipid®) of 3000 mg/dL [33.9 mmol/L] at a gentamicin concentration of 12.03 µg/mL [25.98 µmol/L] tripped a test report message; therefore the magnitude of the interference could not be determined.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Therapeutic Range: Therapeutic gentamicin concentrations vary significantly depending on the individual patient. A range of 4.0 – 10.0 µg/mL [8.64 – 21.6 µmol/L] for peak drug levels indicates effective plasma or serum levels for many patients; however, some individuals are best treated at concentrations outside this range. Although it has not yet been established exactly what specimen concentration is toxic, trough concentration in excess of 2 µg/mL [4.32 µmol/L] for longer than 10 days have been associated with toxicity.^{1,2} The physician must determine the most appropriate therapeutic range for each patient.

Specific Performance Characteristics^a

Precision^d

Material	Mean µg/mL [µmol/L]	Standard Deviation (% CV)	
		Within-run	Total
Serum Pool	2.1 [4.54]	0.02 [0.04] (1.06)	0.05 [0.11] (2.41)
Dade® TDM Control	3.9 [8.42]	0.03 [0.06] (0.86)	0.07 [0.15] (1.75)
Plasma Pool	8.0 [17.28]	0.09 [0.19] (1.19)	0.22 [0.48] (2.74)

h. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).

i. Specimens at each level were analyzed in duplicate, two runs per day for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

Method Comparison

Regression Statistics^j

Comparative Method	Regression Statistics ^j		Correlation Coefficient	n
	Slope	Intercept µg/mL [µmol/L]		
Abbott Gentamicin on Abbott AxSYM®	1.02	-0.15 [-0.32]	0.987	197 ^k
Abbott Gentamicin on Abbott TDx®	1.00	-0.004 [-0.006]	0.997	205 ^l

j. Model equation for regression statistics is: results of Dimension® system = [slope x (comparative method results)] + intercept.

k. Range of samples: µg/dL [µmol/L]: 0.1 - 10.6 [0.03 - 3.4]

l. Range of samples: µg/dL [µmol/L]: 0.1 - 9.9 [0.03 - 3.2]

AxSYM® and TDx® are registered trademarks of Abbott Diagnostics, Abbott Park, IL 60064.

Specificity

HIL Interference

The GENT method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	Gentamicin Concentration µg/mL [µmol/L]	Bias ^m %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L] (monomer)	11.12 [24.02]	<10
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	12.03 [25.98]	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	11.25 [24.3]	<10

m. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances do not interfere with the GENT method when present in serum and plasma at the concentrations indicated. Systemic inaccuracies (bias) due to these substances are less than 10% at gentamicin concentrations of 4 and 10 µg/mL [8.6 and 21.6 µmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1323 µmol/L
Acetylsalicylic Acid	50 mg/dL	1.67 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Amobarbital	10 mg/dL	442 µmol/L
Ampicillin	5 mg/dL	143 µmol/L
Ascorbic Acid	3 mg/dL	170.3 µmol/L
Caffeine	10 mg/dL	515 µmol/L
Carbamazepine	12 mg/dL	508 µmol/L
Carbenicillin	200 mg/dL	528 µmol/L
Cephalothin	100 mg/dL	2522 µmol/L
Clindamycin	30 mg/dL	71 µmol/L
Chloramphenicol	25 mg/dL	774 µmol/L
Chlordiazepoxide	2 mg/dL	67 µmol/L
Chlorpromazine	5 mg/dL	157 µmol/L
Cimetidine	10 mg/dL	396 µmol/L
Codeine	10 mg/dL	334 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Diazepam	2 mg/dL	70 µmol/L
Erythromycin	20 mg/dL	273 µmol/L
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 µmol/L
Furosemide	2 mg/dL	61 µmol/L
Heparin (Porcine)	8000 U/L	8000 U/L
Ibuprofen	40 mg/dL	1939 µmol/L
Kanamycin A	25 mg/dL	52 µmol/L
Kanamycin B	25 mg/dL	52 µmol/L
Lidocaine	6 mg/dL	256 µmol/L
Lithium	3.5 mg/dL	5.07 mmol/L
Neomycin	10 mg/dL	815 µmol/L
Nicotine	2 mg/dL	123 µmol/L
Nortriptyline	14 mg/dL	532 µmol/L
Penicillin V	80 mg/mL	2247 µmol/L
Pentobarbital	10 mg/dL	442 µmol/L
Phenobarbital	15 mg/dL	646 µmol/L
Phenytoin	10 mg/dL	396 µmol/L
Primidone	10 mg/dL	458 µmol/L
Propoxyphene	0.4 mg/dL	12 µmol/L
Protein-Albumin	12 g/dL	120 g/L
Protein - IgG	12 g/dL	120 g/L
Rheumatoid Factor	750 IU/L	750 IU/L
Salicylic Acid	50 mg/dL	3.62 mmol/L
Secobarbital	5 mg/dL	210 µmol/L
Streptomycin	20 mg/dL	343.9 µmol/L
Sulfanilamide	150 mg/dL	8.71 mmol/L
Tetracycline	4 mg/dL	90 µmol/L
Theophylline	25 mg/dL	1388 µmol/L
Tobramycin	10 mg/dL	214 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Cross-reactivity

Aminoglycosides structurally similar to gentamicin (e.g. netilmicin, sagamicin, and sisomicin) may significantly cross-react with this method. Therefore, GENT cannot reliably be used for patients receiving these antibiotics separately or in combination with gentamicin.

Analytical Sensitivity: 0.2 µg/mL [0.3 µmol/L]

This sensitivity represents the lowest concentration of gentamicin that can be distinguished from zero. The sensitivity is defined as the mean value (n = 20) plus two standard deviations of the level 1 (0.0 µg/mL [0.0 µmol/L]) Drug Calibrator II.

Symbols Key: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics

All rights reserved.

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Burtis CA, Ashwood ER, Bruns DE. Tietz Textbook of Clinical Chemistry and Molecular Biology. Fourth Edition, Elsevier Saunders, St. Louis, MO 2006; pp. 1262 (clinical significance), pp. 2305 (reference values).
2. Finn AL, Taylor WJ, Individualizing Drug Therapy, Practical Applications of Drug Monitoring. New York: Gross, Townsend, Frank, Inc., 1981: 114-147.
3. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
4. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
LOT	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
REF	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante
EC REP	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos
IVD	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
CE	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
CONTENTS	Contents / Inhalt / Contenu / Contenido / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
LEVEL	Level / Konzentration / Niveau / Livello / Nivel

2014-10_EPGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare