Dimension® clinical chemistry system

Flex® reagent cartridge

MALB

See shaded sections: Updated information from 2015-03 version.


Issue Date 2019-04-01

Microalbumin

Intended Use: The MALB method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended to quantitatively measure albumin in human urine.

Summary: Diabetic nephropathy occurs in 20 – 40% of patients with diabetes and is the single leading cause of end-stage renal disease. Measurements of albumin in urine are used as an aid in monitoring the development of incipient diabetic nephropathy and for the diagnosis and treatment of a variety of diseases characterized by albuminuria.¹

Principles of Procedure: The Microalbumin (MALB) method is based on a particle-enhanced turbidimetric inhibition immunoassay (PETINIA) adapted to the Dimension® clinical chemistry system which allows direct quantitation of albumin in urine samples. The MALB Flex® reagent cartridge contains a particle reagent (PR) consisting of synthetic particles with human albumin bound to the surface. Aggregates of these particles are formed when a monoclonal antibody (Ab) to human albumin is introduced. Albumin (ALB) present in the sample competes with the particles for the antibody, thereby decreasing the rate of aggregation. Hence, the rate of aggregation is inversely proportional to the concentration of albumin in the sample. The rate of aggregation is measured using bichromatic turbidimetric reading at 340 and 700 nm.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1, 2	Liquid	Particle Reagent Microbial inhibitors	2 mg/mL ^c	
3, 6	Liquid	NaOH ^d	0.5 N	
4, 5	Liquid	Antibody to human albumin Microbial inhibitors	140 µg/mL ^c	Mouse, monoclonal
7, 8	Liquid	Buffer Microbial inhibitors		

- a. Wells are numbered consecutively from the wide end of the cartridge.
b. Nominal value per well in a cartridge.
c. The particle reagent and antibody are a matched pair and may vary from lot to lot.
d. Sodium hydroxide is used as a probe cleaning solution and is not used in the reaction.

Risk and Safety:

H290, H314, H317
P280, P304 + P340 + P310, P301 + P310 + P331,
P303 + P361 + P353 + P310, P305 + P310, P501

Danger!

May be corrosive to metals. Causes severe skin burns and eye damage. May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Immediately call a POISON CENTER or doctor/physician. IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. Do NOT induce vomiting. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Immediately call a POISON CENTER or doctor/physician. IF IN EYES: Immediately call a POISON CENTER or doctor/physician. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: Sodium hydroxide; 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone

Safety data sheets (MSDS/SDS) available on www.siemens.com/diagnostics

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 8

Specimen Collection and Handling:^{1,2,3,4} Normal procedures for collecting urine may be used for samples to be analyzed by this method. The following samples are acceptable:

- 24-hour collection
- Overnight (8 – 12-hour) collection
- 1 to 2-hour collection
- First-morning sample for simultaneous albumin and creatinine measurement

Samples should not be collected after exertion, in the presence of urinary tract infection, during acute illness, immediately after surgery, or after an acute fluid load. Specimens should be collected without preservatives. Specimens visibly contaminated with blood are not suitable for analysis of albumin concentration.

Urine specimens are stable at room temperature for 2 days and may be stored at 8 °C for up to 14 days. The use of frozen urine samples is not recommended.⁵

Procedure

Materials Provided

MALB Flex® reagent cartridge, Cat. No. DF114

Materials Required But Not Provided

Microalbumin Calibrator, Cat. No. DC114
Quality Control Materials

Test Steps

Sampling,⁶ reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

- e. The sample container must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

Test Conditions

Sample Size	17 µL
Particle Reagent Volume	76 µL
Antibody Reagent Volume	76 µL
Assay Buffer Volume	288 µL
Temperature	37 °C
Wavelengths	340 and 700 nm
Type of Measurement	Turbidimetric rate

Calibration

Assay Range	1.3 – 100 mg/L
Calibration Material	Microalbumin Calibrator, Cat. No. DC114
Calibration Scheme	5 levels, n = 2
Units	mg/L
Typical Calibration Levels	0, 12.5, 25.0, 50.0, 110.0 mg/L
Calibration Frequency	Every 30 days for any one lot.
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Starting coefficients	C ₀ 590 C ₁ -605 C ₂ -1.6 C ₃ 37.4 C ₄ 0.5

Quality Control

At least once each day of use, analyze two levels of a quality control material with known microalbumin concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of microalbumin in mg/L using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 1.3 – 100 mg/L

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 100 mg/L should be repeated on dilution.

Manual Dilution: Make appropriate dilution with Reagent grade water or Microalbumin Calibrator (Cat. No. DC114) Level 1, to obtain results within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution: Refer to your Dimension® Operator's Guide. The recommended autodilute volume is 2 µL.

Results less than 1.3 mg/L should be reported as "less than 1.3 mg/L" instead of the numerical value.

Limitations of Procedure

The instrument reporting system contains error messages to warn the user of specific malfunctions. MALB results followed by such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
12.5 mg/L	>1.9 mg/L
50 mg/L	>2.5 mg/L

Expected Values:^{1,2}Less than 20 µg/min^f

Less than 30 mg/24 hour

Excretion rate: less than 30 mg albumin/g creatinine^gf. $\text{MALB (mg/L)} \times [\text{Urine Volume (mL)} / \text{Time (minutes)}] = \mu\text{g MALB/min}$ g. $[\text{MALB (mg/L)} / \text{Urine creatinine (mg/dL)}] \times 100 = \text{mg MALB/g creatinine}$

To minimize intra-individual variation, analysis of three random urine samples collected over the course of a week has also been recommended.¹

Each laboratory should establish its own reference interval for microalbumin as performed on the Dimension® system.

Specific Performance Characteristics^h

Material	Precision ^{h,i}		
	Mean mg/L	Standard Deviation (% CV) Within-run	Total
Bio-Rad Liquichek™ Urine Control			
Level 1	10	0.65 (6.2)	1.06 (10.2)
Level 2	68	0.69 (1.0)	1.17 (1.7)
Urine Pool 1	20	0.64 (3.2)	1.57 (7.8)
Urine Pool 2	50	0.65 (1.3)	2.07 (4.2)

h. All Specific Performance Characteristics testing were run on the Dimension® XL and RxL systems after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).

i. Reproducibility testing was done in accordance with the CLSI/NCCLS Approved Guideline for Evaluation of Precision Performance of Clinical Chemistry Devices (EP5-A, Feb. 1999)

j. Specimens at each level were analyzed in duplicate, twice a day, for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

Bio-Rad Liquichek™ Urine Control is a trademark of Bio-Rad Laboratories Diagnostic Group, Irvine, CA 92618.

Method Comparison

Comparative Method	Regression Statistics ^k			
	Slope	Intercept mg/L	Correlation Coefficient	n
N Antiserum to Human Albumin	0.97	0.45	0.999	50 ^l

k. Model equation for regression statistics is: [Result of Dimension® system] =

[Slope x (result of BNI albumin method)] + Intercept

l. The range of microalbumin values in the correlation study was 4.2 to 98.4 mg/L.

Specificity**Non-Interfering Substances**

The MALB was evaluated for interference according to CLSI/NCCLS EP7-A. The following substances do not interfere with the MALB method when present in urine concentrations indicated. Systemic inaccuracies (bias) due to these substances are less than 10% at a urinary albumin concentration of 18.4 to 22.5 mg/L.

Substance	Test Concentration
Alpha-1 microglobulin	0.3 mg/L
Alpha-2 macroglobulin	7 mg/L
Acetoacetate	750 mg/L
Acetone	150 mg/L
β-hydroxybutyrate	2250 mg/L
β-2 Microglobulin	0.6 mg/L
Bilirubin	20 mg/dL [342 µmol/L] ^m
Captopril	750 mg/L
Creatinine	8540 mg/L
Enalapril	2000 mg/L
Glucose	5000 mg/L
Hemoglobin	10 mg/L
Immunoglobulin L Chain (k and l)	5000 mg/L
Immunoglobulin L Chain (k)	30 mg/L
Immunoglobulin G	100 mg/L
Retinol binding protein	1 mg/L
Sodium chloride	200 mmol/L
Transferrin	40 mg/L
Urea	10000 mg/L [1670 mmol/L]

m. système International d'Unités [SI Units] are in brackets.

Analytical Sensitivity: 1.3 mg/L

The analytical sensitivity represents the lowest concentration of microalbumin that can be distinguished from zero. The sensitivity is defined as the mean value (n = 20) plus two standard deviations of a microalbumin-free sample.

Symbols Key: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics

All rights reserved.

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. American Diabetes Association: Standards of Medical Care for Patients with Diabetes Mellitus (Position Statement). *Diabetes Care* 26 (Supp. 1): S33-S50, 2003.
2. Burtis CA, Ashwood ER. Tietz Fundamentals of Clinical Chemistry, Fifth Edition, W.B. Saunders Company, Philadelphia, PA 2001; pp. 457 (clinical significance), pp. 457-458 (specimen collection and storage).
3. Kaplan LA, Pesce AJ. Clinical Chemistry: Theory, analysis, and correlation, Third Edition, Mosby, Inc., St. Louis, MO 1996; pp. 78-79 (specimen collection).
4. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, Second Edition, AACCPress, Washington, DC 1997; pp. 3-18, 19 (specimen storage).
5. Tietz N. Clinical Guide to Laboratory Tests, Fourth Edition, 2006: p 70.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EPGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare