Dimension® clinical chemistry system

Flex® reagent cartridge

CA

See shaded sections: Updated information from 2017-12 version.

Issue Date 2019-05-21


Calcium

Intended Use: The CA method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of calcium in human serum, plasma and urine.

Summary: In blood, approximately 50% of the plasma calcium is free, 40% is protein bound and 10% is complexed. About 80% of protein-bound calcium is associated with albumin, with the remaining 20% associated with globulins.¹ Calcium measurements are used in the diagnosis and treatment of parathyroid disease, a variety of bone diseases, chronic renal diseases and tetany (intermittent muscular contractions or spasms).

The calcium method is a modification of the calcium o-cresolphthalein complexone (OCPC) reaction originally reported by Schwarzenbach, et al.² Stern and Lewis³ later adapted this reaction to a colorimetric calcium assay. Connerty and Briggs⁴ demonstrated the use of 8-quinolinol to reduce magnesium interference, and Sarkar and Chauhan⁵ reported the use of this procedure for serum calcium without protein precipitation. This method incorporates the use of 8-quinolinol to reduce magnesium interference and glycine buffer for pH control.

Principles of Procedure: Calcium reacts with OCPC to form a purple complex. The amount of complex thus formed is proportional to the calcium concentration and is measured using a bichromatic (577, 540 nm) endpoint technique. Magnesium ions, which also form a colored complex with OCPC, are removed from the reaction by complexation with 8-quinolinol.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1 – 6	Liquid	Glycine buffer	0.22 mol/L
7 – 8	Liquid	OCPC/ 8-Quinolinol	0.39 mmol/L 6.6 mmol/L

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal final value per test at manufacture.

Risk and Safety

Safety data sheets (MSDS/SDS) available on [siemens.com/healthcare](https://www.siemens.com/healthcare)

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 1 days for wells 1 – 6
10 days for wells 7 – 8

Specimen Collection and Handling: Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁶

Follow the instructions provided with your specimen collection device for use and processing.⁷

Complete clot formation should take place before centrifugation.⁸

Serum should be separated from red cells and analyzed promptly.⁹ Specimens are stable for 8 hours at room temperature, 2 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C or colder.⁸

Collection of 24 hour urine specimen should be made in a container with 10 to 20 mL of 6M HCl.

If preservatives are not used during 24 hour urine collection, acidify the sample to below pH 3.0 before analysis.¹⁰

Corvac® and SST® collection tubes do not affect the CA method.

Corvac® is a registered trademark of Monoject, Division of Sherwood Medical, St. Louis, MO.
SST® is a registered trademark of Becton Dickinson, Rutherford, NJ.

Procedure

Materials Provided

CA Flex® reagent cartridge, Cat. No. DF23A

Materials Required But Not Provided

CHEM I Calibrator, Cat. No. DC18B or DC18C

Quality Control Materials

Test Steps

Sampling,⁶ reagent delivery, mixing, processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

c. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus the dead volume; precise container filling is not required.

Test Conditions

Sample Size	5 µL
Reagent 1 Volume	145 µL
Reagent 2 Volume	33 µL
Diluent Volume	258 µL
Temperature	37 °C
Wavelength	577 and 540 nm
Type of Measurement	Bichromatic endpoint

Calibration

Assay Range	5.0 – 15.0 mg/dL [1.25 – 3.75 mmol/L] ^d
Calibration Material	CHEM I Calibrator, Cat. No. DC18B or DC18C
Calibration Scheme	3 levels, n = 3
Units	mg/dL [mmol/L] (mg/dL x 0.25) = [mmol/L]
Typical Calibration Levels	7, 10, 14 mg/dL [1.8, 2.5, 3.5 mmol/L]
Calibration Frequency	Every 3 months for any one lot
A new calibration is required	- For each lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ 1.000 C ₁ 0.090

d. Système International d'Unités [SI units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known calcium concentration. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of calcium in mg/dL [mmol/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 5.0 – 15.0 mg/dL [1.25 – 3.75 mmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 15.0 mg/dL [3.75 mmol/L] should be repeated on dilution.

Manual Dilution: Serum/Plasma: Make appropriate dilutions with Reagent grade water to obtain result within the assay range. Enter dilution factor.

Urine: Dilute 1 part urine: 1 part Reagent grade water. Enter dilution factor of 2.

Serum/Plasma/Urine: Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): Refer to your Dimension® Operator's Guide.
(for serum/plasma/urine)

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
7.0 mg/dL [1.75 mmol/L]	> 0.17 mg/dL [0.04 mmol/L]
13.5 mg/dL [3.37 mmol/L]	> 0.23 mg/dL [0.06 mmol/L]

Interfering Substances

Interference due to magnesium is negligible at magnesium levels normally encountered in human serum. A maximum positive interference of 0.7 mg/dL [0.17 mmol/L] occurs at a magnesium level of 7 mg/dL [2.9 mmol/L].¹¹

Calcium values may be falsely decreased in the presence of gadolinium-containing contrast agents such as Omniscan™. Therefore the manufacturer of this product recommends to avoid drawing samples for serum calcium determination 24 hours after administration of Omniscan™.^{12,13}

Bilirubin (unconjugated) of 80 mg/dL [1368 µmol/L] decreases calcium at 6.4 mg/dL [1.60 mmol/L] by 11%.

Lipemia (Intralipid®) of 600 mg/dL [6.78 mmol/L] and above tripped a test report message; therefore the magnitude of the interference could not be determined.

Omniscan™ is a trademark of Amersham Health AS, Oslo, Norway.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

EDTA when present at 200 mg/dL [2 g/L]^a and potassium oxalate when present at 500 mg/dL [5 g/L]^a depresses the CA result to less than the assay range of the method.

e. Concentrations normally found in blood collection tubes.

Expected Values¹¹

Serum: 8.5 – 10.1 mg/dL [2.12 – 2.52 mmol/L]
Urine:¹⁴ 42 – 353 mg/24 hr [1.0 – 8.8 mmol/24 hr]

This reference population consisted of the following:
75 males, ages 18 – 54
75 females, ages 18 – 55

The reference interval was calculated non-parametrically and represents the central 95% of the population.
Each laboratory should establish its own reference interval for calcium as performed on the Dimension® system.

Specific Performance Characteristics¹

Material	Precision ^a		
	Mean mg/dL [mmol/L]	Standard Deviation (% CV)	
		Within-run	Total
Moni-Trol® Serum			
Level 1	9.2 [2.3]	0.08 [0.02] (0.8)	0.16 [0.04] (1.7)
Level 2	12.9 [3.2]	0.12 [0.03] (0.9)	0.24 [0.06] (1.9)
LYPHOCHEK® Urine Control			
Urine			
Normal	8.4 [2.1]	0.06 [0.02] (0.7)	0.08 [0.02] (0.9)

- f. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (see your Dimension® Operator's Guide).
- g. Specimens at each level were analyzed in duplicate for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

Moni-Trol® is a registered trademark of Medical Analysis Systems, Inc., Camarillo, CA 93012, USA.
LYPHOCHEK® is a registered trademark of Bio-Rad, Anaheim, CA 92806, USA.

Method Comparison

Regression Statistics^b

Comparative Method	Slope	Intercept mg/dL [mmol/L]	Correlation Coefficient	n
Serum CALCM method on the aca® analyzer	0.95	0.22 [0.06]	0.96	106
Urine CALCM method on the aca® analyzer	0.97	0.12 [0.03]	0.987	112

h. Model equation for regression statistics is: [results of Dimension® system] = slope x [comparative method results] + intercept.

Specificity

HIL Interference

The CA method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration [SI Units]	CA Concentration mg/dL [mmol/L]	Bias ⁱ %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L] (monomer)	6.4 [1.60]	<10
Bilirubin	60 mg/dL [1026 µmol/L]	6.4 [1.60]	<10
Lipemia (Intralipid®)	200 mg/dL [2.26 mmol/L]	6.3 [1.58]	<10

i. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances do not interfere with the CA method when present in serum at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at the CA concentration of 10 mg/dL [2.5 mmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	0.025 mg/dL	1.66 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	5 mg/dL	284 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	17.6 µmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	421 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	4 mg/dL	222 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Analytical Sensitivity: 5.0 mg/dL [1.25 mmol/L]

The analytical sensitivity represents the lowest concentration of calcium that can be distinguished from zero. This sensitivity is defined as the low value of the analytical measurement range.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, aca® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía:

1. Burtis C, Ashwood E. Tietz Fundamentals of Clinical Chemistry, 5th ed. Philadelphia: WB Saunders, 2001:797-799.
2. Anderegg G, Flaschka H, Sallmann R, Schwartzbach G. Metallindikatoren VII ein auf Erdalkalitionen ansprechendes Phtalein und seine analytische Verwendung HELV, Chim Acta 1954; 37:113-120.
3. Stern J, Lewis WHP. The colorimetric estimation of calcium in serum with O-cresolphthalein complexone, Clin Chim Acta 1957; 2:576-580.
4. Connerty HV, Briggs AR. Determination of serum calcium by means of O-cresolphthalein complexone, Am J Clin Path 1966; 45:290-296.
5. Sarker B, Ray C, Chauhan UPS. A new method for determining micro-quantities of calcium in biological materials, Anal Biochem 1967; 20:155-166.
6. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
7. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
8. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
9. Henry RJ. Clinical Chemistry, Principles and Technics, Harper and Row, New York, NY, 1974, p 656.
10. Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, PA 1999, p 51 and 1405.
11. Henry JB, and LaFever D. An Evaluation of a Calcium Method for the DuPont aca®, DuPont Company, Wilmington, DE (1974).
12. Norman PT, Froysa A, Svaland M. Interference of gadodiamide injection (Omniscan®) on the colorimetric determination of serum calcium. Scand J. Clin Lab Invest 1995; 55; 421-426.
13. Lin J, Idee JM, Port M, Dai A, Berthommier C, Robert M, Raynal I, Devoldere L, Corot C. Interference of magnetic resonance imaging contrast agents with the calcium measurement technique using colorimetric reagents. Journal of Pharmaceutical and Biomedical analysis 21 (1999) 931-943.
14. Barger JD, Dicht JJ. Reference Intervals for Urine on AMY, CA, MG, Na/K, and UP. DuPont Company, Wilmington, DE, October 1983.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenido / Contenuto

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFSS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare