Dimension® clinical chemistry system

Flex® reagent cartridge

C4

See shaded sections: Updated information from 2016-03 version.

Issue Date 2019-04-01

Complement C4


Intended Use: The C4 method for the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended to quantitatively measure Complement C4 (C4) in human serum as an aid in the diagnosis of immunologic disorders associated with Complement C4 protein.

Summary: The complement system is an integral part of the nonspecific antigen immune defenses. It can be activated via two reaction pathways; the classical pathway, which is triggered primarily by cell-bound immune complexes, and the alternative pathway, which is activated primarily by foreign bodies such as microorganisms. The complement component C4 belongs to the classical pathway of complement activation. Complement activation is associated with consumption of components C3 and/or C4 so that a reduction in their concentrations is diagnostically useful.^{1,2,3,4,5}

Diminished serum concentrations of C4 are observed primarily in active systemic lupus erythematosus (SLE), in forms of membrane proliferative glomerulonephritis and in immune complex diseases (serum sickness). In the case of SLE the serum concentrations of the complement factors reflect the activity of the disease. Isolated diminished levels of C4 can occur in heredity angioneurotic oedema and in cases of cryoglobulinaemia. Complement component C4 reacts as an acute-phase protein and may therefore show elevated serum concentrations in patients with inflammatory diseases. Hereditary deficiency states of complement C4 have been reported.^{1,2,3,4,5}

C4 is a relatively fragile protein, especially in the presence of Ca⁺⁺. Activation of complement and splitting of C4 can occur at room temperature, or colder, in many samples with or without cryoglobulins.⁴

Principles of Procedure: The C4 method is a quantitative, turbidimetric assay using endpoint detection, based on the precipitation of C4 by its polyclonal antibody.³ C4 from serum reacts with its polyclonal antibody to form immune complexes. Addition of polyethylene glycol accelerates formation of these complexes. The resulting turbidity is measured by bichromatic end-point measurements at 340 and 700 nm. The increase in turbidity is proportional to the concentration of C4 in the sample and the results are reported in analyte units (mg/dL or g/L).

a. The antibody is manufactured by Siemens Healthcare Diagnostics, Marburg, Germany.

Reagents

Wells ^b	Form	Ingredient	Concentration ^c	Source
1, 2, 3	Liquid	Polyethylene Glycol, Buffer, and stabilizers	60 mg/mL 60 mM	Rabbit
4, 5, 6	Liquid	Anti C4 Antibody and stabilizers	1.3 mg/mL ^d	

b. Wells are numbered consecutively from the wide end of the cartridge.

c. Nominal value in final reaction mixture.

d. Antibody titer varies from lot to lot.

Risk and Safety

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 4 days for wells 1 – 6

Specimen Collection and Handling: Serum can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁶ Follow the instructions provided with your specimen collection device for use and processing.⁷

Specimens should be free of particulate matter. To prevent the appearance of fibrin in serum samples, complete clot formation should take place before centrifugation.⁸ Visibly lipemic samples or frozen samples which are turbid after thawing must be clarified by centrifugation (10 minutes at approximately 15000 x g) prior to the assay.

Separated samples are stable for 8 hours at ambient temperatures, 2 days at 2 – 8 °C, or frozen at -20 °C for up to three weeks. Repeated freezing and thawing should be avoided.⁹ During prolonged storage the concentration of C4 decreases and the concentration of C4c increases. Therefore, complement protein results for stored samples need to be assessed against reference intervals determined under similar conditions.

Aged samples can be expected to give higher values than fresh samples.

Plasma samples are not recommended.

Procedure

Materials Provided

C4 Flex® reagent cartridge, Cat. No. DF102

Materials Required But Not Provided

Special Protein Calibrator, Cat. No. DC51

Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

Test Conditions

Sample Size	3 µL
Reagent 1 Volume	105 µL
Reagent 2 Volume	50 µL
Temperature	37 °C
Wavelength	340 and 700 nm
Type of Measurement	Bichromatic endpoint

Calibration

Assay Range	5.0 – 140.0 mg/dL [0.050 – 1.400 g/L] ^e
Calibration Material	Special Protein Calibrator, Cat. No. DC51
Calibration Scheme	5 levels, n = 2
Units	mg/dL [g/L] (mg/dL x 0.01) = [g/L]
Typical Calibration Levels	0.0, 10.0, 40.0, 100.0, 180.0 mg/dL [0.000, 0.100, 0.400, 1.000, 1.800 g/L]
Calibration Frequency	Every 60 days for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ -3000.0 C ₁ 7500.0 C ₂ -0.50 C ₃ 5.0 C ₄ 0.50

e. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known C4 concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of C4 in mg/dL [g/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 5.0 – 140.0 mg/dL [0.050 – 1.400 g/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 140.0 mg/dL [1.400 g/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilution with saline to obtain results within assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): None

C4 results less than 5.0 mg/dL [0.050 g/L] should be reported as "less than 5.0 mg/dL [0.050 g/L]" instead of the numerical value.

Limitations of Procedure

The instrument reporting system contains error messages to warn the user of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

C4 Concentration	SD
15.0 mg/dL [0.150 g/L]	>0.8 mg/dL [0.008 g/L]
40.0 mg/dL [0.400 g/L]	>2.0 mg/dL [0.020 g/L]

Interfering Substances

Lipemia (Intralipid®) of 600 mg/dL [6.78 mmol/L] and greater tripped a test report message; therefore the magnitude of the interference could not be determined.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Expected Values: 17.4 – 52.2 mg/dL [0.174 – 0.522 g/L]^f

This reference population applies to serum samples from 155 males and females, ages 18 – 61. This reference interval was calculated non-parametrically and represents the central 95% of the population.⁹ Each laboratory should establish its own reference interval for C4 as performed on the Dimension® system.

f. The values obtained from aged samples may be higher.

Specific Performance Characteristics^a

Precision^{a,i}

Material	Mean mg/dL [g/L]	Standard Deviation (% CV)	
		Within-run	Total
Bio-Rad Liquichek™ Immunology Control			
Level 1	15.8 [0.158]	0.46 [0.005] (2.9)	0.57 [0.006] (3.6)
Level 2	28.5 [0.285]	0.81 [0.008] (2.8)	1.24 [0.012] (4.3)
Level 3	39.4 [0.394]	0.40 [0.004] (1.0)	0.69 [0.007] (1.7)
Serum Pool	28.4 [0.284]	0.58 [0.006] (2.0)	0.96 [0.010] (3.4)

- g. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).
- h. Reproducibility testing was done in accordance with the CLSI/NCCLS Approved Guideline for Evaluation of Precision Performance of Clinical Chemistry Devices (EP5-A, Feb. 1999).
- i. Specimens at each level were analyzed in duplicate, once a day, for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

Liquichek™ is a trademark of Bio-Rad Laboratories, Inc., Irvine, CA 92618.

Method Comparison

Regression Statistics^j

Comparative Method	Slope	Intercept mg/dL [g/L]	Correlation Coefficient	n
Beckman Array 360	1.02	0.37 [0.004]	0.936	125 ^k

- j. Model equation for Passing-Bablok regression statistics is: [Result of Dimension® system] = [Slope x comparable method result] + Intercept.
- k. The range of C4 values in the correlation study was 10.0 – 73.9 mg/dL [0.100 – 0.739 g/L].

Specificity

HIL Interference

The C4 method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	C4 Concentration mg/dL [g/L]	Bias ⁱ %
Hemoglobin (hemolysate)	1000 mg/dL	20 [0.20]	<10
	[0.62 mmol/L] (monomer)		
Bilirubin (unconjugated)	80 mg/dL	21 [0.21]	<10
	[1368 µmol/L]		
Lipemia (Intralipid®)	200 mg/dL	20 [0.20]	<10
	[2.26 mmol/L]		

- i. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances at the concentrations indicated have no significant effect (less than 10%) on the C4 method when added to a serum pool containing 19.0 mg/dL [0.190 g/L] C4.

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1323 µmol/L
Albumin	6 g/dL	60 g/L
Amikacin	15 mg/dL	256 µmol/L
Ascorbic Acid	3 mg/dL	170 µmol/L
Caffeine	10 mg/dL	515 µmol/L
Carbamazepine	12 mg/dL	508 µmol/L
Chloramphenicol	25 mg/dL	774 µmol/L
Chlordiazepoxide	2 mg/dL	67 µmol/L
Chlorpromazine	5 mg/dL	157 µmol/L
Cholesterol Ester	500 mg/dL	12.9 mmol/L
Cimetidine	10 mg/dL	396 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 75	2500 mg/dL	333 µmol/L
Diazepam	2 mg/dL	70 µmol/L
Digoxin	5 ng/mL	6.4 nmol/L
Erythromycin	20 mg/dL	273 µmol/L
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 µmol/L
Furosemide	2 mg/dL	61 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Heparin (Sodium)	8 U/mL	8000 U/L
Hydralazine	1 µg/mL	6.2 µmol/L
Ibuprofen	40 mg/dL	1939 µmol/L
Lidocaine	6 mg/dL	256 µmol/L
Lithium Chloride	3.5 mg/dL	5.07 mmol/L
Nicotine	2 mg/dL	123 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	442 µmol/L
Phenobarbital	15 mg/dL	646 µmol/L
Phenytin	10 mg/dL	396 µmol/L
Primidone	10 mg/dL	458 µmol/L
Procainamide	10 µg/mL	36.8 µmol/L
Propoxyphene	0.4 mg/dL	12 µmol/L
Rheumatoid Factors	1040 IU/mL	1040 IU/mL
Salicylic Acid	50 mg/dL	3.62 mmol/L
Theophylline	25 mg/dL	1388 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Recovery

Recovery of reference material CRM470/RPPHS¹⁰ ranged from 96 – 101%, with a mean recovery of 97%.

Analytical Sensitivity: ≤ 5 mg/dL [0.050 g/L]

The analytical sensitivity represents the low end of the assay range of C4.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

- Herbert LE et al., *Kidney International* 1991; 39: 811.
- West CD. *Complement Inflamm* 1989; 6: 49.
- Dalmasso P. *Critical Reviews in Clinical Laboratory Sciences* 1986; 4: 123.
- Anderson A. *Complement Component C4, Serum Proteins in Clinical Medicine Volume I Laboratory Section, First Edition*, Foundation for Blood Research, Scarborough, ME, 1996.
- Thomas L. *Clinical Laboratory Diagnostics*, T-H Books Verlagsgesellschaft, Frankfurt, 1998, First English Edition, pp 794-806 (The complement system).
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests – Third Edition*. Washington, DC: AACC Press, 2007; p 273.
- Whicher JT et al. *New International Reference Preparation for Proteins in Human Serum (RPPHS)*. *Clin. Chem.* 1994; 40: 934.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EPGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare