Dimension® clinical chemistry system

Flex® reagent cartridge

IGA

See shaded sections: Updated information from 2016-03 version.

Issue Date 2019-04-01

Immunoglobulin A


Intended Use: The IGA method for the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended to quantitatively measure Immunoglobulin A (IgA) in human serum and plasma.

Summary: Immunoglobulins are formed by plasma cells as a humoral response to contact of the immune system with antigens.¹ The initial reaction is production of IgM antibodies, followed later by IgG and IgA antibodies. Quantitative determination of the immunoglobulins can provide important information on the humoral immune status.

IgA is the second most abundant immunoglobulin in normal adult plasma. It is the predominant immunoglobulin in body secretions such as colostrum, saliva, and sweat.² Synthesis of IgA begins during the first few weeks of life; there is essentially no IgA in the plasma of newborns. Decreased IgA concentrations occur in primary immunodeficiency conditions as well as in secondary immune insufficiencies.³ Increased IgA concentrations occur due to polyclonal or oligoclonal immunoglobulin proliferation.⁴ In serum, IgA exists in monomeric and polymeric forms. Monoclonal immunoglobulinemia⁴ requires detailed differential diagnostic investigation in addition to the quantitative determination.

Principles of Procedure: The IGA method is a quantitative, turbidimetric assay using end-point detection, based on the precipitation of IgA by its polyclonal antibody.⁵

IgA from serum or plasma reacts with its polyclonal antibody to form immune complexes. Addition of polyethylene glycol accelerates formation of these complexes. The resulting turbidity is measured by bichromatic end-point measurements at 340 and 700 nm. The increase in turbidity is proportional to the concentration of IgA in the sample and the results are reported in analyte units (mg/dL or g/L).

a. The antibody is manufactured by Siemens Healthcare Diagnostics, Marburg, Germany.

Reagents

Wells ^b	Form	Ingredient	Concentration ^c	Source
1, 2, 3	Liquid	Polyethylene Glycol, Buffer, and stabilizers	120 mg/mL 60 mM	
4, 5, 6	Liquid	Anti IgA Antibody and stabilizers	1.5 mg/mL ^d	Rabbit

b. Wells are numbered consecutively from the wide end of the cartridge.

c. Nominal value in final reaction mixture.

d. Antibody titer varies from lot to lot.

Risk and Safety

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 4 days for wells 1 – 6

Specimen Collection and Handling: Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁶ Follow the instructions provided with your specimen collection device for use and processing.⁷

Specimens should be free of particulate matter. To prevent the appearance of fibrin in serum samples, complete clot formation should take place before centrifugation.⁸ Very lipemic samples or frozen samples which are turbid after thawing must be clarified by centrifugation prior to the assay.

Separated samples are stable for 3 days at ambient temperature, 7 days at 2 – 8 °C or stored frozen. Samples can be stored at below -20 °C for up to six months, if they are frozen within 24 hours after collection and if repeated freeze-thaw cycles are avoided.⁹

Procedure

Materials Provided

IGA Flex® reagent cartridge, Cat. No. DF74

Materials Required But Not Provided

Special Protein Calibrator, Cat. No. DC51
Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

Test Conditions

Cuvette 1	
Sample Size	10 µL
PEG Buffer Volume	135 µL
Diluent Volume	255 µL
Temperature	37 °C
Cuvette 2	
First Transfer Volume	4 µL
PEG Buffer Volume	110 µL
Diluent Volume	280 µL
Antibody Volume	50 µL
Second Transfer Volume	36 µL
Temperature	37 °C
Wavelengths	340 and 700 nm
Type of Measurement	Bichromatic

Calibration

Assay Range	10 – 1000 mg/dL [0.1 – 10.0 g/L]*
Calibration Material	Special Protein Calibrator, Cat. No. DC51
Calibration Scheme	5 levels, n = 2
Units	mg/dL [g/L] (mg/dL x 0.01) = [g/L]
Typical Calibration Levels	0, 90, 300, 600, 1200 mg/dL [0.0, 0.9, 3.0, 6.0, 12.0 g/L]
Calibration Frequency	Every two months for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C_0 -50.0 C_1 400.0 C_2 -1.0 C_3 500.0 C_4 0.5

e. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known IgA concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of IgA in mg/dL [g/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 10 – 1000 mg/dL [0.1 – 10.0 g/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 1000 mg/dL [10.0 g/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilution with saline to obtain results within assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): Refer to your Dimension® Operator's Guide.

The autodilute sample volume is 2 µL.

IgA results less than 10 mg/dL [0.1 g/L] should be reported as "less than 10 mg/dL [0.1 g/L]" instead of the numerical value.

Limitations of Procedure

The instrument reporting system contains error messages to warn the user of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

IGA Concentration	SD
86 mg/dL [0.86 g/L]	> 7.4 mg/dL [0.074 g/L]
298 mg/dL [2.98 g/L]	> 10.9 mg/dL [0.109 g/L]

Interfering Substances

Samples containing a monoclonal immunoglobulin may result in a condition of antigen excess and artificially decreased values. Since the presence of a monoclonal protein can normally be detected using protein electrophoresis, the validity of immunochemical results may be confirmed with an electrophoretic pattern.

Lipemia (Intralipid®) of 600 mg/dL [6.78 g/L] increases an IGA result of 87 mg/dL [0.87 g/L] by 53%.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Expected Values: 87 – 474 mg/dL [0.87 – 4.74 g/L]

This reference population applies to serum samples from 150 healthy adults (75 males ages 19 – 65 and 75 females ages 19 – 59).

This reference interval was calculated non-parametrically and represents the central 95% of the population.¹⁰

Each laboratory should establish its own reference interval for IGA as performed on the Dimension® system.

Specific Performance Characteristics¹

Material	Precision ^{a,b}		Standard Deviation (% CV)	
	Mean mg/dL [g/L]	Within-run	Total	
Serum Pool				
Level 1	261 [2.61]	3.9 [0.04] (1.5)	6.9 [0.07] (2.7)	
Bio-Rad Liquichek™ Immunology Control				
Level 1	110 [1.10]	2.0 [0.02] (1.8)	5.4 [0.05] (4.9)	
Level 2	294 [2.94]	4.7 [0.05] (1.6)	7.0 [0.07] (2.4)	

f. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).

g. Reproducibility testing was done in accordance with the CLSI/NCCLS Tentative Guideline for User Evaluation of Precision Performance of Clinical Chemistry Devices (EP5-T2, 1992).

h. Specimens at each level were analyzed in duplicate, twice a day, for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

Liquichek™ is a trademark of Bio-Rad Laboratories, Inc., Irvine, CA 92618.

Method Comparison

Comparative Method	Regression Statistics ⁱ		Correlation Coefficient	n ^j
	Slope	Intercept mg/dL [g/L]		
Beckman Array 360	0.95	-1.48 [-0.01]	0.993	224

i. Model equation for regression statistics is: [Result of Dimension® analyzer] = [Slope x comparable method result] + Intercept.

j. The range of IgA values in the correlation study was 5 – 1050 mg/dL [0.05 – 10.50 g/L].

Specificity

HIL Interference

The IGA method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	IGA Concentration mg/dL [g/L]	Bias ^b %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L] (monomer)	89 [0.89]	<10
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	91 [0.91]	<10
Lipemia (Intralipid®)	200 mg/dL [2.26 mmol/L]	87 [0.87]	<10

k. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances at the concentrations indicated have no significant effect (less than 10%) on the IGA method when added to a serum pool containing 214 mg/dL [2.14 g/L] of IgA.

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1323 µmol/L
Albumin	6 g/dL	60 g/L
Amikacin	15 mg/dL	256 µmol/L
Ascorbic Acid	3 mg/dL	170 µmol/L
Caffeine	10 mg/dL	515 µmol/L
Carbamazepine	12 mg/dL	508 µmol/L
Chloramphenicol	25 mg/dL	774 µmol/L
Chlordiazepoxide	2 mg/dL	67 µmol/L
Chlorpromazine	5 mg/dL	157 µmol/L
Cholesterol Ester	500 mg/dL	12.9 mmol/L
Cimetidine	10 mg/dL	396 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 75	2500 mg/dL	333 µmol/L
Diazepam	2 mg/dL	70 µmol/L
Digoxin	5 ng/mL	6.4 nmol/L
Erythromycin	20 mg/dL	273 µmol/L
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 µmol/L
Furosemide	2 mg/dL	61 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Heparin (Sodium)	8 U/mL	8000 U/L
Ibuprofen	40 mg/dL	1939 µmol/L
Lidocaine	6 mg/dL	256 µmol/L
Lithium Chloride	3.5 mg/dL	5.07 mmol/L
Nicotine	2 mg/dL	123 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	442 µmol/L
Phenobarbital	15 mg/dL	646 µmol/L
Phenytoin	10 mg/dL	396 µmol/L
Primidone	10 mg/dL	458 µmol/L
Propoxyphene	0.4 mg/dL	12 µmol/L
Salicylic Acid	50 mg/dL	3.6 mmol/L
Theophylline	25 mg/dL	1388 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Recovery

Recovery of reference material CRM470/RPPHS¹¹ ranged from 101.0 – 110.7%, with a mean recovery of 106.1%.

Analytical Sensitivity: ≤ 10 mg/dL [0.1 g/L]

The analytical sensitivity represents the low end of the assay range of the IGA method.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Benjamini E, Sunshine G. Leskowitz. Immunology: A Short Course, 3rd edition, Wiley, New York, NY, 1996.
2. Tomasi TB. The Discovery of Secretory IgA and the Mucosal Membrane System. Immunol. Today 1992;13:416.
3. Hayward AR. Immunodeficiency. In: Lachmann PJ and Peters DK (eds). Clinical Aspects of Immunology, vol II, 4th ed, Blackwell Scientific Publications, Oxford, UK, 1982.
4. Hobbs JR. Immunoglobulins in Clinical Chemistry. Adv. Clin. Chem. 1971;14:220.
5. Kyle RA and Lust JA. The Monoclonal Gammopathies (Paraproteins). Adv. Clin. Chem. 1990;28:146.
6. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
7. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
8. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
9. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, 3rd ed. Washington, DC: AACC Press, 2007; p 529.
10. National Committee for Clinical Laboratory Standards. How to Define, Determine, and Utilize Reference Intervals in the Clinical Laboratory. CLSI/NCCLS Document C 28-P, 1992, Villanova, PA.
11. Whicher JT et al. New International Reference Preparation for Proteins in Human Serum (RPPHS). Clin. Chem. 1994;40:934.

**Symbols Key
Symbolschlüssel
Explication des Symboles
Interpretazione simboli
Clave de los Símbolos**

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EPGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare