## Dimension® clinical chemistry system Heterogeneous Immunoassay Module

### Flex® reagent cartridge

### MMB

See shaded sections: Updated information from 2017-08 version.


Issue Date 2019-04-01

#### Mass Creatine Kinase MB Isoenzyme

**Intended Use:** The MMB method for Dimension® clinical chemistry system with the heterogeneous immunoassay module is an *in vitro* diagnostic test intended to quantitatively measure the MB isoenzyme of creatine kinase (EC 2.7.3.2) in human serum and plasma for confirmation of acute myocardial infarction.

**Summary:** The creatine kinase MB isoenzyme (CKMB) is found primarily in cardiac tissue, with substantially lower concentrations also seen in skeletal muscle.<sup>1</sup> The quantitation of CKMB is routinely ordered as part of the cardiac panel and is useful in the diagnosis of acute myocardial infarction (AMI).<sup>2</sup> Typically, in cases of uncomplicated AMI, serial determinations show a pattern wherein CKMB levels become elevated within 4 – 8 hours after onset of pain, peak between 12 – 24 hours and then drop to normal by 48 hours. Mass CKMB is the biochemical marker of choice for perioperative myocardial infarction during the first 48 hours after the onset of pain.<sup>3</sup> CKMB concentrations have also been used to assess the extent of AMI and subsequent reinfarction.<sup>4,5,6</sup> The diagnostic sensitivity, specificity and efficiency of mass CKMB is superior to that of CK isoenzymes by electrophoresis.<sup>7,8</sup>

**Principles of Procedure:** The MMB method is a one-step enzyme immunoassay based on the "sandwich" principle. The sample is incubated with chromium dioxide particles coated with monoclonal antibodies<sup>1</sup> specific for CKB subunit, and conjugate reagent (β-galactosidase labeled monoclonal antibodies<sup>1</sup> specific for CKMB isoenzyme).<sup>9</sup> A particle/CKMB/conjugate sandwich forms during the incubation period. Unbound conjugate is removed by magnetic separation and washing. The sandwich bound β-galactosidase is combined with a chromogenic substrate chlorophenol red-β-d-galactopyranoside (CPRG). Hydrolysis of CPRG releases a chromophore (CPR). The concentration of CKMB present in the patient sample is directly proportional to the rate of color change due to formation of CPR measured at 577 nm. The amount of CKMB protein is measured immunologically and the results are reported in mass units (ng/mL or µg/L).

#### Reagents

Wells <sup>a</sup>	Form	Ingredient	Concentration <sup>b</sup>	Source
1, 2	Liquid	Conjugate	c	Mouse, monoclonal
3, 4	Tablet <sup>d</sup>	Antibody-CrO <sub>2</sub>	0.9 mg/mL <sup>c</sup>	Mouse, monoclonal
5, 6	Tablet <sup>d</sup>	CPRG	15.7 mM	
7	Liquid	Substrate Diluent		
		Buffer	175 mM	
8	Liquid	Chrome Diluent		
		Buffer	200 mM	

- a. Wells are numbered consecutively from the wide end of the cartridge.  
 b. Nominal value in hydrated cartridge.  
 c. Antibody titer and conjugate activity may vary from lot to lot.  
 d. Tablets contain excipients, buffers, and stabilizers.

#### Risk and Safety:

H317  
 P280, P272, P302 + P352, P333 + P313, P501

#### Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

**Contains:** 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone

Safety data sheets (MSDS/SDS) available on [siemens.com/healthcare](http://siemens.com/healthcare)

**Precautions:** Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

**Reagent Preparation:** Hydrating, diluting and mixing are automatically performed by the instrument.

**Store at:** 2 – 8 °C

**Expiration:** Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge when on the instrument are stable for 30 days.

**Open Well Stability:** 3 days for wells 1 – 6  
 10 days for well 7 – 8

**Specimen Collection and Handling:** Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.<sup>10</sup> Follow the instructions provided with your specimen collection device for use and processing.<sup>11</sup>

EDTA, lithium heparin and sodium heparin anticoagulants do not interfere with the MMB method.

Blood collected in the presence of oxalate potentially can cause clumping of the chrome particles and should not be used.

Complete clot formation should take place before centrifugation.<sup>12</sup>

Specimens should be free of particulate matter. To prevent the appearance of fibrin in serum samples, complete clot formation should take place before centrifugation. If clotting time is increased due to thrombolytic or anticoagulant therapy, the use of plasma specimens will allow for faster sample processing and reduce the risk of particulate matter.

Separated samples may be stored at room temperature for 12 hours, at 4 °C for 3 days and at -20 °C for 1 month.<sup>13</sup>

Frozen plasma samples with insufficient anticoagulants potentially can cause clumping of chrome particles and should not be used.

#### Procedure

#### Materials Provided

MMB Flex® reagent cartridge, Cat. No. RF420

#### Materials Required But Not Provided

MMB Calibrator, Cat. No. RC420  
 Reaction Vessels, Cat. No. RXV1A  
 Chemistry Wash, Cat. No. RD701  
 Sample Diluent, Cat. No. 791092901  
 Quality Control Materials

#### Test Steps

Sampling, reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system with the heterogeneous immunoassay module. For details of this processing, refer to your Dimension® Operator's Guide.

#### Test Conditions

Reaction Vessel	
Sample Size	60 µL
Antibody-CrO <sub>2</sub>	50 µL
MMB Ab-β-galactosidase	50 µL
Incubation Temperature	42 °C*
Incubation Period	8 minutes (sample addition to wash)
Wash Steps	5 minutes (3 washes)

\*Dimension® EXL™ with LOCI® Module: 37 °C

#### Cuvette

Transfer Volume	40 µL
Substrate Reagent Volume (CPRG)	91 µL
Diluent Volume	279 µL
Temperature	37 °C
Incubation Period	4 minutes (sample addition to read)
Wavelength	577 and 700 nm
Type of Measurement	Bichromatic rate

#### Calibration

Assay Range 0.5 – 300 ng/mL [µg/L]<sup>†</sup>  
 Calibration Material MMB Calibrator, Cat. No. RC420  
 Calibration Scheme Levels 1, 2 n = 3  
 Level 3, n = 2  
 Level 4, n = 3  
 Level 5, n = 2

#### Units

Typical Calibration Levels 0.0, 10.0, 50.0, 150.0, 330.0 ng/mL [µg/L]  
 Calibration Frequency Every 60 days for any one lot

- For each new lot of Flex® reagent cartridges
- After major maintenance or service, if indicated by quality control results
- As indicated in laboratory quality control procedures
- When required by government regulations

#### Assigned Coefficients

C<sub>0</sub> -130  
 C<sub>1</sub> 1800  
 C<sub>2</sub> -3.7  
 C<sub>3</sub> 750  
 C<sub>4</sub> 0.5

e. Système International d'Unités [SI Units] are in brackets.

#### Quality Control

At least once each day of use analyze two levels of a Quality Control material with known MMB concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

**Results:** The instrument automatically calculates and prints the concentration of MMB result in ng/mL [µg/L].

**Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.**

**Analytical Measurement Range (AMR): 0.5 – 300 ng/mL [µg/L]**

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 300 ng/mL [µg/L] should be repeated on dilution.

**Manual Dilution:** Results in excess of 300 ng/mL [ $\mu\text{g/L}$ ] should be repeated after diluting the sample with Sample Diluent to obtain results within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

**Autodilution (AD):** Refer to your Dimension® Operator's Guide.

Results less than 0.5 ng/mL [ $\mu\text{g/L}$ ] should be reported as "less than 0.5 ng/mL [ $\mu\text{g/L}$ ]" instead of the numerical value.

#### Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

Patient samples may contain heterophilic antibodies that could react with immunoassays to give falsely elevated or depressed results. This assay has been designed to minimize interference from heterophilic antibodies.<sup>5</sup>

Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed. A test result that is inconsistent with the clinical picture and patient history should be interpreted with caution.<sup>14</sup>

A system malfunction may exist if the following 5-test precision is observed:

MMB Concentration	SD
0 ng/mL [ $\mu\text{g/L}$ ]	> 1.0 ng/mL [ $\mu\text{g/L}$ ]
10 ng/mL [ $\mu\text{g/L}$ ]	> 1.1 ng/mL [ $\mu\text{g/L}$ ]
300 ng/mL [ $\mu\text{g/L}$ ]	> 20.0 ng/mL [ $\mu\text{g/L}$ ]

**Expected Values:** Levels of CKMB in normal individuals are low to undetectable. The reference interval was determined from a population of healthy adults. A total of 257 serum samples, comprised of 142 females and 115 males, ages 19 – 76 were tested. A mean of 1.03 ng/mL [ $\mu\text{g/L}$ ] was obtained with an SD of 1.06 ng/mL [ $\mu\text{g/L}$ ]. The reference interval for MMB was calculated non-parametrically and represents the central 95% of the population.

The calculated reference interval for the Dimension® immunoassay system is 0 – 3.6 ng/mL [ $\mu\text{g/L}$ ].

Each laboratory should establish its own reference intervals for MMB as performed on the Dimension® system.

#### Diagnostic Criteria

	MB isoenzyme of creatine	Relative Index (RI)*
Non-AMI	$\leq 5$	NA
Gray Zone	$> 5$	$\leq 4$
AMI	$> 5$	$> 4$

\*RI = (mass CKMB/CK) x 100.

The current recommendation for diagnosis of AMI from the World Health Organization (WHO) is based on the presence of two of the following criteria: chest pain, electrocardiographic (ECG) changes, and cardiac enzyme abnormalities.<sup>15,16</sup> Of the cardiac enzymes, CKMB is currently considered the definitive serum marker.

Following an AMI without reperfusion, CKMB is typically detected in the serum 4 – 8 hours after the onset of pain. CKMB peaks in the serum 12 – 24 hours and returns to normal levels 48 – 72 hours after the MI. Since elevated CKMB results can be obtained in disorders or disease states other than AMI, no single result can be considered diagnostic.<sup>17,18</sup> A sequential sampling protocol to demonstrate the classical rise and fall of CKMB should be followed.<sup>17,19,20</sup> Sampling intervals should include time of admission, 8 – 12 hours and 24 hours after admission.<sup>21,22</sup> For very early detection of myocardial infarction an accelerated (two-hour) sequential sampling protocol has been suggested.<sup>23</sup> To minimize myocardial tissue damage, thrombolytic therapy is commonly administered within 4 hours of the onset of pain.<sup>24</sup>

For the reperfused patient the CK and CKMB isoenzyme peaks are earlier and the enzyme activity is greater.<sup>25</sup> Sampling protocol for the reperfused patient is suggested at the time of admission, the time of thrombolytic therapy, and then at 3, 6, and 12 hours after the therapy to document the early rise and early peak associated with successful reperfusion.<sup>26</sup>

For the clinical diagnosis of acute myocardial infarction, a sequential sampling protocol is recommended. A cutoff of  $> 5$  ng/mL [ $\mu\text{g/L}$ ] and  $> 4$  RI is suggested as being consistent with AMI. Since the diagnostic sensitivity and specificity as well as predictive value are population dependent, this cutoff should be confirmed in your laboratory. Because of the overlap of results between the non-AMI and AMI populations, the discriminating power and predictive value of the result will vary with the cutoff selected.

For elevated total CK and CKMB results, a relative index (RI) can be used to differentiate CKMB elevations due to skeletal muscle injury from AMI. A cutoff of  $< 1.5$  RI is suggested as being consistent with skeletal muscle injury.<sup>27</sup>

#### Specific Performance Characteristics<sup>1</sup>

Material	Mean ng/mL [ $\mu\text{g/L}$ ]	Precision <sup>a,b</sup>	
		Standard Deviation (% CV)	
		Within-run	Between-day
Serum Pool			
Level 1	1.27	0.29 (23.0)	0.29 (23.0)
Level 2	7.76	0.27 (3.5)	0.69 (8.8)
Level 3	13.07	0.32 (2.4)	0.57 (4.4)
Dade® CK-MB/Myoglobin Immunoassay Control			
Level 1	3.58	0.24 (6.7)	0.35 (9.7)
Level 2	16.69	0.43 (2.6)	1.10 (6.6)
Level 3	46.05	0.64 (1.4)	1.86 (4.0)

f. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).

g. Reproducibility testing was done in accordance with the CLSI/NCCLS Tentative Guideline for User Evaluation of Precision Performance of Clinical Chemistry Devices (EP5-T2, 1992).

h. Specimens at each level were analyzed in duplicate for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

#### Method Comparison

Comparative Method	Regression Statistics <sup>1</sup>			n
	Slope	Intercept ng/mL [ $\mu\text{g/L}$ ]	Correlation Coefficient	
Stratus®	1.28	-1.85	0.994	88 <sup>f</sup>
IMx® Stat CKMB	0.85	0.52	0.998	132 <sup>k</sup>
AxSYM®	0.80	-0.60	0.999	102 <sup>l</sup>

i. Model equation for regression statistics is: Result of Dimension® analyzer = [Slope x (comparative method result)] + Intercept.

j. The range of MMB values in the correlation study was 0.5 to 118.2 ng/mL [ $\mu\text{g/L}$ ].

k. The range of MMB values in the correlation study was 0.5 to 283.4 ng/mL [ $\mu\text{g/L}$ ].

l. The range of MMB values in the correlation study was 0.5 to 190.8 ng/mL [ $\mu\text{g/L}$ ].

IMx® and AxSYM® are registered trademarks of Abbott Laboratories, Abbott Park, IL 60064.

#### Specificity

##### Non-Interfering Substances

No clinically significant interference was observed from icterus (bilirubin 60 mg/dL [1026 mmol/L]), hemolysis (hemoglobin 1000 mg/dL [0.62 mmol/L]) or lipemia (triglyceride 1500 mg/dL [16.93 mmol/L]).

The following substances have no significant effect (less than 10%) on the MMB method when added to a 5 ng/mL MMB serum pool at the concentrations indicated.

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1323 $\mu\text{mol/L}$
Amikacin	15 mg/dL	256 $\mu\text{mol/L}$
Ampicillin	5 mg/dL	143 $\mu\text{mol/L}$
Ascorbic Acid	3 mg/dL	170 $\mu\text{mol/L}$
Caffeine	10 mg/dL	515 $\mu\text{mol/L}$
Carbamazepine	12 mg/dL	508 $\mu\text{mol/L}$
Chloramphenicol	25 mg/dL	774 $\mu\text{mol/L}$
Chlordiazepoxide	2 mg/dL	67 $\mu\text{mol/L}$
Chlorpromazine	5 mg/dL	157 $\mu\text{mol/L}$
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	10 mg/dL	396 $\mu\text{mol/L}$
Creatinine	30 mg/dL	2652 $\mu\text{mol/L}$
Dextran 75	1250 mg/dL	167 $\mu\text{mol/L}$
Diazepam	2 mg/dL	70 $\mu\text{mol/L}$
Digoxin	5 ng/mL	6.4 nmol/L
Erythromycin	20 mg/dL	273 $\mu\text{mol/L}$
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 $\mu\text{mol/L}$
Furosemide	2 mg/dL	61 $\mu\text{mol/L}$
Gentamicin	12 mg/dL	251 $\mu\text{mol/L}$
Heparin (sodium)	8 U/mL	8000 U/L
Ibuprofen	40 mg/dL	1939 $\mu\text{mol/L}$
Lidocaine	6 mg/dL	256 $\mu\text{mol/L}$
Lithium	3.5 mg/dL	5.04 mmol/L
Nicotine	2 mg/dL	123 $\mu\text{mol/L}$
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	442 $\mu\text{mol/L}$
Phenobarbital	15 mg/dL	646 $\mu\text{mol/L}$
Phenytoin	10 mg/dL	396 $\mu\text{mol/L}$
Primidone	10 mg/dL	458 $\mu\text{mol/L}$
Protein, Human IgG	3 g/dL	30 g/L
Propoxyphene	0.4 mg/dL	12 $\mu\text{mol/L}$
Salicylic Acid	50 mg/dL	3.6 mmol/L
Theophylline	25 mg/dL	1388 $\mu\text{mol/L}$
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic acid	50 mg/dL	3467 $\mu\text{mol/L}$

#### Hook Effect

One-step sandwich immunometric assays are susceptible to a high-dose "hook effect," where an excess of antigen prevents simultaneous binding of the capture and detection antibodies to a single analyte molecule.<sup>28</sup> Such samples must be diluted and reassayed prior to reporting the results (see dilution information). The MMB method shows no hook effect up to at least 2000 ng/mL [ $\mu\text{g/L}$ ].

#### Recovery

Samples containing known amounts of CKMB were mixed with human serum samples with baseline MMB values of 0.1 – 0.5 ng/mL. The sample concentrations were measured and the percent recovery calculated. Recovery ranged from 94.2 to 104% with a mean recovery of 98%.

$$\% \text{ Recovery} = \frac{\text{Value obtained}}{(\text{high sample} + \text{low sample}) / 2} \times 100$$

#### Analytical Sensitivity: 0.5 ng/mL [ $\mu\text{g/L}$ ]

The analytical sensitivity represents the lowest concentration of MMB that can be distinguished from zero. This sensitivity is defined as the concentration at two standard deviations above the 0 ng/mL [ $\mu\text{g/L}$ ] MMB Calibrator (n = 20).

**Symbols Key:** See adjacent panel.

**Bibliography:** See adjacent panel.

Dimension®, Flex™, EXL™, LOCI® and Stratus® are trademarks of Siemens Healthcare Diagnostics

©2009 Siemens Healthcare Diagnostics

All rights reserved.

Siemens Healthcare Diagnostics Inc.  
500 GBC Drive  
Newark, DE 19714 USA

## Dimension® clinical chemistry system Heterogeneous Immunoassay Module

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Siemens Healthcare Diagnostics Inc.  
500 GBC Drive  
Newark, DE 19714 USA

**Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía**

1. Tsung SH. Creatine Kinase Isoenzyme Patterns in Human Tissue Obtained at Surgery. *Clin Chem* 1976;22:173-175.
2. Lott JA. Serum Enzyme Determinations in the Diagnosis of Acute Myocardial Infarction: An Update. *Hum Pathol* 1984;15:706-716.
3. Gulbis B. Mass Concentration of Creatine Kinase MB Isoenzyme and Lactate Dehydrogenase Isoenzyme 1 in Diagnosis of Perioperative Myocardial Infarction after Coronary Bypass Surgery. *Clin Chem* 1990;36:1784-1788.
4. Galen RS, Gambino SR. Isoenzymes of CPK and LDH in Myocardial Infarction and Certain Other Diseases. *Pathobiol Ann* 1975;5:283-315.
5. Vaidya HC, Maynard Y, Dietzler DN, Ladenson JH. Direct Measurement of Creatine Kinase-MB Activity in Serum After Extraction with a Monoclonal Antibody Specific to the MB Isoenzyme. *Clin Chem* 1986;32:657-663.
6. Vaidya HC, Beatty BG. Eliminating Interference from Heterophilic Antibodies in a Two- Site Immunoassay for Creatine Kinase MB by using F(ab)<sub>2</sub> Conjugate and Polyclonal Mouse IgG. *Clin Chem* 1992;38:1737-1742.
7. Ryan F, Scott D, Gregory M, Turner M. Measurement of Creatine Kinase MB by the Baxter Stratus® Fluorometric Analyzer. *Clin Chem* 1989;35:1123.
8. Lee K, Kypridakis G. Evaluation of the Stratus® CK-MB Assay in a Community Hospital. *Clin Chem* 1990;36:1-133.
9. Birkmeyer R C, Diaco R, Hutson DK., et al. Applications of Novel Chromium dioxide Magnetic Particles to Immunoassay Development. *Clin Chem* 1987;33:1543-1547.
10. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
11. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
12. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
13. Tietz NW. *Textbook of Clinical Chemistry*, 3<sup>rd</sup> ed. Carl A. Burtis, Edward R. Ashwood, W.B. Saunders Co., Philadelphia, PA 1999:662.
14. Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. *Clin Chem* 1999; 45:7:942-956.
15. World Health Organization criteria for the diagnosis of acute myocardial infarction. Proposal for the multinational monitoring of trends and determinants in cardiovascular disease. Geneva Cardiovascular Disease Unit of World Health Organization; 1981.
16. Wyngaarden JB, Smith LH. *Cecil Textbook of Medicine*. W. B. Saunders Co., Philadelphia, PA 1988:330-332.
17. Henry JB. *Clinical Diagnosis and Management by Laboratory Methods*. W. B. Saunders Co., Philadelphia, PA 1991:273-274.
18. Sasse EA, Madiedo G, Kopinski W. Evaluation of Abbott IMx CKMB Immunoassay. *Clin Chem* 1990;36:1858-1859.
19. Chan KM, Ladenson JH, Pierce GF. Increased Creatine Kinase MB in the Absence of Acute Myocardial Infarction. *Clin Chem* 1986;32:2044-2051.
20. Lee TH, et al. Evaluation of Creatine Kinase and Creatine Kinase-MB for Diagnosing Myocardial Infarction. Clinical Impact in the Emergency Room. *Arch Intern Med* 1987;147:115-120.
21. Koumanhkis G, Radcliff FJ. Diagnostic sensitivity and specificity of serum enzymes and isoenzymes: optimization for the diagnosis of myocardial infarction. *Austra J of Med Lab Sci* 1986;7:114-121.
22. Moerman J, Nellens P, Beckers F, Sierens L. Diagnostic value of an immunoenzymetric assay compared with an enzymic immunoinhibition procedure for creatine kinase isoenzyme MB. *Clin Chim Acta* 1986;160:313-317.
23. Marin MM, Teichman S L. Use of rapid serial sampling of creatine kinase MB for very early detection of myocardial infarction in patients with acute chest pain. *Am Heart J* 1992;123:354-361.
24. Sharkey SW, Brunette DD, Ruiz E, Hession WT et al. An Analysis of Time Delays Preceding Thrombolysis for Acute Myocardial Infarction. *J Am Med Assoc* 1989;262:3171-3174.
25. Apple FS, Sharkey SW, Werdick M, Elsparger KJ, Tilbury RT. Analysis of Creatine Kinase Isoenzymes and Isoforms in Serum to Detect Reperfusion after Acute Myocardial Infarction. *Clin Chem* 1987;33:507-511.
26. Apple FS. Creatine Kinase-MB. *Lab Med* 1992;23:298-302.
27. Skogen W, Osborne D, Price G. Guidelines for Interpreting the CK-MB Relative Index Value. *Clin Chem* 1991;37:916.
28. Ryall RG, Story CJ, Turner DR. Reappraisal of the causes of the “hook effect” in two-site immunoradiometric assays. *Anal Biochem* 1982;127:308.

† U.S. Patent No. Re. 33,405

**Symbols Key  
Symbol Schlüssel  
Explication des Symboles  
Interpretazione simboli  
Clave de los Símbolos**

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
<b>LOT</b>	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
<b>REF</b>	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante
<b>EC REP</b>	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos
<b>IVD</b>	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
<b>CE</b>	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
<b>CONTENTS</b>	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
<b>LEVEL</b>	Level / Konzentration / Niveau / Livello / Nivel

2014-10\_EFR6S

Siemens Healthcare Diagnostics Inc.  
500 GBC Drive  
Newark, DE 19714 USA

Global Siemens  
Headquarters  
Siemens AG  
Wittelsbacherplatz 2  
80333 Muenchen  
Germany

Global Siemens  
Healthcare Headquarters  
Siemens AG  
Healthcare Sector  
Henkestrasse 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0  
siemens.com/healthcare

Global Division  
Siemens Healthcare  
Diagnostics Inc.  
511 Benedict Avenue  
Tarrytown, NY 10591  
USA  
siemens.com/healthcare