Dimension® clinical chemistry system

Flex® reagent cartridge

GGT

See shaded sections: Updated information from 2018-09 version.


Issue Date 2019-04-01

γ-Glutamyl Transferase

Intended Use: The GGT method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of γ-glutamyl transferase activity in human serum and plasma.

Summary: The γ-glutamyl transferase method is an adaptation of the methodology recommended by the International Federation of Clinical Chemistry (IFCC).¹ The method uses the substrate L-gamma-glutamyl-3-carboxy-4-nitranilide with glycylglycine.

Principles of Procedure: γ-glutamyl transferase catalyzes the transfer of the glutamyl moiety from γ-glutamyl-3-carboxy-4-nitranilide (GCNA) to glycylglycine thereby releasing 5-amino-2-nitrobenzoate which absorbs at 405 nm. This change is proportional to the γ-glutamyl transferase activity and is measured using a bichromatic (405, 600 nm) rate technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1 – 6	Tablet ^c	GCNA	0.69 mg
7 – 8	Liquid	Glycylglycine	90 mM/L

- a. Wells are numbered consecutively from the wide end of the cartridge.
b. Nominal value per test at manufacture.
c. Tablet contains excipients.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501

Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 2-chloroacetamide

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: Hydrating, diluting and mixing are automatically performed by the instrument.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6
30 days for wells 7 – 8

Specimen Collection and Handling: Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method.²

Follow the instructions provided with your specimen collection device for use and processing.³

Clot formation should take place before centrifugation. Separated serum/plasma specimens are stable for 8 hours at room temperature, 2 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C or colder.⁴

Blood collection tubes containing EDTA, lithium heparin, potassium oxalate or sodium fluoride do not interfere in this method.

Blood collection tubes containing a combination of potassium oxalate and sodium fluoride will decrease GGT activity by 15%.

Corvac® and SST® collection tubes do not affect the GGT method.

Corvac® is a registered trademark of Monoject, Division of Sherwood Medical, St. Louis, MO.
SST® is a registered trademark of Becton-Dickinson, Rutherford, NJ.

Procedure

Materials Provided

GGT Flex® reagent cartridge, Cat. No. DF45A

Materials Required But Not Provided

Enzyme Verifier, Cat. No. DC19
Quality Control Materials

Test Steps

Sampling,^{4a} reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

- d. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.
e. An alternate sample size of 15 µL can be programmed; refer to the Dimension® Operator's Guide for the use of alternate sample size

Test Conditions

Sample Size	32 µL, (15 µL) ^e
Reagent 1 Volume	75 µL
Diluent Volume	313 µL
Temperature	37° C
Wavelength	405 and 600 nm
Type of Measurement	Bichromatic rate

Verification

Assay Range (@ 37° C)	0 – 800 U/L
Verification Material	Enzyme Verifier, Cat. No. DC19
Verification Scheme	3 levels, n = 3
Units	U/L
Typical Verification Levels	70, 400, 700 U/L
Verification Frequency	Every 3 months for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	Standard sample size = 32 µL C ₀ -1.000 C ₁ 3.360 Alternate sample size = 15 µL C ₀ -8.000 C ₁ 7.230

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known γ-glutamyl transferase activity.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of γ-glutamyl transferase in U/L using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0 – 800 U/L

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range. Samples with results in excess of 800 U/L should be repeated on dilution.

Manual Dilution: Results in excess of 800 U/L should be repeated after diluting the sample with Enzyme Diluent (Cat. No. 790035901) or equivalent to produce a sample within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): Refer to your Dimension® Operator's Guide.

Samples with results less than 7 U/L should be reported as "less than 7 U/L".

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed at the standard sample size (32 µL):

GGT Activity	SD
80 U/L	>3 U/L
400 U/L	>7 U/L

Interfering Substances

Hemolysis (hemolysate) of 500 mg/dL [0.31 mmol/L]^l (monomer) and above tripped a test report message; therefore the magnitude of the interference could not be determined.

Bilirubin (unconjugated) of 60 mg/dL [1024 µmol/L] and above tripped a test report message; therefore the magnitude of the interference could not be determined.

Triglycerides as (Intralipid®) of 200 mg/dL [2.299 mmol/L] decreases a GGT at 128 U/L by 16%.

f. Système International d'Unités [SI Units] are in brackets.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Expected Values

Females 5 – 55 U/L

Males 15 – 85 U/L

This reference population consisted of the following:

112 males, ages 18 – 60

121 females, ages 18 – 60

The reference interval was calculated non-parametrically and represents the central 95% of the population.

Heavy drinkers have been reported to have elevated γ-glutamyl transferase levels.^{5,6}

Each laboratory should establish its own reference interval for γ-glutamyl transferase as performed on the Dimension® system.

Specific Performance Characteristics^a

Material	Precision ^a		
	Mean U/L	Standard Deviation (% CV)	
		Within-run	Between-day
Multiqua [®]			
Nor Unasy	31	1.0 (3.3)	1.3 (4.2)
Abn Unasy	112	0.9 (0.8)	1.4 (1.2)
Moni-Trol [®] Control ^l			
Level 1	26	1.9 [7.4]	2.7 [10.5]
Level 2	106	2.4 [2.3]	2.3 [2.2]

g. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (Refer to your Dimension® Operator's Guide).

h. Specimens at each level were analyzed in duplicate for 20 days. The within-run and between-day standard deviations were calculated by the analysis of variance method.

i. Using reduced sample size (15 µL).

Multiqua[®] is a registered trademark of Bio-Rad., Irvine, CA 92714, USA.

Moni-Trol[®] is a registered trademark of Medical Analysis Systems Inc., Camarillo, CA 93012-8058, USA.

Comparative Method	Method Comparison		
	Slope	Regression Statistics	
		Intercept U/L	Correlation Coefficient
GGT method on the aca® discrete clinical analyzer ^l	1.114	-5.6	0.999
Reduced vs. Standard sample size ^k	1.015	-1.0	0.999
			n
			104
			56 ^l

j. Model equation for regression statistics is:

Results of Dimension® GGT = [slope x comparative method results] + intercept.

k. Model equation for regression statistics is: Results of Dimension® GGT using reduced sample size (15 µL) = [slope x results of Dimension® GGT system using standard sample size (32 µL)] + intercept.

l. Range of samples: 12 – 545 U/L.

Specificity

HIL Interference

The GGT method (using the standard sample size of 32 µL) was evaluated for interference from hemolysis, icterus, and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	GGT Conc. U/L	Bias ^m %
Hemoglobin	300 mg/dL [0.19 mmol/L] (monomer)	129	<10
(hemolysate)	500 mg/dL [0.31 mmol/L] (monomer)	129	n
Bilirubin	40 mg/dL [684 µmol/L]	129	<10
(unconjugated)			
Lipemia	50 mg/dL [0.57 mmol/L]	128	<10
(Intralipid®)	200 mg/dL [2.29 mmol/L]	128	-16

m. Analyte results should not be corrected based on this bias.

n. Hemolysis (hemolysate) of 500 mg/dL [0.31 mmol/L] (monomer) and above tripped an error flag, so the magnitude of the interference is not available.

The GGT method (using an alternate sample size of 15 µL) was evaluated for interference from hemolysis, icterus, and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI units	GGT Activity U/L	Bias %
Hemoglobin	300 mg/dL [0.19 mmol/L] (monomer)	87	<10
(hemolysate)	500 mg/dL [0.31 mmol/L] (monomer)	87	-14
Bilirubin	60 mg/dL [1026 µmol/L]	89	<10
(unconjugated)	80 mg/dL [1368 µmol/L]	89	13
Lipemia	50 mg/dL [0.57 mmol/L]	86	<10
(Intralipid®)	200 mg/dL [2.29 mmol/L]	86	-25

Non-Interfering Substances

The following substances have no measurable effect on the GGT method at the concentrations indicated:

Substance	Test Concentration	SI Units
Acetaminophen	200 µg/mL	1323 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	20 µg/mL	57 µmol/L
Ascorbic Acid	5 mg/dL	284 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	20 µg/mL	70 µmol/L
Digoxin	20 ng/mL	25.6 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	800 mg/dL	174 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	16 µg/mL	29.4 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Nortriptyline	1000 ng/mL	3797 nmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	80 µg/mL	344 µmol/L
Phenytoin	30 µg/mL	119 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Salicylate	100 mg/dL	7.24 mmol/L
Theophylline	100 µg/mL	555 µmol/L
Triglycerides (endogenous)	3000 mg/dL	33.9 mmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Analytical Sensitivity: 7 U/L

The analytical sensitivity represents the lowest activity of GGT that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations of Reagent grade water.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, aca® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics

All rights reserved.

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Shaw LM, Stromme JH, London JL, Theodorsen L. IFCC methods for the determination of enzymes part 4. IFCC method for gamma-glutamyl transferase [(gamma-glutamyl)-peptide; amino acid gamma-glutamyl transferase, EC 2.3.2.2]. Clinica Chimica Acta 1983; 15F-338F.
2. Tietz NW. Textbook of Clinical Chemistry, W. B. Saunders Co., Philadelphia, PA 1986, pp 478-497 (specimen collection and storage recommendations).
3. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
4. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
5. Rollason JG, Pincherle G, and Robinson D. Serum – gamma glutamyl transpeptidase in relation to alcohol consumption, Clinica chimica acta; international journal of clinical chemistry (Netherlands), Jun 1972, 39, (1), p 75-80.
6. Rosalki SB, and Rau D. Serum g -glutamyl transpeptidase activity in alcoholism, Clinica chimica acta; international journal of clinical chemistry (Netherlands), Jun 1972, 39, (1), p 41-47.

Symbols Key Symbol Schlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

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Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare