## Dimension® clinical chemistry system

### Flex® reagent cartridge

MBI

See shaded sections: Updated information from 2017-08 version.

Issue Date 2019-04-08

#### Creatine Kinase MB

**Intended Use:** The creatine kinase MB (MBI) method is an *in vitro* diagnostic test for the quantitative measurement of creatine kinase MB isoenzyme activity in human serum and plasma on the Dimension® clinical chemistry system.

**Summary:** The Dimension MBI method is a modification of the International Federation of Clinical Chemistry (IFCC) creatine kinase (CK) primary reference 37 °C procedure, adapted for use on the Dimension® clinical chemistry system.<sup>1</sup>


The MB isoenzyme of creatine kinase is used as a marker of cardiac muscle cell damage, and may be especially useful to detect re-infarction after a prior event. Normally an intracellular enzyme, CK-MB begins to appear in blood about 4 to 6 hours after the onset of myocardial damage, peaking at about 24 hours after a single event. CK-MB activity declines and returns to normal within about three days.<sup>2,3</sup>

Calculation of the relative percentage of CK-MB isoenzyme activity in the total CK activity can aid in distinguishing acute coronary syndromes from non-cardiac events.<sup>4</sup>

**Principles of Procedure:** The activity of the CK-MM isoenzyme is inhibited by an antibody specific for the CK-M subunit. The activity of the B subunit of creatine kinase MB isoenzyme is not inhibited, and it is on this basis that CK-MB can be measured.

In an enzyme coupled reaction, creatine kinase in patient samples catalyzes the transphosphorylation of creatine phosphate to adenosine-diphosphate (ADP), producing adenosine-triphosphate (ATP). Hexokinase (HK) uses the ATP to phosphorylate glucose. The resulting glucose-6-phosphate is oxidized by glucose-6-phosphate dehydrogenase (G-6-PDH) with the simultaneous reduction of nicotinamide adenine dinucleotide phosphate (NADP) to NADPH.

The rate of formation of NADPH is measured bichromatically at 340, 540 nm and is directly proportional to CK-B activity in the sample.

#### Reagents

Wells*	Form	Ingredient	Concentration <sup>b</sup>	Source
1 – 3 <sup>c</sup>	Liquid	NADP	2.46 mmol/L	
Reagent 1		ADP	2.46 mmol/L	
		AMP	6.14 mmol/L	
		Diadenosine Pentaphosphate (AP5A)	19 µmol/L	
		N-acetyl cysteine	24.6 mmol/L	
		Hexokinase	≥ 36.7 µkat/L	Yeast
		G-6-PDH	≥ 23.4 µkat/L	Bacterial
		EDTA	2.46 mmol/L	
		Mg acetate	12.3 mmol/L	
	Imidazole buffer	123 mmol/L		
6 <sup>c</sup>	Liquid	Glucose	120 mmol/L	
Reagent 2		Creatine phosphate	184 mmol/L	
		Anti CK-M antibody	0.24 mg/mL	Mouse
		Mg Acetate	11 mmol/L	
		3-(cyclohexylamino)-2-hydroxyl-1-propane sulfonic acid (CAPSO Buffer)	20 mmol/L	
		EDTA	2.46 mmol/L	

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

c. Wells 1 – 3 and 6 contain stabilizers and preservatives.

#### Risk and Safety

H360d  
P201, P202, P280, P281, P308 + P313, P501

**Danger!**  
May damage unborn child.

Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Wear protective gloves/protective clothing/eye protection/face protection. Use personal protective equipment as required. If exposed or concerned: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

**Contains:** Imidazole

Safety data sheets (MSDS/SDS) available on [siemens.com/healthcare](http://siemens.com/healthcare)

**Precautions:** Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

**Reagent Preparation:** All reagents are liquid and ready to use.

**Store at:** 2 – 8 °C

**Expiration:** Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

**Open Well Stability:** 5 days for wells 1 – 3  
15 days for well 6

**Specimen Collection and Handling:** Recommended specimen types: serum, sodium and/or lithium heparin plasma, and EDTA.

Grossly hemolyzed samples should not be used with the MBI method.

Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.<sup>5</sup>

Follow the instructions provided with your specimen collection device for use and processing.<sup>6</sup>

For serum, complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.<sup>7</sup> Separated serum/plasma samples should be stored at 2 – 8 °C and analyzed within 24 hours. For longer storage, samples should be frozen at or below -20 °C.<sup>7</sup>

#### Procedure

##### Materials Provided

MBI Flex® reagent cartridge, Cat. No. DF32

##### Materials Required But Not Provided

Dimension® CKI/MBI Calibrator, Cat. No. DC32  
Quality control materials

##### Test Steps

Sampling<sup>d</sup>, reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® clinical chemistry system. For details of this processing, refer to your Dimension® Operator's Guide.

d. The sample container must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

##### Test Conditions

Sample Volume	20 µL
Reagent 1 Volume	288 µL
Reagent 2 Volume	58 µL
Temperature	37.0 °C ± 0.1 °C
Reaction Time	8.7 minutes <sup>e</sup>
Wavelengths	340 and 540 nm
Type of Measurement	Bichromatic rate

e. Calculated as time from test initiation to final result

##### Calibration

Assay Range	3 –125 U/L [0.05 – 2.08 µkat/L] <sup>f</sup>
Calibration Material	Dimension® CKI/MBI Calibrator, Cat. No. DC32
Calibration Scheme	3 levels in triplicate
Units	U/L [µkat/L] (U/L ÷ 60) = [µkat/L]

Typical Calibration Levels	Level 1: 0 U/L [0.00 µkat/L] <sup>g</sup>
	Level 2: 66 U/L [1.10 µkat/L]
	Level 3: 138 U/L [2.30 µkat/L]

Calibration Frequency: Every 90 days for any one lot

- A new calibration is required
- For each new lot of Flex® reagent cartridges
  - After major maintenance or service, if indicated by quality control results
  - As indicated in laboratory quality control procedures
  - When required by government regulations

Assigned Coefficients	C <sub>0</sub>	-1.000
	C <sub>1</sub>	14.400

f. Système International d'Unités [SI Units] are in brackets.

g. Level 1 calibrator for CKI/MBI is not included in the CKI/MBI carton. Purified Water Diluent (Cat No. 710615901) or Reagent grade water should be used as the level 1 calibrator for the MBI method.

##### Quality Control

Follow government regulations or accreditation requirements for quality control frequency.

At least once each day of use, analyze two levels of a Quality Control (QC) material with known creatine kinase concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

**Results:** The instrument calculates the activity concentration of creatine kinase in U/L [µkat/L] using the calculation scheme described in your Dimension® Operator's Guide.

**Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.**

**Analytical Measurement Range (AMR): 3 – 125 U/L [0.05 – 2.08 µkat/L]**

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- **Serum/plasma samples** with results in excess of 125 U/L [2.08 µkat/L] are reported as “Above Assay Range” and should be repeated on dilution.

**Manual Dilution:** Dilute with Reagent grade water to obtain results within the analytical measurement range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.

**Autodilution (AD):** The autodilute sample volume is 5 µL (dilution factor = 4) for serum/plasma. Refer to your Dimension® Operator’s Guide.

- **Serum/plasma** with results less than 3 U/L [0.05 µkat/L] should be reported as “less than 3 U/L [0.05 µkat/L]”.

**Limitations of Procedure**

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in creatine kinase results. Refer to your Dimension® Operator’s Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory’s procedure manual and not reported.

Venipuncture should occur prior to sulfapyridine administration due to the potential for falsely depressed results.

A system malfunction may exist if the following 5-test precision is observed:

MBI Activity Concentration	SD
15 U/L [0.25 µkat/L]	>3.0 U/L [0.05 µkat/L]
50 U/L [0.83 µkat/L]	>3.0 U/L [0.05 µkat/L]

**Interfering Substances**

The MBI was evaluated for interference according to CLSI/NCCLS EP7-A2.<sup>8</sup> Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Hemoglobin (hemolysate) at 30 mg/dL [0.019 mmol/L] increases MBI results by 26% at a creatine kinase activity of 23 U/L [0.38 µkat/L].

Lipemia (Intralipid®) at 1000 mg/dL [11.3 mmol/L] increases MBI results by 18% at a creatine kinase activity of 28 U/L [0.47 µkat/L].

Methotrexate at 100 mg/dL [2201 µmol/L] increases MBI results by 35% at a creatine kinase activity of 23 U/L [0.38 µkat/L].

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

**Expected Values: 7 – 25 U/L [0.12 – 0.42 µkat/L]<sup>9,10</sup>**

Each laboratory should establish its own expected values for MBI as performed on the Dimension® clinical chemistry system.

**Specific Performance Characteristics**

The following data represent typical performance for the Dimension® clinical chemistry system.

Material	Precision <sup>11,h</sup>		
	Mean U/L [µkat/L]	Standard Deviation (% CV)	
		Repeatability	Within-Lab
<b>Serum</b>			
Serum Pool 1	26 [0.434]	1.0 [0.017] (3.9)	1.3 [0.022] (5.0)
<b>BioRad® Liquechek™ Cardiac Marker Control LT</b>			
Level 2	15 [0.251]	1.5 [0.025] (9.8)	1.7 [0.028] (11.5)
Level 3	92 [1.536]	3.1 [0.052] (3.4)	3.4 [0.057] (3.7)

h. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Liquechek™ is a trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

Comparative Method	Method Comparison <sup>12</sup> Regression Statistics <sup>i</sup>			
	Slope	Intercept U/L [µkat/L]	Correlation Coefficient	n
Roche CKMB Hitachi® 917	1.01	-1.3 [-0.022]	0.986	98 <sup>j</sup>

i. CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was orthogonal regression fit.

j. The range of plasma creatine kinase – MB values in the correlation study was 5 to 122 U/L.

Roche/Hitachi CKMB is a product of F. Hoffmann-La Roche Ltd, Diagnostics Division, Grenzacherstrasse 124, CH-4070 Basel, Switzerland.

**Specificity****HIL Interference**

The MBI was evaluated for interference according to CLSI/NCCLS EP7-A2.<sup>8</sup> Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration SI Units	Creatine Kinase MB U/L [µkat/L]	Bias* %
Hemoglobin (hemolysate)	20 mg/dL [0.012 mmol/L]	30 U/L [0.50] 53 U/L [0.88]	< 10
Bilirubin (unconjugated)	20 mg/dL [342 µmol/L]	22 U/L [0.37] 53 U/L [0.88]	< 10
Bilirubin (conjugated)	40 mg/dL [684 µmol/L]	20 U/L [0.33] 53 U/L [0.88]	< 10
Lipemia (Intralipid®)	600 mg/dL [6.8 mmol/L]	20 U/L [0.33] 53 U/L [0.88]	< 10

\*Analyte results should not be corrected based on this bias.

**Non-Interfering Substances – serum/plasma**

The following substances do not interfere with the MBI method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at MBI method concentrations of 20 U/L [0.33 µkat/L].

Substance	Test Concentration	SI Units
Acetaminophen	20.0 mg/dL	1324 µmol/L
Allopurinol	2.5 mg/dL	184 µmol/L
Amikacin	8.0 mg/dL	137 µmol/L
Amiodarone	2.5 mg/dL	3.9 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	6.0 mg/dL	342 µmol/L
Atenolol	1.0 mg/dL	37.6 µmol/L
Caffeine	6.0 mg/dL	308 µmol/L
Captopril	5.0 mg/dL	230 µmol/L
Carbamazepine	3.0 mg/dL	127 µmol/L
Chloramphenicol	5.0 mg/dL	155 µmol/L
Chlordiazepoxide	1.0 mg/dL	33.3 µmol/L
Chlorpromazine	0.20 mg/dL	6.27 µmol/L
Cholesterol	503 mg/dL	13 mmol/L
Cimetidine	2.0 mg/dL	79.2 µmol/L
Cinnarizine	3.0 mg/dL	81.4 µmol/L
Creatine	30.0 mg/dL	2.7 mmol/L
Cyclosporin A	4000 ng/mL	3.3 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.51 mg/dL	18 µmol/L
Digitoxin	350 ng/mL	458 nmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
Diltiazem	120 µg/mL	266.1 µmol/L
Disopyramide	4.0 mg/dL	118 µmol/L
Dopamine	65 mg/dL	3428 µmol/L
Erythromycin	6.0 mg/dL	81.6 mmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25.0 mg/dL	1770 µmol/L
Furosemide	6.0 mg/dL	181 µmol/L
Gentamicin	1.0 mg/dL	21 µmol/L
Heparin	3.0 U/mL	3000 U/L
Ibuprofen	50.0 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5000 mg/dL	50 g/L
Immunoglobulin M (IgM)	5000 mg/dL	50 g/L
Isosorbide dinitrate	6.0 mg/dL	254 µmol/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lisinopryl	16 µg/mL	36.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
Lovastatin	16 µg/mL	39.6 µmol/L
L-thyroxine	60 µg/dL	0.77 µmol/L
Methotrexate	60 mg/dL	1321 µmol/L
Methylolopa	2.5 mg/dL	118 µmol/L
Methylprednisolone	4.0 mg/dL	107 µmol/L
Mexiletine	24 mg/dL	1339 µmol/L
N-acetyl-procainamide	30 mg/dL	1082 µmol/L
Nicotine	0.10 mg/dL	6.2 µmol/L
Nifedipine	6.0 mg/dL	173 µmol/L
Nitroglycerin	0.16 µg/mL	0.58 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8.0 mg/dL	354 µmol/L
Phenobarbital	10.0 mg/dL	431 µmol/L
Phenytoin	5.0 mg/dL	198 µmol/L
Primidone	4.0 mg/dL	183 µmol/L
Procainamide	10 mg/dL	443 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Propranolol	0.5 mg/dL	19 µmol/L
Protein: Albumin	6000 mg/dL	60 g/L
Protein: Total	12000 mg/dL	120 g/L
Quinidine	20 µg/mL	61.6 µmol/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Simvastatin	32 µg/mL	76.5 µmol/L
Streptokinase	300 IU/mL	300,000 IU/L
Theophylline	4.0 mg/dL	222 µmol/L
Tocainide	10 mg/dL	437 µmol/L
Urea	500 mg/dL	83 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L
Verapamil	16 mg/dL	0.33 µmol/L
Warfarin	10 mg/dL	324 µmol/L

**Limit of Detection:<sup>13</sup> 3 U/L [0.05 µkat/L]**

The limit of detection represents the lowest concentration of creatine kinase MB that can be detected with at least 95% probability. CLSI/NCCLS EP17-A, Protocols for Determination of Limits of Detection and Limits of Quantitation, was followed using 288 replicate determinations at two concentration levels to determine Limit of Blank and Limit of Detection.

**Symbols Key:** See adjacent panel.

**Bibliography:** See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics  
All rights reserved.

Siemens Healthcare Diagnostics Inc.  
500 GBC Drive  
Newark, DE 19714 USA

## Dimension® clinical chemistry system

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Siemens Healthcare Diagnostics Inc.  
500 GBC Drive  
Newark, DE 19714 USA

**Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía**

1. Siekmann L, Ceriotti F, Ferard G, Kanno T, Schumann G, et. al. IFCC Primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 2. Reference procedure for the measurement of catalytic concentrations of creatine kinase [ATP: creatine N-phosphotransferase (CK), EC 2.7.3.2] Clin Chem Lab Med. 2002; 40(6):635-642.
2. McPherson RA, Pincus MR, eds. Henry's clinical diagnosis and management by laboratory methods 21st ed. St. Louis: Saunders Elsevier, 2007: 221, 260.
3. Anderson JL, Adams CD, Antman EM, Bridges Cr, Califf RM, Casey DE, et. al. ACC/ACC 2007 Guidelines for the management of patients with unstable angina/non-ST-elevation myocardial infarction-executive summary. JACC. 50(7). 2007:652-726.
4. Burtis CA, Ashwood ER, eds. Tietz fundamentals of clinical chemistry, 5th ed. Philadelphia: WB Saunders Co., 2001:358-62.
5. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA, 2003.
6. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
7. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimen; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1- 56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
8. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
9. Klein G, Berger A, Bertholf R et al. Abstract: Multicenter Evaluation of Liquid Reagents for CK, CK-MB and LDH with Determination of Reference Intervals on Hitachi Systems. Clin Chem 2001; 47:Suppl. A30.
10. Thomas L, Müller M, Schumann G, Weidemann G et al. Consensus of DGKL and VDGH for interim reference intervals on enzymes in serum. J Lab Med 2005;29:301–308.
11. Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition*. CLSI/NCCLS Document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
12. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
13. Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

<b>Symbols Key</b> <b>Symbolschlüssel</b> <b>Explication des Symboles</b> <b>Interpretazione simboli</b> <b>Clave de los Símbolos</b>	
	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabricante / Fabricante
	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contents sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenido suficiente per "n" saggi / Contenido suficiente para <n> ensayos
	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
	Level / Konzentration / Niveau / Livello / Nivel

2014-10\_EFH5

Siemens Healthcare Diagnostics Inc.  
500 GBC Drive  
Newark, DE 19714 USA

Global Siemens  
Headquarters  
Siemens AG  
Wittelsbacherplatz 2  
80333 Muenchen  
Germany

Global Siemens  
Healthcare Headquarters  
Siemens AG  
Healthcare Sector  
Henkestrasse 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0  
siemens.com/healthcare

Global Division  
Siemens Healthcare  
Diagnostics Inc.  
511 Benedict Avenue  
Tarrytown, NY 10591  
USA  
siemens.com/healthcare