Dimension® clinical chemistry system

Flex® reagent cartridge

DBI

See shaded sections: Updated information from 2018-12 version.

Issue Date 2019-04-30

Direct Bilirubin

Intended Use: The DBI method for the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended to quantitatively measure direct (conjugated) bilirubin in human serum and plasma. Measurements of direct bilirubin are used in the diagnosis and treatment of liver, hemolytic, hematological, and metabolic disorders, including hepatitis and gall bladder disease.

Summary: There are at least four distinct bilirubin fractions that make up total bilirubin in serum. The direct reacting fractions are mono- and diconjugated bilirubin (β - and γ -bilirubin) and the delta fraction (δ -bilirubin), which is tightly bound to albumin. Unconjugated bilirubin (α -bilirubin) is water-insoluble and reacts only after addition of an accelerator such as caffeine.¹ The DBI method is a modification of the Doumas reference method,² which is a modification of the diazo method described by Jendrassik and Grof in 1938.³

Principles of Procedure: Diazotized sulfanilic acid is formed by combining sodium nitrite and sulfanilic acid at low pH. The sample is diluted in 0.5M HCl. A blank reading is taken to eliminate interference from non-bilirubin pigments. Upon addition of the diazotized sulfanilic acid, the conjugated bilirubin is converted to diazo-bilirubin, a red chromophore which absorbs at 540 nm and is measured using a bichromatic (540, 700 nm) endpoint technique.

Conjugated bilirubin + Diazotized sulfanilic acid $\longrightarrow$ Red chromophore (absorbs at 540 nm)

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1 – 4		Empty ^c	
5	Liquid	Sodium nitrite	72.5 mM
6	Liquid	Hydrochloric acid	500 mM
7 – 8	Liquid	Sulfanilic acid	25.89 mM
		Hydrochloric acid	132 mM

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

c. Diazotized sulfanilic acid will be prepared in these wells automatically by the instrument.

Risk and Safety

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes or vessels contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use.

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 2 days for wells 1 – 4 (once the diazotized sulfanilic acid has been prepared)
30 days for wells 5 – 8

Specimen Collection and Handling: Serum, lithium heparin plasma, and EDTA plasma can be collected by normal procedures.^{4,5}

Follow the instructions provided with your specimen collection device for use and processing.⁶

Serum and plasma specimens should be separated from cells within 2 hours after venipuncture.⁴

Specimens should be free of particulate matter. To prevent the appearance of fibrin in serum samples, complete clot formation should take place before centrifugation.⁸ Clotting time may be increased due to thrombolytic or anticoagulant therapy.

Bilirubin is photosensitive. Care should be taken to protect sample from both daylight and fluorescent light to avoid photodegradation. Separated specimens are stable for 8 hours at room temperature, 7 days at 2 – 8 °C or 6 months frozen at -20 °C or colder. Protection from light is required when specimens are stored for more than 8 hours.⁷

Hemolysis may depress DBI results. Follow your laboratory's procedures for reporting results when the sample is hemolyzed.

Procedure

Materials Provided

DBI Flex® reagent cartridge, Cat. No. DF125

Materials Required But Not Provided

TBI/DBI Calibrator, Cat. No. DC167

Quality Control Materials

Purified Water Diluent (Cat. No. 710615901) or Reagent grade water

Test Steps

Sampling,⁴ reagent delivery, mixing, processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

d. The sample container must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

Test Conditions

	Cuvette 1
Sample Volume	10 μ L
Reagent 1 Volume	25 μ L
Reagent 2 Volume	50 μ L
Temperature	37 °C
Wavelength	540, 700 nm
Type of Measurement	Bichromatic endpoint

Calibration

Assay Range	0.05 – 16.0 mg/dL [0.86 – 274 μ mol/L] ^e
Calibration Material	TBI/DBI Calibrator, Cat. No. DC167
Calibration Scheme	3 levels, n = 3
Units	mg/dL [μ mol/L] (mg/dL x 17.1) = [μ mol/L]
Typical Calibration Levels	0.0, 7.0, 17.5 mg/dL [0, 120, 299 μ mol/L]
Calibration Frequency	Every 90 days for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ 0.06 C ₁ 0.075

e. Système International d'Unités [SI Units] are in brackets.

Note: Level 1 calibrator for DBI is not included in the TBI/DBI calibrator carton. Purified Water Diluent (Cat. No. 710615901) or Reagent grade water should be used as the Level 1 calibrator for the DBI method.

Quality Control

At least once each day of use, analyze two levels of a Quality Control material (QC) with known direct bilirubin concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of direct bilirubin in mg/dL [μ mol/L] using the calculation scheme described in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.05 – 16.0 mg/dL [0.86 – 274 μ mol/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 16.0 mg/dL [274 μ mol/L] should be repeated on dilution.

Manual Dilution: Dilute with Reagent grade water 1:2 (1 part of sample with 1 part of Reagent grade water) to obtain results within assay range. Enter dilution factor (2). Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): The recommended autodilute sample volume is 5 μ L for serum and plasma. Refer to your Dimension® Operator's Guide.

Results less than 0.05 mg/dL [0.86 μ mol/L] should be reported as "less than 0.05 mg/dL [0.86 μ mol/L]".

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

DBI Concentration	SD
0.6 mg/dL [10.3 μ mol/L]	> 0.06 mg/dL [1.0 μ mol/L]
16.0 mg/dL [287 μ mol/L]	> 0.34 mg/dL [5.8 μ mol/L]

Interfering Substances

The DBI method was evaluated for interference according to CLSI/NCCLS EP7-A.⁹ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in mg/dL [$\mu\text{mol/L}$]. Bias exceeding 10% is considered interference.

Interferent	Interferent Concentration SI Units	Direct Bilirubin SI Units	Bias SI units	Bias' %
Albumin	6 g/dL [60 g/L]	0.4 mg/dL [7 $\mu\text{mol/L}$]	-0.04 mg/dL [-0.7 $\mu\text{mol/L}$]	-11
Albumin	6 g/dL [60 g/L]	5.8 mg/dL [99 $\mu\text{mol/L}$]	-0.8 mg/dL [-14 $\mu\text{mol/L}$]	-13
Ascorbic acid	5 mg/dL [227 $\mu\text{mol/L}$]	0.3 mg/dL [5 $\mu\text{mol/L}$]	+0.03 mg/dL [0.5 $\mu\text{mol/L}$]	+11
Carbenicillin	3 mg/dL [7.1 mmol/L]	0.4 mg/dL [7 $\mu\text{mol/L}$]	+0.04 mg/dL [0.7 $\mu\text{mol/L}$]	+11
Cholesterol	500 mg/dL [12.9 mmol/L]	0.3 mg/dL [5 $\mu\text{mol/L}$]	-0.03 mg/dL [-0.5 $\mu\text{mol/L}$]	-11
Hemoglobin	20 mg/dL [0.013 mmol/L] (monomer)	0.4 mg/dL [7 $\mu\text{mol/L}$]	-0.2 mg/dL [-3 $\mu\text{mol/L}$]	-44
Hemoglobin	50 mg/dL [0.03 mmol/L] (monomer)	0.3 mg/dL [5 $\mu\text{mol/L}$]	-0.25 mg/dL [-4 $\mu\text{mol/L}$]	-83
Hemoglobin	50 mg/dL [0.03 mmol/L] (monomer)	3.0 mg/dL [51 $\mu\text{mol/L}$]	-0.7 mg/dL [-12 $\mu\text{mol/L}$]	-23
Hemoglobin	50 mg/dL [0.03 mmol/L] (monomer)	5.0 mg/dL [86 $\mu\text{mol/L}$]	-0.8 mg/dL [-14 $\mu\text{mol/L}$]	-16
Hemoglobin	100 mg/dL [0.06 mmol/L] (monomer)	0.3 mg/dL [5 $\mu\text{mol/L}$]	-0.3 mg/dL [-5 $\mu\text{mol/L}$]	-100
Hemoglobin	100 mg/dL [0.06 mmol/L] (monomer)	3.0 mg/dL [51 $\mu\text{mol/L}$]	-1.0 mg/dL [-17 $\mu\text{mol/L}$]	-33
Hemoglobin	100 mg/dL [0.06 mmol/L] (monomer)	5.0 mg/dL [86 $\mu\text{mol/L}$]	-1.3 mg/dL [-22 $\mu\text{mol/L}$]	-26
Hemoglobin	100 mg/dL [0.06 mmol/L] (monomer)	14.1 mg/dL [241 $\mu\text{mol/L}$]	-1.6 mg/dL [-27 $\mu\text{mol/L}$]	-11
Immunoglobulin G (IgG)	5 g/dL [50 g/L]	0.4 mg/dL [7 $\mu\text{mol/L}$]	-0.08 mg/dL [-1 $\mu\text{mol/L}$]	-20
Levodopa	300 $\mu\text{g/mL}$ [1.52 mmol/L]	0.4 mg/dL [7 $\mu\text{mol/L}$]	+0.3 mg/dL [5 $\mu\text{mol/L}$]	+76
Lipemia (Intralipid®)	200 mg/dL [2.26 mmol/L]	0.4 mg/dL [7 $\mu\text{mol/L}$]	-0.09 mg/dL [-2 $\mu\text{mol/L}$]	-22
Oxytetracycline	50 mg/dL [0.10 mmol/L]	0.4 mg/dL [7 $\mu\text{mol/L}$]	+0.05 mg/dL [0.9 $\mu\text{mol/L}$]	+13
Phloroglucinol	1500 ng/mL [11.9 $\mu\text{mol/L}$]	0.2 mg/dL [2.4 $\mu\text{mol/L}$]	+0.04 mg/dL [0.5 $\mu\text{mol/L}$]	+17
Rheumatoid factor (RF)	510 IU/mL [510 IU/mL]	0.4 mg/dL [7 $\mu\text{mol/L}$]	-0.06 mg/dL [-1 $\mu\text{mol/L}$]	-14
Total protein	12 g/dL [120 g/L]	0.4 mg/dL [7 $\mu\text{mol/L}$]	-0.04 mg/dL [-0.7 $\mu\text{mol/L}$]	-11
Total protein	12 g/dL [120 g/L]	5.8 mg/dL [99 $\mu\text{mol/L}$]	-1.4 mg/dL [-24 $\mu\text{mol/L}$]	-13

f. Analyte results should not be corrected based on this bias.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Expected Values: 0.0 – 0.2 mg/dL [0 – 3 $\mu\text{mol/L}$]¹⁰

The reference interval was validated in a confirmatory study using 30 serum samples.

Each laboratory should establish its own reference interval for conjugated bilirubin as performed on the Dimension® system.

Specific Performance Characteristics

All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).

Material	Precision ^{11g, h}		Standard Deviation (% CV)	
	Mean mg/dL [$\mu\text{mol/L}$]	Repeatability	Within-Lab	
Serum pool	18.2 [311]	0.27 [6] (1.5)	0.47 [9] (2.6)	
Bio-Rad Lyphochek® Control Level 1	0.3 [5]	0.01 [0.2] (1.8)	0.01 [0.2] (2.9)	
Bio-Rad Lyphochek® Control Level 2	1.5 [25]	0.03 [0.5] (2.3)	0.04 [0.7] (2.8)	
MAS® Liquid Bilirubin QC Level 3	6.3 [108]	0.1 [2] (1.5)	0.2 [3] (2.4)	

g. Reproducibility testing was done in accordance with the CLSI/NCCLS Approved Guideline for Evaluation of Precision Performance of Quantitative Measurement Methods (EP5-A2, 2004).

h. Specimens at each level were analyzed in duplicate, twice a day, for 20 days. The repeatability (within-run) and within-lab (total) standard deviations were calculated by analysis of variance method.

Lyphochek® is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618.

MAS® is a registered trademark of Medical Analysis Systems, Inc., Camarillo, CA 93012.

Comparative Method	Method Comparison Regression Statistics ⁱ		Correlation Coefficient	r^{1k}
	Slope	Intercept mg/dL [$\mu\text{mol/L}$]		
Dimension® DBIL	1.02	-0.02 [-0.3]	1.000	156
Vitros™ BuBc (conjugated)	0.911	-0.096 [-1.6]	0.989	156

i. Model equation for regression statistics is: results for Dimension® system = [slope x comparative method results] + intercept.

j. The range of DBI values in the Dimension® correlation study was: 0.04 to 15.7 mg/dL [0.7 to 269 $\mu\text{mol/L}$].

k. The range of DBI values in the Vitros™ BuBc correlation study was: 0.02 to 15.7 mg/dL [0.3 to 269 $\mu\text{mol/L}$].

Vitros™ is a trademark of Ortho Clinical Diagnostics, Inc., Rochester, NY 14626.

Recommended Sample Types

Matched human serum and plasma samples were analyzed for direct bilirubin. As shown in the table below, no clinically significant difference was observed between different sample types using ordinary least squares analysis to fit the regression line.

Sample Type	n	Slope (m)	Intercept (b)	Correlation Coefficient (r)
Lithium heparin plasma versus serum	25	0.98	0.23	0.989
EDTA plasma versus serum	25	0.95	0.08	0.991
EDTA plasma versus lithium heparin plasma	25	0.97	-0.13	0.999

Hemolysis and Lipemia Information

The DBI method was evaluated for interference according to CLSI/NCCLS EP7-A. Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Test Concentration SI Units	DBI Concentration SI Units %	Bias'
Hemoglobin ^m	20 mg/dL [0.13 mmol/L] (monomer)	0.4 mg/dL [7 $\mu\text{mol/L}$]	-44
Lipemia (Intralipid®) ⁿ	50 mg/dL [0.57 mmol/L]	0.4 mg/dL [7 $\mu\text{mol/L}$]	<10
	200 mg/dL [2.26 mmol/L]	0.4 mg/dL [7 $\mu\text{mol/L}$]	-22

i. Analyte results should not be corrected for this bias.

m. Samples with hemolysis greater than 50 mg/dL [0.03 mmol/L] of hemoglobin will be flagged with a "Hemoglobin" error message. Refer to your Dimension® Operator's Guide.

n. Samples with Intralipid® concentrations greater than 600 mg/dL [6.78 mmol/L] may be flagged with an "Abnormal Reaction" error message. Refer to your Dimension® Operator's Guide.

Non-Interfering Substances

The following substances do not interfere with the DBI method when present in serum in the amounts indicated. Systematic inaccuracies (bias) due to these substances are less than 10% at a direct bilirubin concentration of 0.31 – 0.37 mg/dL [5.1 – 6.3 $\mu\text{mol/L}$] and 4.5 – 5.8 mg/dL [71 – 98 $\mu\text{mol/L}$].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1328 $\mu\text{mol/L}$
Amikacin	15 mg/dL	256 $\mu\text{mol/L}$
4-aminosalicylic acid	80 $\mu\text{g/dL}$	5 $\mu\text{mol/L}$
Amphotericin B	4 $\mu\text{g/mL}$	4 $\mu\text{mol/L}$
Ampicillin	5.3 mg/dL	152 $\mu\text{mol/L}$
Ascorbic acid ^o	5 mg/dL	227 $\mu\text{mol/L}$
Biliverdin	4 mg/dL	0.061 mmol/L
Caffeine	6 mg/dL	308 $\mu\text{mol/L}$
Carbamazepine	3 mg/dL	127 $\mu\text{mol/L}$
Carbenicillin ^p	3 mg/dL	7.1 mmol/L
Cefotiam	840 $\mu\text{g/mL}$	1.4 mmol/L
Chloramphenicol	5 mg/dL	155 $\mu\text{mol/L}$
Chlordiazepoxide	1 mg/dL	33.3 $\mu\text{mol/L}$
Chlorpromazine	1 mg/dL	31.4 $\mu\text{mol/L}$
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 $\mu\text{mol/L}$
Creatinine	30 mg/dL	2652 $\mu\text{mol/L}$
Dextran 40 ^q	6000 mg/dL	1500 $\mu\text{mol/L}$
Diazepam	0.5 mg/dL	17.6 $\mu\text{mol/L}$
Digoxin	5 ng/mL	6.15 nmol/L
Erythromycin	6 mg/dL	81.6 $\mu\text{mol/L}$
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 $\mu\text{mol/L}$
Furosemide	6 mg/dL	181 $\mu\text{mol/L}$
Gentamicin	12 mg/dL	251 $\mu\text{mol/L}$
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 $\mu\text{mol/L}$
Immunoglobulin G ^p	5 g/dL	50 g/L
Levodopa ^p	300 $\mu\text{g/mL}$	1.52 mmol/L
Lidocaine	1.2 mg/dL	51.2 mmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Methotrexate	50 $\mu\text{g/mL}$	110 $\mu\text{mol/L}$
Nicotine	0.1 mg/dL	6.2 $\mu\text{mol/L}$
Nitrofurantoin	20 $\mu\text{g/mL}$	0.08 mmol/L
Oxytetracycline ^p	50 mg/dL	0.1 mmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	443 $\mu\text{mol/L}$
Phenazopyridine	80 $\mu\text{g/mL}$	320 $\mu\text{mol/L}$
Phenobarbital	15 mg/dL	647 $\mu\text{mol/L}$
Phenytoin	5 mg/dL	198 $\mu\text{mol/L}$
Pindolol	5 mg/L	0.02 mmol/L
Piroxicam	10 $\mu\text{g/mL}$	30 $\mu\text{mol/L}$
Primidone	4 mg/dL	183 $\mu\text{mol/L}$
Propoxyphene	1 mg/dL	25 $\mu\text{mol/L}$
Reserpine	60 mg/L	0.10 mmol/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Sulfasalazine	40 $\mu\text{g/mL}$	100 $\mu\text{mol/L}$
Theophylline	4 mg/dL	222 $\mu\text{mol/L}$
Triamterene	60 $\mu\text{g/mL}$	237 $\mu\text{mol/L}$
Triglycerides ^q	724 mg/dL	8.18 mmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 $\mu\text{mol/L}$
Valproic Acid	50 mg/dL	3467 $\mu\text{mol/L}$

o. Testing performed at a DBI concentration of 4.5 – 5.8 mg/dL [71 – 98 $\mu\text{mol/L}$].

p. Interference testing for Dextran 40 was performed only at a DBI concentration of 0.4 mg/dL [7 $\mu\text{mol/L}$].

q. Interference testing for Triglyceride was performed only at a DBI concentration of 0.3 mg/dL [5 $\mu\text{mol/L}$].

Analytical Sensitivity: 0.05 mg/dL [0.86 $\mu\text{mol/L}$]

The analytical sensitivity represents the lower limit of detection (lowest concentration of direct bilirubin that can be distinguished from zero). This sensitivity is defined as the mean value (n = 20) plus two standard deviations of the low level (0 mg/dL [0 $\mu\text{mol/L}$]) TBI/DBI Calibrator. Reagent grade water should be used as the Level 1 calibrator for the DBI method.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Dimension® clinical chemistry system

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Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia: W.B. Saunders Company, 1999: p. 1136.
2. Dumas BT, Perry BW, Sasse EA, et al., Standardization in bilirubin assays: Evaluations of selected methods and stability of bilirubin solutions, Clin Chem 1973; 19:984-993.
3. Jendrassik L, Grof P. Vereinfachte photometrische Methode zur Bestimmung des Blutbilirubin, Biochem Z 938; 297-81.
4. Burtis CA, Ashwood ER. Tietz Fundamentals of Clinical Chemistry, 5th ed. Philadelphia: W.B. Saunders Company, 2001: pp. 31-41.
5. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
6. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
7. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, 3rd ed. Washington, DC: AACC Press, 2007: p 175.
8. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline–Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
9. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline*. CLSI/NCCLS document EP7-A [ISBN 1-56238-480-5]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2002.
10. Kaplan LA, Pesce AJ. Clinical Chemistry: Theory, Analysis, and Correlation, 3rd ed. St. Louis, MO: Mosby, Inc, 1996: p. 712.
11. Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2004.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFIS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare