Dimension® clinical chemistry system

Flex® reagent cartridge

AMY

See shaded sections: Updated information from 2017-11 version.


Issue Date 2019-04-01

Amylase

Intended Use: The AMY method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of amylase activity in human serum, plasma and urine.

Summary: The AMY method on the Dimension® system utilizes a chromogenic substrate, 2-chloro-4-nitrophenol linked with maltotriose.¹ The direct reaction of α-amylase with the substrate results in the formation of 2-chloro-4-nitrophenol, which is monitored spectrophotometrically. Amylase measurements are used primarily for the diagnosis and treatment of pancreatitis. The AMY method responds to both pancreatic and salivary amylase isoenzymes.

Principles of Procedure: α-amylase (α-1, 4-glucan, 4-glucanohydrolase; EC 3.2.1.1) catalyzes the hydrolysis of a defined synthetic substrate, 2-chloro-4-nitrophenyl-α-D-maltotriose (CNP3), to yield 2-chloro-4-nitrophenol (CNP), 2-chloro-4-nitrophenyl-α-D-maltoside (CNP2), maltotriose (G3) and glucose. After an incubation of 70 seconds at 37 °C, the absorbance due to the formation of 2-chloro-4-nitrophenol (CNP) is measured using a bichromatic (405, 577 nm) rate technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1–6	Liquid	CNP3	1.24 mmol/L

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per test at manufacture.

Risk and Safety

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Contains sodium azide (< 0.1%) as a preservative. Sodium azide can react with copper or lead pipes in drain lines to form explosive compounds. Dispose of properly in accordance with local regulations.

Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6

Specimen Collection and Handling: Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.²

Blood collection tubes containing EDTA, citrate and oxalate have been reported to inhibit α-amylase activity and should not be used.³

Corvac® and SST® collection tubes and tubes containing sodium fluoride do not affect the AMY method.

Follow the instructions provided with your specimen collection device for use and processing.⁴

Complete clot formation should take place before centrifugation.⁵

Specimens should be free of particulate matter.

Separated samples are stable for 7 days at room temperature and six months at 2 to 8 °C. For longer storage specimens may be frozen at -20 °C or colder.⁶ Urine amylase is unstable in acid urine. Adjust urine to a pH of 7.0 then store refrigerated.^{7,8}

The purpose of specimen storage information is to provide guidance to users; however, users may validate their own procedures for storing patient samples.

Albumin must be added to all urine specimens to maximize amylase activity.^{9,10} The final albumin concentration in urine should be at least 3.0 g/dL [30 g/L].^c Refer to Manual Dilution.

c. Système International d'Unités [SI Units] are in brackets.

Corvac® is a registered trademark of Monoject, Division of Sherwood Medical, St. Louis, MO. SST® is a registered trademark of Becton-Dickinson, Rutherford, NJ.

Procedure

Materials Provided

AMY Flex® reagent cartridge, Cat. No. DF17A

Materials Required But Not Provided

Enzyme Verifier, Cat. No. DC19
Quality Control Materials

Test Steps

Sampling,^d reagent delivery, mixing, and processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

d. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus the dead volume; precise container filling is not required.

Test Conditions

Sample Size	14 µL (10 µL) ^e
Reagent Volume	220 µL
Diluent Volume	166 µL
Temperature	37 °C
Wavelength	405 and 577 nm
Type of Measurement	Bichromatic rate
e. An alternate sample size of 10 µL can be programmed; refer to the Operator's Guide for the use of alternate sample size.	

Verification

Assay Range (@37 °C)	0 – 650 U/L
Verification Material	Enzyme Verifier, Cat. No. DC19
Verification Scheme	3 levels, n = 3
Units	U/L
Typical Verification Levels	60, 400, 725 U/L
Verification Slope Range	0.90 – 1.10
Verification Frequency	Every 3 months for any one lot
A new verification is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	Standard Sample size = 14 µL C ₀ 0.000 C ₁ 5.400 Alternate Sample size = 10 µL C ₀ 0.000 C ₁ 7.560

Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once daily analyze two levels of a Quality Control (QC) material with known amylase activity. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

To run urine quality control materials, treat material as a urine specimen described under specimen collection.

Results: The instrument automatically calculates and prints the activity of amylase in U/L using the calculation scheme illustrated in your Dimension® Operator's Guide. A change of 0.2 milli-absorbance units (mAU) per minute corresponds to an α-amylase activity of 1 U/L at 37 °C.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0 – 650 U/L

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 650 U/L should be repeated on dilution.

Manual Dilution:	Serum/Plasma: Make appropriate dilution with Enzyme Diluent (Cat. No. 790035901) or equivalent to obtain result within the assay range. Enter dilution factor.
	Urine: Dilute 1 part urine: 1 part Enzyme Diluent or equivalent. Enter dilution factor of 2. If readout is not within assay range, make appropriate dilution to obtain result within the assay range. Enter dilution factor.
	Serum/Plasma/Urine: Reassay. Resulting readout is corrected for dilution.

Autodilution (AD)

(for serum, plasma):	Refer to your Dimension® Operator's Guide.
(for urine):	Not recommended for urines.

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed at the standard sample size (14 µL):

Activity	SD
50 U/L	>4 U/L
600 U/L	>10 U/L

Interfering Substances

The AMY method was evaluated for interference according to CLSI/NCCLS EP7-A2.¹¹ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Hemoglobin at 1000 mg/dL [0.62 mmol/L] (monomer) decreased AMY results of 141 U/L by 19% (at the standard sample size of 14 µL) and decreased AMY results of 108 U/L by 12% (at an alternate sample size of 10 µL).

Immunoglobulin G at 5 g/dL [50 g/L] increased AMY results of 134 U/L by 32%.

Lipemia (Intralipid®) at 3000 mg/dL [33.9 mmol/L] and above tripped an error flag on this method, so the magnitude of the interference is not available.

Total protein at 12 g/dL [120 g/L] increased AMY results of 134 U/L by 95%.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Expected Values (@37 °C):

Serum: ¹²	25 – 115 U/L
Urine: ¹¹³	59 – 401 U/24 hr
AMY/CREA Clearance Ratio: ¹⁴	1.3 – 4.3%

The reference interval was calculated non-parametrically and represents the central 95% of the population. The reference population consisted of the following:

serum study	118 adults
urine study	164 adults
AMY/CREA clearance ratio study	107 adult random serum and urines

A urine reference interval in U/L is of limited usefulness because it does not account for specimen volume or kidney efficiency. Timed urine collections can be used to express the urine reference interval in U/24 hr. The urine reference interval can also be expressed as the urine AMY/CREA ratio.^{15,16}

Each laboratory should establish its own reference interval for amylase as performed on the Dimension® system.

$$f. \text{ UAMY/24 hr} = \frac{\text{urine AMY (U/L)}}{1000} \times \text{mL urine/24 hr}$$

$$g. \text{ AMY/CREA Clearance Ratio} = \frac{\text{urine AMY (U/L)}}{\text{serum AMY (U/L)}} \times \frac{\text{serum CREA (mg/dL)}}{\text{urine CREA (mg/dL)}} \times 100$$

Specific Performance Characteristics^h

Material	Mean U/L	Precision ⁱ	
		Within-run	Standard Deviation (% CV) Total
Serum Pool			
Low	50	0.4 (0.81)	0.68 (1.36)
High	408	1.1 (0.27)	3.25 (0.80)
Urine Pool			
Low	41	0.4 (1.1)	0.8 (2.0)
High	205	1.8 (0.9)	6.3 (3.1)
UrichemTRAK Control ^j			
Level 1	54	0.7 (1.3)	1.2 (2.2)
Level 2	178	1.3 (0.7)	4.0 (2.3)
Moni-Trol® Control			
Low	57	0.58 (1.01)	3.05 (5.32)
High	338	0.99 (0.29)	5.52 (1.63)
Moni-Trol® Control ^j			
Level 1	46	0.9 (1.9)	1.0 (2.2)
Level 2	209	1.1 (0.5)	1.7 (0.8)

h. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (see your Dimension® Operator's Guide).

i. Specimens at each level were analyzed in duplicate, twice a day for 20 days. The within-run and total standard deviations were calculated by EP5-T2: Evaluation of Precision Performance of Clinical Chemistry Devices, 2nd Edition, CLSI/NCCLS 1992; 12(4):146.

j. Using reduced sample size (10 µL).

Moni-Trol® is a registered trademark of Medical Analysis Systems Inc., Camarillo, CA 93012-8058, USA.

Comparative Method	Method Comparison			
	Regression Statistics		Correlation Coefficient	n
	Slope	Intercept U/L		
Original Dimension® ^k (Serum)	0.999	-2.6	0.995	254 ^l
Original Dimension® (Urine)	1.108	-2.9	0.999	93 ^m
Reduced vs. Standard Sample Size ⁿ (Serum)	1.02	-1.0	0.999	59 ^o

k. Model equation for regression statistics is: Results of Revised Dimension® AMY system = [slope x results of original Dimension® AMY system] + intercept.

Range of Samples (U/L)

l. 10 – 579.

m. 1 – 601.

n. Model equation for regression statistics is: Results of Dimension® AMY system using reduced sample size (10 µL) = [slope x results of Dimension® AMY system using standard sample size (14 µL)] + intercept.

o. Range of Samples 18 – 653 U/L.

Specificity

Hemolysis, Icterus, Lipemia (HIL) Interference

The AMY method (using the standard sample size of 14 µL) was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	AMY Activity U/L	Bias % ^p
Hemoglobin (hemolysate)	500 mg/dL [0.31 mmol/L] (monomer)	141	<10
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	141	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	140	<10
	3000 mg/dL [33.9 mmol/L]	140	q

The AMY method (using an alternate sample size of 10 µL) was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	AMY Activity U/L	Bias % ^p
Hemoglobin (hemolysate)	500 mg/dL [0.31 mmol/L] (monomer)	108	<10
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	110	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	106	<10

p. Analyte results should not be corrected based on this bias.

q. The interference testing at this level tripped a test report message; therefore the magnitude of the interference could not be determined.

Non-Interfering Substances

The following substances do not interfere with the AMY method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at the AMY activity of 134 U/L.

Substance	Test Concentration	SI Units
Acetaminophen	0.025 mg/dL	1.66 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	5 mg/dL	284 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	17.6 µmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	421 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	4 mg/dL	222 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Analytical Sensitivity: 2 U/L

The analytical sensitivity represents the lowest activity of amylase that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations of Reagent grade water.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Chavez RG. U.S. Patent 4,963,479.
2. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
3. Zakowski JJ, Bruns DE. Biochemistry of Human Alpha Amylase Isoenzymes, Crit Rev Clin Lab Sci, 1985; 21:283-322.
4. Clinical and Laboratory Standards Institute /NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
6. Henry JB. Clinical Diagnosis and Management by Laboratory Methods. 20th ed., 2001, Sanders WB Co, Philadelphia, PA, pp 464.
7. Wu AHB, Tietz NW. Clinical Guide to Laboratory Tests. 4th ed., 2006, Saunders WB Co, Philadelphia, PA., pp 102.
8. Kaplan LA. Clinical Chemistry Theory, Analysis, Correlation. Third ed., 1996, Mosby Inc., St. Louis.
9. Irie A, Hunaki M B, Anda K, Kawai K. Activation of α -Amylase in Urine. Clin Chem Acta 1974; 5:241-245.
10. Berger GMB, Liesching W. Interbatch Variability in Measurement of Urinary Amylase with a Commercial Kit, Clin Chem 1976; 22:556.
11. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
12. Dicht JJ. Reference Intervals for Serum Amylase and Urine Creatinine on the aca® discrete clinical analyzer, DuPont Company, Wilmington, DE (November, 1984).
13. Barger JD, Dicht JJ. Reference Intervals for Urine on AML, CALCM, MG, Na/K and UP, DuPont Company, Wilmington, DE (October 1983).
14. Gadsden RH, Phelps CA. A Normal Range Study of Amylase in Urine and Serum on the aca®, DuPont Company, Wilmington, DE (March 1978).
15. McCroskey R, Chang T, David H, Winn E. p-Nitrophenylglycosides as Substrates for Measurement of Amylase in Serum and Urine, Clin Chem 1982; 28:1787-1791.
16. Salt WB, Schenker S. Amylase-Its Clinical Significance: A Review of the Literature, Medicine, 1976; 55:269.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EPGS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare