Dimension® clinical chemistry system

Flex® reagent cartridge

LDI

See shaded sections: Updated information from 2016-02 version.

Issue Date 2019-04-08

Lactate Dehydrogenase

Intended Use: The lactate dehydrogenase (LDI) method is an *in vitro* diagnostic test for the quantitative measurement of total lactate dehydrogenase in human serum and plasma on the Dimension® clinical chemistry system.


Summary: The Dimension® LDI method has been standardized against the International Federation of Clinical Chemistry (IFCC) lactate dehydrogenase (LD) primary reference method procedure at 37 °C, adapted to the Dimension® clinical chemistry system.¹

Lactate dehydrogenase (LDI) is present in the cytoplasm of all cells in the body. The concentration of lactate dehydrogenase in tissues is several hundred-fold higher than in serum or plasma and even a small amount of tissue damage can lead to an elevation in lactate dehydrogenase activity.^{2,3} This makes lactate dehydrogenase especially useful in the diagnosis and monitoring of disease states where tissue turnover is accelerated such as the liver, cardiac muscle, skeletal muscle, kidneys and erythrocytes.^{2,3}

Lactate dehydrogenase is elevated in myocardial or pulmonary infarction, leukemias, hemolytic anemias, non-viral hepatitis, sickle cell disease, lymphoma, renal infarction, acute pancreatitis and any condition that results in the leaking of cytoplasm. It is moderately elevated in cirrhosis, obstructive jaundice, renal disease, skeletal muscle diseases, neoplastic diseases and congestive heart failure. Lactate dehydrogenase is markedly elevated in megaloblastic and pernicious anemia, metastatic carcinoma, viral hepatitis, shock, hypoxia and extreme hyperthermia.

Principles of Procedure: The LDI method uses as a substrate L-lactate buffered at a pH of 9.4. Lactate dehydrogenase oxidizes the substrate in the presence of NAD⁺ to yield pyruvate and NADH which absorbs at 340 nm. Lactate dehydrogenase activity concentration is measured as a rate reaction at 340/700 nm, proportional to the amount of lactate dehydrogenase in the sample.

The reversible reaction catalyzed by lactate dehydrogenase is as follows:

Reagents

Wells*	Form	Ingredient	Concentration*
1 – 4	Liquid	N-Methyl-D-glucamine	1091 mmol/L
Reagent 1		L(+)-lactate	168 mmol/L
		NaCl	513 mmol/L
5 – 6 ^b	Liquid	β-NAD ⁺	16.5 mmol/L
Reagent 2		β-Nicotinamide Adenine Dinucleotide Lithium Salt (NAD)	36.0 mmol/L
		Preservative	
		Stabilizer	

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

c. Wells 5 – 6 contain preservative and stabilizer.

Risk and Safety

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 4
6 days for wells 5 – 6

Specimen Collection and Handling: Recommended specimen types: serum, and sodium and/or lithium heparin plasma.

Phlebotomies to obtain samples for lactate dehydrogenase activity concentration should be carefully done to avoid hemolysis since the abundant lactate dehydrogenase in red blood cells can contaminate such samples.⁴

Hemolyzed samples should not be used with the LDI method.

Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁵

Follow the instructions provided with your specimen collection device for use and processing.⁶

For serum, complete clot formation should take place before centrifugation.

Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.⁷

For plasma, make sure separated specimen is free of platelet contamination. Platelets contain high activity concentrations of lactate dehydrogenase.⁸

Separated serum/plasma samples are stable for 3 days at 20 – 25 °C.^{8,9} Do not refrigerate or freeze.^{8,9}

Note: Activity of LD-4 and LD-5 are lost if the samples are stored below ambient temperatures.⁸

Procedure

LDI Flex® reagent cartridge, Cat. No. DF54

Materials Provided

Materials Required But Not Provided

Dimension® Enzyme I Calibrator, Cat. No. DC35

Quality control materials

Enzyme Diluent, Cat. No. 790035.901

Test Steps

Sampling^d, reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® clinical chemistry system. For details of this processing, refer to your Dimension® Operator's Guide.

d. The sample container must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

Test Conditions

Sample Volume	8 µL
Reagent 1 Volume	106 µL
Reagent 2 Volume	50 µL
Temperature	37.0 °C ± 0.1 °C
Reaction Time	7.5 minutes ^d
Wavelength	340/700 nm
Type of Measurement	Bichromatic rate

e. calculated as time from test initiation to final result

Calibration

Assay Range	6 – 1000 U/L [0.10 – 16.70 µkat/L] ^f
Calibration Material	Dimension® ENZ I CAL (Cat. No. DC35)
Calibration Scheme	Three levels in triplicate
Units	U/L [µkat/L] (U/L x 0.0167) = [µkat/L]
Typical Calibration Levels	Level 1: 0 U/L [0.00 µkat/L] ^g Level 2: 545 U/L [9.10 µkat/L] Level 3: 1075 U/L [17.95 µkat/L]
Calibration Frequency	Every 90 days for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ 0.000 C ₁ 1.000

f. Système International d'Unités [SI Units] are in brackets.

g. Level 1 calibrator for LDI is not included in the ENZ I CAL carton. Purified Water Diluent (Cat. No. 710615901) or Reagent grade water should be used as the Level 1 calibrator for the LDI method.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of a Quality Control (QC) material with known lactate dehydrogenase activity concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the activity concentration of lactate dehydrogenase in U/L [µkat/L] using the calculation scheme described in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 6 – 1000 U/L [0.10 – 16.70 µkat/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

• **Serum/plasma samples** with results in excess of 1000 U/L [16.70 µkat/L] are reported as "Above Assay Range" and should be repeated on dilution.

Manual Dilution: Dilute with Enzyme Diluent Cat. No. 790035.901 to obtain results within the analytical measurement range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): The recommended autodilute sample volume is 2 µL (dilution factor = 4) for serum/plasma. Refer to your Dimension® Operator's Guide.

• **Serum/plasma samples** with results less than 6 U/L [0.10 µkat/L] should be reported as "less than 6 U/L [0.10 µkat/L]".

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in lactate dehydrogenase results. Refer to your Dimension® Operator’s Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory’s procedure manual and not reported if indicated.

A system malfunction may exist if the following 5-test precision is observed:

LDI Activity Concentration	SD
100 U/L [1.67 µkat/L]	>4 U/L [> 0.07 µkat/L]
374 U/L [6.25 µkat/L]	>7 U/L [> 0.12 µkat/L]

Interfering Substances

The LDI method was evaluated for interference according to CLSI/NCCLS EP7-A2.¹⁰ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Hemoglobin (hemolysate) at 50 mg/dL [0.03 mmol/L] increases LDI results by 13.7% at a lactate dehydrogenase activity of 300 U/L [5.01 µkat/L] and 500 U/L [8.35 µkat/L].

Note: Lactate dehydrogenase released from red blood cells especially intravascularly may vary by disease state. Therefore bias from hemolysis in any individual sample may vary from patient to patient.¹¹

Dopamine at 65 mg/dL [4243 µmol/L] increases LDI results by 113% at a lactate dehydrogenase activity of 300 U/L [5.01 µkat/L].

Samples from patients with unexpectedly high levels of lactate dehydrogenase should be investigated for the presence of immunoglobulin-complexed enzyme as a cause for such results.¹²

Expected Values

Males:	85 – 227 U/L [1.42 – 3.79 µkat/L] ^h
Females:	81 – 234 U/L [1.35 – 3.91 µkat/L] ⁱ

Lactate dehydrogenase activity concentration was determined in 270 apparently healthy adult individuals. The reference interval was calculated nonparametrically and represents the central 95% of results.
h. The age range of males in years in the reference interval study was 19 to 61.
i. The age range of females in years in the reference interval study was 18 to 57.

Each laboratory should establish its own expected values for LDI as performed on the Dimension® clinical chemistry system.

Specific Performance Characteristics

The following data represent typical performance for the Dimension® clinical chemistry system.

Precision ^{13,j}			
Material	Mean	Standard Deviation (% CV)	
	U/L [µkat/L]	Repeatability	Within-Lab
Plasma Pool	157 [2.62]	1.4 [0.02] (1.0)	3.0 [0.05] (1.9)
Serum Pool	470 [7.85]	4.8 [0.08] (1.0)	21.1 [0.35] (4.5)
BioRad® Multiqual®			
Level 1	113 [1.89]	1.9 [0.03] (1.7)	3.1 [0.05] (2.8)
Level 3	376 [6.28]	3.1 [0.05] (0.8)	4.2 [0.07] (1.1)
MAS® chemTRAK® H			
Level 2	273 [4.56]	1.7 [0.03] (0.6)	4.1 [0.07] (1.5)

j. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Multiqual® is a registered trademark of Bio-Rad® Laboratories, Irvine, CA 92618, USA.
chemTRAK® is a registered trademark of Microgenics Corporation, Fremont, CA 94538, USA.

Method Comparison ¹⁴ Regression Statistics ^k				
Comparative Method	Slope	Intercept U/L [µkat/L]	Correlation Coefficient	n
IFCC LD at 37 °C reference method	0.985	6.5 [0.11]	0.996	52 ^l
Roche LDH IFCC on Hitachi® 917	0.985	+15.5 [0.26]	0.994	178 ^m

k. CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was ordinary least squares.
l. The range of lactate dehydrogenase values in the correlation study was 66 – 566 U/L [1.10 – 9.45 µkat/L].
m. The range of lactate dehydrogenase values in the correlation study was 33 – 957 U/L [0.55 – 15.98 µkat/L].
International Federation of Clinical Chemistry and Laboratory Medicine (IFCC) primary reference procedure for the measurement of catalytic concentration of LD at 37 °C.¹
Roche LDH IFCC/ Hitachi® is a registered trademark of Roche Diagnostics Division, GmbH, D68298 Mannheim.

Specificity

HIL interference

The LDI method was evaluated for interference according to CLSI/NCCLS EP7-A2.¹⁰ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration SI Units	Lactate Dehydrogenase U/L [µkat/L]	Bias* %
Hemoglobin (hemolysate)	50 mg/dL [0.03 mmol/L]	300 [5.01]	13.7
		500 [8.35]	13.7
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	300 [5.01]	<10
		500 [8.35]	<10
Bilirubin (conjugated)	80 mg/dL [1368 µmol/L]	300 [5.01]	<10
		500 [8.35]	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	300 [5.01]	<10
		500 [8.35]	<10
		300 [5.01]	n

n. Lipemia (Intralipid®) of 3000 mg/dL [33.9 mmol/L] generated an error flag on this method, so the magnitude of the interference is not available.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

***Analyte results should not be corrected based on this bias.**

Non-Interfering Substances

The following substances do not interfere with the LDI method when present at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at LDI method activity concentrations of 300 U/L [5.01 µkat/L].

Substance	Test Concentration	SI Units
Acetaminophen	20.0 mg/dL	1324 µmol/L
Allopurinol	2.5 mg/dL	184 µmol/L
Amikacin	8.0 mg/dL	137 µmol/L
Amiodarone	2.5 mg/dL	39 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	6.0 mg/dL	342 µmol/L
Atenolol	1.0 mg/dL	37.5 µmol/L
Caffeine	6.0 mg/dL	308 µmol/L
Captopril	5.0 mg/dL	230 µmol/L
Carbamazepine	3.0 mg/dL	127 µmol/L
Chloramphenicol	5.0 mg/dL	155 µmol/L
Chlordiazepoxide	1.0 mg/dL	33.3 µmol/L
Chlorpromazine	0.20 mg/dL	6.27 µmol/L
Cholesterol	503 mg/dL	13 mmol/L
Cimetidine	2.0 mg/dL	79.2 µmol/L
Cinnarizine	3.0 mg/dL	81.4 µmol/L
Creatinine	30 mg/dL	2.65 mmol/L
Cyclosporine A	4000 ng/mL	3.3 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.51 mg/dL	18 µmol/L
Digitoxin	350 ng/mL	458 nmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
Diltiazem	120 µg/mL	289 µmol/L
Disopyramide	4.0 mg/dL	118 µmol/L
Dopamine	32.5 mg/dL	2122 µmol/L
Erythromycin	6.0 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25.0 mg/dL	1770 µmol/L
Furosemide	6.0 mg/dL	181 µmol/L
Gentamicin	1.0 mg/dL	21 µmol/L
Heparin	3.0 U/mL	3000 U/L
Ibuprofen	50.0 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5000 mg/dL	50 g/L
Immunoglobulin M (IgM)	1041 mg/dL	10.41 g/L
Isosorbide Dinitrate	6.0 mg/dL	254 µmol/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
L-Thyroxine	60 µg/dL	0.77 µmol/L
Lisinopryl	16 µg/mL	36.2 µmol/L
Lovastatin	16 µg/mL	39.6 µmol/L
Methotrexate	225 mg/dL	4951 µmol/L
Methyldopa	2.5 mg/dL	105 µmol/L
Methylprednisolone	4.0 mg/dL	107 µmol/L
Mexiletine	24 mg/dL	1113 µmol/L
N-acetyl-procainamide	30 mg/dL	1082 µmol/L
Nicotine	0.10 mg/dL	6.2 µmol/L
Nifedipine	6.0 mg/dL	173 µmol/L
Nitroglycerine	0.16 µg/mL	0.577 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8.0 mg/dL	354 µmol/L
Phenobarbital	10.0 mg/dL	431 µmol/L
Phenytoin	5.0 mg/dL	198 µmol/L
Primidone	4.0 mg/dL	183 µmol/L
Procainamide	10 mg/dL	425 µmol/L
Propranolol	0.5 mg/dL	19 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Quinidine	20 µg/mL	61.6 µmol/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Simvastatin	32 µg/mL	76.4 µmol/L
Streptokinase	300 IU/mL	300 IU/mL
Theophylline	4.0 mg/dL	222 µmol/L
Tocainide	10 mg/dL	437 µmol/L
Triglyceride	1000 mg/dL	10 g/L
Urea	500 mg/dL	83 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L
Vancomycin	10 mg/dL	69 µmol/L
Verapamil	16 mg/dL	326 µmol/L
Warfarin	10 mg/dL	324 µmol/L

Limit of Detection¹⁵: 6 U/L [0.10 µkat/L]

The Limit of Detection (LoD) for the Dimension® LDI method was determined to be 6 U/L [0.10 µkat/L], determined consistent with CLSI/NCCLS guideline EP17-A, with false positives (α) and false negatives (β) less than 5%, and a Limit of Blank (LoB) of 3 U/L [0.05 µkat/L].

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía

1. Siekmann L, Ceriotti F, Ferard G, Kanno T, Schumann G, et. al. IFCC Primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37o C. Part 3. Reference procedure for the Measurement of Catalytic Concentration of Lactate Dehydrogenase. Clin Chem Lab Med. 2002; 40(6):643-648.
2. Henry JB. Clinical Diagnosis and Management by Laboratory Methods, 21st ED. WB Saunders Co.; 2007: 258-259.
3. Burtis CA, Ashwood ER, eds. Tietz Fundamentals of Clinical Chemistry, 5th ed. Philadelphia, PA: W.B. Saunders Company, 2001: pp. 362-366.
4. Thomas L, Haemolysis as influence and interference factor, eIJFCC vol 13 no 4: <http://www.ifcc.org/eijfcc/vol13no4/130401002.htm>.
5. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA, 2003.
6. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
7. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimen; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1- 56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
8. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed., W.B. Saunders Co. Philadelphia, PA 2006: pp 602.
9. Wu A, ed. Tietz Clinical Guide to Laboratory Tests, 4th ed. Philadelphia; W.B. Saunders Company, 2006: p. 648.
10. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
11. Henry JB. Clinical Diagnosis and Management by Laboratory Methods, 21st ED. WB Saunders Co.; 2007: 516.
12. Klonoff DC. Macroamylasemia and other immunoglobulin-complexed enzyme disorders. West J Med. 1980 November; 133(5): 392-407.
13. Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition*. CLSI/NCCLS Document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
14. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline- Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
15. Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

**Symbols Key
Symbol Schlüssel
Explication des Symboles
Interpretazione simboli
Clave de los Símbolos**

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante
	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenido suficiente per "n" saggi / Contenido suficiente para <n> ensayos
	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
	Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFIS

 Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare