Dimension® clinical chemistry system

Flex® reagent cartridge

CKI

See shaded sections: Updated information from 2016-12 version.

Issue Date 2019-04-08

Creatine Kinase


Intended Use: The CKI method is an *in vitro* diagnostic test for the quantitative measurement of creatine kinase in human serum and plasma on the Dimension® clinical chemistry system.

Summary: The Dimension® creatine kinase (CKI) method is standardized to the International Federation of Clinical Chemistry (IFCC) primary reference procedure at 37 °C, adapted to the Dimension® clinical chemistry system.¹

Measurements of creatine kinase are used in the diagnosis and treatment of myocardial infarction and muscle diseases. Creatine kinase may also be elevated following muscle injury or strenuous exercise.^{2,3}

Principles of Procedure: In a coupled enzyme reaction, the creatine kinase in patient samples catalyzes the transphosphorylation of phosphate from creatine phosphate to adenosine-diphosphate (ADP) producing adenosine-triphosphate (ATP). Hexokinase (HK) phosphorylates glucose from the ATP to phosphorylate glucose. The resulting glucose-6-phosphate is oxidized by glucose-6-phosphate dehydrogenase (G-6-PDH) with the simultaneous reduction of nicotinamide adenine dinucleotide phosphate (NADP).

The rate of formation of NADPH is directly proportional to the CK activity in the sample and is measured bichromatically at 340 and 540 nm.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1 – 4 ^c	Liquid	Hexokinase	13700 U/L	Yeast
		Glucose -6-phosphate dehydrogenase	14900 U/L	Bacterial
		ADP	7.05 mM	
		AMP	16.15 mM	
		EDTA	6.47 mM	
		Mg acetate	32.3 mM	
		Diadenosine Pentaphosphate (AP5A)	0.05 mM	
		NADP	6.28 mM	
		N-acetylcysteine	64.7 mM	
		Imidazole buffer	123 mM	
5, 6 ^c	Liquid	Creatine phosphate	184 mM	
		Glucose	120 mM	
		EDTA	2.46 mM	
		3-(cyclohexylamino)-2-hydroxyl-1-propane sulfonic acid (CAPSO) buffer	20 mM	

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

c. Wells 1 – 6 contain preservative and stabilizer.

Risk and Safety

H360d
P201, P202, P280, P281, P308 + P313, P501

Danger!

May damage unborn child.

Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Wear protective gloves/protective clothing/eye protection/face protection. Use personal protective equipment as required. If exposed or concerned: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: Imidazole

Safety data sheets (MSDS/SDS) available on [siemens.com/healthcare](https://www.siemens.com/healthcare)

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 5 days for wells 1 – 4, 10 days for wells 5 – 6.

Specimen Collection and Handling: Recommended specimen types: serum, sodium and/or lithium heparin plasma, and EDTA.

Grossly hemolyzed samples should not be used with the CKI method.³

Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁴

Follow the instructions provided with your specimen collection device for use and processing.⁵

For serum, complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.⁶ Separated serum/plasma samples should be stored at 2 – 8 °C and analyzed within 7 days. For longer storage, samples should be frozen at or below -20 °C.⁷

Procedure

CKI Flex® reagent cartridge, Cat. No. DF38

Materials Provided

Materials Required But Not Provided

Dimension® CKI/MBI Calibrator, Cat. No. DC32
Quality Control Materials

Test Steps

Sampling^d, reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® System. For details of this processing, refer to your Dimension® Operator's Guide.

d. The sample container must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

Test Conditions

Sample Volume	14 µL
Reagent 1 Volume	112 µL
Reagent 2 Volume	55 µL
Temperature	37.0 °C ± 0.1 °C
Reaction Time	8.7 minutes ^e
Wavelength	340 and 540 nm
Type of Measurement	Bichromatic rate

Calibration

Assay Range	7 – 1000 U/L [0.12 – 16.67 µkat/L] ^f
Calibration Material	Dimension® CKI/MBI Calibrator, Cat. No. DC32
Calibration Scheme	3 levels in triplicate
Units	U/L [µkat/L] (U/L ÷ 60) = [µkat/L]
Typical Calibration Levels	Level 1: 0 U/L [0.00 µkat/L] ^g Level 2: 525 U/L [8.75 µkat/L] Level 3: 1100 U/L [18.33 µkat/L]
Calibration Frequency	Every 90 days for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ 1.708 C ₁ 8.044

e. Calculated as time from test initiation to final result.

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g. Level 1 calibrator for CKI/MBI is not included in the CKI/MBI carton. Purified Water Diluent (Cat No. 710615901) or Reagent grade water should be used as the level 1 calibrator for the CKI method.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency.

At least once each day of use, analyze two levels of a Quality Control (QC) material with known creatine kinase activity. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the enzymatic activity of creatine kinase in U/L using the calculation scheme described in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 7 – 1000 U/L [0.12 – 16.67 µkat/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

• **Serum/plasma samples** with results in excess of 1000 U/L [16.67 µkat/L] are reported as "Above Assay Range" and should be repeated on dilution.

Manual Dilution: Dilute with Reagent grade water to obtain results within the analytical measurement range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): The autodilute sample volume is 2 µL (dilution factor = 7) for serum/plasma. Refer to your Dimension® Operator's Guide.

• **Serum/plasma** with results less than 7 U/L [0.12 µkat/L] should be reported as "less than 7 U/L [0.12 µkat/L]"

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in creatine kinase results. Refer to your Dimension® Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual and not reported.

A system malfunction may exist if the following 5-test precision is observed:

CKI Concentration	SD
150 U/L [2.50 µkat/L]	> 7 U/L [0.12 µkat/L]
800 U/L [13.33 µkat/L]	> 16 U/L [0.27 µkat/L]

Interfering Substances

The CKI was evaluated for interference according to CLSI/NCCLS EP7-A2.⁸ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Expected Values⁹

Serum/Plasma

Males: 39 – 308 U/L [0.65 – 5.11 µkat/L]
Females: 26 – 192 U/L [0.43 – 3.19 µkat/L]

Each laboratory should establish its own expected values for CKI as performed on the Dimension® clinical chemistry system.

Specific Performance Characteristics

The following data represent typical performance for the Dimension® clinical chemistry system.

Precision^{10,h}

Material	Mean	Standard Deviation (% CV)	
	U/L [µkat/L]	Repeatability	Within-Lab
Serum			
Serum Pool 1	46 [0.768]	0.7 [0.012] (1.4)	1.3 [0.022] (2.8)
Serum Pool 2	166 [2.772]	1.4 [0.023] (0.8)	3.5 [0.058] (2.1)
Multiqual®			
Level 1	108 [1.804]	1.1 [0.018] (1.0)	3.5 [0.058] (3.3)
Level 3	788 [13.160]	3.3 [0.055] (0.4)	13.4 [0.224] (1.7)

h. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Multiqua[®] is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

Method Comparison¹¹
Regression Statisticsⁱ

Comparative Method	Slope	Intercept U/L [µkat/L]	Correlation Coefficient	n
Serum/Plasma				
IFCC reference method	1.06	0.35 [0.01]	0.999	54 ^j
Plasma				
Roche CK Hitachi [®] 917	1.04	6.07 [0.10]	1.000	90 ^k

i. CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was orthogonal regression fit.

j. The range of serum and plasma creatine kinase values in the correlation study was 10 to 985 U/L.

k. The range of plasma creatine kinase values in the correlation study was 38 to 953 U/L.

International Federation of Clinical Chemistry and Laboratory Medicine (IFCC) primary reference procedure for the measurement of catalytic concentration of CK.¹

Roche/Hitachi CK is a product of F. Hoffmann-La Roche Ltd, Diagnostics Division, Grenzachstrasse 124, CH-4070 Basel, Switzerland.

Specificity

HIL Interference

The CKI was evaluated for interference according to CLSI/NCCLS EP7-A2.⁸ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Hemoglobin (hemolysate) at 300 mg/dL [0.186 mmol/L] increases CKI results by 22% at a creatine kinase activity of 200 U/L [3.34 µkat/L].

Substance Tested	Substance Concentration SI Units	Creatine Kinase U/L	Bias* %
Hemoglobin (hemolysate)	Hemoglobin	200 U/L [3.33 µkat/L]	< 10
	100 mg/dL [0.06 mmol/L]	611 U/L [10.18 µkat/L]	
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	200 U/L [3.33 µkat/L]	< 10
		611 U/L [10.18 µkat/L]	
Bilirubin (conjugated)	80 mg/dL [1368 µmol/L]	200 U/L [3.33 µkat/L]	< 10
		611 U/L [10.18 µkat/L]	
Lipemia (Intralipid [®])	3000 mg/dL [33.9 mmol/L]	200 U/L [3.33 µkat/L]	< 10
		611 U/L [10.18 µkat/L]	

Intralipid[®] is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

*Analyte results should not be corrected based on this bias.

Non-Interfering Substances – serum/plasma

The following substances do not interfere with the CKI method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at CKI method concentrations of 200 U/L [3.33 µkat/L].

Substance	Test Concentration	SI Units
Acetaminophen	20.0 mg/dL	1324 µmol/L
Allopurinol	2.5 mg/dL	184 µmol/L
Amikacin	8.0 mg/dL	137 µmol/L
Amiodarone	2.5 mg/dL	3.9 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	6.0 mg/dL	342 µmol/L
Atenolol	1.0 mg/dL	37.6 µmol/L
Caffeine	6.0 mg/dL	308 µmol/L
Captopril	5.0 mg/dL	230 µmol/L
Carbamazepine	3.0 mg/dL	127 µmol/L
Chloramphenicol	5.0 mg/dL	155 mmol/L
Chlordiazepoxide	1.0 mg/dL	33.3 µmol/L
Chlorpromazine	0.20 mg/dL	6.27 µmol/L
Cholesterol	544 mg/dL	13 mmol/L
Cimetidine	2.0 mg/dL	79.2 µmol/L
Cinnarizine	3.0 mg/dL	81.4 µmol/L
Creatine	30.0 mg/dL	2.7 mmol/L
Cyclosporin A	4000 ng/mL	3.3 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.51 mg/dL	18 µmol/L
Digitoxin	350 ng/mL	458 nmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
Diltiazem	120 µg/mL	266.1 µmol/L
Disopyramide	4.0 mg/dL	1017 µmol/L
Dopamine	65 mg/dL	118 µmol/L
Erythromycin	6.0 mg/dL	81.6 mmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25.0 mg/dL	1770 µmol/L
Furosemide	6.0 mg/dL	181 µmol/L
Gentamicin	1.0 mg/dL	21 µmol/L
Heparin	3.0 U/mL	3000 U/L
Ibuprofen	50.0 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5000 mg/dL	50 g/L
Immunoglobulin M (IgM)	5000 mg/dL	50 g/L
Isosorbide dinitrate	6.0 mg/dL	254 µmol/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lisinopril	16 µg/mL	36.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
Lovastatin	16 µg/mL	39.6 µmol/L
L-thyroxine	60 µg/dL	0.77 µmol/L
Methotrexate	100 mg/dL	2201 µmol/L
Methylidopa	2.5 mg/dL	118 µmol/L
Methylprednisolone	4.0 mg/dL	107 µmol/L
Mexiletine	24 mg/dL	1339 µmol/L
N-acetyl-procainamide	30 mg/dL	1082 µmol/L
Nicotine	0.10 mg/dL	6.2 µmol/L
Nifedipine	6.0 mg/dL	173 µmol/L
Nitroglycerin	0.16 µg/mL	0.58 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8.0 mg/dL	354 µmol/L
Phenobarbital	10.0 mg/dL	431 µmol/L
Phenytoin	5.0 mg/dL	198 µmol/L
Primidone	4.0 mg/dL	183 µmol/L
Procainamide	10 mg/dL	443 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Propranolol	0.5 mg/dL	19 µmol/L
Protein: Albumin	6000 mg/dL	60 g/L
Protein: Total	12000 mg/dL	120 g/L
Quinidine	20 µg/mL	61.6 µmol/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Simvastatin	32 µg/mL	76.5 µmol/L
Streptokinase	300 IU/mL	300,000 IU/L
Theophylline	4.0 mg/dL	222 mmol/L
Tocainide	10 mg/dL	437 µmol/L
Urea	500 mg/dL	83 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L
Verapamil	16 mg/dL	0.33 µmol/L
Warfarin	10 mg/dL	324 µmol/L

Limit of Detection¹²: 7 U/L [0.12 µkat/L]

The limit of detection represents the lowest concentration of creatine kinase that can be detected with at least 95% probability. CLSI/NCCLS EP17-A, Protocols for Determination of Limits of Detection and Limits of Quantitation, was followed using 288 replicate determinations at two concentration levels to determine Limit of Blank and Limit of Detection.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension[®] and Flex[®] are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Dimension® clinical chemistry system

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Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía

1. Siekmann L, Ceriotti F, Ferard G, Kanno T, Schumann G, et. al. IFCC Primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 2. Reference procedure for the measurement of catalytic concentrations of creatine kinase [ATP: creatine N-phosphotransferase (CK), EC 2.7.3.2] Clin Chem Lab Med. 2002; 40(6):635-642.
2. McPherson RA, Pincus MR, eds. Henry's clinical diagnosis and management by laboratory methods 21st ed. St. Louis: Saunders Elsevier, 2007: 221, 260.
3. Burtis CA, Ashwood ER, eds. Tietz fundamentals of clinical chemistry, 5th ed. Philadelphia: WB Saunders Co., 2001:358-62.
4. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
6. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimen; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1- 56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
7. Wu AHB, ed. Tietz Clinical Guide to Laboratory Tests 4th Ed. St. Louis MO. Elsevier Saunders. 2006: pp.306 – 7.
8. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
9. Klein G, Berger A, Bertholf R et al. Abstract: Multicenter Evaluation of Liquid Reagents for CK, CK-MB and LDH with Determination of Reference Intervals on Hitachi Systems. Clin Chem 2001; 47:Suppl. A30.
10. Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition*. CLSI/NCCLS Document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
11. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
12. Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "<n"> tests / Contenuto sufficiente per "<n"> saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic *in vitro* / Dispositivo medico-diagnostico *in vitro* / Producto sanitario para diagnóstico *in vitro*

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFIS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare