Dimension® clinical chemistry system

Flex® reagent cartridge

ACTM

See shaded sections: Updated information from 2017-12 version.


Issue Date 2019-04-01

Acetaminophen

Intended Use: The ACTM method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended to measure acetaminophen, an analgesic and antipyretic, in human serum and plasma. ACTM test results are used in the diagnosis and treatment of acetaminophen overdose.

Summary: Acetaminophen is an analgesic found in many "over-the-counter" pain remedies. It is rapidly and completely absorbed from the gastrointestinal tract. After oral administration, peak plasma concentrations are reached in less than an hour. Approximately 90% of a therapeutic dose is eliminated by conjugation with glucuronic acid (and to a slight extent, sulfuric acid) in the liver. Another 3 – 5% is catabolized by the P-450 mixed function oxidase enzyme system to the acid and cysteine conjugates. All of these metabolites are excreted in the urine. Only a slight amount of the drug is excreted unchanged. Intermediate metabolites of uncertain structure formed during the biotransformation in the liver are believed to be responsible for the hepatotoxicity. After a therapeutic dose of acetaminophen, the biologic half-life in normal adults is 2 – 3 hours. Metabolism is more rapid in children (except newborns). Because the hepatic conjugation is the rate-limiting step in the catabolic pathway, the half-life is prolonged in patients with liver disease, alcoholics, or in the presence of other drugs which compete for the hepatic conjugation mechanism. Acetaminophen does not have anti-inflammatory activity and it does not affect blood clotting (hemostasis). It is preferred over aspirin when the hemostatic side effects of aspirin must be avoided. Severe liver damage in adults is generally associated with ingestion of 15 grams or more. Since the drug is catabolized in the liver, hepatotoxicity will result in elevated plasma drug levels and prolonged half-life. The availability of a rapid accurate plasma acetaminophen assay is of extreme importance in cases of suspected intoxication because effective antidotes are available. Therapy with N-acetylcysteine (NAC) must be started within eight hours after ingestion to prevent hepatic injury as signified by elevations in AST and ALT.^{1,2}

Principles of Procedure: The methodology for ACTM is based on the enzymatic hydrolysis of acetaminophen producing acetate and p-aminophenol. The p-aminophenol is determined colorimetrically by reaction with o-cresol and ammoniacal copper sulfate.³ The amount of aminophenol produced is proportional to the acetaminophen concentration and is measured using a bichromatic endpoint technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1, 2	Tablet	Enzyme	0.04 U/mL
3, 4	Liquid	Copper Sulfate	0.7 mM
5, 6	Liquid	o-Cresol	17.2 mM

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per test at manufacture.

Risk and Safety:

H319, H315, H411
P280, P264, P273, P305 + P351 + P338, P337 + P313, P391, P501

Warning!

Causes serious eye irritation. Causes skin irritation. Toxic to aquatic life with long lasting effects.

Wear protective gloves/protective clothing/eye protection/face protection. Wash hands thoroughly after handling. Avoid release to the environment.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. Collect spillage. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: Ammonium chloride, Copper sulphate, o-Cresol, Sodium carbonate

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: Hydrating, diluting and mixing are automatically performed by the instrument.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6

Specimen Collection and Handling: Recommended specimen types: serum or plasma (lithium heparin). Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁴

Follow the instructions provided with your specimen collection device for use and processing.⁵

Heparin, citrate, potassium oxalate, sodium fluoride and EDTA anticoagulants do not interfere with the ACTM method at concentrations normally found in blood collection tubes.

Complete clot formation should take place before centrifugation.⁶ Specimens should be free of particulate matter.

Each laboratory must determine its appropriate sample timing/collection protocol.

Separated specimens are stable for 8 hours at 20 – 25 °C, 2 weeks at 2 – 8 °C and 45 days at -20 °C or colder.⁷

Procedure

Materials Provided

ACTM Flex® reagent cartridge, Cat. No. DF88

Materials Required But Not Provided

Drug Calibrator II, Cat. No. DC49D
Quality Control Materials

Test Steps

Sampling,^c reagent delivery, mixing, processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

c. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus dead volume; precise container filling is not required.

Test Conditions

	Blank	Test
Sample Size	4 µL	4 µL
Enzyme Volume	–	80 µL
Copper Sulfate Volume	80 µL	80 µL
O-Cresol Volume	80 µL	80 µL
Temperature	37 °C	
Wavelength	600 and 700 nm	
Type of Measurement	Bichromatic endpoint	

Calibration

Assay Range	0.0 – 300.0 µg/mL [0.0 – 1985.4 µmol/L] ^d
Calibration Material	Drug Calibrator II, Cat. No. DC49D
Calibration Scheme	3 levels, n = 3
Units	µg/mL [µmol/L] (µg/mL x 6.62) = [µmol/L]
Typical Calibration Levels	0.0, 75.0, 300.0 µg/mL [0.0, 496.4, 1985.4 µmol/L]
Calibration Frequency	Every 90 days for any one lot
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, as indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations
Assigned Coefficients	C ₀ 0.000 C ₁ 1.990
d. Système International d'Unités [SI Units] are in brackets.	

Quality Control

At least once each day of use, analyze two levels of a quality control material with known ACTM concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of acetaminophen in µg/mL [µmol/L] using the calculation scheme outlined in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.0 – 300.0 µg/mL [0.0 – 1985.4 µmol/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 300.0 µg/mL [1985.4 µmol/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilution with Level 1 (0 µg/mL) Drug Calibrator II, or drug-free serum to obtain results within assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution: Refer to your Dimension® Operator's Guide.

For samples with acetaminophen concentrations extremely beyond the method's assay range, test results will be accompanied by a message flag, however, the autodilute feature may not be initiated.

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
16.2 µg/mL [107.2 µmol/L]	> 1.93 [12.78 µmol/L]
75.0 µg/mL [496.4 µmol/L]	> 2.35 [15.56 µmol/L]

Interfering Substances

The ACTM method was evaluated for interference according to CLSI/NCCLS EP7-A2.⁸ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Bilirubin (unconjugated) of 20 mg/dL [342 µmol/L] decreases an ACTM result of 152.5 µg/mL [1009.1 µmol/L] by 21%.

Lipemia (Intralipid®) of 3000 mg/dL [33.9 mmol/L] tripped an error flag on this method, so the magnitude of the interference is not available.

Albumin of 6 g/dL [60 g/L] decreases an ACTM result of 30 µg/mL [19.9 µmol/L] by 13%.

Immunoglobulin G (IgG) of 5 g/dL [50 g/L] decreases an ACTM result of 30 µg/mL [19.9 µmol/L] by 26%.

Total protein of 12 g/dL [120 g/L] decreases an ACTM result of 30 µg/mL [19.9 µmol/L] by 31%.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Therapeutic Range: Therapeutic acetaminophen concentrations vary significantly depending on the individual patient. A range of 10.0 – 30.0 µg/mL [66.2 – 198.5 µmol/L] may be an effective sample concentration for many patients;⁹ however, some individuals are best treated at concentrations outside this range. The physician must determine the appropriate therapeutic range for each patient.

Acetaminophen concentrations > 150.0 µg/mL [992.7 µmol/L] at four hours after ingestion and 50.0 µg/mL [330.9 µmol/L] at twelve hours after ingestion are often associated with toxic reactions.¹⁰

Specific Performance Characteristics^a

Precision^d

Material	Mean µg/mL [µmol/L]	Standard Deviation (% CV)		
		Within-run	Total	
Liquichek™ QC	42.9 [284]	0.91 [6.02] (2.12)	1.09 [7.22] (2.55)	
Serum Pool	146.5 [970]	1.07 [7.08] (0.73)	2.13 [14.1] (1.45)	
Plasma Pool	13.4 [88.7]	0.39 [2.58] (2.89)	0.59 [3.91] (4.44)	

e. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).

f. Specimens at each level were analyzed in duplicate, two runs/day for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

Liquichek™ is a trademark of Bio-Rad Laboratories, Anaheim, CA 92806.

Method Comparison

Regression Statistics^a

ACTM Flex® cartridge on Dimension® System Compared to	Slope	Intercept µg/mL [µmol/L]	Correlation Coefficient	n ^b
aca® ACTMN test pack on aca® clinical analyzer	1.04	-3.27 [21.67]	0.998	154

g. Model equation for regression statistics is: [results of Dimension® system] = slope x [comparative method results] + intercept.

h. Range of samples (µg/mL): 10.2 – 295.0.

Specificity

HIL Interference

The ACTM method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	ACTM Concentration µg/mL [µmol/L]	Bias ⁱ %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L] (monomer)	150.2 [996.6]	<10
Bilirubin (unconjugated)	5 mg/dL [86 µmol/L]	152.5 [1009.1]	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	149.1 [986.6]	<10

i. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances do not interfere with the ACTM method when present in serum and plasma at the concentrations indicated. Systemic inaccuracies (bias) due to these substances are less than 10% at acetaminophen concentration of approximately 50 µg/mL [330.9 µmol/L].

Substance	Test Concentration	SI Units
Acetylsalicylic Acid	50 mg/dL	1.67 mmol/L
Amikacin	15 mg/dL	256 µmol/L
Amitriptyline	20 mg/dL	721 µmol/L
Amobarbital	10 mg/dL	442 µmol/L
Ampicillin	5 mg/dL	143 µmol/L
Ascorbic Acid	3 mg/dL	170 µmol/L
Caffeine	10 mg/dL	515 µmol/L
Carbamazepine	12 mg/dL	508 µmol/L
Chloramphenicol	25 mg/dL	774 µmol/L
Chlordiazepoxide	2 mg/dL	67 µmol/L
Chlorpromazine	5 mg/dL	157 µmol/L
Cimetidine	10 mg/dL	396 µmol/L
Codeine	10 mg/dL	334 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 75	2500 mg/dL	333 µmol/L
Diazepam	2 mg/dL	70 µmol/L
Digoxin	20 ng/mL	26 nmol/L
Erythromycin	20 mg/dL	273 µmol/L
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 µmol/L
Furosemide	2 mg/dL	61 µmol/L
Gentamicin	12 mg/dL	221 µmol/L
Heparin (Porcine)	8000 U/L	8000 U/L
Ibuprofen	40 mg/dL	1939 µmol/L
Imipramine	20 mg/dL	713 µmol/L
Lidocaine	6 mg/dL	256 µmol/L
Methadone HCL	20 mg/dL	578 µmol/L
Methionine	20 mg/dL	1340 µmol/L
Nicotine	2 mg/dL	123 µmol/L
Nortriptyline	1000 ng/mL	3797 nmol/L
Penicillin V	80 mg/dL	2247 µmol/L
Pentobarbital	10 mg/dL	442 µmol/L
Phenobarbital	15 mg/dL	646 µmol/L
Phenytoin	10 mg/dL	395 µmol/L
Phenylephrine	10 mg/dL	491 µmol/L
Primidone	10 mg/dL	458 µmol/L
Procainamide	20 mg/dL	736 µmol/L
Propoxyphene	0.4 mg/dL	12 µmol/L
Protein - Albumin	4 g/dL	40 g/L
Quinidine	20 mg/dL	616 µmol/L
Rheumatoid Factor	375 IU/L	375 IU/L
Salicylic Acid	50 mg/dL	3.62 mmol/L
Secobarbital	5 mg/dL	210 µmol/L
Theophylline	25 mg/dL	1389 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3472 µmol/L

Cross-reactivity

The following method specific compounds (metabolites/cross-reactants) have no measurable effect on the ACTM method at the concentrations indicated:

Substance	Test Concentration	SI Units
Acetaminophen Cysteine	20 mg/dL	740 µmol/L
Acetaminophen Glucuronide	20 mg/dL	573 µmol/L
Acetaminophen Mercapturate	20 mg/dL	640 µmol/L
Acetaminophen Sulfate	20 mg/dL	743 µmol/L
Acetanilide	10 mg/dL	740 µmol/L
Acetophenetidin	10 mg/dL	558 µmol/L
N-Acetyl/cysteine	20 mg/dL	1225 µmol/L
N-Acetylprocainamide	20 mg/dL	721 µmol/L
Salicylamide	10 mg/dL	729 µmol/L

Analytical Sensitivity: 2.0 µg/mL [13.2 µmol/L]

The analytical sensitivity represents the lowest concentration of acetaminophen that can be distinguished from zero. The sensitivity is defined as the mean value (n = 20) plus two standard deviations of the level 1 (0.0 µg/mL) [0.0 µmol/L] Drug Calibrator II.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, aca® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Burtis CA, Ashwood ER, Bruns DE. Tietz Textbook of Clinical Chemistry and Molecular Biology, Fourth Edition, Elsevier Saunders, St. Louis, MO 2006; pp. 1253 (clinical significance), pp. 2315 (reference values).
2. Rumack BH, Matthew H. Acetaminophen poisoning and toxicity. Pediatrics, 1975;55, 871-876.
3. Price CP, Hammond PM, Scawen MD. Evaluation of an Enzymatic Procedure for the Measurement of Acetaminophen. Clin Chem 29/2; pp 358-61 (1983).
4. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
6. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
7. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, 3rd ed. Washington, DC: AACC Press, 2007; p 46.
8. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
9. Tietz NW. Clinical Guide to Laboratory Tests. Philadelphia: W.B. Saunders Co., 1990; pp684-685.
10. Rumack BH, Matthew H. Acetaminophen poisoning and toxicity, Pediatrics, 1975; p 55; pp871-876.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFIS

Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare