Dimension® clinical chemistry system

Flex® reagent cartridge

CSA

See shaded sections: Updated information from 2015-02 version.

Issue Date 2019-04-08

Cyclosporine A


Intended Use: The CSA method is an *in vitro* diagnostic test intended to quantitatively measure cyclosporine A (CSA) in human whole blood for the Dimension® clinical chemistry system. Measurements of CSA are used as an aid in the management of heart, liver and kidney transplant patients.

Summary: Cyclosporine is a cyclic undecapeptide of fungal origin and a potent immunosuppressive agent. Inadequate cyclosporine levels may result in rejection of the transplanted organ. Toxic levels of cyclosporine are associated with many serious side effects, including nephrotoxicity and hepatotoxicity. Monitoring parent drug cyclosporine concentrations in whole blood and interpreting these values in conjunction with other laboratory data and clinical considerations is the most effective means of ensuring adequate immunosuppressant therapy for recipients of solid-organ transplants.^{1,2,3}

Cyclosporine is extensively metabolized by the liver. Circulating cyclosporine levels may be influenced by drugs that affect hepatic microsomal enzymes, particularly the cytochrome P-450 system. Substances known to inhibit these enzymes, such as cimetidine, ketoconazole and methylprednisolone, will decrease hepatic metabolism and increase cyclosporine levels. Substances that are inducers of cytochrome P-450 activity, including rifampin, phenytoin, phenobarbital, and carbamazepine, will increase hepatic metabolism and decrease cyclosporine levels.⁴

Principles of Procedure: The automated Dimension® CSA method uses an immunoassay technique in which free and CSA-bound antibody-enzyme species are separated using magnetic particles. The assay is performed using a method specific Flex® reagent cartridge. The Flex® cartridge contains a pretreatment reagent, β -galactosidase-CSA antibody conjugate, CSA immobilized on chromium dioxide particles, chlorophenol red β -D-galactopyranoside (CPRG) substrate, and diluent to hydrate the tablets.

To perform the CSA assay, a sample cup containing the whole blood sample to be analyzed and a CSA Flex® reagent cartridge are placed appropriately on the Dimension® system. The Dimension® system mixes and lyses the whole blood sample. The lysed sample is then mixed with the antibody conjugate reagent. The CSA present in the sample is bound by the CSA antibody conjugate reagent. Magnetic particles coated with cyclosporine A are added to bind free (unbound) antibody-enzyme conjugate. The reaction mixture is then separated magnetically. Following separation, the supernatant containing the CSA-antibody-enzyme complex is transferred to another cuvette and mixed with the substrate. β -galactosidase catalyzes the hydrolysis of CPRG (chlorophenol red β -D-galactopyranoside) to produce CPR (chlorophenol red) that absorbs light maximally at 577 nm. The change in absorbance at 577 nm due to the formation of CPR is directly proportional to the amount of CSA in the patient's sample and is measured using a bichromatic (577, 700 nm) rate technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1,2	Liquid	Ab- β -galactosidase	c	Mouse, monoclonal
3,4	Tablet ^d	CSA-CrO ₂	3.4 mg/tab	
5,6	Tablet ^d	CPRG	9.3 mg/tab	
7	Liquid	Substrate Diluent		
8	Liquid	Pretreatment reagent		

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per tablet at manufacture.

c. Antibody concentration and conjugate activity varies from lot to lot.

d. Tablets contain protease-free BSA, bovine gamma globulin, 30% BSA, excipients, buffers and stabilizers.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501

Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 5-chloro-2-methyl-3(2H)-isothiazolone mixture with 2-methyl-3(2H)-isothiazolone; Chromium oxide.

Safety data sheets (MSDS/SDS) available on www.siemens.com/diagnostics

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: Hydrating, diluting, and mixing are automatically performed by the instrument.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6
10 days for wells 7 – 8

Specimen Collection and Handling: Recommended specimen types: Use only whole blood in EDTA collected by normal procedures.⁵

Blood can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁵

Follow the instructions provided with your specimen collection device for use and processing.⁶

A trough sample is recommended for measurement of cyclosporine. The time when blood is drawn relative to the last dose of cyclosporine should be noted to ensure proper interpretation of results.

Mix freshly collected, refrigerated and thawed samples thoroughly prior to use to allow equilibration and uniform distribution of cyclosporine in whole blood.⁷

Specimens that contain particulate matter and clots should not be used.

Whole blood specimens are stable for 7 days at ambient temperatures.⁸ Whole blood specimens may be refrigerated at 2 – 8 °C for up to 7 days. For longer periods of storage, samples should be frozen at -20 °C.⁹

Avoid repeated freezing and thawing.

Procedure

Materials Provided

CSA Flex® reagent cartridge, Cat. No. DF89A

Materials Required But Not Provided

Sample Probe Cleaner, Cat. No. RD703

Reagent Probe Cleaner, Cat. No. RD702

CSA Calibrator, Cat. No. DC89

Reaction Vessels (required for HM only), Cat. No. RXV1A

Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, separation, processing and printing of results are automatically performed by the Dimension® System with the heterogeneous immunoassay module.

- Mix all calibrators, controls and samples gently but thoroughly just before use. Mix by hand or on an inverter or rocker. **DO NOT VORTEX.**

- Transfer **200 μ L** to a sample cup for processing on the Dimension® system. The specimen must be sampled within 30 minutes of being placed on the sample wheel.

- Do not process primary tubes directly on the Dimension® clinical chemistry system for this method.**

For details of this processing, refer to your Dimension® Operator's Guide.

Test Conditions

Reaction Vessel

Pretreatment Reagent	70 μ L
Sample Size	5 μ L
Antibody- β -galactosidase	50 μ L
CSA - CrO ₂	50 μ L
Incubation Temperature	42 °C*

*Dimension® EXL™ with LOCi® Module: 37°C

Reaction Cuvette

Transfer Volume	54 μ L
CPRG Reagent Volume	150 μ L
Diluent Volume	191 μ L
Temperature	37.0 °C
Wavelength	577 and 700 nm
Type of Measurement	Bichromatic rate

Calibration

Assay Range	25 – 500 ng/mL [20.8 – 416.3 nmol/L] ^e										
Calibration Material	CSA Calibrator, Cat. No. DC89										
Calibration Scheme	5 levels, Four @ level 1, three @ level 2, two @ level 3, and three @ levels 4 and 5.										
Units	ng/mL [nmol/L] (ng/mL x 0.832) = [nmol/L]										
Typical Calibration Levels	0, 100, 200, 350, 500 ng/mL [0.0, 83.2, 166.4, 291.2, 416.0 nmol/L]										
Calibration Frequency	Every 30 days for any one lot										
A new calibration is required	- For each new lot of Flex® reagent cartridges - After major maintenance or service, if indicated by quality control results - As indicated in laboratory quality control procedures - When required by government regulations										
Assigned Coefficients	<table> <tr> <td>C₀</td><td>41.9</td></tr> <tr> <td>C₁</td><td>746</td></tr> <tr> <td>C₂</td><td>-2.25</td></tr> <tr> <td>C₃</td><td>562</td></tr> <tr> <td>C₄</td><td>0.5</td></tr> </table>	C ₀	41.9	C ₁	746	C ₂	-2.25	C ₃	562	C ₄	0.5
C ₀	41.9										
C ₁	746										
C ₂	-2.25										
C ₃	562										
C ₄	0.5										

e. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known cyclosporine concentration. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of CSA in ng/mL [nmol/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 25 – 500 ng/mL [20.8 – 416.3 nmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 500 ng/mL [416.3 nmol/L] are reported as "Above Assay Range" and should be run on CSAE method (DF108), which provides quantitative values for CSA from 350 – 2000 ng/mL [291.7 – 1666.7], or repeated on dilution.

Manual Dilution: Dilute 1:4 [1 part of whole blood and 3 parts of Dimension® Level 1 (0 ng/mL), CSA Calibrator] to obtain an estimate of the cyclosporine concentration. Enter dilution factor (4). Reassay. Resulting readout is corrected for dilution.

If the CSA result of the 1:4 manually diluted sample is still flagged "Above Assay Range or Assay Range", the CSA result should be reported as "greater than 2000 ng/mL [1664.0 nmol/L]."

Autodilution (AD): Not available for whole blood methods.

Samples with results less than 25 ng/mL [20.8 nmol/L] should be reported as "less than 25 ng/mL [20.8 nmol/L]."

Limitations of Procedure

Patient samples may contain heterophilic antibodies that could react in immunoassays to give falsely elevated or depressed results. This assay has been designed to minimize interference from heterophilic antibodies. Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed. A test result that is inconsistent with the clinical picture and patient history should be interpreted with caution.^{10,11}

The assay has also been designed to minimize interference from antibodies to β-galactosidase. These antibodies can be encountered in samples as a consequence of bacterial infection and may produce falsely elevated results which may not be consistent with clinical evaluation. In very rare instances, immunoassays may produce falsely elevated or decreased results due to other patient-specific interferents.¹²

For patients with impaired liver function and patients receiving other drugs which may induce or inhibit microsomal enzyme activity, the routine use of the Dimension® Cyclosporine Whole Blood Assay analysis may be supported by HPLC data to assess possible changes in biotransformation and elimination.

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
200 ng/mL [166.4 nmol/L]	> 20.4 ng/mL [17.0 nmol/L]
350 ng/mL [291.2 nmol/L]	> 35.2 ng/mL [29.3 nmol/L]

NOTE: To do a 5-test precision check, pipette 300 µL of CSA Calibrators Level 3 and 4 in separate sample cups. Request five replicates of CSA measurements for each calibrator level. Evaluate the SD obtained for each level using the above guidelines.

Expected Values

CAUTION: There is no universally established therapeutic range for cyclosporine in whole blood. Values obtained with different assay methods cannot be used interchangeably due to differences in assay methods and reagent specificity.

Individual cyclosporine values cannot be used as the sole indicator for making changes in the treatment regimen. Each patient should be thoroughly evaluated clinically before treatment adjustments are made and each user must establish his or her own ranges based on these clinical experiences. Because of varying patterns of cross-reactivity with metabolites, the consistent use of one assay for individual patients is recommended. Conversion factors obtained from assay method comparisons should not be used to predict values for individual patients.

The complexity of the clinical state, individual differences in sensitivity to immunosuppressive and nephrotoxic effects of cyclosporine, co-administration of other immunosuppressants, type of transplant, time post-transplant, and a number of other factors contribute to different requirements for optimal blood levels of cyclosporine.¹³

Specific Performance Characteristics¹

The following data represent typical performance for the Dimension® system.

Precision^{a,b}

Material	Standard Deviation (% CV)		
	Mean ng/mL [nmol/L]	Within-Run ng/mL [nmol/L]	Total ng/mL [nmol/L]
QC Level 1*	65 [54.1]	3.4 [2.8] (5.2)	5.3 [4.4] (8.2)
QC Level 2	151 [125.6]	6.5 [5.4] (4.3)	8.1 [6.7] (5.4)
QC Level 3	400 [332.8]	17.3 [14.4] (4.3)	19.0 [15.8] (4.8)
Whole Blood Pool-1**	83 [69.1]	4.9 [4.1] (5.9)	6.0 [5.0] (7.3)
Whole Blood Pool-2	132 [109.8]	5.5 [4.6] (4.2)	8.1 [6.7] (6.2)
Whole Blood Pool-3	219 [182.2]	7.7 [6.4] (3.5)	9.5 [7.9] (4.3)
Whole Blood Pool-4	369 [307.0]	13.5 [11.2] (3.6)	14.4 [12.0] (3.9)

*More Diagnostics Tac/CsA Immunosuppressant Controls, More Diagnostics Inc, Los Osos, CA 93402.

**Whole Blood Pool (1,2,3) prepared from EDTA whole blood and gravimetrically supplemented with cyclosporine at above levels.

f. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).

g. Reproducibility testing was done in accordance with the NCCLS Guideline for User Evaluation of Precision Performance of Clinical Chemistry Devices (EP5-A2, 2004).

h. Specimens at each level were analyzed in duplicate, once a day, for 20 days. The within-run and total standard deviations were calculated by analysis of variance method.

Method Comparison Regression Statistics¹



Comparative Method	Slope	Intercept ng/mL [nmol/L]	Correlation Coefficient	n ¹
Abbott TDx®				
All	0.83	-21.3 [-17.7]	0.957	667
Heart	0.83	-20.9 [-17.4]	0.927	111
Liver	0.82	-26.3 [-21.9]	0.952	201
Kidney	0.83	-18.8 [-15.6]	0.966	279
HPLC				
All	1.19	-13.7 [-11.4]	0.935	661
Heart	1.28	-16.0 [-13.3]	0.893	111
Liver	1.18	-13.3 [-11.1]	0.961	201
Kidney	1.13	-8.9 [-7.4]	0.920	273




i. Passing-Bablok¹⁴ method was used to fit the linear regression line. The model equation for regression statistics is: Result of Dimension® system = (Slope x comparative method result) + Intercept.

Note: Regression of Dimension® CSA method (y) vs. HPLC (x) was also performed using the ordinary least squares regression method. The regression error of estimate (sy/x) was 32.3 ng/mL [26.9 nmol/L] (ALL); 40.5 ng/mL [33.7 nmol/L] (HEART); 24.7 ng/mL [20.6 nmol/L] (LIVER); 25.4 ng/mL [21.1 nmol/L] (Kidney).

j. The range of CSA values in the correlation study was 27.3 to 494.3 ng/mL [22.7 – 411.6 nmol/L]. The number of individual patients (adults, 23 to 76 years old), represented in the above method comparisons for each organ transplant group is: 36 (Heart); 44 (Liver); 68 (Kidney); 32 (Bone Marrow and Lung).

TDx® is the registered trademark of Abbott Diagnostics, Abbott Park, IL 60064.

Specificity

Non-Interfering Substances

The following substances do not interfere with the CSA method when present in whole blood at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at a CSA level of 200 ng/mL [166.5 nmol/L]:

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1322 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Amikacin Sulfate	100 µg/mL	127.9 µmol/L
Amphotericin B	100 µg/mL	108.2 µmol/L
Ampicillin	5 mg/dL	143 µmol/L
Apresoline	100 µg/mL	508.7 µmol/L
Ascorbic Acid	3 mg/dL	170.3 µmol/L
Azathioprine	100 µg/mL	360.6 µmol/L
Bilirubin	60 mg/dL	1026 µmol/L
Caffeine	10 mg/dL	515 µmol/L
Carbamazepine	12 mg/dL	508 µmol/L
Cephalexin	100 µg/mL	208.9 µmol/L
Chloramphenicol	25 mg/dL	774 µmol/L
Chlorpromazine	5 mg/dL	157 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	10 mg/dL	396 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 75	2500 mg/dL	333 µmol/L
Diazepam	2 mg/dL	70 µmol/L
Digoxin	5 ng/mL	6.4 nmol/L
Digitoxin	100 µg/mL	130.7 µmol/L
Dipyridamole	100 µg/mL	198.2 µmol/L
Disopyramide	100 µg/mL	294.6 µmol/L
Erythromycin	20 mg/dL	273 µmol/L
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 µmol/L
FK506 (Tacrolimus)	0.1 µg/mL	0.12 µmol/L
Furosemide	2 mg/dL	61 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Hemoglobin	1000 mg/dL	0.62 mmol/L (monomer)
Heparin	8000 U/L	8000 U/L
Ibuprofen	40 mg/dL	1939 µmol/L
Immunoglobulin G	6 g/dL	60 g/L
Kanamycin A	100 µg/mL	171.6 µmol/L
Kanamycin sulfate B	100 µg/mL	206.8 µmol/L
Ketoconazole	100 µg/mL	188.2 µmol/L
Lidocaine	6 mg/dL	256 µmol/L
Lincomycin	100 µg/mL	225.7 µmol/L
Lipemia	1500 mg/dL	16.95 mmol/L (triglyceride)
Lithium	3.5 mg/dL	5.07 mmol/L (lithium chloride)
Methotrexate	100 µg/mL	220 µmol/L
Methylprednisolone	100 µg/mL	267 µmol/L
Mycophenolic Acid	100 µg/mL	312.1 µmol/L
N-acetyl procainamide	100 µg/mL	360.5 µmol/L
Neomycin sulfate	100 µg/mL	110 µmol/L
Nicotine	2 mg/dL	123 µmol/L
Oxytocin	100 µg/mL	99.3 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	442 µmol/L
Phenobarbital	15 mg/dL	646 µmol/L
Phenytoin	10 mg/dL	396 µmol/L
Prazosin	100 µg/mL	238.2 µmol/L
Prednisolone	100 µg/mL	277.5 µmol/L
Prednisone	100 µg/mL	279 µmol/L
Primidone	10 mg/dL	458 µmol/L
Procainamide	100 µg/mL	367.9 µmol/L
Propoxyphene	0.4 mg/dL	12 µmol/L
Propranolol	100 µg/mL	338 µmol/L
Protein (Albumin)	6 g/dL	60 g/L
Protein (Total)	4 g/dL	40 g/L
Protein (Total)	12 g/dL	120 g/L
Quinidine	100 µg/mL	308.3 µmol/L
Rapamycin	0.1 µg/mL	0.11 µmol/L
Rheumatoid Factor	500 IU/mL	500 IU/mL
Rifampin	100 µg/mL	121.5 µmol/L
Salicylic Acid	50 mg/dL	3.62 mmol/L
Spectinomycin	100 µg/mL	246.7 µmol/L
Streptomycin	100 µg/mL	68.6 µmol/L
Theophylline	25 mg/dL	1388 µmol/L
Tobramycin	100 µg/mL	213.9 µmol/L
Tocainide	100 µg/mL	438.6 µmol/L
Triamterene	100 µg/mL	394.8 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L
Vancomycin	100 µg/mL	67.3 µmol/L
Verapamil	100 µg/mL	203.6 µmol/L

Cross-reactivity

Six major metabolites were evaluated for cross-reactivity in the presence of 200 ng/mL [166.5 nmol/L] cyclosporine. The percent cross-reactivity was calculated as follows:

% cross-reactivity =
$$\frac{\text{measured CSA} - \text{control CSA (200 ng/mL [166.5 nmol/L])}}{\text{metabolite added}} \times 100$$

Cyclosporine Metabolite	Metabolite Level Tested ng/mL [nmol/L]	Cross-Reactivity%
AM1 (M17)	1000 [821]	1.8
AM1c (M18)	1000 [821]	1.1
AM1c9 (M26)	1000 [810]	1.7
AM4N (M21)	1000 [842]	6.0
AM9 (M1)	1000 [821]	2.1
AM19 (M8)	1000 [810]	1.8

Recovery

Known amounts of CSA were added to human whole blood samples with concentrations of 100, 200, 300 and 400 ng/mL [83.3, 166.5, 249.8, 333.0 nmol/L]. CSA concentrations on these samples were then measured and the calculated percent recovery ranged from 93.0% to 101.4% with a mean recovery of 96.6%.

Percent Recovery =
$$\frac{\text{Value obtained} - \text{Baseline}}{\text{Amount Added}} \times 100$$

Linearity

The analysis of cyclosporine on the Dimension® RxL by the Dimension® CSA method shows a deviation from the estimated line of less than 10 ng/mL [8.3 nmol/L] with 95% confidence at the analytical concentration of 200 ng/mL [166.5 nmol/L]. The lower limit of the linearity claim is 25 ng/mL [20.8 nmol/L].

Limit of Detection: 25 ng/mL [20.8 nmol/L]

The limit of detection for CSA is 25 ng/mL, determined consistent with CLSI guideline EP-17A and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 120 determinations, with 5 blank samples and 5 low level samples.¹⁵

Functional Sensitivity: 30 ng/mL [25.0 nmol/L]

The Functional Sensitivity of CSA was determined by assaying multiple specimens (10 – 70 ng/mL) and assessing the total precision for these specimens. The study was conducted over a six week period, with 12 days of testing, 1 run per day with samples run in singlet on two instruments with two reagent lots. The testing included data from two calibrations. The inter-assay coefficient of variation was plotted by nonlinear regression versus concentration. The functional sensitivity (at 20% CV) was 30 ng/mL [25.0 nmol/L].

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, Flex®, EXL™ and LOCI® are trademarks of Siemens Healthcare Diagnostics.

©2009 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía:

- McMillan MA. Clinical pharmacokinetics of cyclosporin. *Pharmac Ther* 1989;42:135-156.
- Kahan BD. Cyclosporine. *N Engl J Med* 1989;321:1725-1738.
- Borel JF. Cyclosporin A - present experimental status. *Transplant Proc* 1981;13:344-8.
- Physicians Desk Reference, ed 49 Montvale, NJ Medical Economics Data Production Company, 1995, pp 2184-2185.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne PA 19087-1898 USA, 2003.
- Lemaire M, Tillement JP. Role of lipoproteins and erythrocytes in the *in vitro* binding and distribution of cyclosporin A in the blood. *J Pharm Pharmacol* 1982;34:715.
- Warner A, Annesley T. *Guidelines for Therapeutic Drug Monitoring Services*, NACB/AACC Press, 2029 K Street, N.W., Washington, D.C., 1999, p 46.
- Kaplan LA, Pesce, AJ. *Clinical Chemistry-Theory, analysis and correlation*, Third Edition, Mosby Co., 1996, p. 1105.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem* 1992;38:1737-1742.
- Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. *Clin Chem* 1999; 45:7,942-956.
- Ismail AA, Walker, PL, Cawood ML and Arth JH. Interference in immunoassay is an underestimated problem. *Ann. Clin. Biochem* 2002; 39:1737-1742.
- Cyclosporine Task Force. Critical issues in cyclosporine monitoring: Report of the task force on cyclosporine monitoring. *Clin Chem* 1987;33:1269-88.
- Passing HW. Bablok (1983): A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. *J. Clin. Chem. Biochem.* 21: 709-720.
- Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar

Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad

Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote

Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo

Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso

Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante

Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea

Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos

In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic *in vitro* / Dispositivo medico-diagnostico *in vitro* / Producto sanitario para diagnóstico *in vitro*

Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura

Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso

Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril

CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE

Contents / Inhalt / Contenu / Contenuto / Contenido

Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución

Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFIS

 Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA

Global Siemens
Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens
Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591
USA
siemens.com/healthcare